



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

OptimAI Line, Aluminum 16 Bin Cassette Medication Cart #AL16CS-5

Specifications

Construction:

- Painted aluminum cabinet and drawers

Features:

- 16 each 5" wide (12.7 cm) patient bin drawers in 2 removable cassettes; bin dividers and labels are included
- Cassettes are located in pull-out drawers for ergonomic removal
- BEST® lock key locking system with two keys
- Pull-out shelf on right side of cart
- Pre-threaded accessory mounting holes
- Self-closing, ball-bearing full extension drawer slides
- Removable plastic top with integrated push handles
- Wrap-around bumper on separate stabilizing frame
- 5" ball-bearing casters, all full swivel - one directional, one with brake
- Cart is delivered fully assembled
- Latex Free
- Twelve year warranty
- Hard-baked powder coat finish

Suggested Cart Options:

- Cassette exchange package: 16CS5-PKG (includes 16 each 5" wide (12.7 cm) patient bin drawers and 2 double-tier aluminum cassettes)
- Narcotics drawer: CS-4FLN (separately locking narcotics drawer for 4" multi-purpose storage drawer)

Internal Drawer Dimensions:

- 20" W x 17" D (50.8 cm W x 43.2 cm D)

Bin Dimensions:

- 5" W x 9.5" L x 3" H
(12.7 cm W x 24.1 cm L x 7.6 cm H)

Assembled Cart Dimensions

- 45.75" H x 30.5" W x 23.875" D
(116.2 cm H x 77.5 cm W x 60.6 cm D)
- Cart net weight: 132 pounds (59.9 kg)
- Approximate shipping weight: 165 pounds (78.8 kg)



Drawer Configuration:

- (1) - 4" (10.2 cm) multi-purpose storage drawer
- (2) - 8" (20.3 cm) double tier cassettes, with eight 5.0" wide patient bin drawers each
- (1) - 4" (10.2 cm) multi-purpose storage drawer
- (1) - 8" (20.3 cm) multi-purpose storage drawer

Sixteen available color choices:

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

OptimAI Line, Aluminum 24 Bin Cassette Medication Cart #AL24CS-5

Specifications

Construction:

- Painted aluminum cabinet and drawers

Features:

- 24 each 5" wide (12.7 cm) patient bin drawers in 3 removable cassettes; bin dividers and labels are included
- Cassettes are located in pull-out drawers for ergonomic removal
- BEST® lock key locking system with two keys
- Pull-out shelf on right side of cart
- Pre-threaded accessory mounting holes
- Self-closing, ball-bearing full extension drawer slides
- Removable plastic top with integrated push handles
- Wrap-around bumper on separate stabilizing frame
- 5" ball-bearing casters, all full swivel - one directional, one with brake
- Cart is delivered fully assembled
- Latex Free
- Twelve year warranty
- Hard-baked powder coat finish

Suggested Cart Options:

- Cassette exchange package: 24CS5-PKG (includes 24 each 5" wide (12.7 cm) patient bin drawers and 3 double-tier aluminum cassettes)
- Narcotics drawer: CS-4FLN (separately locking narcotics drawer for 4" multi-purpose storage drawer)

Internal Drawer Dimensions:

- 20" W x 17" D (50.8 cm W x 43.2 cm D)

Bin Dimensions:

- 4.2" W x 9.5" L x 3" H
(10.7 cm W x 24.1 cm L x 7.6 cm H)

Assembled Cart Dimensions

- 45.75" H x 30.5" W x 23.875" D
(116.2 cm H x 77.5 cm W x 60.6 cm D)
- Cart net weight: 132 pounds (59.9 kg)
- Approximate shipping weight: 165 pounds (78.8 kg)



Drawer Configuration:

- (1) - 4" (10.2 cm) multi-purpose storage drawer
- (3) - 8" (20.3 cm) double tier cassettes, with eight 5.0" wide patient bin drawers each
- (1) - 4" (10.2 cm) multi-purpose storage drawer

Sixteen available color choices:

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

OptimAI Line, Aluminum 36 Bin Cassette Medication Cart #AL36CS-5

Specifications

Construction:

- Painted aluminum cabinet and drawers

Features:

- 36 each 3.5" wide (8.9cm) patient bin drawers in 3 removable cassettes; bin dividers and labels are included
- Cassettes are located in pull-out drawers for ergonomic removal
- BEST® lock key locking system with two keys
- Pull-out shelf on right side of cart
- Pre-threaded accessory mounting holes
- Self-closing, ball-bearing full extension drawer slides
- Removable plastic top with integrated push handles
- Wrap-around bumper on separate stabilizing frame
- 5" ball-bearing casters, all full swivel - one directional, one with brake
- Cart is delivered fully assembled
- Latex Free
- Twelve year warranty
- Hard-baked powder coat finish

Suggested Cart Options:

- Cassette exchange package: 36CS3-PKG (includes 36 each 3.5" wide (8.9 cm) patient bin drawers and 3 double-tier aluminum cassettes)
- Narcotics drawer: CS-4FLN (separately locking narcotics drawer for 4" multi-purpose storage drawer)

Internal Drawer Dimensions:

- 20" W x 17" D (50.8 cm W x 43.2 cm D)

Bin Dimensions:

- 3.5" W x 9.5" L x 3" H
(8.9cm W x 24.1cm L x 7.6 cm H)

Assembled Cart Dimensions

- 45.75" H x 30.5" W x 23.875" D
(116.2 cm H x 77.5 cm W x 60.6 cm D)
- Cart net weight: 132 pounds (59.9 kg)
- Approximate shipping weight: 165 pounds (78.8 kg)



Drawer Configuration:

- (1) - 4" (10.2 cm) multi-purpose storage drawer
- (3) - 8" (20.3 cm) double tier cassettes, with 12 3.5" wide patient bin drawers each
- (1) - 4" (10.2 cm) multi-purpose storage drawer

Sixteen available color choices:

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black