

Model DS5100

Digital Physician's Scale

Today's busy physician's office needs Doran's DS5100 Digital Physician's Scale that provides a faster and easier way to weigh and measure patients accurately. With Doran's Touchless Start feature, you can acquire a patient's weight in as little as two seconds without any buttons to press. Simply have the patient step onto the large DS5100 platform and an accurate weight up to 500 pounds will be held on the display. Once the patient leaves the platform, the DS5100 is immediately ready for another weighing. Imagine the time saved over a year in weighing patients alone!

The sturdy aluminum column houses the height measurement rod. Simply extend the measurement arm and pull upwards for adult heights or push down to measure a child's height. To calculate BMI, enter the height using the scroll buttons, press the BMI button and the patient's Body Mass Index is instantly displayed in the right-hand LCD. When not in use, the BMI display functions as a clock in 24 or 12 hour format.

MAXIMUM CAPACITY

500 x 0.1 lb (200 x 0.05kg)

ACCURACY

+/- 0.1 lb (0.05kg)

PLATFORM SIZE

15" W x 13" D x 2" H

DISPLAYED UNITS

lb, kg

DISPLAY

Dual 1" High LCD

CONTROLS

Zero, LB, KG, Height Input/BMI, FT/IN, CM, Clock

HEIGHT ROD ADJUSTMENT

Minimum 2' - 7" (81 cm)

Maximum 6' - 10" (210 cm)

CLOCK

12 or 24 Hour Display

BMI CALCULATOR

Input data in FT/IN or CM

POWER

Four C Batteries (Included)

BATTERY LIFE

Three Years

CONSTRUCTION

ABS and Metal Constructio

WARRANTY

3 Three Years



Doran Scales, Inc. 1-800-264-4106
1315 Paramount Parkway Fax: 630-879-0073
Batavia, IL 60510 www.doranmedical.com