

ENOCHS ENCORE 4200

High-Low Examination Table

PATIENT LOAD: 500 LBS.

evenly distributed on the table's main couch

PAPER ROLL HOLDER is concealed at headend of table and accommodates 18" or 21" paper

PAPER CUTTER positioned at the headend of table, encourages staff to remove the entire length of table paper contaminated by patient

FULL POWER CONTROLLED BACK smoothly adjusts to any angle from horizontal to an 80° chair position

FIELD REVERSIBLE SIDE DRAWER

Durable 6" deep **MOULDED LARGE SIDE DRAWER** with rounded internal corners and edges for easy and efficient cleaning and maintenance

Extra strength **STEEL REINFORCED SIDE GUIDE MOUNTING SYSTEM** for discipline non-binding movement

The 100% seamless **VACUUM-FORMED TOP** is designed without decorative embossing, cross-seams and fabric folds for easy and efficient cleaning. Top is padded with 1.5" of high-density foam with an extra .5" along the edges for patient comfort.

100% PERINEAL RECESS seamless for infection control maintenance

PELVIC TILT - 7° and 11° positions with release handles on either side of table

4" deep STAINLESS STEEL IRRIGATION BASIN moves independent of leg extension

POSITIVE LOCK MULTI-POSITION STIRRUPS with length and lateral adjustment

PADDED PULLED LEG EXTENSION works independent of Irrigation Basin and allows full access to Basin when unit is extended

4" Deep DRAWER with rounded internal corners and edges for quick and easy cleaning

MAINS POWER LIGHT

'EASY TOUCH' FOOT CONTROL adjusts table's height and back with reliable hands-free action. Pedals are covered with non-slip tread and utilize a natural left to right rocking motion of operator's foot.

TABLE FEATURES:

- 2-YEAR Warranty
- Hand Control available as an accessory
- Fully welded All Steel Base Construction - NO PLASTIC
- Full Power adjustment of back
- Easy to clean 100% Seamless Vacuum Formed Top with additional padding on edges for patient comfort
- Pelvic Tilt 7° and 11°
- Stainless Steel Irrigation Basin operates independent of Leg Extension
- Easy to adjust Stirrups with length and lateral adjustment
- 1.4 cubic feet of Storage
- Certified by UL to Standard 60601-1 CAN/CSA C22.3 No. 601-1
- ENOCHS Encore 4200 replaces the Power 4000-01



Securely attached to the table base, the ENOCHS vacuum-formed top is the strongest top on the market with an all steel frame and back. Top is padded for patient comfort with 1.5" of high density foam with an additional .5" of foam contouring the edges. The sleek Euro-look top features a seamless surface, moulded perineal recess and a pelvic tilt with positions of 7° and 11°. The top offers proven performance of strength, durability and long service life.



The working end of the table has the exclusive ENOCHS design of the pullout padded leg extension moving independent of the 4" deep stainless steel irrigation basin. This design reduces surface area to allow fast easy clean up after irrigation procedure. Patients with limited rotation can be placed in a GYN position utilizing the pad board instead of stirrups while physician maintains unrestricted access to the irrigation basin.

Specifications for the ENOCHS ENCORE 4200 High-Low Examination Table

Warranty:	2-YEAR
Base Construction:	Fully welded all steel base
Power System:	Electro-Mechanical provides smooth and quiet movement
Power Adjustments:	Height: 24" - 40" Back: horizontal - 80 °
Foot Control:	"Hands free", standard on every table, SELV
Hand Control:	"Easy to read", optional accessory, add at anytime, SELV
Couch Top:	100% Seamless vacuum formed style, all steel frame construction
Perineal Recess:	Moulded, no seams, 10" W x 2" D at center point
Pelvic Tilt:	7°, 11°
Leg Extension:	Extends out 20" from end of seat section
Irrigation Basin:	4" Deep stainless steel basin
Electrical Outlets: with Circuit Breaker	Not available on this model
Instrument Warmer: with ON/OFF Switch, Circuit Breaker & Reset	Not available on this model
Stirrups:	Multi-positional, lateral and length adjustment
Storage:	Side Drawer: 15.5" W x 15" D x 6" H Footend Drawer: 15.5" W x 15" D x 4" H Total Storage Space: 1.4 cubic feet
Paper Roll Holder:	Accommodates 18" or 21" paper widths
Accessories:	hand control, knee crutches, padded arm board, drawer locks, restraint straps, IV pole, urology basin with screen and drain, headend accessory rails, footend accessory rails, caster base
Mains Power Entry: with Circuit Breaker	IEC connection, internal circuit breaker with reset
Duty Cycle:	Table: 1 minute ON / 9 minutes OFF Receptacles: n/a Instrument Warming Drawer: n/a
Ratings:	Table: 115V 60 Hz 10 amps Receptacles: n/a Total Receptacle Load: n/a
Table Dimensions:	Width: 27", Height: 24" - 40", Length: 52" (extends to 72")
Table Weight:	416 lbs.
Patient Load:	500 lbs. evenly distributed on main couch
Certifications:	By UL to Standard 60601-1, CAN/CSA C22.3 No. 601-1