

Secure-A-Scope™ Storage System

- Revolutionary new design sets the gold standard for scope storage and accessibility
- Rotating scope holders make loading and removing scopes easy, safe, and secure
- Scope holders engineered to protect scope heads and support umbilicals
- Back wall of cabinet padded vinyl to prevent damage to distal tips
- Removable plastic drip tray
- Filtered ambient air vents
- Specialized holders designed for EUS, Enteroscopy, TEE



Ask to see a demonstration in your hospital today.

Call MASS™ or find your local rep on our website.



M1200LHG
SECURE-A-SCOPE™ STORAGE SYSTEM
Locking Hinged Door with Glass Inserts.
Holds 9 Flexible Endoscopes.
13.62"D x 40.87"W x 88.62"H



M1800
SECURE-A-SCOPE™ STORAGE SYSTEM
with FlexLok™ Door.
Holds 18 Flexible Endoscopes.
18"D x 40.87"W x 95.25"H

Secure-A-Scope™ Exclusively From MASS™
Unique, Rotating Scope Holders
The Storage System To Protect All Your Scopes

- Cabinets to store 4, 8, 9, or 18 flexible Endoscopes
- Door Options - all models
 - FlexLok™ Door
 - Hinged Doors with Glass Inserts
 - Solid Hinged Doors



M1200
SECURE-A-SCOPE™
STORAGE SYSTEM
with FlexLok™ Door. Holds 9 Flexible Endoscopes.
12.87"D x 40.87"W x 94.12"H



M1202
SECURE-A-SCOPE™
STORAGE SYSTEM
with FlexLok™ Door. Holds 4 flexible Endoscopes.
12.87"D x 21.62"W x 94.12"H



M1802
SECURE-A-SCOPE™
STORAGE SYSTEM
with FlexLok™ Door. Holds 8 Flexible Endoscopes. 18"D x 21.62"W x 95.25"H



M2000LHG
Pass-through SECURE-A-SCOPE™
Scope Cabinet
with Locking Hinged Glass Doors. Holds 4 scopes on each side wall for capacity of 8 scopes. 23.12"D x 40.87"W x 93"H

- Cabinets to store 4/8 or 8/16 small diameter scopes such as Bronchoscopes, Cystoscopes, or Ureteroscopes



M1208
SECURE-A-SCOPE™
STORAGE SYSTEM
with FlexLok™ Door. Holds 9 small diameter scopes. 12.87"D x 40.87"W x 77"H



M1204
SECURE-A-SCOPE™
STORAGE SYSTEM
with FlexLok™ Door. Holds 4 small diameter scopes. 12.87"D x 21.62"W x 77"H



- Vertically and independently adjustable upper and lower scope supports
- Holds 4-8 small diameter scopes
- Hinged, 3-point locking, clear hinged door
- Wall mounted for loading and removing scopes easily, safely and securely
- Silicone inserts gently cushion and support any scopes. The inserts are easily, quickly cleaned or replaced
- Available for immediate shipment

ENDOSCOPY PROCEDURE CARTS



M3103
Full Size ERCP Cart
with FlexLok™ Door on upper
and lower sections.
26"D x 45.12"W x 71"H



M3104
Mid Size ERCP Cart
with FlexLok™ Door.
26"D x 45.12"W x 42"H



M3102ED
Dilator Storage Cart
with FlexLok™ Door.
26"D x 45.12"W x 42"H



M1201HG
**Small Diameter Wall
Mounted SECURE-A-SCOPE™
Scope Cabinet**
with hinged, 3-point locking, clear hinged
door. Holds 4-8 small scopes.
11.25"D x 19.125"W x 48.125"H, 51 lbs

**SCOPE
TRANSPORT
CARTS**

- Carts available with four or nine autoclavable trays (lids included)
- Heavy-duty handles for easier mobility

**This Cart is Steam or
Machine Washable!**



M3124
HDPE Scope Transport Cart
22.62"D x 27"W x 44"H



M6204SS
**Stainless Steel Scope
Transport Cart**
21"D x 24.4"W x 35.75"H