

ISOLATION STATION CARTS

Reduce the Risk of Pathogen Transmission



ISO6550

Isolation Station with Antimicrobial Paint,
Open Well Design

FEATURES & BENEFITS

- Isolation Station™ (patent pending) designed to counter the alarming rise in drug-resistant bacteria in hospitals.
- Makes following hygiene practices easier and less time consuming.
- Provides convenient access to supplies that assist infection control practitioners to follow their isolation protocols. As a result, the spread of infection decreases, hospitals avoid infection cases, and organization costs are reduced.
- Open well design eliminates searching in drawers and touching handles.
- Incorporates silver ion technology in the paint as an antimicrobial agent, which resists nosocomial pathogens when compared to current carts.
- Improves adherence to infection control and stops the spread of bacteria on contaminated hands and clothes by using a storage and supply system that makes it easier for healthcare workers to properly disinfect between patient visits.
- Two each, folding side shelves with quick-release latch (**ISO6550**)
- Well lid can be latched open
- Built-in glove box holder and mask box holder
- Area for hand sanitizer/soap dispenser
- Adjustable coat rack
- Powder coat paint finish
- Patent pending on design
- Ball-bearing, full extension drawer slides (**ISO6550**)
- 12.70 cm (5 inch) full swivel casters, one with brake, one directional
- Cart is delivered fully assembled



ISO6548

Isolation Station (small footprint)
with Antimicrobial Paint,
Open Well Design



Open well showing supplies that are accessed without having to open a drawer



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Classic Line, Four Drawer Infection Control Cart with Key Lock, Standard Package #6520



Color Shown is Yellow
Additional Colors Available

Isolation Cart Upgrade:

- Lemon drop yellow finish available (antimicrobial)

Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- Key lock with two keys
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Ball-bearing full extension drawer slides
- 5" ball-bearing casters, all full swivel - one directional, one with brake
- Removable plastic top
- Pre-drilled threaded holes for side mounted accessories
- Aluminum side rail / push handle
- Stabilizing pontoon bumpers
- Cart is delivered fully assembled
- Twelve year warranty on manufactured parts
- Latex free
- Durable, hard-baked powder coat finish

Drawer Configuration:

- (1) – 3" (7.6 cm) multi-purpose storage drawer
- (3) – 9" (22.9 cm) multi-purpose storage drawer

Internal Drawer Dimensions:

23" W x 17" D (58.4 cm W x 43.2 cm D)

Assembled Cart Dimensions:

- 42.5" H x 32" W x 22" D
(108 cm H x 81.3 cm W x 55.9 cm D)
- Cart net weight: 113 pounds (51.3 kg)
- Approximate shipping weight: 156 pounds (70.8 kg)
- Approximate shipping dims: 53" H x 33" W x 26" D
(134.69 cm H x 83.8 cm W x 66 cm D)



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Classic Line, Three Drawer Infection Control Cart with Key Lock, Specialty Package

#6511



Color Shown is Yellow
Additional Colors Available

Isolation Cart Upgrade:

- Lemon drop yellow finish available (antimicrobial)

Assembled Cart Dimensions:

- 45.75" H x 32" W x 22" D
(116.2 cm H x 81.3 cm W x 55.9 cm D)
- With coat and hat rack: 63.5" H (161.3 cm H)
- Cart net weight: 114 pounds (63.6 kg)
- Approximate shipping weight: 150 pounds (68 kg)
- Approximate shipping dims: 53" H x 33" W x 26" D
(134.6 cm H x 83.8 cm W x 68 cm D)

Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- Key lock with two keys
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Stabilizing pontoon bumpers
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Ball-bearing full extension drawer slides
- Removable plastic top
- Painted steel top rail / push handle (68160)
- Coat and hat rack (680412)
- Utility hooks (4 ea) (680402)
- Coat and hat rack / customer installed
- Cart is delivered mostly assembled
- Pre-drilled threaded holes for side mounted accessories
- Latex free
- Durable, hard-baked powder coat finish
- Twelve year warranty on manufactured parts

Quick Ship:

QuickShip applies to 6511Q
Yellow only, in quantities of 5 or less-
Ships within 5 business days after receipt of order

Drawer Configuration:

- (1) - 6" (15.2 cm) multi-purpose storage drawer
- (2) - 12" (30.5 cm) multi-purpose storage drawer

Internal Drawer Dimensions:

23" W x 17" D (58.4 cm W x 43.2 cm D)



650 Ford Street
Colorado Springs, CO 80915
Toll: (800)433-4064 * Ph: (719)637-0300
Fax: (719)597-8273 * www.harloff.com

Classic Line, Three Drawer Infection Control Cart with Key Lock, Standard Package

#6510

Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- Key lock with two keys
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Ball-bearing full extension drawer slides
- 5" ball-bearing casters, all full swivel - one directional, one with brake
- Removable plastic top
- Pre-drilled threaded holes for side mounted accessories
- Aluminum side rail / push handle / Customer mount
- Stabilizing pontoon bumpers
- Cart is delivered fully assembled
- Twelve year warranty on manufactured parts
- Latex free
- Durable, hard-baked powder coat finish

Quick Ship:

QuickShip applies to 6510Q
Yellow only, in quantities of 5 or less-
Ships within 5 business days after receipt of order

Drawer Configuration:

- (1) - 6" (15.2 cm) multi-purpose storage drawer
- (2) - 12" (30.5 cm) multi-purpose storage drawer

Internal Drawer Dimensions:

23" W x 17" D (58.4 cm W x 43.2 cm D)

Assembled Cart Dimensions:

- 42.5" H x 32" W x 22" D
(108 cm H x 81.3 cm W x 55.9 cm D)
- Cart net weight: 112 pounds (50.8 kg)
- Approximate shipping weight: 148 pounds (67.2 kg)
- Approximate shipping dims: 53" H x 33" W x 26" D
(134.6 cm H x 83.8 cm W x 66 cm D)



Color Shown is Yellow
Additional Colors Available

Isolation Cart Upgrade:

- Lemon drop yellow finish available (antimicrobial)



650 Ford Street
Colorado Springs, CO 80915
Toll: (800)433-4064 * Ph: (719)637-0300
Fax: (719)597-8273 * www.harloff.com

Classic Line, Three Drawer Infection Control Cart, with Basic Electronic Pushbutton Lock, Standard Package #7510E

Specifications



**Color Shown is Yellow
Additional Colors Available**

Construction:

- Sturdy, uni-body steel construction

Features:

- Basic electronic pushbutton lock
- Key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Stabilizing pontoon bumpers
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Ball-bearing full extension drawer slides
- Removable plastic top
- Aluminum side rail / push handle
- Cart is delivered fully assembled
- Pre-drilled threaded holes for side mounted accessories
- Latex free
- Durable, hard-baked powder coat finish
- Twelve year warranty on manufactured parts

Drawer Configuration:

- (1) - 6" (15.2 cm) multi-purpose storage drawer
- (2) - 12" (30.5 cm) multi-purpose storage drawer

Internal Drawer Dimensions:

- 23" W x 17" D (58.42 cm W x 43.18 cm D)

Assembled Cart Dimensions:

- 42.5" H x 32" W x 22" D
(108 cm H x 81.3 cm W x 55.9 cm D)
- Cart net weight: 141 pounds (63.9 kg)
- Approximate shipping weight: 183 pounds (83 kg)
- Approximate shipping dims: 53" H x 33" W x 26" D
(134.6 cm H x 83.8 cm W x 66 cm D)

Isolation Cart Upgrade:

- Lemon drop yellow finish available (antimicrobial)



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

V-Series, Four Drawer Procedure Cart with Key Lock #V30-4K

Harloff's V-Series Carts are value-priced for use in a variety of hospital departments. This drawer configuration is often a popular choice for treatment, procedure and isolation applications.

Specifications

Construction:

- Steel construction

Features:

- Key lock with two keys
- Aluminum side rail / push handle
- 30" of vertical drawer space
- Stabilizing bumper frame with extended front casters and (4) corner bumpers
- 5" (12.7 cm) ball-bearing casters, all full swivel - two with brake
- Pre-drilled threaded holes for V-Series side mounted accessories
- Ball-bearing full extension drawer slides
- Twelve year manufacturer's warranty
- Ships with wheels installed and ready to use
- V-Series accessories are customer installed
- Latex free
- Durable, hard-baked powder coat finish

Internal Drawer Dimensions:

- 23" W x 16.75" D (58.4 cm W x 42.5 cm D)

Assembled Cart Dimensions

- 40" H x 29.5" W x 22" D (101.6cm H x 74.9 W x 55.9 D)
- Cart weight: 105 lbs (47.6 kg)
- Approx shipping weight: 128 pounds (58 kg)

Drawer Configuration:

- (1) – 3" (15.2 cm)
 - (2) – 12" (30.5 cm)
- Integrated Pull Handle on all drawers



Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

V-Series, Three Drawer Procedure Cart with Key Lock #V30-3K

Harloff's V-Series Carts are value-priced for use in a variety of hospital departments. This drawer configuration is often a popular choice for treatment, procedure and isolation applications.

Specifications



Construction:

- Steel construction

Features:

- Key lock with two keys
- Aluminum side rail / push handle
- 30" of vertical drawer space
- Stabilizing bumper frame with extended front casters and (4) corner bumpers
- 5" (12.7 cm) ball-bearing casters, all full swivel - two with brake
- Pre-drilled threaded holes for V-Series side mounted accessories
- Ball-bearing full extension drawer slides
- Twelve year manufacturer's warranty
- Ships with wheels installed and ready to use
- V-Series accessories are customer installed
- Latex free
- Durable, hard-baked powder coat finish

Internal Drawer Dimensions:

- 23" W x 16.75" D (58.4 cm W x 42.5 cm D)

Assembled Cart Dimensions

- 40" H x 29.5" W x 22" D (101.6cm H x 74.9 W x 55.9 D)
- Cart weight: 102 lbs (46.2 kg)
- Approx shipping weight: 125 pounds (56.7 kg)

Drawer Configuration:

- (1) – 6" (15.2 cm)
- (2) – 12" (30.5 cm)

Integrated Pull Handle on all drawers

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

V-Series, Three Drawer Procedure Cart with Key Lock #V24-3K

Harloff's V-Series Carts are value-priced for use in a variety of hospital departments. This drawer configuration is often a popular choice for treatment, procedure and isolation applications.

Specifications

Construction:

- Steel construction

Features:

- Key lock with two keys
- Aluminum side rail / push handle
- 24" of vertical drawer space
- Stabilizing bumper frame with extended front casters and (4) corner bumpers
- 5" (12.7 cm) ball-bearing casters, all full swivel - two with brake
- Pre-drilled threaded holes for V-Series side mounted accessories
- Ball-bearing full extension drawer slides
- Twelve year manufacturer's warranty
- Ships with wheels installed and ready to use
- V-Series accessories are customer installed
- Latex free
- Durable, hard-baked powder coat finish

Internal Drawer Dimensions:

- 23" W x 16.75" D (58.4 cm W x 42.5 cm D)

Assembled Cart Dimensions

- 34" H x 29.5" W x 22" D
(86.4 cm H x 74.9 cm W x 55.9 cm D)
- Cart weight: 88 pounds (39.9 kg)
- Approximate shipping weight: 110 pounds (49.9 kg)

Drawer Configuration:

- (1) – 6" (15.2 cm)
 - (2) – 9" (22.9 cm)
- Integrated Pull Handle on all drawers



Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

ETC Line, Four Drawer, Treatment Cart #ETC-4



Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- Low height make these carts ideal for bedside or under a counter
- Key lock with two keys
- 3" stem casters, two with brake
- Roller drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Hard-baked powder coat finish

Drawer Configuration:

- (1) – 2.75" (6.98 cm)
- (3) – 5.75" (14.6 cm)

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Assembled Cart Dimensions

- 27.25" H x 23.38" W x 15.13" D
(69.2 cm H x 59.4 W x 38.4 D)
- Cart weight: 48 lbs (21.8 kg)
- Approx. Shipping Weight: 63 pounds (28.7 kg)

** PLEASE NOTE THIS CART IS SHORTER THAN STANDARD COUNTER HEIGHTS



650 Ford Street
Colorado Springs, CO 80915
Toll: (800) 433-4064 * Ph: (719) 637-0300
Fax: (719) 597-8273 * www.harloff.com

Four Drawer Mini Line Isolation Cart Standard Package #3154K

Specifications



Construction

- Painted steel cabinet and drawers

Standard Features

- Cabinet with 30 vertical inches (76.2 cm) drawer space
- Keylocking system with two keys
- 3-inch (7.6 cm) single-wheel stem caster set, full swivel – two with brake
- Pre-threaded accessory mounting holes
- Ball-Bearing full extension drawer slides
- Removable plastic top
- Certified “Latex Free”
- Five Year Warranty
- Hard-baked powder coat finish

Twelve available color choices:

- Beige (-BG), Burgundy (-BU), Light Blue (-LB), Navy (-NY), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammer Tone Blue (-HTB), Hammer Tone Gray (-HTGY) and Hammer Tone Green (-HTGN)

Drawer configuration:

- 1 x 3-inch (7.6 cm)
- 3 x 9-inch (22.9 cm)

Assembled Cart Dimensions

- 38.25" H x 18" W x 18" D (97.2 cm H x 45.7cm W x 45.7 cm D)
- Cart Net Weight: 90 pounds (41 kg)
- Approx. Shipping Weight: 105 pounds (47.8 kg)
- Cart requires caster installation and accessory assembly after delivery



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Three Drawer Mini Line Infection Control Cart (Short Cabinet), Key-locking, Standard Package #3143K

Specifications



Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Cabinet with 24 vertical inches (61 cm) drawer space
- Key-locking system with 2 keys
- 3" (7.6 cm) single-wheel stem caster set, full swivel – two with brakes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Pre-threaded accessory mounting holes
- Latex free
- Twelve year warranty
- Cart requires caster installation and accessory assembly after delivery
- Durable hard-baked powder coat finish

Drawer Configuration:

- (1) - 6-inch (15.2 cm)
- (2) - 9-inch (22.9 cm)

Internal Drawer Dimensions:

- 12.6" W x 16.75" D (32 cm x 42.5cm)

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Assembled Cart Dimensions

- 34.4" H x 18" W x 18"D
(87.4 cm H x 45.7 cm W x 45.7 cm D)
- Cart net weight: 71 pounds (32.2 kg)
- Approximate shipping weight: 90 pounds (41 kg)



650 Ford Street
Colorado Springs, CO 80915
Toll: (800)433-4064 * Ph: (719)637-0300
Fax: (719)597-8273 * www.harloff.com

Four Drawer Mini Line Isolation Cart with Extended Stabilizing Bumper and 5" Wheels, Key Lock, Standard Package #3254K

Specifications



Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Cabinet with 30 vertical inches (76.2 cm) drawer space
- Key lock with two keys
- Bumper system extends wheel base and stabilizes cart
- 5" ball-bearing full swivel casters – one brake, one directional
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Cart requires caster installation and accessory assembly after delivery
- Durable hard-baked powder coat finish

Drawer Configuration:

- (1) - 3-inch (7.6 cm)
- (3) - 9-inch (22.9 cm)

Internal Drawer Dimensions:

- 12.6" W x 16.75" D (32 cm x 42.5cm)

Assembled Cart Dimensions

- 43.75" H x 19.5" W x 25" D
(110.5 cm H x 49.5 cm W x 63.5 cm D)
- Cart net weight: 87 pounds (39.5 kg)
- Approximate shipping weight: 107 pounds (48.5 kg)

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black



650 Ford Street
Colorado Springs, CO 80915
Toll: (800)433-4064 * Ph: (719)637-0300
Fax: (719)597-8273 * www.harloff.com

Three Drawer Mini Line Infection Control Cart with Extended Stabilizing Bumper and 5" Wheels, with Key Lock #3243K

Harloff's Mini Line of carts offer a narrow design with reduced footprint. Suitable for hospitals, clinics and doctor offices with both space and budget constraints

Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Cabinet with 24 vertical inches (61 cm) drawer space
- Key lock with two keys
- Bumper system extends wheel base and stabilizes cart
- 5" ball-bearing full swivel casters – one brake, one directional
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish

Drawer Configuration:

- (1) - 6-inch (15.2 cm)
- (2) - 9-inch (22.9 cm)

Internal Drawer Dimensions:

- 12.6" W x 16.75" D (32cm x 42.5cm)

Assembled Cart Dimensions:

- 37" H x 19.5" W x 25" D
(94.0 cm H x 49.5 cm W x 63.5 cm D)
- Cart net weight: 94 pounds (42.6 kg)
- Approximate shipping weight: 114 pounds (51.7 kg)



Sixteen available color choices:

*Antimicrobial Finish Available

Sand, Gloss Beige, Burgundy, Light Blue, Navy, Light Gray, Mauve, Red, Teal, Yellow, Hammer Tone Blue, Hammer Tone Gray, Hammer Tone Green, Hammer Tone Red, Hammer Tone Black and White



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Isolation Station, Wall or Door Mount, Standard Package #ISO2800

Specifications



Construction:

- Painted aluminum

Features:

- Saves floor space, no more carts getting in the way of patient flow
- Allows storage of PPE (personal protective equipment) in areas with space limitations
- Pre-drilled holes for wall mount (mounting hardware included)
- Door-hanging bracket included to hang on 2" door
- Designed to hold 2 boxes of gloves, bulk pack of disposable gowns, face masks, shoe covers and caps, or face masks with splash guards
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish

Other Wall Mount / Door Hanging Isolation Station Options:

- Silver ion antimicrobial paint upgrade available
- Custom configurations are available based on your requirements

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Overall Dimensions:

- Wall mount: 23" H x 22.75" W x 5"D
(58.4 cm H x 57.8 cm W x 12.7 cm D)
- Door-hanging: 35" H x 22.75" W x 7.5" D
(88.9 cm H x 57.8 cm W x 19.1 cm D)
- Net weight: 12 pounds (5.5 kg)
- Approx. shipping weight: 17 pounds (7.7 kg)



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Isolation Station, Wall or Door Mount, Standard Package #ISO2825

Specifications

Construction:

- Painted aluminum

Features:

- Saves floor space, no more carts getting in the way of patient flow
- Allows storage of PPE (personal protective equipment) in areas with space limitations
- Pre-drilled holes for wall mount (mounting hardware included)
- Door-hanging bracket included to hang on 2" door
- Designed to hold 3 boxes of gloves, bulk pack of disposable gowns, face masks or shoe covers and caps
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish

Other Wall Mount / Door Hanging Isolation Station Options:

- Silver ion antimicrobial paint upgrade available
- Custom configurations are available based on your requirements

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Overall Dimensions:

- Wall mount: 22" H x 24.9" W x 3.7"D (55.9 cm H x 63.2 cm W x 9.4 cm D)
- Net weight: 15 pounds (6.8 kg)





650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Isolation Station, Wall or Door Mount, Standard Package #ISO2831

Specifications



Color Shown is SAND
Additional Colors Available

Construction:

- Painted aluminum

Features:

- Saves floor space, no more carts getting in the way of patient flow
- Allows storage of PPE (personal protective equipment) in areas with space limitations
- Pre-drilled holes for wall mount (mounting hardware not included)
- Door-hanging bracket included to hang on 2" door
- Designed to hold bulk packs of disposable gowns, face masks or shoe covers, caps and gloves
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish

Other Wall Mount / Door Hanging Isolation Station Options:

- Custom configurations are available based on your requirements

Overall Dimensions:

- Wall mount: 30.5" H x 26" W x 4" D
(77.5 cm H x 66 cm W x 10.2 cm D)
- Net weight: 15 pounds (6.8 kg)



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Isolation Station, Wall or Door Mount, Standard Package #ISO2835

Specifications

Construction:

- Painted aluminum

Features:

- Saves floor space, no more carts getting in the way of patient flow
- Allows storage of PPE (personal protective equipment) in areas with space limitations
- Pre-drilled holes for wall mount (mounting hardware included)
- Door-hanging bracket included to hang on 2" door
- Designed to hold 2 boxes of gloves, 2 bulk packs of disposable gowns, face masks with splash guards or shoe covers and caps
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish

Other Wall Mount / Door Hanging Isolation Station Options:

- Silver ion antimicrobial paint upgrade available
- Custom configurations are available based on your requirements



Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Overall Dimensions:

- Wall mount: 33" H x 28" W x 5"D
(83.82 cm H x 71.52 cm W x 12.8 cm D)
- Net weight: 24 pounds (10.9 kg)