

SS8343

SS8300 Series Stands and Receptacles



Stainless Steel I.V. Stand, Two Hook, Hand Wheel



Standard Features

- 1" (2.53cm) diameter tubular stainless steel construction
- 2 Hooks
- 24" (61cm) non-tip aluminum base (5 legs)
- 2" (5.1cm) swivel casters
- Positive locking hand wheel
- 3 year limited warranty

Specifications

- **Extended height:** 92.5" (235cm)
- **Non-extended height:** 50.5" (128.27cm)
- **Base diameter:** 24" (61cm)
- **Caster diameter:** 2" (5.1cm)
- **Tube diameter:** 1" (2.54cm)



SS8342

SS8300 Series Stands and Receptacles



Stainless Steel I.V. Stand, Four Hook, Hand Wheel



Standard Features

- 1" (2.53cm) diameter tubular stainless steel construction
- 4 Hooks
- 24" (61cm) non-tip aluminum base (5 legs)
- 2" (5.1cm) swivel casters
- Positive locking hand wheel
- 3 year limited warranty

Specifications

- **Extended height:** 92.5" (235cm)
- **Non-extended height:** 50.5" (128.27cm)
- **Base diameter:** 24" (61cm)
- **Caster diameter:** 2" (5.1cm)
- **Tube diameter:** 1" (2.54cm)



SS8310

SS8300 Series Stands and Receptacles



Stainless Steel Mayo Stand 13" x 19" (33cm x 48cm) with Foot Pedal



Standard Features

- All stainless steel construction
- Foot pedal height control
- Removable stainless steel tray 13"x19" (33cm x 48.3cm)
- Tip-resistant, weighted flat base
- 2" (5.1cm) ball bearing swivel casters
- 3 year limited warranty

Specifications

- **Min. Height:** 38.25" (97.2cm)
- **Max. Height:** 38" (96.5cm)
- **Tray Dimensions:** 13"x19" (33cm x 48.3cm)



SS8311

SS8300 Series Stands and Receptacles



Stainless Steel Mayo Stand 16" x 21" (41cm x 53cm) with Foot Pedal



Standard Features

- All stainless steel construction
- Foot pedal height control
- Removable stainless steel tray 16" x 21" (40.6cm x 53cm)
- Tip-resistant, weighted flat base
- 2" (5.1cm) ball bearing swivel casters
- 3 year limited warranty

Specifications

- **Min. Height:** 38.25" (97.2cm)
- **Max. Height:** 63.875" (162.2cm)
- **Tray Dimensions:** 16" x 21" (40.6cm x 53cm)



SS8525

SS8500 Series Bassinets



Bassinet, One Shelf, One Drawer with Basket and Mattress



Standard Features

- Heavy duty tubular stainless steel framework
- Full width stainless steel drawer
- Convenient storage shelf
- 3" (7.6cm) ball bearing swivel casters
- 3 year limited warranty
- Includes basket and mattress pad
- Name card holder

Optional Accessories

- 132 – Swivel casters locking 3" (set of 4)

Specifications

- **Height:** 35.375" (89.9cm)
- **Width:** 31.25" (79.4cm)
- **Depth:** 17" (43.2cm)



SS8528

SS8500 Series Bassinets



Bassinet, Two Door with Basket and Mattress



Standard Features

- Heavy duty tubular stainless steel framework
- Full width stainless steel drawer
- Convenient storage shelf
- 3" (7.6cm) ball bearing swivel casters
- 3 year limited warranty
- Includes basket and mattress pad
- Name card holder

Optional Accessories

- 132 – Swivel casters locking 3" (set of 4)

Specifications

- **Height:** 35.375" (89.9cm)
- **Width:** 31.25" (79.4cm)
- **Depth:** 17" (43.2cm)



SS8538

SS8500 Series Bassinets



Bassinet, One Shelf with Basket and Mattress



Standard Features

- Heavy duty tubular stainless steel framework
- Convenient storage shelf
- 3" (7.6cm) ball bearing swivel casters
- 3 year limited warranty
- Includes basket and mattress pad
- Name card holder

Optional Accessories

- 132 – Swivel casters locking 3" (set of 4)

Specifications

- **Height:** 35.375" (89.9cm)
- **Width:** 31.25" (79.4cm)
- **Depth:** 17" (43.2cm)



SS8350

SS8300 Series Stands and Receptacles



Stainless Steel 13 Quart (12.3L) Kick Bucket



Standard Features

- Heavy duty stainless steel mobile base
- 13 quart (12.3L) stainless steel pail
- Rubber cushions inside bucket support ring
- Rubber bumper on outside of mobile frame
- 2" (5.1cm) ball bearing swivel casters
- 3 year limited warranty

Specifications

- **Height:** 14" (35.6cm)
- **Base Diameter:** 14.5" (36.8cm)
- **Pail capacity:** 13 Qt. (12.3L)

