

MEDICINE-ON TIME® & MULTI-DOSE

MD30-4K (pictured)
Multi-Dose Cart with 4" drawers, key lock on cabinet, 30 bin capacity, internal narcotics box in 1/3 of the bottom drawer, drawer dividers included. Standard Package

MD50-6K (pictured)
Multi-Dose Cart with 4" drawers, key lock on cabinet, 50 bin capacity, drop-shelf, drawer dividers included. Standard Package



680540 - Multi-Dose Divider Set MD30-4K and MD50-6K (1 per drawer) holds 10 bins
Dimensions: 23"W x 17"D

680545 - Multi-Dose Divider Set SL72MD and SL60MD (1 per drawer) holds 12 bins

*Bins shown in photo are not included with carts

AL64MDBiN (pictured)
OptimAL Line Multi-Dose Cart with Best® Lock in cabinet. Includes internal narcotics box and pull-out shelf. 4" drawers designed to hold multi-dose bins. (bins not included)
Capacity of 64 bins.

AL48MDBiN (pictured)
OptimAL Line Multi-Dose Cart with Best® Lock in cabinet. Includes internal narcotics box and pull-out shelf. 4" drawers designed to hold multi-dose bins. (bins not included) Capacity of 48 bins.



The carts on this page have 4" drawers designed to accommodate bin styles used with automated packaging. Drawer dividers available to mark bin locations in drawers.


NeW!!

These carts provide a low maintenance pontoon bumper system and an ergonomic height. Writing surface on these carts is 40 inches.

Accessory packages available for all cart styles

SL72MD (pictured)
Standard Line Multi-Dose Cart with BeST® lock on cabinet, narcotics drawer, pull-out shelf, 4" drawers with dividers. Capacity of 72 bins, Standard Package

SL60MD (pictured)
Standard Line Multi-Dose Cart with BeST® lock on cabinet, internal narcotics box, pull-out shelf, 4" drawers with dividers. Capacity of 60 bins. Standard Package



Multi-Dose Bins

40425 - Multi-Dose Bin (includes label, 2 dividers)

40426 - Bin Labels

40427 - Bin Divider (ea)

*Bins shown in photo are not included with carts



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

Standard Line Multi-Dose Bin Medication Cart, Standard Package #SL72MD

Specifications

Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Dual column medication cart with 72 multi-dose bin capacity, single width narcotics drawer
- BEST® brand lock on cabinet; Illinois lock on narcotics drawer
- Ball-bearing, full extension drawer slides
- Painted steel pull-out side shelf
- Plastic replaceable top with integrated push handles
- Full wrap around vinyl bumper, protecting the cart
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish
- Offered in specialty package (#SL72MDSP) Included are: #68640-P1 drawer tray with adjustable dividers, #680551 eye drop tray, #AL680435 aluminum waste container and #684810 locking sharps container



Capacity of 12 bins per drawer.
 (Bins not included)

Internal Drawer Dimensions:

Left Column Drawers:

- 6 ea – 4" drawers with internal dividers, capacity 12 bins per drawer (Bins are not included)
- 1 ea – 8" multipurpose storage drawer

Right Column Drawers:

- 2 ea – 4" multipurpose storage drawer
- 1 ea – 8" locking narcotics drawer
- 2 ea – 8" multipurpose storage drawer

Thirteen available color choices:

Sand, Gloss Beige, Burgundy,
 Light Blue, Navy, Light Gray,
 Mauve, Red, Teal, Yellow,
 Hammer Tone Blue, Hammer Tone Gray
 and Hammer Tone Green

Assembled Cart Dimensions:

- 45.75" H x 43" W x 23.5" D
 (116.2 cm x 109.2 cm x 59.7 cm)
- Writing surface: 44" H (111.8 cm)
- Cart weight: 210 lbs (95.5 kg)
- Approx shipping weight: 255 lbs (116 kg)



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

**Standard Line Multi-Dose Bin Medication Cart, Single Column,
 Standard Package
 #SL60MD**

Specifications

Construction:

- Painted steel cabinet and aluminum drawers

Features:


- Single column medication cart with 60 multi-dose bin capacity
- Internal narcotics box with Illinois Lock
- BEST® brand lock on cabinet
- Ball-bearing, full extension drawer slides
- Painted steel pull-out side shelf
- Plastic replaceable top with integrated push handles
- Full wrap around vinyl bumper, protecting the cart
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish
- Offered in specialty package (#SL60MDSP) Included are:
 #68640-P1 drawer tray with adjustable dividers, #AL680435 aluminum waste container and #684810 locking sharps container

Internal Drawer Dimensions:

- 1 ea – 4" multi-purpose drawer
- 5 ea – 4" drawers with internal dividers, capacity 12 bins per drawer (Bins are not included)
- 1 ea – 9.75" multipurpose storage drawer with locking narcotics box and adjustable divider

Assembled Cart Dimensions:

- 45.75" H x 36" W x 23.5" D
 (116.2 cm x 91.5 cm x 59.7 cm)
- Writing surface: 44" H (111.8 cm)
- Cart weight: 210 lbs (95.5 kg)
- Approx shipping weight: 240 lbs (109.1 kg)



Capacity of 12 bins per drawer.
 (Bins not included)

Thirteen available color choices:

Sand, Gloss Beige, Burgundy,
 Light Blue, Navy, Light Gray,
 Mauve, Red, Teal, Yellow,
 Hammer Tone Blue, Hammer Tone Gray
 and Hammer Tone Green



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

OptimAL Line 64 Bin Capacity Medication Cart, Standard Package #AL64MDBIN

Specifications

Construction:

- Painted aluminum cabinet and drawers

Features:

- Internal narcotics box with Illinois Lock
- BEST® brand lock on cabinet
- Painted pull-out side shelf with recessed work area
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with durable work surface
- Ergonomic top rail/push handle
- Easy clean protective pontoon bumpers
- 5" ball-bearing premium rubber casters, all full swivel, 1 brake and 1 directional
- Hard-baked powder coat finish (Customer chooses color for cabinet)



Drawer Configuration:

- 4 ea – 4" drawers
- 1 ea – 11" multipurpose storage drawer with locking narcotics box

Assembled Cart Dimensions:

- 42" H x 36" W x 24" D (includes top rail)
(106.6 cm x 91.4 cm x 61 cm)
- Work surface: 38" H (96.5 cm)

Thirteen available color choices:

Sand, Gloss Beige, Burgundy,
Light Blue, Navy, Light Gray,
Mauve, Red, Teal, Yellow,
Hammer Tone Blue, Hammer Tone Gray
and Hammer Tone Green



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

OptimAL Line Multi-Dose Bin Medication Cart, Standard Package #AL48MDBIN



Construction:

- Painted lightweight aluminum cabinet and drawers

Standard Accessories:

- Internal narcotics box with Illinois Lock
- BEST® brand lock on cabinet
- Painted pull-out side shelf with recessed work area
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with durable work surface
- Ergonomic top rail/push handle
- Easy clean protective pontoon bumpers
- 5" ball-bearing premium rubber casters, all full swivel, 1 brake and 1 directional
- Hard-baked powder coat finish (Customer chooses color for cabinet)

Drawer Configuration:

- 4 ea – 4" drawers
- 1 ea – 11" multipurpose storage drawer with locking narcotics box

Thirteen available color choices:

Beige, Burgundy, Light Blue, Red, Teal, Yellow,
Dark Blue, Light Gray, Mauve, Hammertone Grey
Hammertone Blue, Hammertone Green and Sand

Colors may be viewed at www.harloff.com

Assembled Cart Dimensions:

- 42" H x 31" W x 24" D (includes top rail)
- Work surface: 38" H (96.5 cm)



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Classic Line Medication Cart for Multi-Dose Packaging with Key Lock #MD30-4K

Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Key lock with two keys
- Stabilizing pontoon bumpers
- Five-inch sealed ball-bearing premium casters, all full swivel – one with brake, one directional
- Internal locking narcotics box in 1/3 of bottom drawer
- Pull-out side shelf
- Aluminum push handle/side rail
- Removable plastic top
- Ball-bearing full extension drawer slides
- Cart is delivered fully assembled
- Certified "Latex Free"
- Twelve year warranty on manufactured parts
- Compatible with most accessories
- Pre-threaded accessory mounting holes
- Durable hard-baked powder coat finish

- Available in Specialty Package. (MD30-4KSP)
"SP" package includes sharps and waste container with brackets.

Drawer Configuration:

- 3 – 4" Drawers with plastic drawer dividers (10 spaces per drawer)
- 1 – 12" Drawer with 1 adjustable divider and locking narcotics box

Assembled Cart Dimensions:

- 39.25" H x 32" W x 22" D
(100.3 cm H x 81.3 cm W x 55.9 cm D)
- Cart weight: 139 pounds (63.1 kg)
- Approximate shipping dims 47" H x 33" W x 26" D
(119.4 cm H x 83.8 cm W x 66 cm D)
- Approximate shipping weight 180 lbs (81.7 kg)



Thirteen available color choices:

Sand, Gloss Beige, Burgundy,
Light Blue, Navy, Light Gray,
Mauve, Red, Teal, Yellow,
Hammer Tone Blue, Hammer Tone Gray
and Hammer Tone Green



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Classic Line Medication Cart for Multi-Dose Packaging with Key Lock #MD50-6K

Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Key lock with two keys
 - Stabilizing pontoon bumpers
 - Five-inch sealed ball-bearing premium casters, all full swivel – one with brake, one directional
 - Internal locking narcotics box in 1/3 of bottom drawer
 - Pull-out side shelf
 - Aluminum push handle/side rail
 - Removable plastic top
 - Ball-bearing full extension drawer slides
 - Cart is delivered fully assembled
 - Certified "Latex Free"
 - Twelve year warranty on manufactured parts
 - Compatible with most accessories
 - Pre-threaded accessory mounting holes
 - Durable hard-baked powder coat finish
- Available in Specialty Package. (MD30-4KSP)
"SP" package includes sharps and waste container with brackets.

Drawer Configuration:

- 5 – 4" Drawers with plastic drawer dividers (10 spaces per drawer)
- 1 – 10" Drawer with 1 adjustable divider and locking narcotics box

Assembled Cart Dimensions:

- 42.5" H x 32" W x 22" D (108 cm H x 81.3 cm W x 55.9 cm D)
- Cart weight: 149 pounds (67.6 kg)
- Approximate shipping dims: 53" H x 33" W x 26" D (134.62 cm H x 83.8 cm W x 66 cm D)
- Approximate shipping weight: 189 lbs (85.7 kg)



Thirteen available color choices:

Sand, Gloss Beige, Burgundy,
Light Blue, Navy, Light Gray,
Mauve, Red, Teal, Yellow,
Hammer Tone Blue, Hammer Tone Gray
and Hammer Tone Green