

Eliminate Clutter & Save Time!

BIG CASE *back tables*



Part No. 429
Big Case Back Table



Part No. 429-ADJ
Adjustable Back Table

OR Specific is a leading provider of stainless steel operating room back tables and surgical instrument tables. The Big Case Back tables are specifically designed to accommodate large cases such as orthopedics, spinal fusions, neuro, endoscopy, open heart and craniotomies. The additional space created by the upper tier shelf promotes improved visibility, organization and arrangement of instrument trays without the need for stacking.



Suggested Drape

Part No.	Description	Dimensions	418-S	419-S	420-S	418-HD-S	419-HD-S	420-HD-S
427	4' Back Table with removable 21½"D second tier shelf	48"L x 30"D x 33"H	•			•		
428	5' Back Table with removable 21½"D second tier shelf	60"L x 30"D x 33"H		•			•	
428-ADJ	5' Back Table with removable 21½"D second tier shelf (height adjustment from 29" to 38")	60"L x 30"D x 33"H (height adjustment from 29" to 38")		•			•	
429	6' Back Table with removable 21½"D second tier shelf	72"L x 30"D x 33"H			•			•
429-ADJ	6' Back Table with removable 21½"D second tier shelf (height adjustment from 29" to 38")	72"L x 30"D x 33"H (height adjustment from 29" to 38")			•			•



Store the shelf and use as a standard back table

BIG CASE 6' adjustable

Shelf

Upper tier may be set at 12" or 15" above main table surface.


Table Height

Main table surface may be adjusted for a maximum of 38" or minimum of 29".

- May be adjusted intraoperatively

Drape

Unique one piece drape easily covers table. Clear plastic backdrop allows light into all areas of the sterile field.



Safety locking pin for securing upper tier shelf



Casters are 4" in diameter with the front two locking.



Upper tier weight capacity of 150 lbs.



Quickly raise or lower table using the table crank.



Post hole for use in storing the upper tier shelf

Two-Tier Table Drapes

Single piece disposable drapes are available in standard or heavy-duty weight to cover both surfaces of the sterile field. Use the reference chart to the right to match the table with suggested drape.

Reinforced Single Level Table Drapes

This *very heavy* flat drape is intended for use with our Specialty Back Tables as a regular single level drape.



STERILE

For use with Table	Two-Tier Standard	Two-Tier Heavy-duty	Per Case
427	418-S	418-HD-S	10
428	419-S	419-HD-S	10
428-ADJ	419-S	419-HD-S	10
429	420-S	420-HD-S	10
429-ADJ	420-S	420-HD-S	10

For use with Table	Single Level Heavy-duty	Drape Dimensions	Per Case
428	315-S	80" x 90"	10
429	320-S	80" x 110"	10