



23001

Injection molded
onboard regulator &
canister holder



23002



22000

Stainless steel
onboard regulator &
canister holder

22002

SSSCOR, Inc. provides portable and onboard vacuum solutions to meet your emergency suction requirements. SSCOR onboard vacuum systems include options for different pre-hospital applications, conditions and ambulance configurations.

All SSCOR suction products utilize disposable collection for ease of use and operator safety, meet standards for suction efficiency and provide variable vacuum regulation without restricting flow. The new 23000 series rugged plastic regulator and canister bracket provide rounded edges conducive to the interior of the patient care area.

Onboard Suction Solutions from SSCOR, Inc.



SSCOR/board™ CUSTOM CONFIGURATIONS

- 1 Choose a regulator.
- 2 Choose a canister holder.
- 3 Include a pump for a complete onboard vacuum system.

- 1 Panel mount regulators

Model #23001
Rugged plastic



Model #22000
Stainless steel



Standard items included: Variable regulator & gauge assembly, 10' vinyl tubing (5/16" ID x 9/16" OD) [304.8cm tubing (.79cm ID x 1.428cm OD)], mounting screws, straight hose barbs, and hose clamps.

- 2 Canister holders

Model #23002
Rugged plastic



Model #22002
Stainless steel



Standard items included: Canister holder, patient tubing, HI-D®, "The Big Stick®" suction tip, 1200cc/ml Bemis disposable canister, and mounting screws.

- 3 Vacuum pump

Part #SB90120-RT

Description: The SSCOR/board pump package includes a small quiet 12V DC, diaphragm pump capable of generating 30 LPM clinical airflow and over 525 mmHg negative pressure.

Call or look online at www.sscor.com for ordering information and for an illustrated parts and accessory list.

* shown with charging/retention bracket