

# ENOCHS

EXAMINING ROOM FURNITURE

## PEDIATRIC TABLE



# Specifications for *Ped Table*™ Pediatric Table



	50 PEDIATRIC TABLE	51 SIDE CABINET
<b>Warranty:</b>	4 year	4 year
<b>Base Construction:</b>	all steel, fully welded	all steel, fully welded
<b>Base Finish:</b>	durable powder coat finish	durable powder coat finish
<b>Top Design:</b>	seamless Formica™ work surface with 4" back splash	seamless Formica™ work surface with 4" back splash
<b>Height:</b>	36"	32"
<b>Width:</b>	24"	18"
<b>Length:</b>	50"	25"
<b>Work Surface Area:</b>	18"W x 50" L	25"W x 16" D
<b>Storage:</b>	7 - drawers 1- double door storage area with half shelf	2- drawers 1- double door storage area with half shelf
<b>Drawer Type:</b>	high impact moulded drawers, easy to clean rounded corners and sides	high impact moulded drawers, easy to clean rounded corners and sides
<b>Drawer System:</b>	steel reinforced side-guide suspension, engineered for smooth non-binding drawer movement	steel reinforced side-guide suspension, engineered for smooth non-binding drawer movement
<b>Drawer Dimensions:</b>	21.5" W x 15" D x 4" H	21.5" W x 15" D x 4" H
<b>Bulk Storage Dimensions:</b>	25" W x 18" D x 17" H	25" W x 18" D x 17" H
<b>Total Storage Area:</b>	8.57 cubic feet	6 cubic feet
<b>Paper Roll Holder:</b>	concealed, accommodates 18" or 21" paper widths	N/A
<b>Patient Load (evenly distributed):</b>	infant to toddler	when used with infant scale - follow scale capacity
<b>Product Weight:</b>	184 lbs.	88 lbs.
<b>Shipping:</b>	some assembly required	shipped fully assembled
<b>Color Selection:</b>	Top: Cambridge Blue, Earthen Brown, Emerald, Ocean Jade, Fawn, Ninja, Smoke or Windsor Wine Base: Cream White or Pearl Grey Drawers/doors: Beige, Grey or Natural Oak	Base: Cream White or Pearl Grey Drawers/doors: Beige, Grey or Natural Oak