

# ENOCHS RESOLUTE 6220

## General Procedural Table

**PATIENT LOAD: 500 LBS.**  
evenly distributed on the table's main couch



### TABLE FEATURES:

- **4-YEAR WARRANTY**
- "Easy Touch" Foot Control for hands free adjustment of table
- Electrical Receptacles, 2 per side. Overload Protection: quick reset Circuit Breaker — no fuses
- Safety Lockout, accessible from either side of table, prevents accidental movement of table during exam or procedure, but maintains power to electrical receptacles
- Hand Control Connection Port, standard on both sides of table, allows hand control accessory to be added at any time
- Receptacles, Safety Lockout and Hand Control Ports located on table's moveable column for convenience and easy access
- Items plugged into side of table travel with the height of the table, eliminating over-stretched cords
- Adjustable Padded Headrest
- Multi-Position Positive Lock Stirrups
- 2" deep Stainless Steel Irrigation Basin
- Stainless Steel Headend Accessory Rails



Power adjustment of height (24"–40"), back (flat – 80°), tilt (flat – 30°) and foot sections (vertical – horizontal) allows table to be positioned for various exams and procedures including chair, contoured chair, Trendelenburg, Lithomy, OBGYN and flat.



Quick access to electrical receptacles and the safety lockout switch on the side of the table saves time. The innovative Lockout System disables the hand and foot control to eliminate accidental table movement during exams or procedures but retains power to table's receptacles.



# Specifications for the RESOLUTE 6220

<b>Warranty:</b>	4-YEAR
<b>Base Construction:</b>	Fully welded all steel base
<b>Power System:</b>	Electro-mechanical
<b>Attainable Positions:</b>	Chair, Contoured Chair, or OB/GYN, Lithotomy, Urology, Flat
<b>Power Adjustments:</b>	<b>Height:</b> 24"to 40", <b>Back:</b> horizontal 80° <b>Tilt:</b> horizontal to 30°, <b>Leg:</b> vertical to horizontal
<b>Programmable:</b>	Programmable with optional hand control accessory
<b>Hand Control:</b>	Optional accessory that may be added later
<b>Hand Control Connection Port:</b>	Standard feature, located on both sides of table's moveable column for easy access and elimination of over-stretched cords.
<b>Foot Control:</b>	"Easy touch" hands-free foot control Safety Extra-Low Voltage (SELV) and IPX1
<b>Upholstered Couch Top:</b>	One-piece style, no cross seams, embossing or folded fabric, securely fastened to all steel frame
<b>Headend Accessory Rails:</b>	Stainless steel rails, mounted on both sides of the table
<b>Irrigation Basin:</b>	2" deep stainless steel basin, integrated safety switch prohibits operation of foot section when basin is extended
<b>Electrical Receptacles:</b>	Four hospital grade electrical receptacles (2 per side) Overload Protection: Dedicated circuit breaker with quick reset
<b>Safety Lockout:</b>	Shuts power off to all articulating sections of table, but retains power to receptacles
<b>Stirrups:</b>	Positive lock, multi-position stirrups, length and lateral adjustment
<b>Paper Roll Holder:</b>	Accommodates 18" or 21" paper widths
<b>Accessories:</b>	Knee crutches, padded arm board, restraint straps, IV pole, urology basin with screen and drain, footend accessory rails, caster base, universal clamp, multi-positional clamp
<b>Mains Power Entry:</b>	Master ON/OFF, IEC connection, Overload Protection: Integrated circuit breaker with quick reset
<b>Duty Cycle:</b>	1 minute ON / 9 minutes OFF
<b>Ratings:</b>	<b>Table:</b> 115V 60 Hz 10 amps <b>Receptacles:</b> 115V 60 Hz <b>Total Receptacle Load:</b> 3.9 amps
<b>Table Dimensions:</b>	<b>Length:</b> 71", extends with headrest to 83" <b>Width top:</b> 30" at seat, tapers to 25" at headend <b>Headrest:</b> 19" wide, tapers to 15" <b>Legrest:</b> 19" wide, tapers to 18"
<b>Table Weight:</b>	485 lbs.
<b>Patient Load:</b>	500 lbs. evenly distributed on main couch
<b>Certifications:</b>	UL-60601-1, Can/CSA C22.3 No. 601-1