

ENOCHS RESOLUTE 6250

Programmable Procedural Table

PATIENT LOAD: 500 LBS.
evenly distributed on the table's main couch



TABLE FEATURES:

- 4-YEAR WARRANTY
- Programmable Hand Control with 4 memory positions
- "Easy Touch" foot control for hands free table positioning
- Power Adjustment of Back, Tilt, Foot and Height
- Safety Lockout to prevent accidental movement
- Circuit Breakers for overload protection-no fuses
- Electrical Receptacles, Safety Lockouts and Hand Control Ports located on table's movable column for convenience and easy access
- Cords plugged into table's receptacles, move up and down with the table, eliminating over-stretched cords
- Master ON/OFF
- Multi-Position Positive Lock Stirrups
- Stainless Steel Irrigation Basin
- Stainless Steel Headend Accessory Rails



Positioning includes chair, contoured chair, Trendelenburg, Lithomy, OBGYN and flat. A safety anti-pinch system stops movement of table if a pinch situation exists between headrest and floor.



Power adjustment of height (24"-40"), back (flat to 80°), tilt (flat to 30°) and foot sections (vertical to horizontal) allows the 6250 to be positioned for numerous exams and procedures.

Specifications for the RESOLUTE 6250

Warranty:	4-YEAR
Base Construction:	Fully welded all steel base
Power System:	Electro-mechanical
Attainable Positions:	Chair, Contoured Chair, Trendelenburg, OB/GYN, Lithotomy, Urology, Flat
Power Adjustments:	Height: 24" to 40", Back: horizontal to 80° Tilt: horizontal to 30°, Legrest: vertical to horizontal
Programmable:	Quick, easy set-up of 4 memory positions
Hand Control:	Easy to operate, simple to program table positions
Hand Control Connection Port:	Located on both sides of table's moveable column for easy access and elimination of over-stretched cords
Foot Control:	"Easy touch" hands-free foot control
Upholstered Couch Top:	One-piece style, no cross seams, embossing or folded fabric, securely fastened to all steel frame
Headend Accessory Rails:	Stainless steel rails, mounted on both sides of the table
Irrigation Basin:	2" deep stainless steel basin
Electrical Receptacles:	Four hospital grade electrical receptacles (2 per side) Overload Protection: Dedicated circuit breaker and quick reset
Safety Lockout:	Shuts power off to all articulating sections of tables, but retains power to receptacles
Stirrups:	Positive lock multi-position stirrups, length and lateral adjustment
Paper Roll Holder:	Accommodates 18" or 21" paper widths
Accessories:	Knee crutches, padded arm board, restraint straps, IV pole, urology basin with screen and drain, footend accessory rails, caster base, universal clamp, multi-positional clamp
Mains Power Entry:	Master ON/OFF, IEC connection Overload Protection: Integrated circuit breaker with quick reset
Duty Cycle:	1 minute ON / 9 minutes OFF
Ratings:	Table: 115V 60 Hz 10 amps Receptacles: 115V 60Hz Total Receptacle Load: 3.9 amps
Table Dimensions:	Length: 71", extends to 83" with headrest Width (top): 30" at seat, tapers to 25" at headend Headrest: 19" wide, tapers to 15" Legrest: 19" wide, tapers to 18"
Table Weight:	485 lbs.
Patient Load:	500 lbs. evenly distributed on main couch
Certifications:	UL Standard 60601-1 CAN/CSA C22.3 No. 601-1