

SPECIALTY PACKAGING CARTS

VLT150MOT (pictured)

Value Line Med Cart, key lock, internal narcotics box designed to hold MOT cards, pull-out shelf. Capacity of 150 Medicine-On-Time® cards, Standard Package



VLT360PC (not shown)

Value Line Medication Cart with key lock, adjustable dividers. Capacity of 150 MOT cards and internal narcotics box designed to hold punch cards



SL270MOT (pictured)

Standard Line cart with BeST® Lock on cabinet. Narcotics box with Illinois lock in top 9" drawer. Capacity of 270 Medicine-on-Time® cards.



MOT Narcotics Box



SL270MOT Drawer

SL100DSP (pictured)

Standard Line Cart with Best® Lock on cabinet. Includes pull-out shelf and internal drawer dividers designed to hold DISPILL® packaging. Capacity of over 100 cards.

SL100DSPN (not shown)

Adds internal narcotics box to SL100DSP

SL70DSP (pictured)

Standard Line Cart with Best® Lock on cabinet. Includes pull-out shelf and internal drawer dividers designed to hold DISPILL® packaging. Capacity of over 70 cards.

SL70DSPN (not shown)

Adds internal narcotics box to SL70DSP





650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Standard Line Medication Cart with Key Lock, Standard Package (MEDICINE ON TIME® Compatible) #SL270MOT

Specifications

Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Single column medication delivery cart
- Capacity of 270 Medicine On Time® cards
- Key lock system with Best ® brand lock on cabinet
- Illinois lock on narcotics drawer located in top 9" drawer, designed to hold MOT® cards
- Painted pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Full wrap-a-round vinyl bumper
- 5" ball-bearing premium casters, all full swivel, 2 with toe locking brake & 1 directional,
- Also available in a specialty package (#SL270MOTSP) which includes #68640-P1 top drawer tray with adjustable dividers, #AL689435 aluminum waste container and #684810 locking sharps container

Internal Drawer Dimensions:

- 29.5" W x 17" D (75cm x 43cm)

Drawer Configuration:

- 1 ea – 4" multipurpose storage drawer
- 1 ea – 9.75" drawer with one aluminum riveted divider and locking narcotics box.
- 2 ea – 9.75" drawer with two aluminum riveted dividers.

Assembled Cart Dimensions:

- 45.75" H x 36" W x 23.5" D
(116.2 H x 91.4 W x 59.7 D)
- Writing surface: 44" (111.8 cm)
- Cart weight: 175 lbs (100kg)
- Approx. shipping weight: 218 lbs.



Interior view of drawer



Internal narcotics box

Thirteen available color choices:

Sand, Gloss Beige, Burgundy,
Light Blue, Navy, Light Gray,
Mauve, Red, Teal, Yellow,
Hammer Tone Blue, Hammer Tone Gray
and Hammer Tone Green



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Value Line Medication Cart with Key Lock (MEDICINE ON TIME®) # VLT150MOT

Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Tall cabinet with key locking system
- Capacity of up to 150 Medicine On Time® cards
- Internal locking narcotics box made for M.O.T. packaging in bottom 9" drawer
- Pull out side shelf
- Bumper frame with projected front casters and (4) corner bumpers
- 5" ball-bearing premium casters, all full swivel – one with brake, one directional
- Painted steel top rail/push handle
- Removable plastic top
- Full extension ball-bearing drawer slides
- Cart is delivered fully assembled
- Pre-threaded accessory mounting holes
- Hard-baked powder coat finish
- Certified "Latex Free"
- Five year warranty

Thirteen available color choices:

(may be viewed at www.harloff.com)

- Sand, Gloss Beige, Burgundy, Light Blue, Yellow, Navy, Light Gray, Mauve, Red, Teal, Hammer Tone Blue, Hammer Tone Gray and Hammer Tone Green

Drawer Configuration:

- 1 - ea 3" multipurpose storage drawer
- 2 - ea 9.75" drawers with 1 adjustable divider
- 1 - ea 9.75" drawer with locking narcotics compartment made for M.O.T. packaging.

Internal Drawer Dimensions:

- 23" W x 17" D (58.5 cm W x 43.2 cm D)

Assembled Cart Dimensions:

- 46.75" H x 29.5" W x 22.5" D
(121.3 cm H x 75 cm W x 57.2 cm D)
- Writing surface: 43.25" (109.9 cm)
- Cart weight: 140 pounds (63.7 kg)
- Approx. shipping weight: 175 pounds (79.6 kg)



Value Line Drawers Have the Ability to Easily Convert in the Field:

- If your medication packaging changes, the drawers in your Value Line Cart can change
- Replace 1 large drawer with 3 small drawers



Unit Dose/Multi Dose Drawer



Bin Drawer



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

**Standard Line Medication Cart with Key Lock, Standard Package
(DISPILL USA® Compatible)
#SL100DSPN**

Specifications



Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Single column medication delivery cart
- Capacity of 100 plus DISPILL® cards
- Key lock system with Best® brand lock on cabinet
- Illinois lock on narcotics drawer located in top 13" drawer, designed to hold DISPILL® cards
- Painted pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Full wrap-a-round vinyl bumper
- 5" ball-bearing premium casters, all full swivel, 2 with toe locking brake & 1 directional
- Pre-threaded accessory mounting holes

Internal Drawer Dimensions:

- 29.5" W x 18" D (75cm x 45.7cm)

Drawer Configuration:

- 1 ea – 4" multipurpose storage drawer
- 1 ea – 3.25" multipurpose storage drawer
- 1 ea – 13" drawer with aluminum riveted dividers and locking narcotics box
- 1 ea – 13" drawer with aluminum riveted dividers

Thirteen available color choices:

Sand, Gloss Beige, Burgundy,
Light Blue, Navy, Light Gray,
Mauve, Red, Teal, Yellow,
Hammer Tone Blue, Hammer Tone Gray
and Hammer Tone Green

Assembled Cart Dimensions:

- 45.75" H x 36" W x 23.5" D
(116.2 H x 91.4 W x 59.7 D)
- Writing surface: 44" (111.8 cm)
- Cart weight: 175 lbs (100kg)
- Approx. shipping weight: 205 lbs. (93.2 kg)



650 Ford Street
Colorado Springs, CO 80915
Toll: (800)433-4064 * Ph: (719)637-0300
Fax: (719)597-8273 * www.harloff.com

Standard Line Medication Cart, DISPILL® Compatible, Standard Package #SL70DSPN

Specifications

Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Single column medication cart with single width internal narcotics box
- Capacity of 70 plus DISPILL® cards
- Internal dividers designed to hold DISPILL® packaging included
- BEST® brand lock on cabinet. Illinois lock on narcotics box
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Painted steel pull-out shelf
- Full wrap around vinyl bumper, protecting the cart
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish

Internal Drawer Dimensions:

- 20.5" W x 18" D (52.1cm x 45.7cm)

Drawer Configuration:

- 1 ea – 4" multipurpose storage drawer
- 1 ea – 3.25" multipurpose storage drawer
- 1 ea - 13" drawer with locking narcotics box
- 1 ea – 13" drawers with internal dividers

Assembled Cart Dimensions:

- 45.75" H x 26.5" W x 23.5" D
(116.2 cm x 67.3 cm x 59.7 cm)
- Writing Surface: 44" (111.8)
- Cart weight: 137 lbs (62.2 kg)
- Approx shipping weight: 167 lbs (75.9 kg)

Thirteen available color choices:

Sand, Gloss Beige, Burgundy,
Light Blue, Navy, Light Gray,
Mauve, Red, Teal, Yellow,
Hammer Tone Blue, Hammer Tone Gray
and Hammer Tone Green

