

ENOCHS

EXAMINING ROOM FURNITURE

TREATMENT
TABLES



Specifications for Treatment Tables



	49 EXTRA-WIDE TABLE	49D EXTRA-WIDE TABLE WITH DRAWERS	48 ECONOMY TABLE
Warranty:	4 year	4 year	4 year
Base Construction:	H-brace with steel support angles	H-brace with steel support angles	fully welded, all steel
Base Finish:	laminate finish	laminate finish	durable powder coat finish
Top Design:	extra-wide, padded, one-piece upholstered	extra-wide, padded, one-piece upholstered top	padded one-piece upholstered top
Height:	32"	32"	33"
Width:	32"	32"	24"
Length:	72"	72"	72"
Drawers:	not available on this model	two, 4" deep drawers with laminate fronts	not available on this model
Drawer Type:	N/A	high impact moulded drawers, rounded corners and sides, easy to clean	N/A
Drawer System:	N/A	steel reinforced side-guide suspension, engineered for smooth non-binding drawer movement	N/A
Drawer Dimensions:	N/A	21.5" W x 15" D x 4" H	N/A
Paper Roll Holder:	concealed, accommodates 18" or 21" paper widths	concealed, accommodates 18" or 21" paper widths	concealed, accommodates 18" or 21" paper widths
Patient Load (evenly distributed):	300 lbs.	300 lbs.	250 lbs.
Table Weight:	119 lbs.	139 lbs.	78 lbs.
Shipping:	shipped fully assembled	shipped fully assembled	shipped knocked-down-assembly required
Color Selection:	Top: Cambridge Blue, Earthen Brown, Emerald, Fawn, Ninja, Smoke, Windsor Wine, Ocean Jade Base finish: Grey, Beige, Natural Oak Drawer fronts: Grey, Beige, Natural Oak		Top: Cambridge Blue, Earthen Brown, Emerald, Fawn, Ninja, Smoke, Windsor Wine, Ocean Jade Base finish: Cream White, Pearl Grey