

TITAN SERIES

**DURABLE, RELIABLE
TEMPERATURE CONTROL**



EC250▲

- Capacity: 4-6 blankets
- Interior size: 2.5 ft³
- Exterior size: 21.9" H; 18.5" W; 23.4" D
- Countertop model
- Stackable

EC350▼

- Capacity: 10-12 blankets
- Interior size: 3.5 ft³
- Exterior size: 27.9" H; 18.5" W; 23.5" D
- Floor or countertop model
- Stackable



Shown with optional casters



Shown with optional left-hinge window door, lock and casters

EC750▲

- Capacity: 20-25 blankets
- Interior size: 7.5 ft³
- Exterior size: 35.5" H; 24.0" W; 25.6" D
- Stackable

EC1540▼

- Capacity: 40-50 blankets
- Interior size: 15.4 ft³; 7.7 ft³ per chamber
- Exterior size: 72.7" H; 23.9" W; 34.5" D



Shown with optional bumper

EC2060

- Capacity: 60-70 blankets
- Interior size: 20.6 ft³
- Exterior size: 74.6" H; 31.9" W; 29.1" D



**COMBINATION BLANKET &
FLUID WARMERS**

EC1540BL

- Total capacity:
 - Blanket chamber: 20-25 blankets; 7.7 ft³ interior size
 - Fluid chamber: 54 1-liter bottles or 42 1-liter bags
- Exterior size: 72.7" H; 25.2" W; 35.1" D



COMBINATION BLANKET & FLUID WARMERS

EC1260BL

- Total capacity:
 - Blanket chamber: 25-30 blankets; 8.7 ft³ interior size
 - Fluid chamber: 18 1-liter bottles or 28 1-liter bags
- Exterior size: 71.0" H; 31.0" W; 21.7" D
- 18" chamber depth for built-in applications



EC1730BL

- Total capacity:
 - Blanket chamber: 30-35 blankets; 11.8 ft³ interior size
 - Fluid chamber: 30 1-liter bottles or 24 1-liter bags
- Exterior size: 74.6" H; 31.9" W; 28.8" D



BLANKET WARMERS

YOU CAN AFFORD THE BEST



ARE YOUR BLANKETS
WARM ENOUGH?



DESIGNER SERIES

WARMSAFE™ ZONE HEAT – WARMS SAFELY, QUIETLY & EVENLY



DC150 ▲

- Capacity: 3-4 blankets
- Interior size: 1.5 ft³
- Exterior size: 16.9" H; 18.5" W; 21.8" D
- Stackable

DC250 ▼

- Capacity: 4-6 blankets
- Interior size: 2.5 ft³
- Exterior size: 21.9" H; 18.5" W; 24.8" D
- Stackable



DC350 ▲

- Capacity: 10-12 blankets
- Interior size: 3.5 ft³
- Exterior size: 27.9" H; 18.5" W; 24.8" D
- Stackable

DC400 ▼

- Capacity: 14-18 blankets
- Interior size: 4.0 ft³
- Exterior size: 21.9" H; 24.0" W; 27.0" D
- Stackable



DC750

- Capacity: 20-25 blankets
- Interior size: 7.5 ft³
- Exterior size: 34.9" H; 24.0" W; 27.0" D
- Stackable



STACKED BLANKET & FLUID WARMERS

DC400L over DC750

- Capacity: 14-18 blankets and 20 1-liter bottles or 28 1-liter bags
- Interior size: 4.0 ft³
- Exterior size: 56.9" H; 24.0" W; 27.0" D
- Numerous combinations available



MIX & MATCH

Blankets, irrigation fluids and injection fluids need to be warmed to different temperatures

Match the right temperature to the right product in the right place at the right time.

Only Enthermics offers such a full line of fluid and blanket warmers to stack to create the warming configuration to fit your needs.



DESIGNER & TITAN SERIES WARMERS

WarmSafe™ zone heat – warms safely, quietly & evenly

Provide the warmest possible blankets and fluids with Designer Series blanket warmers. These convenient, patient-friendly warmers are perfect for point-of-care placement.


Both Designer Series and durable Titan Series blanket warmers feature WarmSafe zone heating that ensure every blanket is safely and evenly warmed. With no fans or moving parts, they are engineered for years of trouble-free operation.



DC150



EC250 over EC750



Many more combinations available