

ilex₃₂

Portable image display when you just need a bigger monitor

Designed for use as a large table side monitor, the Ilex 32 combines practical features with easy use.

- 32" HD monitor with dual DVI, SDI and SOG inputs
- Counter weighted stand with minimal footprint
- Single, side-by-side or picture-in-picture imaging modes
- Typical video input sources; C Arm, Hemodynamics, and ultrasound
- Accessories Included:
 - Accessory basket
 - Cable hook
 - Transport handle

Stand

- Rugged counter-balanced base
- Four casters, two lockable and two non-lockable
- Adaptable column with vertical rails for accessory attachment
- Heavy duty VESA mount with tilt
- Cable sets: three DVI/HDMI, three SDI, one S Video, one VGA

Monitor

- FSN 32" (81cm) high-bright LED monitor
- Aspect Ratio - 16 x 9 " (41cm x 23cm)
- Width: 30.31" (77cm) x Height: 27.64" (70cm)
- Contrast Ratio – 1300:to 1
- Pixel Pitch - .27 (W) x .27 (H)
- Nits - 450
- Dual DVI, SDI, SOG, S-Video and composite video inputs
- Medical device bios software



A veteran-owned company