

# ENOCHS

EXAMINING ROOM FURNITURE

PHYSICIAN  
SEATING



# Specifications for Physician Seating



	92 OPERATOR STOOL	96 OPERATOR STOOL WITH BACK	101 PNUEMATIC STOOL	106 PNUEMATIC STOOL WITH BACK
<b>Warranty:</b>	4 year	4 year	4 year	4 year
<b>Base Construction:</b>	fully welded, tubular steel base	fully welded, tubular steel base	cast aluminum	cast aluminum
<b>Base Finish:</b>	durable powder coated finish	durable powder coated finish	polished aluminum	polished aluminum
<b>Number of Legs:</b>	5	5	5	5
<b>Foot Ring:</b>	yes	yes	no*	no*
<b>Caster Style:</b>	2" ball bearing, hard surface caster	2" ball bearing, hard surface caster	2" carpet caster	2" carpet caster
<b>Adjustment System:</b>	threaded stem	threaded stem	pneumatic gas cylinder	pneumatic gas cylinder
<b>Adjustment By:</b>	rotation of seat	rotation of seat	hand lever under seat	hand lever under seat
<b>Height:</b>	18.5" to 25"	18.5" to 25"	20" to 30"	20" to 30"
<b>Back Support:</b>	no	yes	no	yes
<b>Back Adjustment:</b>	N/A	3" vertical	N/A	3" vertical
<b>Seat Dimensions:</b>	14"W x 2.5" H	14"W x 2.5" H	14"W x 2.5" H	14"W x 2.5" H
<b>Weight Capacity:</b>	350 lbs.	225 lbs.	350 lbs.	225 lbs.
<b>Stool Weight:</b>	16 lbs.	21 lbs.	13 lbs.	18 lbs.
<b>Shipping:</b>	assembly required	assembly required	assembly required	assembly required
<b>Color Selection:</b>	Seat: Black, Cambridge Blue, Desert Sand, Earthen Brown, English Brown, Emerald, Fawn, Ninja, Ocean Jade, Smoke, Windsor Wine			

\*available as an add on accessory