

# TITAN SERIES

**DURABLE, RELIABLE  
TEMPERATURE CONTROL**



Shown with optional casters

## EC2060

- Capacity: 60-70 blankets
- Interior size: 20.6 ft<sup>3</sup>
- Exterior size: 74.6" H; 31.9" W; 29.1" D



**COMBINATION BLANKET &  
FLUID WARMERS**

## EC1540BL

- Total capacity:
  - Blanket chamber: 20-25 blankets; 7.7 ft<sup>3</sup> interior size
  - Fluid chamber: 54 1-liter bottles or 42 1-liter bags
- Exterior size: 72.7" H; 25.2" W; 35.1" D



# BLANKET WARMERS

YOU CAN AFFORD THE BEST



## EC250▲

- Capacity: 4-6 blankets
- Interior size: 2.5 ft<sup>3</sup>
- Exterior size: 21.9" H; 18.5" W; 23.4" D
- Countertop model
- Stackable

## EC350▼

- Capacity: 10-12 blankets
- Interior size: 3.5 ft<sup>3</sup>
- Exterior size: 27.9" H; 18.5" W; 23.5" D
- Floor or countertop model
- Stackable



## EC750▲

- Capacity: 20-25 blankets
- Interior size: 7.5 ft<sup>3</sup>
- Exterior size: 35.5" H; 24.0" W; 25.6" D
- Stackable

## EC1540▼

- Capacity: 40-50 blankets
- Interior size: 15.4 ft<sup>3</sup>; 7.7 ft<sup>3</sup> per chamber
- Exterior size: 72.7" H; 23.9" W; 34.5" D



Shown with optional bumper

**COMBINATION BLANKET & FLUID WARMERS**

## EC1260BL

- Total capacity:
  - Blanket chamber: 25-30 blankets; 8.7 ft<sup>3</sup> interior size
  - Fluid chamber: 18 1-liter bottles or 28 1-liter bags
- Exterior size: 71.0" H; 31.0" W; 21.7" D
- 18" chamber depth for built-in applications



## EC1850BL

- Total capacity:
  - Blanket chamber: 30-35 blankets; 13ft<sup>3</sup> interior size
  - Fluid chamber: 30 1-liter bottles or 24 1-liter bags
- Exterior size: 74.8" H; 30.0" W; 28.0" D



**ENTHERMICS MEDICAL SYSTEMS**

An ISO 13485:2003 certified company  
Tel 262-251-8356 | 800-TO-B-WARM  
generalinfo@enthermics.com  
www.enthermics.com

ARE YOUR BLANKETS  
WARM ENOUGH?



## DESIGNER SERIES

WARMSAFE ZONE HEAT – WARMS SAFELY, QUIETLY & EVENLY



### DC150 ▲

- Capacity: 3-4 blankets
- Interior size: 1.5 ft<sup>3</sup>
- Exterior size: 16.9" H; 18.5" W; 21.8" D
- Stackable

### DC250 ▼

- Capacity: 4-6 blankets
- Interior size: 2.5 ft<sup>3</sup>
- Exterior size: 21.9" H; 18.5" W; 24.8" D
- Stackable



### DC350 ▲

- Capacity: 10-12 blankets
- Interior size: 3.5 ft<sup>3</sup>
- Exterior size: 27.9" H; 18.5" W; 24.8" D
- Stackable

### DC400 ▼

- Capacity: 14-18 blankets
- Interior size: 4.0 ft<sup>3</sup>
- Exterior size: 21.9" H; 24.0" W; 27.0" D
- Stackable



### DC750

- Capacity: 20-25 blankets
- Interior size: 7.5 ft<sup>3</sup>
- Exterior size: 34.9" H; 24.0" W; 27.0" D
- Stackable

## MIX & MATCH

Blankets, irrigation fluids and injection fluids need to be warmed to different temperatures

Match the right temperature to the right product in the right place at the right time.

Only Enthermics offers such a full line of stackable fluid and blanket warmers to create the warming configuration to fit your needs.



## STACKED BLANKET & FLUID WARMERS



### DC400 L over DC750

- Capacity: 14-18 blankets and 20 1-liter bottles or 28 1-liter bags
- Interior size: 4.0 ft<sup>3</sup>
- Exterior size: 56.9" H; 24.0" W; 27.0" D
- Numerous combinations available

## DESIGNER & TITAN SERIES WARMERS

WarmSafe™ zone heat – warms safely, quietly & evenly

Provide the warmest possible blankets and fluids with Designer Series blanket warmers. These convenient, patient-friendly warmers are perfect for point-of-care placement.


Both Designer Series and durable Titan Series blanket warmers feature WarmSafe zone heating that ensure every blanket is safely and evenly warmed. With no fans or moving parts, they are engineered for years of trouble-free operation.



DC150



EC750 with optional casters



Additional combinations available