

 enthermics®

An ISO 13485:2003 certified company



PATIENT WARMING  
TECHNOLOGIES

[www.enthermics.com](http://www.enthermics.com)



# DESIGNER SERIES

## **WarmSafe™ zone heat – warms safely, quietly & evenly**

Provide the warmest possible blankets and fluids with these convenient, patient-friendly warmers for point-of-care placement.

Designer Series blanket warmers with WarmSafe zone heating system ensure every blanket is safely and evenly warmed. With no fans or moving parts, they are engineered for years of trouble-free operation.

Designer Series fluid warmers maintain precise temperature for IV and irrigation fluids.



DC350



DC400L over DC750

# DESIGNER SERIES

WarmSafe™ zone heat – warms safely, quietly & evenly

## CABINET WARMERS



- **Patient-friendly** combination of warmth and style
- **Improved safety** with WarmSafe™ zone heat
- **Energy efficient** temperature control for low energy usage and minimal heat loss
- **Easy to clean** epoxy-coated interior surfaces
- **Simple loading** with hands-free push button door
- **Stackable design** allows for customized configurations in a small footprint
- **No fans or moving parts** means less maintenance\*

**Most commonly used in...** Designer Series blanket warmers have no fans and operate virtually silently, making them perfect for patient rooms and point-of-care in departments like L&D and NICU where patients need warmth and quiet. Fluid warmers are great for point-of-care use in departments like OR, L&D, Dialysis, Oncology, ICU and NICU.

\*Blanket warmers only

## Options & Capacities

	Fluid Warmer?	Blanket Warmer?	1L bottle capacity	1L bag capacity	Blanket capacity	Countertop or floor standing?	Stackable?
DC150	✗	✓	✗	✗	3-4	Counter	✓
DC250	✓	✓	12	16	4-6	Counter	✓
DC350	✗	✓	✗	✗	10-12	Floor	✓
DC400	✓	✓	20	28	14-18	Both	✓
DC750	✗	✓	✗	✗	20-25	Floor	✓



DC350

Blanket warmer also available in this size



DC250L



DC150



DC750

Blanket warmer also available in this size



DC400L

# IMPROVE CLINICAL OUTCOMES & BOOST HCAHPS SCORES



## TITAN SERIES

### **Durable, reliable temperature control**

Experience the benefits of patient warming during medical procedures.

Titan Series blanket warming cabinets provide safe, efficient and even heating. The durable, affordable stainless steel cabinets are built for hard work and years of dependable service.

Designed for busy nurses and aides, the warmers come in various sizes and options to best meet your warming and space needs.



EC1850BL



EC750



EC250L

# TITAN SERIES

Durable, reliable temperature control

## CABINET WARMERS



EC250L over EC350

- **Affordable, rugged stainless steel cabinets** are built for years of worry free service
- **Improved safety** with WarmSafe™ zone heat
- **Easy-grip handles** are designed for heavy use
- **Heavy-duty casters** for mobile usage
- **Optional bumpers**
- **Optional locking doors**
- **Insulated glass** increases energy efficiency
- **Over-temp alarm** alerts of unsafe operating conditions
- **Pass-through** options available
- **No fans or moving parts** means less maintenance\*

**Most commonly used in...** Small, countertop fluid warmers work well in the ED, ICU, Dialysis and NICU. Larger fluid and blanket warmers are perfect for high-volume areas including pre-OP, PACU, L&D, Oncology, Cath Lab and other areas.

\*Blanket warmers only



## Options & Capacities

	Fluid Warmer?	Blanket Warmer?	1L bottle capacity	1L bag capacity	Blanket capacity	Countertop or floor standing?	Stackable?
EC250	✓	✓	8	16	4-6	Counter	✓
EC350	✓	✓	24	24	10-12	Both	✓
EC750	✗	✓	✗	✗	20-25	Floor	✓
EC770L	✓	✗	54	42	✗	Floor	✓
EC1260BL	✓	✓	18	28	25-30	Floor	✗
EC1540	✗	✓	✗	✗	40-50	Floor	✗
EC1540BL	✓	✓	54	42	20-25	Floor	✗
EC1850BL	✓	✓	30	24	30-35	Floor	✗
EC2060	✗	✓	✗	✗	60-70	Floor	✗



EC350L



EC2060

WARM FLUIDS  
WHEN & WHERE  
YOU NEED THEM  
without disposables





## ivNow®

### **Modular Fluid Warmers**

Improve processes and recovery time by efficiently warming patients from admission to discharge. Convenient, easy-to-use ivNow Fluid Warmers quickly heat and maintain safe temperatures of intravenous fluids while saving space and staff time.

ivNow fluid warmers record each bag's warming time to increase safety, comply with healthcare regulations and reduce waste. Save space in patient rooms, surgery centers, and labor and deliver units with these compact units. Simplify procedures by keeping fluids warm and ready to go at the point of care.

# ivNow®

Records the shelf life and temperature of every IV bag



ivNow-1



ivNow-2



ivNow-3

## OPTIONS

Warm fluids easily and cost effectively in every department. Flexible mounting options make ivNow well suited for L & D ED, OR, PACU, ICU, NICU — and everywhere in between.

ivNows can be placed on the counter top, mounted to the wall, attached to heavy-duty mobile equipment pole stands and mast arms, and even anesthesia carts.




- Know the time and temperature of every IV bag
- Warms 1 liter in 30 minutes or less
- Warms 20-140 liters in 12 hours
- Set the expiration date of each bag
- Point-of-care use saves space, time & footsteps
- Holds 1L -3L IV bags
- Requires no disposables
- Modular set-up



- Increase regulatory compliance**
- Save money**
- Improve processes**
- Increase reimbursements**



# PERITONEAL DIALYSIS



ivNow works well as a PD fluid warmer too! Avoid patient discomfort such as cramping and chills with simple, precise, cost-effective fluid warming.

- Install almost anywhere: on counters or walls to keep them adjacent to your exchange area
- Warms PD fluid with 40 minutes
- Pre-set, tamper-proof heating element quickly warms fluids and keeps them at the perfect temperature

# TECHNOLOGY

## WARMSAFE™

### SAFE AND RELIABLE HEAT

WarmSafe™ heating technology, exclusive to Entermics warming cabinets incorporates a multiple zone warming technology that heats where and when it is needed.


- **Improved safety** features
- **Even heat** throughout the cavity
- **Silent operation** without fans in blanket warmers


## STACKING OPTIONS

### Blankets, irrigation fluids and injection fluids need to be warmed to different temperatures




Match the right temperature to the right product in the right place at the right time.

Only Entermics offers such a full line of fluid and blanket warmers to stack to create the warming configuration to fit your needs.





**Blanket and fluid warmers** can be stacked in a wide variety of configurations. However, it is recommended that fluid warming units be placed at, or near, the bottom of the stack whenever possible for maximum stability.





Enthermics is committed to providing the best value and return on investment for our customers. Our never ending drive for innovation and quality has made Enthermics a leader in patient warming. As an ISO 13485:2003 certified company, we provide sales and technical support to our customers worldwide. Each warming unit is proudly manufactured and inspected in the U.S.A.



## ENTHERMICS MEDICAL SYSTEMS

An ISO 13485:2003 certified company

PO Box 443, Menomonee Falls WI 53052-0443  
W164 N9221 Water St, Menomonee Falls WI 53051  
Phone 262-251-8356 | 800-862-9276 | Fax 262-251-7067  
generalinfo@enthermics.com | www.enthermics.com

Made in the USA