

CHAMPION™

Quality Healthcare Seating Products

Clinical Stools

CAREGIVER SEATING

RATED UP TO 300 LBS.

Part #507

Part #509




STANDARD FEATURES

SERIES #	507	509
Adjustable seat – hand lever	•	
Adjustable seat – foot ring		•
Mounts on a 22" pedestal	•	•
Composite base	•	
Polished aluminum base		•

SPECIFICATIONS

SERIES #	507	509
Weight capacity	300 lbs.	300 lbs.
Overall width – wheel base	22"	22"
Overall height – with backrest	34"-43"	34"-43"
Seat width	15"	15"
Seat height – adjustable	21"-28"	21"-28"

Champion's Caregiver Clinic stools are all-purpose seats to meet your exact needs. All stools offer a pneumatic cylinder for easily adjustable height options along with either a composite or polished aluminum base. Three and a half inch foam padding, a 15-inch seat and sealed seams provide caregivers with comfort and durability. Optional backrests are available on all clinic stools. Available in a wide range of standard and up-graded vinyls, they can be designed to complement any interior.




OPTIONS

- Backrest (B)
- 5 - Brake casters (5BR)
- Glides (G) - available on 507 model only
- Vinyl upgrades
 - Champion, graded vinyl selections*
 - Customer's own material (COM)
 - PVC-free material

*For available patterns and colors please contact your Champion sales representative.




GLIDES



BACKREST
(shown on 509)

STANDARD VINYL SELECTION



Printed color swatches may not be an exact match to actual fabrics and are enlarged to show detail.

CHAMPION™

Champion Manufacturing, Inc.
Elkhart, IN 46516

1.800.998.5018
www.championchair.com

Information contained in this literature is based on the latest product data available at the time of printing. Photography, artwork, text and specifications are subject to change. Photography may show optional items. Due to photographic and printing processes, product colors in this brochure are not necessarily accurate. Champion Manufacturing, Inc. reserves the right to make changes to colors, materials, prices, standards, options, equipment and any specification at any time without notice or obligation. Contact Champion direct for additional information and the latest specifications.