

# Apron Options

- Fabrics
- Embroidery
- Tu-Tone Covering (specify Inside and Outside fabric choice)
- Additional "QR" Quick Release Buckle (specify 1" or 2")
- Built-in Vest Belt
- Skirt Suspenders
- Shoulder Strap
- No-Slip Skirt Backing
- Extension Panel

Custom Binding Colors			
 Silver	 Black	 Brown	 Purple
 Maroon	 Gold	 Neon Orange	 Red
 White	 Navy	 Slate	 Royal
 Aqua	 Green	 Orange	 Pink

Bindings



Optional Sleeve (model 812)



Tethered Collar (model T)




Attached Collar (model 600A)



Pockets (model PKT)



Removable Name Tag (model NTG)



No-Slip Skirt Backing



Extension Panel for Vest or Skirt



# Ordering an apron is easy!

Select your Model No.

Use loose measurements to find the correct apron width and length from the charts below.

Choose your protection from charts on facing page.

## Example Part # for Frontal and TQW Aprons Standard Coat Apron with Ties - Reg. Lead

Model - Width X Length - Protection  
ST - 24 X 36 - RL

## Example Part # for Wrap Around Apron Wrap Around - Small Female - TL Lite Ply

Model - Size X Length - Protection  
WA - SF X 36 - TL

### Frontal Aprons

Apron Width	Chest
18" 45.72 cm	30" - 34" 76.2 cm - 86.4 cm
20" 50.8 cm	34" - 38" 86.4 cm - 96.5 cm
22" 55.88 cm	38" - 42" 96.5 cm - 106.7 cm
24" 60.96 cm	42" - 46" 106.7 cm - 116.8 cm
27" 68.58 cm	47" - 52" 119.4 cm - 132.1 cm

### Three Quarter Wrap Aprons

Apron Width	Chest	Wrap Width
18" 45.72 cm	30" - 34" 76.2 cm - 86.4 cm	4" 10.2 cm
20" 50.8 cm	34" - 38" 86.4 cm - 96.5 cm	4" 10.2 cm
22" 55.88 cm	38" - 42" 96.5 cm - 106.7 cm	5" 12.7 cm
24" 60.96 cm	42" - 46" 106.7 cm - 116.8 cm	5" 12.7 cm
27" 68.58 cm	47" - 52" 119.4 cm - 132.1 cm	6" 15.2 cm

### Unisex Wrap Around Aprons

### Male Wrap Around Aprons

### Female Wrap Around Aprons

Size	Chest	Waist	Hips up to	Size	Chest	Waist	Hips up to	Size	Chest	Waist	Hips up to
Small (S)	34" - 37" 86.4 cm - 94 cm	27" - 31" 68.6 cm - 78.7 cm	39" 99.1 cm	Small (SM)	35" - 38" 88.9 cm - 96.5 cm	28" - 32" 71.1 cm - 81.3 cm	40" 101.6 cm	Small (SF)	34" - 35" 86.4 cm - 88.9 cm	26" - 30" 66.4 cm - 76.2 cm	38" 96.5 cm
Medium (M)	38" - 41" 96.5 cm - 104.1 cm	31" - 35" 78.7 cm - 88.9 cm	43" 109.2 cm	Medium (MM)	39" - 42" 99.1 cm - 106.7 cm	32" - 36" 81.3 cm - 91.4 cm	44" 111.8 cm	Medium (MF)	36" - 37" 91.4 cm - 94 cm	30" - 34" 76.2 cm - 86.36 cm	40" 101.6 cm
Large (L)	42" - 45" 106.7 cm - 114.3 cm	35" - 39" 88.9 cm - 99.1 cm	47" 119.4 cm	Large (LM)	43" - 46" 109.2 cm - 116.8 cm	36" - 40" 91.4 cm - 101.6 cm	48" 129.5 cm	Large (LF)	38" - 40" 96.5 cm - 101.8 cm	34" - 38" 86.36 cm - 96.52 cm	43" 109.2 cm
X-Large (XL)	46" - 49" 116.8 cm - 124.5 cm	39" - 43" 99.1 cm - 109.2 cm	51" 129.5 cm	X-Large (XLM)	47" - 50" 119.4 cm - 127 cm	40" - 44" 101.6 cm - 111.8 cm	52" 132.1 cm	X-Large (XLF)	41" - 44" 104.1 cm - 111.8 cm	38" - 42" 96.52 cm - 106.68 cm	47" 119.4 cm
XX-Large (XXL)	50" - 53" 127 cm - 134.6 cm	43" - 47" 109.2 cm - 119.4 cm	55" 139.7 cm	XX-Large (XXLM)	51" - 54" 129.5 cm - 137.2 cm	44" - 48" 111.8 cm - 121.9 cm	56" 142.2 cm	XX-Large (XXLF)	45" - 48" 114.3 cm - 121.9 cm	42" - 46" 106.68 cm - 116.84 cm	51" 129.5 cm

## Sizing

Your Approximate Height	4'11" - 5'3" 149.9 cm - 160 cm	5'3" - 5'5" 160 cm - 165.1 cm	5'5" - 5'8" 165.1 cm - 172.7 cm	5'8" - 5'11" 172.7 cm - 180.3 cm	5'11" - 6'2" 180.3 cm - 187.9 cm
Choose Length	30" 76.2 cm	32" 81.3 cm	34" 86.4 cm	36" 91.4 cm	38" 96.5 cm

Custom made (non-stock) aprons are non-returnable.

### Unisex Vest Skirt

Size	Chest	Vest Length	Waist	Hips	Skirt Length
Small (S)	34" - 37" 86.4 cm - 94 cm	18" 45.7 cm	27" - 31" 68.6 cm - 78.7 cm	39" 99.1 cm	21" 53.3 cm
Medium (M)	38" - 41" 96.5 cm - 104.1 cm	19" 48.3 cm	31" - 35" 78.7 cm - 88.9 cm	43" 109.2 cm	22" 55.9 cm
Large (L)	42" - 45" 106.7 cm - 114.3 cm	20" 50.8 cm	35" - 39" 88.9 cm - 99.1 cm	47" 119.4 cm	22" 55.9 cm
X-Large (XL)	46" - 49" 116.8 cm - 124.5 cm	21" 53.3 cm	39" - 43" 99.1 cm - 109.2 cm	51" 129.5 cm	22" 55.9 cm
XX-Large (XXL)	50" - 53" 127 cm - 134.6 cm	22" 55.9 cm	43" - 47" 109.2 cm - 119.4 cm	55" 139.7 cm	23" 58.4 cm

### Example Part # for Vest Skirt Combo - TL Lite Ply

Vest Skirt - Small Male - TL Lite Ply  
Model - Size - Protection  
VS - SM - TL


### Male Vest Skirt

Size	Chest	Vest Length	Waist	Hips up to	Skirt Length
Small (SM)	35" - 38" 88.9 cm - 96.5 cm	18" 45.7 cm	28" - 32" 71.1 cm - 81.3 cm	40" 101.6 cm	21" 53.3 cm
Medium (MM)	39" - 42" 99.1 cm - 106.7 cm	19" 48.3 cm	32" - 36" 81.3 cm - 91.4 cm	44" 111.8 cm	22" 55.9 cm
Large (LM)	43" - 46" 109.2 cm - 116.8 cm	20" 50.8 cm	36" - 40" 91.4 cm - 101.6 cm	48" 129.5 cm	22" 55.9 cm
X-Large (XLM)	47" - 50" 119.4 cm - 127 cm	21" 53.3 cm	40" - 44" 101.6 cm - 111.8 cm	52" 132.1 cm	22" 55.9 cm
XX-Large (XXLM)	51" - 54" 129.5 cm - 137.2 cm	22" 55.9 cm	44" - 48" 111.8 cm - 121.9 cm	56" 142.2 cm	23" 58.4 cm

### Female Vest Skirt

Size	Chest	Vest Length	Waist	Hips up to	Skirt Length
Small (S)	34" - 35" 86.4 cm - 88.9 cm	17" 43.2 cm	26" - 30" 66.04 cm - 76.2 cm	38" 96.5 cm	20" 50.8 cm
Medium (M)	36" - 37" 91.4 cm - 94 cm	18" 45.7 cm	30" - 34" 76.2 cm - 86.36 cm	40" 101.6 cm	20" 50.8 cm
Large (L)	38" - 40" 96.5 cm - 101.8 cm	19" 48.3 cm	34" - 38" 86.36 cm - 96.52 cm	43" 109.2 cm	21" 53.3 cm
X-Large (XL)	41" - 44" 104.1 cm - 111.8 cm	20" 50.8 cm	38" - 42" 96.52 cm - 106.68 cm	47" 119.4 cm	21" 53.3 cm
XX-Large (XXL)	45" - 48" 114.3 cm - 121.9 cm	21" 53.3 cm	42" - 46" 106.68 cm - 116.84 cm	51" 129.5 cm	22" 55.9 cm

# Finding Your Fit and Protection




(C) CHEST: Measure widest part 1" below armhole  
 (W) WAIST: Measure narrowest part 1" above navel  
 (H) HIPS: Measure widest part of hips  
 LENGTH: We recommend 1" above the knee

Can't find your size?  
**Customize!**

Additional charges may apply for size customization. Contact us. All custom apparel is non returnable/non refundable

Model No. \_\_\_\_\_  
 Chest \_\_\_\_\_  
 Waist \_\_\_\_\_  
 Hips \_\_\_\_\_  
 Length \_\_\_\_\_  
 Length Vest \_\_\_\_\_ Skirt \_\_\_\_\_  
 Protection RL \_\_\_\_\_ TL \_\_\_\_\_ NL \_\_\_\_\_  
 RLL \_\_\_\_\_ TLL \_\_\_\_\_ NLL \_\_\_\_\_



## Protection Options

<b>FRONTAL &amp; HALF APRONS</b> (RLL, TL, TLL, NL, NLL) Add Proper Initials for Optional Protection when ordering, otherwise RL Regular Lead is standard.	<b>RL – REGULAR LEAD<sup>1</sup></b> 0.5mm Pb Eq.	<b>TL – TL Lite Ply<sup>1</sup></b> 0.5mm Pb Eq.	<b>NL – Non Lead<sup>2</sup></b> 0.5mm Pb Eq.
	<b>RLL<sup>1</sup></b> 0.35mm Pb Eq.	<b>TLL<sup>1</sup></b> 0.35mm Pb Eq.	<b>NLL<sup>2</sup></b> 0.35mm Pb Eq.

<b>DENTAL APRONS</b> (TL, NL) Add Proper Initials for Optional Protection when ordering, otherwise RL Regular Lead is standard.	<b>RL – REGULAR LEAD<sup>1</sup></b> 0.3mm Pb Eq.	<b>TL – TL Lite Ply<sup>1</sup></b> 0.25mm Pb Eq.	<b>NL – Non Lead<sup>2</sup></b> 0.25mm Pb Eq.
	<i>Alternate protection also available in 0.5mm and 0.35mm Pb equivalencies.</i>		

<b>VESTS &amp; SKIRTS</b> (RLL, TL, TLL, NL, NLL) Add Proper Initials for Optional Protection when ordering, otherwise RL Regular Lead is standard.  <i>*Overlap does not apply to Flex Back Vest</i>	<b>RL – REGULAR LEAD<sup>1</sup></b>			<b>TL – TL Lite Ply<sup>1</sup></b>		<b>NL – Non Lead<sup>2</sup></b>	
	Vest:	Front	Front Overlap*	Back	Front	Front Overlap*	Back
	Skirt:	0.5mm Pb Eq. 0.3mm Pb Eq.	1.0mm Pb Eq. 0.6mm Pb Eq.	0.3mm Pb Eq. 0.3mm Pb Eq.	0.5mm Pb Eq. 0.25mm Pb Eq.	1.0mm Pb Eq. 0.5mm Pb Eq.	0.25mm Pb Eq. 0.25mm Pb Eq.
	Vest:	<b>RLL<sup>1</sup></b>			<b>TLL<sup>1</sup></b>		<b>NLL<sup>2</sup></b>
Skirt:	Front	Front Overlap*	Back	Front	Front Overlap*	Back	
	0.35mm Pb Eq. 0.25mm Pb Eq.	0.7mm Pb Eq. 0.5mm Pb Eq.	0.25mm Pb Eq. 0.25mm Pb Eq.	0.35mm Pb Eq. 0.25mm Pb Eq.	0.7mm Pb Eq. 0.5mm Pb Eq.	0.25mm Pb Eq. 0.25mm Pb Eq.	

<b>WRAP AROUND APRONS</b> (RLL, TL, TLL, NL, NLL) Add Proper Initials for Optional Protection when ordering, otherwise RL Regular Lead is standard.  <i>*Overlap does not apply to Flex Back Wrap Around</i>	<b>RL – REGULAR LEAD<sup>1</sup></b>			<b>TL – TL Lite Ply<sup>1</sup></b>		<b>NL – Non Lead<sup>2</sup></b>	
	Front	Front Overlap*	Back	Front	Front Overlap*	Back	
	0.5mm Pb Eq.	1.0mm Pb Eq.	0.3mm Pb Eq.	0.5mm Pb Eq.	1.0mm Pb Eq.	0.25mm Pb Eq.	
	<b>RLL<sup>1</sup></b>			<b>TLL<sup>1</sup></b>		<b>NLL<sup>2</sup></b>	
Front	Front Overlap*	Back	Front	Front Overlap*	Back		
0.35mm Pb Eq.	0.7mm Pb Eq.	0.25mm Pb Eq.	0.35mm Pb Eq.	0.7mm Pb Eq.	0.25mm Pb Eq.		

<sup>1</sup> Denotes Lead Equivalency of RL & RLL Regular Lead and TL & TLL Lite Ply @ 100kVp.

<sup>2</sup> Denotes Lead Equivalency of NL & NLL Non Lead @ 90kVp.

NOTE: For Proper protection in your working environment please contact a radiation physicist regarding adequate apron Pb equivalency.

Contact us regarding prices on alternate protection options. Discounts may apply when reducing Pb equivalency.





# Prints Perfect for your Personality and Profession

## DESIGNER FABRICS (DF)



Orange Floral (DF)



Bermuda (DF)



Flower Power (DF)



Hearts of Color (DF)



Pastel Butterfly (DF)



Farm Animals (DF)



Chocolate Floral (DF)



Grey Paisley (DF)



Green Surf (DF)



Tiki Blue (DF)



Blue Metallic (DF)



Denim (DF)



Stars (DF)



Red and White Stars (DF)



Patriotic (DF)



White Lil' Cherries (DF)



Lollipop (DF)



Red Flames (DF)



Splatter (DF)



Pink Camo (DF)



Desert Camo (DF)



Storm Camo (DF)



Digi Camo (DF)



Grey Digi Camo (DF)



Royal Paws (DF)



Dog Paws (DF)



Leopard (DF)



Wild Cat Pink (DF)



Zebra (DF)



Blue Squares (DF)



Royal Plaid (DF)



Blue Plaid (DF)



Green Plaid (DF)



Hot Pink Plaid (DF)



Red Plaid (DF)



Grey Squares (DF)

## DIAMOND RIPSTOP (DRS)



Silver (DRS)



Charcoal (DRS)



Lavender (DRS)



Purple (DRS)



Fuchsia (DRS)



Cinnamon (DRS)



Red (DRS)



Slate (DRS)



Marine (DRS)



Royal Blue (DRS)



Kelly Green (DRS)



Tan (DRS)



Brown (DRS)



Black (DRS)

## FABRIC AND BINDING EXCEPTIONS

Fabric and binding colors are approximate representations of actual colors/patterns available. Economy aprons (**Model E**) are available in Vinyl (V), Nylon (NY) or Diamond Ripstop (DRS) only. Diamond Ripstop (DRS) and Designer Fabrics (DF) feature a black nylon inside cover and black binding. Weblon (W) aprons feature a black or corresponding inside nylon cover. Black binding is standard on most apparel, except for royal binding on royal nylon and light blue vinyl, unless otherwise indicated from binding selection at time of order. Two tone coverings subject to up charge. To avoid snagging or tearing of fabric, all hook and loop fasteners must be secured properly.

## Fabrics & Colors

### WEBLON (W)

\* Anti Microbial Healthcare Fabric




Royal Blue (W) Lt. Blue (W) \* Tan (W) \* Gray (W) \* Green (W) \*

### NYLON (NY)



### VINYL (VI)



Lt. Blue (VI) Royal Blue (VI) Mauve (VI) Beige (VI) Brooks Gray (VI)

### RIPSTOP (RS)



Raspberry (RS) Royal Blue (RS) Kelly Green (RS) Gold (RS) Maroon (RS)

### BINDINGS



Silver Black Brown Purple Maroon Gold Neon Orange Red White Navy  
Slate Royal Aqua Green Orange Pink



# Stock Designs

For additional designs (under 10,000 stitches) see Dakota Collectibles website: [www.dakotacollectibles.com](http://www.dakotacollectibles.com)



Boots



Bone & Paw



Butterfly



Kid's Hands



Balloons



Smiley



Football Helmet



Golf Crest



Sports



Sportsman



Frog



Flip Flops



Swirl Sun



Tropical Fish



Palm Trees



Cherries



Chocolate Flower



Light House



Heart Bear



Heart



Bandaged Heart



EKG Heart



Heart Scope



Ribbon



Skull and Bones



X-Ray Man



Skeleton



Radiology



Neuro



Neurology



Cardiology



Orthopedic



Bone & Joint



Spine



Caducus



Dental Tools



RN Logo



X-Ray



Cath Lab



IR



# Custom Options

## Pocket Styles



Patch Pocket



Slit Pocket



ID Pocket

## Other Embroidery Options



Removable Name Tag



Bottom Corner



Back Shoulder

## Apron Options

- Attached Thyroid Collar (as shown on page 19)
- Tethered Thyroid Collar (as shown on page 19)
- Pocket (Specify style: Patch, Slit or ID)
- Tu-Tone Covering (Specify Inside & Outside Fabric Choice)
- Removable Name Tag
- Additional 'QR' Quick-Release Buckle (Specify 1" or 2")
- Lead Sleeve (Specify Sew On, Hook & Loop or Snap Fastener - as shown on page 19)
- Built-in Vest Belt
- Skirt Suspenders
- **A)** No-Slip Skirt Backing
- **B)** Extension Panel (Specify Vest or Skirt - as shown this page)



\*Example of custom logo option

## Embroidery Options

- Personalization - Two lines per location for one price!
- Stock Designs - Select from a wide variety of designs on opposite page! (Also see [www.dakotacollectibles.com](http://www.dakotacollectibles.com) for additional designs under 10,000 stitches)
- Custom Logos - Subject to one time set up fee.

## Embroidery Personalization - Letter Styles

**Shielding Intl.**

Garamond

**Shielding Intl.**

Block

*Shielding Intl*

Script

**Shielding Intl.**

Micro Scan

**Shielding Intl.**

Impress

*Can't find your letter style? Call for assistance.*

## A) Non-Slip Skirt Option



## B) Extension Panel

