



Armstrong Medical

INDUSTRIES INC.

<http://www.armstrongmedical.com>



Distributed by:

J.D. Honigberg International, Inc.

155 Pfingsten Road, Suite 150

Deerfield, IL 60015 USA.

Medical@jdhintl.com

www.jdhmedical.com





Armstrong Medical
INDUSTRIES INC.

- Started operating in 1957
- Made in Illinois, USA.
- Dedicated to quality
- Present in 65% of Emergency rooms in USA.
- Certified ISO 9001-2000





Armstrong Medical
INDUSTRIES INC.

2 choices of materials: Aluminum versus Steel

- Aluminum does not corrode whereas steel will rust if exposed to humidity and liquids. Without corrosion you can expect a longer life span and personnel is not at risk when in contact with rust.
- Aluminum is much lighter than steel making it easier to move and handle the cart, especially in emergency situations and personnel is less likely to suffer back pain from pushing heavy carts.
- An aluminum cart 76cm high weighs **41.7kg**.
The same size cart in steel weighs **74.8 kg**.





Armstrong Medical
INDUSTRIES INC.

Differences between Standard and Premier carts

- Premier carts are all made of Powder coated Aluminum.
- Most Standard carts are made of Powder coated Steel.
- Premier carts come with a slide thru side shelf.
- Standard carts can receive an optional side shelf on one side.
- Accessories on a Premier cart are affixed to side panels which moves up and down the cart and can be removed or installed in seconds.
- Accessories on Standard carts are affixed with screws directly into the side walls of the cart.



Overview of a Standard Cart



Armstrong Medical
INDUSTRIES INC.

Breakaway Security ABL Series

Soft-grip Handles

Double Side-wall Construction

Ball Bearing Drawer Slides

Stabilizing Frame with Bumper prevents injuries. Casters are moved from under the cart to outside its perimeter. This extra-wide base provides your cart with added stability.

- (1) Swivel
- (2) Swivel w/Brake
- (1) Tracking Caster


Alternative locking systems: Key lock or Push button lock



Overview of a Premier Cart



Armstrong Medical
INDUSTRIES INC.



Stabilizing Frame with Bumper prevents injuries. Casters are moved from under the cart to outside its perimeter. This extra-wide base provides your cart with added stability (*casters and bumpers are moved further away from the base on the drawer side to balance the weight of opened drawers*).

Alternative locking systems: Key lock or Push button lock



Armstrong Medical
INDUSTRIES INC.

Characteristics of Armstrong Carts

- ❖ Built with double wall all around to protect inside drawers from shocks.
- ❖ Wall panels are made from one metal sheet all around instead of being united with soldering which is susceptible to rust.
- ❖ High quality Tente[®] swivel Casters (2 locking & 1 tracking).
- ❖ Casters are positioned beyond the cart's base for better stability.
- ❖ Drawers system on ball bearings for easy and smooth movement.
- ❖ Drawers are interchangeable and can be moved up or down.
- ❖ Bumpers all around to absorb any impact.
- ❖ 5 years warranty.
- ❖ Choice of 18 colors (powder coated).





HUNDREDS OF ACCESSORIES AVAILABLE

- Large choice of accessories to customize your carts.
- Accessories are mounted to the walls of Standard carts or affixed to the accessory bars of aluminum carts and can be changed from one side to the other in a few seconds to accommodate each nurse's preferences.
- Tilt bins and shelves can be mounted above the carts .
- Drawer trays with dividers can be inserted to organize items stored in the carts.





Armstrong Medical
INDUSTRIES INC.

Choice of locking systems

- Plastic seal for emergency carts.
- Key lock for all or individual drawers or push button lock on anesthesia, procedure, isolation carts.





Armstrong Medical
INDUSTRIES INC.

Armstrong carts are present in world renown Hospitals

- ✓ John Hopkins Hospital
- ✓ The Cleveland Clinic
- ✓ The Mayo Clinic
- ✓ Duke University Medical Center
- ✓ New York Presbyterian University Hospital
- ✓ UCLA
- ✓ Stanford Hospital Clinics





Armstrong Medical
INDUSTRIES INC.

- Armstrong Medical carts are used by nurses in Hospitals in over 50 countries on all continents.
- Your Hospital staff also deserves to work with Armstrong Medical carts to help provide great healthcare services.

