



650 Ford Street
Colorado Springs, CO 80915
Toll: (800)433-4064 * Ph: (719)637-0300
Fax: (719)597-8273 * www.harloff.com

Stainless Steel Eight Drawer Cast Cart with Top Compartment, Deluxe Package #6025-TC

Specifications

Construction:

- Stainless steel cabinet and drawers
- Painted aluminum top compartment

Features:

- Top compartment for 3M® and BSN® splinting supplies; includes plexiglass front
- Stainless steel folding shelf
- Stainless steel side storage cabinet
- Scissors cable (scissors not included)
- Stainless steel bucket holder with vinyl mat
- Swivel tourniquet / pump pole
- 24" stainless steel ruler mounted on top of cart
- Aluminum towel bar
- Full wrap-around vinyl bumper on separate steel frame
- 5" sealed ball-bearing casters, all full swivel, one directional, one with brake
- Ball-bearing, full extension drawer slides
- Certified "Latex Free"
- Five year warranty



Casting supplies not included

Assembled Cart Dimensions

Cart only: 43" H x 44" W x 21" D

Cart w/ top compartment (TC): 66" H x 44" W x 21" D

Cart weight w/ TC: 239 lbs

NOTE: This product will ship in two packages

Approx shipping weight of cart: 250 lbs (113.7 kg)

Approx shipping weight of top compartment: 27 lb (12.7 kg)

Drawer configuration:

- 1 x 3" with black vinyl mat (7.6 cm)
- 5 x 3" with 4 adjustable dividers per drawer (15.2 cm)
- 2 x 6" with 4 adjustable dividers per drawer (22.9 cm)

Internal Drawer Dimensions:

- 23.75" W x 17" D (60.3 cm x 43.2 cm)



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Painted Steel Eight Drawer Cast Cart with Top Compartment, Deluxe Package #6035-TC

Specifications

Construction:

- Painted steel cabinet and aluminum drawers
- Painted aluminum top compartment

Features:

- Top compartment for 3M® and BSN® splinting supplies; includes plexi-glass front
- Stainless steel folding side shelf
- Stainless steel side storage cabinet
- Scissors cable (scissors not included)
- Stainless steel bucket holder with vinyl mat
- Swivel tourniquet / pump pole
- 24" stainless steel ruler mounted on top of cart
- Aluminum towel bar
- Internal drawer dividers
- Full wrap-around vinyl bumper on separate stabilizing frame
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Ball-bearing, full extension drawer slides
- Cart is delivered fully assembled
- Latex free
- Twelve year warranty
- Designed in conjunction with orthopedic technicians to maximize functionality

Other Cast Cart Options:

- All cast carts are available with three different lock types: breakaway, key-locking, and mechanical push-button combination locks. Call Harloff for details.

Internal Drawer Dimensions:

- 23.75" W x 17" D (60.3 cm x 43.2 cm)

Assembled Cart Dimensions

- Cart only: 43" H x 44" W x 21" D
- Cart w/ top compartment (TC): 66" H x 44" W x 21" D
- Cart weight: 222 lbs
- Cart weight w/ TC: 239 lbs
- Approximate shipping weight of cart: 250 lbs (113.7 kg)
- Approximate shipping weight of top compartment: 27 lb (12.7 kg)

** NOTE: This product will ship in two packages



Drawer Configuration:

- (1) - 3" with black vinyl mat (7.6 cm)
- (5) - 3" with 4 adjustable dividers per drawer (15.2 cm)
- (2) - 6" with 4 adjustable dividers per drawer (22.9 cm)

Sixteen available color choices:

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black



650 Ford Street
Colorado Springs, CO 80915
Toll: (800)433-4064 * Ph: (719)637-0300
Fax: (719)597-8273 * www.harloff.com

Stainless Steel Eight Drawer Cast Cart, Deluxe Package

#6025

Specifications

Construction:

- Stainless steel cabinet and drawers

Features:

- Stainless steel folding side shelf
- Stainless steel side storage cabinet
- Scissors cable (scissors not included)
- Stainless steel bucket holder with vinyl mat
- Swivel tourniquet / pump pole
- 24" stainless steel ruler mounted on top of cart
- Aluminum towel bar
- Internal drawer dividers
- Full wrap-around bumper on stabilizing frame
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Ball-bearing, full extension drawer slides
- Cart is delivered fully assembled
- Latex Free
- Twelve year warranty
- Designed in conjunction with orthopedic technicians

Other Cast Cart Options:

- All cast carts are available with three different lock types:
 - breakaway
 - key-locking
 - mechanical push-button combination locks

Internal Drawer Dimensions:

- 20" W x 17" D (50.8 cm W x 43.2 cm D)

Assembled Cart Dimensions

- Cart only: 43" H x 44" W x 21" D (109 cm x 111.8 cm x 53.3 cm)
- Cart w/ tourniquet / pump pole: 76" H x 44" W x 21" D (193 cm x 111.8 cm x 53.3 cm)
- Cart weight: 209 lbs (95 kg)
- Approximate shipping weight: 250 lbs (113.7 kg)

Drawer Configuration:

- (1) - 3" with black vinyl mat (7.6 cm)
- (5) - 3" with 4 adjustable dividers per drawer (15.2 cm)
- (2) - 6" with 4 adjustable dividers per drawer (22.9 cm)

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Stainless Steel, Eight Drawer Cast Cart, Standard Package #6020

Specifications

Construction:

- Stainless steel cabinet and drawers

Features:

- Stainless steel bucket holder with vinyl mat
- 24" stainless steel ruler mounted on top of cart
- Aluminum towel bar
- Internal drawer dividers
- Full wrap-around bumper on stabilizing frame
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Ball-bearing, full extension drawer slides
- Cart is delivered fully assembled
- Latex free
- Twelve year warranty
- Designed in conjunction with orthopedic technicians

Other Cast Cart Options:

- All cast carts are available with three different lock types:
 - breakaway
 - key-locking
 - mechanical push-button combination locks

Drawer Configuration:

- (1) - 3" with black vinyl mat (7.6 cm)
- (5) - 3" with 4 adjustable dividers per drawer (15.2 cm)
- (2) - 6" with 4 adjustable dividers per drawer (22.9 cm)

Internal Drawer Dimensions:

- 20" W x 17" D (50.8 cm W x 43.2 cm D)

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Assembled Cart Dimensions

- 43" H x 44" W x 21" D (109 cm x 111.8 cm x 53.3 cm)
- Cart weight: 182 lbs (82.8 kg)
- Approximate shipping weight: 222 lbs (101 kg)





650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Painted Steel Eight Drawer Cast Cart, Deluxe Package #6035



Specifications

Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Stainless steel folding side shelf
- Stainless steel side storage cabinet
- Scissors cable (scissors not included)
- Stainless steel bucket holder with vinyl mat
- Swivel tourniquet / pump pole
- 24" stainless steel ruler mounted on top of cart
- Aluminum towel bar
- Internal drawer dividers
- Full wrap-around bumper on separate stabilizing frame
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Ball-bearing, full extension drawer slides
- Cart is delivered fully assembled
- Latex free
- Twelve year warranty
- Designed in conjunction with orthopedic technicians to maximize functionality

Drawer Configuration:

- (1) - 3" with black vinyl mat (7.6 cm)
- (5) - 3" with 4 adjustable dividers per drawer (15.2 cm)
- (2) - 6" with 4 adjustable dividers per drawer (22.9 cm)

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Other Cast Cart Options:

- All cast carts are available with three different lock types: breakaway, key-locking, and mechanical push-button combination locks. Call Harloff for details.

Internal Drawer Dimensions:

- 23.75" W x 17" D (60.3 cm x 43.2 cm)

Assembled Cart Dimensions

- Cart only: 43" H x 44" W x 21" D (109.2cm x 111.8 cm x 53.3 cm)
- Cart w/ tourniquet / pump pole: 76" H x 44" W x 21" D (193 cm x 111.8 cm x 53.3 cm)
- Cart weight: 182 lbs
- Approximate shipping weight of cart: 227 lbs (102.2 cm)



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Painted Steel Eight Drawer Cast Cart, Standard Package #6030



Specifications

Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Stainless steel bucket holder with vinyl mat
- 24" stainless steel ruler mounted on top of cart
- Aluminum towel bar
- Internal drawer dividers
- Full wrap-around bumper on separate stabilizing frame
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Ball-bearing, full extension drawer slides
- Cart is delivered fully assembled
- Latex free
- Twelve year warranty
- Designed in conjunction with orthopedic technicians to maximize functionality

Drawer Configuration:

- (1) - 3" with black vinyl mat (7.6 cm)
- (5) - 3" with 4 adjustable dividers per drawer (15.2 cm)
- (2) - 6" with 4 adjustable dividers per drawer (22.9 cm)

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Other Cast Cart Options:

- All cast carts are available with three different lock types: breakaway, key-locking, and mechanical push-button combination locks. Call Harloff for details.

Internal Drawer Dimensions:

- 23.75" W x 17" D (60.3 cm x 43.2 cm)

Assembled Cart Dimensions

- Cart: 43" H x 44" W x 21" D (109.2cm x 111.8 cm x 53.3 cm)
- Cart weight: 152 lbs (69.1 kg)
- Approximate shipping weight: 197 lbs (89.5 kg)