

Two-Tier Table Drapes (Protected by U.S. Patent No. 6,019,102)

- Unique single piece disposable drapes:
 - Standard:** Two layer system including absorbent layer - 6 mil thick.
 - Heavy Duty:** Three layer system including extra absorbent center layer - 28 mil thick.
- One piece drape easily covers both upper and main table surfaces.
- The clear plastic backdrop protects the sterile field while allowing light through to the main table surface.



| TABLE & DRAPE CROSS REFERENCE CHART | | | | Guaranteed True-Fit Drapes for OR Specific's Big Case Back Tables | | | | | |
|-------------------------------------|---|--|----------|--|------------|----------|----------|----------|----------|
| Part No. | Description | Dimensions | Standard | | Heavy-Duty | | | Per Case | |
| | | | 418-S | 419-S | 420-S | 418-HD-S | 419-HD-S | | 420-HD-S |
| 427 | 4' Back Table with removable 21½" D second tier shelf | 48" L x 30" D x 33" H | ● | | | ● | | | 10 |
| 427-ADJ | 4' Back Table with removable 21½" D second tier shelf | 48" L x 30" D x 29" H (height adjustment from 29" to 38") | ● | | | ● | | | 10 |
| 428 | 5' Back Table with removable 21½" D second tier shelf | 60" L x 30" D x 33" H | | ● | | | ● | | 10 |
| 428-ADJ | 5' Back Table with removable 21½" D second tier shelf | 60" L x 30" D x 29" H (height adjustment from 29" to 38") | | ● | | | ● | | 10 |
| 429 | 6' Back Table with removable 21½" D second tier shelf | 72" L x 30" D x 33" H | | | ● | | | ● | 10 |
| 429-ADJ | 6' Back Table with removable 21½" D second tier shelf | 72" L x 30" D x 29" H (height adjustment from 29" to 38") | | | ● | | | ● | 10 |

Reinforced Single Level Table Drapes

Utilize these single level drapes when the upper tier of the Big Case Back Table is in the storage position and the table intended as a single surface work area.



| SINGLE LEVEL TABLE DRAPE | | | |
|--------------------------|-------------------------|------------------|----------|
| For use with Table | Single Level Heavy-duty | Drape Dimensions | Per Case |
| 428 | 315-S | 80" x 90" | 10 |
| 429 | 320-S | 80" x 110" | 10 |



Draping Video DVD Demo

Request a copy of the features, benefits and draping instructions for the Big Case Back Table. CD-ROM or DVD. Also available as a download at www.orspecific.com.



Efficiency Report

Request a copy of the time saving benefits for the Big Case Back Table from your local rep. Also available as a PDF download from our website at www.orspecific.com.



30 Day Evaluation

To learn more about our table and drape purchasing program, please contact OR Specific's customer service department at 800.937.7949.

BIG CASE back tables

Specifically designed for a variety cases



Maximize Space

Upper shelf creates additional space without compromising visibility



Visibility

Easily arrange & organize instrument trays without the need for stacking



Time

Improve setup efficiency & workflow during large surgical cases



Drape

Create a sterile field with our patented one-piece drape.

Standard or heavy-duty
Protected by U.S. Patent No. 6,019,102



standard BIG CASE Back Tables

The OR Specific standard Big Case Back Table is ideal in supporting any surgical procedure that is implant and/or instrument intensive, especially when space is limited.

STANDARD TABLE SPECIFICATIONS

The standard Big Case Back Table has a fixed table height of 33" (floor to main surface). When upper tier is in the storage position, the table can be used as a single tier back table.

- **Upper Shelf Set Positions**
12" above main surface
15" above main surface
- **Weight Capacity**
150 lbs. upper shelf
200 lbs. main surface
200 lbs. storage shelf
— weight uniformly distributed
- **18 gauge stainless steel**
- **Four 3½" casters**
- **Two front locking casters**



427 Big Case 4' Standard

Table Dimensions: 48"L x 30"D x 33"H

Shelf Dimensions: 48"L x 21½"D

428 Big Case 5' Standard

Table Dimensions: 60"L x 30"D x 33"H

Shelf Dimensions: 60"L x 21½"D

429 Big Case 6' Standard

Table Dimensions: 72"L x 30"D x 33"H

Shelf Dimensions: 72"L x 21½"D

BIG CASE back tables

427-ADJ Big Case 4' Adjustable

Table Dimensions: 48"L x 30"D x 29"H

Adjustable: 29" LOW to 38" HIGH

Shelf Dimensions: 48"L x 21½"D

428-ADJ Big Case 5' Adjustable

Table Dimensions: 60"L x 30"D x 29"H

Adjustable: 29" LOW to 38" HIGH

Shelf Dimensions: 60"L x 21½"D

429-ADJ Big Case 6' Adjustable

Table Dimensions: 72"L x 30"D x 29"H

Adjustable: 29" LOW to 38" HIGH

Shelf Dimensions: 72"L x 21½"D

adjustable BIG CASE Back Tables

OR Specific's adjustable back tables accommodate variable staff heights ranging from 4'11" to 6'6". The main surface of the table may be adjusted intraoperatively from a minimum height of 29" to a maximum height of 38".

ADJUSTABLE TABLE SPECIFICATIONS

Adjustable Big Case Back Table models may be raised or lowered intraoperatively through the use of a hand crank.

- **Adjustable Main Surface**
29" minimum height
38" maximum height
- **Upper Shelf Set Positions**
12" above main surface
15" above main surface
- **Weight Capacity**
150 lbs. upper shelf
200 lbs. main surface
200 lbs. storage shelf
— weight uniformly distributed
- **18 gauge stainless steel**
- **Four 3½" casters**
- **Two front locking casters**



Table shown with upper tier in storage position



Two locking casters on front legs of table hold table stationary.




Secure the upper tier with T-bolts once raised and positioned.




Post hole for storing upper tier when not in use.



Raise or lower the table by using the hand crank at back of table.



Safety locking pin used for securing upper tier shelf.



Set at 15" above main surface

Set at 12" above main surface

38" max. height

29" min. height




Table has a 9" vertical travel to reach a maximum height of 38".

standard
BIG CASE BACK TABLES

The OR Specific standard Big Case Back Table is ideal in supporting any surgical procedure that is implant and/or instrument intensive, especially when space is limited.

427 Big Case 4' Standard

Table Dimensions: 48" L x 30" D x 33" H

Shelf Dimensions: 48" L x 21½" D



I love this table, the Big Case Back Table. I always ask for it in my room. It makes it easier, more organized - don't need 3 or 4 tables, only this one.
- Staff Surgical Technologist

428 Big Case 5' Standard

Table Dimensions: 60" L x 30" D x 33" H


Shelf Dimensions: 60" L x 21½" D



429 Big Case 6' Standard

Table Dimensions: 72" L x 30" D x 33" H

Shelf Dimensions: 72" L x 21½" D




adjustable
BIG CASE BACK TABLES

OR Specific's adjustable back tables accommodate variable staff heights ranging from 4' 11" to 6' 6". The main surface of the table may be adjusted intraoperatively to a maximum height of 38" or a minimum of 29".

427-ADJ Big Case 4' Adjustable

Table Dimensions: 48" L x 30" D x 33" H

Shelf Dimensions: 48" L x 21½" D



428-ADJ Big Case 5' Adjustable

Table Dimensions: 60" L x 30" D x 33" H

Shelf Dimensions: 60" L x 21½" D



429-ADJ Big Case 6' Adjustable

Table Dimensions: 72" L x 30" D x 33" H

Shelf Dimensions: 72" L x 21½" D



OR SPECIFIC

STANDARD TABLE SPECIFICATIONS

The standard Big Case Back Table has a fixed table height of 33" (floor to main surface). When upper tier is in the storage position, the table can be used as a single tier back table.

- Upper Shelf Set Positions
12" above main surface
15" above main surface
- Weight Capacity
150 lbs. upper shelf
200 lbs. main surface
200 lbs. storage shelf
— weight uniformly distributed
- 18 gauge stainless steel
- Four 3½" casters
- Two front locking casters

Secure the upper tier with T-bolts once raised and positioned.

Post hole for storing upper tier when not in use.

Table shown with upper tier in storage position

Two locking casters on front legs of table hold table stationary.

This is great. Perfect for A/P spine cases, revisions - love the table. Minimizes the space you need. Can take the top tier off and use as a single tier table.
- Staff Surgical Technologist

OR SPECIFIC

ADJUSTABLE TABLE SPECIFICATIONS

Adjustable Big Case Back Table models may be raised or lowered intraoperatively through the use of a hand crank.

- Adjustable Main Surface
29" minimum height
38" maximum height
- Upper Shelf Set Positions
12" above main surface
15" above main surface
- Weight Capacity
150 lbs. upper shelf
200 lbs. main surface
200 lbs. storage shelf
— weight uniformly distributed
- 18 gauge stainless steel
- Four 3½" casters
- Two front locking casters

Set at 15" above main surface

Set at 12" above main surface

Raise or lower the table by using the hand crank at back of table.

Safety locking pin used for securing upper tier shelf.

Table has a 9" vertical travel to reach a maximum height of 38".

Two-Tier Table Drapes (Protected by U.S. Patent No. 6,019,102)

- Unique single piece disposable drapes:
 - Standard: 0.006 thick - two layer system including absorbent layer.
 - Heavy Duty: 0.028 thick - three layer system including extra absorbent center layer.
- One piece drape easily covers both upper and main table surfaces.
- The clear plastic backdrop protects the sterile field while allowing light through to the main table surface.



TABLE & DRAPE CROSS REFERENCE CHART

SUGGESTED DRAPE

| Part No. | Description | Dimensions | 418-S | 419-S | 420-S | 418-HD-S | 419-HD-S | 420-HD-S | Per Case |
|----------|--|---|-------|-------|-------|----------|----------|----------|----------|
| 427 | 4' Back Table with removable 21½"D second tier shelf | 48"L x 30"D x 33"H | • | | | • | | | 10 |
| 427-ADJ | 4' Back Table with removable 21½"D second tier shelf | 48"L x 30"D x 33"H (height adjustment from 29" to 38") | • | | | • | | | 10 |
| 428 | 5' Back Table with removable 21½"D second tier shelf | 60"L x 30"D x 33"H | | • | | | • | | 10 |
| 428-ADJ | 5' Back Table with removable 21½"D second tier shelf | 60"L x 30"D x 33"H (height adjustment from 29" to 38") | | • | | | • | | 10 |
| 429 | 6' Back Table with removable 21½"D second tier shelf | 72"L x 30"D x 33"H | | | • | | | • | 10 |
| 429-ADJ | 6' Back Table with removable 21½"D second tier shelf | 72"L x 30"D x 33"H (height adjustment from 29" to 38") | | | • | | | • | 10 |

Reinforced Single Level Table Drapes

This very heavy flat drape is intended for use with our Specialty Back Tables as a regular single level drape.



SINGLE LEVEL TABLE DRAPE

| For use with Table | Single Level Heavy-duty | Drape Dimensions | Per Case |
|--------------------|-------------------------|------------------|----------|
| 428 | 315-S | 80" x 90" | 10 |
| 429 | 320-S | 80" x 110" | 10 |

Draping Video DVD Demo



Request a copy of the features, benefits and draping instructions for the Big Case Back Table. CD-ROM or DVD. Also available as a download at www.orspecific.com.

Time Commission Report




Request a copy of the time saving benefits for the Big Case Back Table from your local rep. Also available as a PDF download from our website at www.orspecific.com.

Special Purchasing Program



To learn more about our table and drape purchasing program, please contact OR Specific's customer service department at 800.937.7949.



BIG CASE BACK TABLES

Specifically designed for large cases



for better utilization of...

Space
Ideal for Smaller ORs

The upper shelf creates additional space without compromising visibility.

Surface Area
No Tray Stacking

Large surface area to organize and arrange instrument trays.

Time
Streamline Workflow

Improve setup efficiency and workflow during large surgical cases.