

CHAMPION™

Quality Healthcare Seating Products

Accessories +

PAPER ROLL DISPENSER

Part# PAPER-STD (STANDARD*)
Part# PAPER-XL (BARIATRIC RECLINERS)
Part# PAPER-CND (CONCORD RECLINER)



Shown on the Ascent recliner.

Minimize the risk of cross contamination quickly and conveniently with Champion's Paper Roll Dispenser. Our dispenser attaches securely to the rear push handle of your chair and allows caregivers to quickly dispense antimicrobial paper. A vinyl strap, that is color-coordinated with your chair, secures the paper to your chair and positions it for maximum patient comfort.

Add an additional layer of protection and increased patient confidence with the Paper Roll Dispenser from Champion.



Functions in both upright and recline positions.

Protects against mold, mildew, fungi and odor-causing bacteria.



SIMPLE AND QUICK TO USE



Acts as a shield against the spread of unwanted bacteria. Additional protection to patient as well as the recliner.

PRODUCT INFORMATION

- Sturdy, steel powder-coated holder.
- Vinyl band with velcro closure in back, secures paper in position (specify color).
- First roll of InhibX™ paper included with initial purchase of paper dispenser (while supplies last).
- Dispenser can be factory-installed or retrofit.

*Standard size recliners - Classic, Ascent, Passage

CHAMPION™

Champion Manufacturing, Inc.
2601 Industrial Parkway, Elkhart, IN 46516

1.800.998.5018
www.championchair.com

www.championchair.com