

Model DS4050

# Digital Baby Scale



The Model DS4050 digital baby scale is exactly what you've been looking for in a professional infant scale. The DS4050 combines accuracy, high quality construction, easy operation - at a price that's right for even the tightest budgets.

Providing precise measurements of babies up to 44 pounds, the DS4050 displays the patient's weight in 0.5 ounce increments. The built-in measuring tape can be used to measure the length of infants up to 22.5 inches. The rugged ABS weighing surface and housing are precision molded with smooth edges and a tight fit throughout.

The weight is displayed on an easy-to-read LCD display with large 1" tall characters. With just three buttons, this scale is easy to use.

The Doran DS4050 Baby Scale...amazing quality, amazingly affordable.

## MAXIMUM CAPACITY

44 lb x 0.5 oz (20 kg x 10 g)

## ACCURACY

+/- 0.5 oz (10 g)

## DISPLAYED UNITS

lb:oz or kg

## POWER

AC Adaptor (Included)  
Four AA Batteries (Not Included)

## CONSTRUCTION

Rugged ABS Weighing Platform and Base

## DIMENSIONS

Weighing Surface 11" W x 23.5" L  
Scale 15" W x 23.6" L x 3.9" H

## OPTIONS

Carrying Case

## WARRANTY

 Two Years



**Doran Scales, Inc.** 1-800-264-4106  
1315 Paramount Parkway Fax: 630-879-0073  
Batavia, IL 60510 www.doranmedical.com