

Models DS4100 & DS4200

# Infant/Pediatric Scales



**DS4100 Infant Scale**



**DS4200 Infant Scale with Seat**

Designed for portability and accuracy, the Model DS4100 is the perfect choice for a traveling nurse or a pediatric office. Light enough for visiting nurses, yet stable and reliable enough for busy pediatric or family practices.

The Model DS4200 offers all the features of the DS4100 with the additional feature of a molded plastic infant seat. For maximum infant safety, the seat is permanently mounted to the scale. The seat can be adjusted to face either right or left so that it can allow easy access to the infant chair in any office configuration.

Although very affordable, the Doran Infant and Pediatric Scale contains features you only see in far more expensive scales. With 0.2 oz (5 g) accuracy, your medical staff will now be able to accurately track even the smallest weight changes. To operate, simply turn on the scale and place the patient onto the platform.

## MAXIMUM CAPACITY

45 lb x 0.2 oz (20 kg x 5 g)

## ACCURACY

+/- 0.2 oz (5 g)

## DISPLAYED UNITS

lb:oz, kg

## DISPLAY

0.66" High LCD

## POWER

Four AA Batteries (Included)

## PLATFORM DIMENSIONS

DS4100 - 25" L x 11" W

DS4200 - 25" L x 11" W

## OVERALL DIMENSIONS

DS4100 - 25" L x 13.75" W x 3.75" H

DS4200 - 25" L x 13.75" W x 18" H

## CONTROLS

On/Off, Reweigh, LB, KG, Zero

## BATTERY LIFE

Two year battery life

## CONSTRUCTION

ABS Plastic Housing

## OPTIONS

Carrying Case (P/N DS4100-C)

## WARRANTY



Three Years



**Doran Scales, Inc.** 1-800-264-4106

1315 Paramount Parkway Fax: 630-879-0073

Batavia, IL 60510 [www.doranmedical.com](http://www.doranmedical.com)