

DORAN SCALES

DS500 Flat Digital Medical Scale



Features

The “Auto-On” feature of the DS500 switches on the scale when the patient taps on the heavy duty glass platform. The weight is shown in a 2” high LCD display. The tempered glass weighing surface is 12” x 12” and the height of the scale is just 3/4” high. The DS500 can display weight in lbs or kg. The accurate results of the DS500 are easy to read, thanks to the easy to read display and 2” tall numbers.

Specifications

Maximum Capacities:	181 Kgs
Display:	5 cms High LCD
Units:	Lb / Kg
Resolution Range:	0.09 Kg
Construction:	Heavy Duty Tempered Glass and Steel
Power:	3V Lithium Battery
Dimensions:	29.8cms L x 31.1cms W x 1.9cm H
Warranty:	2 Years

JD Honigberg International, Inc

155 Pfingsten Road, Ste 150, Deerfield, IL 60015 USA

Tel: (1) 847-412-0200 / Fax: (1) 847-412-0034 / E-Mail: medical@jdhintl.com