

Model DS8150

Digital Chair Scale

Doran's Model 8150 Chair Scale has been carefully designed to provide a digital chair scale that is easy to move and operate. The strong steel frame provides a very durable chair scale that can weigh patients up to 550 pounds and display the weight on an easy to read LCD display. The highly maneuverable chair scale can be rolled to the patient, and easily and safely locked in place, to weigh the patient.

The four full swivel wheels feature a convenient locking mechanism. To secure the wheels, simply press down on the locking arms. To release the wheel, press down on the release arm. Most chair scales on the market today have only two locking swivel casters, we've put them on all four corners to ensure the safety of the patient and ease of movement in tight spaces.

The anti-tip bar at the bottom of the foot rest arm, prevents the patient from accidentally tipping the chair scale forward while seated and leaning too far forward or when they enter or exit the scale. The large foot rests can be folded up and out of the way when the patient enters and exits the scale, and folded down when the scale is used to weigh the patient. The padded arm rests fold up and out of the way to allow patients to slide on to the scale from the side from a bed or wheelchair. The arm rests can then be folded down for the patient's safety.



The DS8150 is EMR Ready and data can be transmitted via either the included USB or RS-232 ports.

MAXIMUM CAPACITY

550 lb x 0.2 lb (250 kg x 0.1 kg)

ACCURACY

+/- 0.2 lb (0.1 kg)

SEAT DIMENSIONS

22" W x 18" H x 15" D

OVERALL SCALE DIMENSIONS

28" W x 44" H x 45" D

DISPLAYED UNITS

lb, kg

DISPLAY

1" High LCD

BMI CALCULATOR

Input data in FT/IN or CM

APPROVALS

CE, UL (AC Adapter)

POWER

AC Adapter (Included)
or Six AA Batteries

BATTERY LIFE

85,000 weighments

CONSTRUCTION

Powder Coated Steel
and Heavy Duty ABS Plastic

WARRANTY

2 Two Years



Doran Scales, Inc. 1-800-264-4106
1315 Paramount Parkway Fax: 630-879-0073
Batavia, IL 60510 www.doranmedical.com