

In-Bed Scale

Complete patient weighing system



Monitoring a patient's weight is critical in a hospital or extended care facility, but moving the patient from the bed to a scale can be dangerous for the patient as well as the staff. Between 700,000 and 1 million patients fall in hospitals each year, according to the Agency for Healthcare Research and Quality. Since 2008, patient falls with injury have been included in the Centers for Medicare & Medicaid Services list of hospital-acquired conditions for which providers will not be reimbursed, so that risk of a fall now becomes a financial risk.

Partner with Doran Scales and offer your customers a bed that can provide a quick, safe, and accurate method to weigh patients.

Put our over 35 years of experience designing and building digital weighing systems to work for you. Our engineering team will design a turn-key system specifically suited for your beds. Stand out from the competition and win sales by providing a in-bed scale option.

MAXIMUM CAPACITY

1000 x 0.2 lb (500 x 0.1 kg)

DISPLAYED UNITS

lb, kg

DIMENSIONS

10" W x 1" D x 2.5" H

DISPLAY

0.66" High LCD

CONTROLS

ZERO, WEIGH, LB KG

CONSTRUCTION

ABS Plastic Housing

POWER

5V DC



Doran Scales, Inc.
Batavia, IL 60510
800-262-6844
doranmedical.com