



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Standard Line 600 Capacity Punch Card Medication Cart with CompX® Electronic Lock, Standard Package #SL600PC-EKC

Specifications

Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Dual column medication cart with 600 punch card capacity, single width narcotics drawer
- CompX® electronic lock on cabinet with BEST® brand key lock override; Illinois lock on narcotics drawer
- Access right may be assigned for up to 250 users per cart
- Drawers automatically re-lock after access (pre-set or programmed time)
- Painted pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Full wrap around vinyl bumper, protecting the cart
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish

Drawer Dimensions:

Left Column Drawers:

- 1 ea – 4" multipurpose storage drawer
- 3 ea - 9.75" punch card drawers

Right Column Drawers:

- 1 ea – 9.75" locking narcotics drawer
- 2 ea – 9.75" multipurpose storage drawer



lockview software available: #41052 Includes
audit software and dongle

Thirteen available color choices:

Sand, Gloss Beige, Burgundy,
Light Blue, Navy, Light Gray,
Mauve, Red, Teal, Yellow,
Hammer Tone Blue, Hammer Tone Gray
and Hammer Tone Green

Assembled Cart Dimensions:

- 45.75" H x 43" W x 23.5" D
(116.2 cm x 109.2 cm x 59.7 cm)
- Writing surface: 44" (111.8" cm)
- Cart weight: 185 lbs (93.2 kg)
- Approx shipping weight: 230 lbs (113.7 kg)



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Standard Line Med-Bin Cart with CompX® Electronic Lock

54 – 3” Bins Standard Package

#SL54BIN3-EKC

Specifications

Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Dual column medication cart with 54 – 3” wide medication bins. Dividers and labels are included, single width narcotics drawer
- CompX® electronic lock on cabinet with key lock override; Illinois lock on narcotics drawer
- Access rights may be assigned for up to 250 users per cart
- Drawers automatically re-lock after access (pre-set or programmed time)
- Painted steel pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Full wrap around vinyl bumper, protecting the cart
- 5” ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish
- Offered in specialty package (#SL54BIN3SP) Included are: #68640-P1 drawer tray with adjustable dividers, #680551 eye drop tray, #AL680435 aluminum waste container and #684810 locking sharps container

Internal Drawer Dimensions:

Left Column Drawers:

- 1 ea – 4” multipurpose storage drawer
- 9 ea – 3.25” med bin drawers with 6 – 3” Bins per drawer

Right Column Drawers:

- 1 ea – 4” Multipurpose Storage Drawer
- 1 ea – 9.75” Locking Narcotics Drawer
- 2 ea – 9.75” Multipurpose Storage Drawer

Assembled Cart Dimensions:

- 45.75” H x 43” W x 23.5” D
(116.2 cm x 109.2 cm x 59.7 cm)
- Writing Surface: 44” (111.8 cm)
- Cart Weight: 220 lbs (100 kg)
- Approx Shipping Weight: 265 lbs (120.5 kg)



lockview software available: #41052 Includes audit software and dongle

Thirteen available color choices:

Sand, Gloss Beige, Burgundy,
Light Blue, Navy, Light Gray,
Mauve, Red, Teal, Yellow,
Hammer Tone Blue, Hammer Tone Gray
and Hammer Tone Green



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

OptimAL Line 64 Bin Capacity Medication Cart with CompX® Electronic Lock #AL64MDBIN-EKC

Specifications

Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Internal narcotics box with Illinois Lock
- CompX® brand electronic keypad lock on cabinet
- Access rights may be assigned for up to 250 users per cart
- Drawers automatically re-lock after access (pre-set or programmed time)
- Painted pull-out side shelf with recessed work area
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with durable work surface
- Ergonomic top rail/push handle
- Easy clean protective pontoon bumpers
- 5" ball-bearing premium rubber casters, all full swivel, 1 brake and 1 directional
- Hard-baked powder coat finish (Customer chooses color for cabinet)

Drawer Dimensions:

- 4 ea – 4" drawers
- 1 ea – 11" multipurpose storage drawer with locking narcotics box

Assembled Cart Dimensions:

- 42" H x 36" W x 24" D (includes top rail)
(106.6 cm x 91.4 cm x 61 cm)
- Work surface: 38" H (96.5 cm)



lockview software available: #41052 Includes
audit software and dongle

Thirteen available color choices:

Sand, Gloss Beige, Burgundy,
Light Blue, Navy, Light Gray,
Mauve, Red, Teal, Yellow,
Hammer Tone Blue, Hammer Tone Gray
and Hammer Tone Green



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Aluminum Universal Line Tall 6 Drawer Procedure/Nurse Supply Cart with Basic Electronic Pushbutton Lock #AL809E6

Harloff's 20 years of innovation in hospital cart design and construction has resulted in the Universal Line. Its modern appearance incorporates industry-leading standard features with flexible configuration options.

Specifications

Construction:

- Aluminum construction

Features:

- Internal Pullout shelf on right side
- Basic Electronic Pushbutton Lock
- Key Lock Override
- Factory Set User and Supervisor Code
- Cabinet may be programmed in the field
- 29.25" of vertical drawer space
- Cabinet design allows for customized drawer configurations
- Painted steel top rail / push handle
- Removable plastic top
- Protective pontoon bumpers
- Extended caster footprint for stability when drawers are extended
- 5" (12.7 cm) ball-bearing premium casters, all full swivel- one tracking, one braking
- Ball-bearing full extension drawer slides
- Cart is delivered wheels-on and mostly assembled
- Pre-drilled threaded holes for standard accessories
- Latex free
- Durable, hard-baked powder coat finish available in 16 colors
- Twelve year warranty on manufactured parts

Drawer Configuration:

- (4) – 3.25" (8.3 cm)
- (1) – 6.5" (16.5 cm)
- (1) – 9.75" (24.8 cm)

Internal Drawer Dimensions:

- 23" W x 17" D (58.4 cm W x 43.2 cm D)

Assembled Cart Dimensions:

- 44.5" H x 31" W x 24.75" D
(113 cm H x 78.7 cm W x 62.9 cm D)
- Cart weight: 110 pounds (50 kg)
- Approximate shipping weight: 140 pounds (64 kg)



**Color Shown is Beige.
Additional Colors Available**