

ENOCHS ENCORE 4250

High-Low Examination Table

PATIENT LOAD: 500 LBS.
evenly distributed on the table's main couch



'EASY TOUCH' FOOT CONTROL adjusts table's height and back with reliable hands-free action. Pedals are covered with non-slip tread and utilize a natural left to right rocking motion of operator's foot.

TABLE FEATURES:

- 4-YEAR Warranty
- Hand Control available as an accessory
- Fully welded All Steel Base Construction -NO PLASTIC
- Full Power adjustment of back
- Easy to clean 100% Seamless Vacuum Formed Top with additional padding on edges for patient comfort
- Pelvic Tilt 7° and 11°
- Electrical Outlets with circuit breaker - NO FUSES
- Instrument Warming Drawer efficiently warms drawer contents to approximate body temperature
- Instrument Warmer ON/OFF Switch with Integrated Circuit Breaker - NO FUSES
- Stainless Steel Irrigation Basin operates independent of Leg Extension
- Easy to adjust Stirrups with length and lateral adjustment
- 1.4 cubic feet of Storage
- Certified by UL to standard 60601-1 CAN/CSA C22-2 No. 601-1
- ENOCHS Encore 4250 replaces the Power 4000



The hospital grade dual electrical outlets are conveniently located at the foot end of the table. The outlets feature a dedicated circuit breaker that will trip if outlets are overloaded. A quick reset of the circuit breaker re-establishes power to the outlets, eliminating costly down time. The Power 4250 also features a SELV (safety extra low voltage) hand control, an accessory that can be ordered at time of purchase or added at a later date.



The footend of the table features a padded pullout leg extension, 4" deep stainless steel irrigation basin and the next generation instrument warming drawer system. The exclusive design of separate operation between the leg extension and basin allows for fast transition from one to the other. The pull-out leg extension enables patients with limited hip rotation to be placed in a GYN position, while still giving the physician unrestricted access to the irrigation basin. The instrument warming drawer is designed to warm speculums and instruments to approximate body temperature.

Specifications for the ENOCHS ENCORE 4250 High-Low Examination Table

Warranty:	4-YEAR
Base Construction:	Fully welded all steel base
Power System:	Electro-Mechanical provides smooth and quiet movement
Power Adjustments:	Height: 24" - 40" Back: horizontal - 80 °
Foot Control:	"Hands free", standard on every table, SELV
Hand Control:	"Easy to read", optional accessory, add at anytime, SELV
Couch Top:	100% Seamless vacuum formed style, all steel frame construction
Perineal Recess:	Moulded, no seams, 10" W x 2" D at center point
Pelvic Tilt:	7°, 11°
Leg Extension:	Extends out 20" from end of seat section
Irrigation Basin:	4" Deep stainless steel basin
Electrical Outlets: with Circuit Breaker	Standard feature
Instrument Warmer: with ON/OFF Switch, Circuit Breaker & Reset	Standard feature
Stirrups:	Multi-positional, lateral and length adjustment
Storage:	Side Drawer: 15.5" W x 15" D x 6" H Footend Drawer: 15.5" W x 15" D x 4" H Total Storage Space: 1.4 cubic feet
Paper Roll Holder:	Accommodates 18" or 21" paper widths
Accessories:	hand control, knee crutches, padded arm board, drawer locks, restraint straps, IV pole, urology basin with screen and drain, headend accessory rails, footend accessory rails, caster base
Mains Power Entry: with Circuit Breaker	IEC connection, internal circuit breaker with reset
Duty Cycle:	Table: 1 minute ON / 9 minutes OFF Receptacles: Continuous operation Instrument Warming Drawer: Continuous operation
Ratings:	Table: 115V 60 Hz 10 amps Receptacles: 115V 60 Hz Total Receptacle Load: 3.9 amps
Table Dimensions:	Width: 27", Height: 24" - 40", Length: 52" (extends to 72")
Table Weight:	416 lbs.
Patient Load:	500 lbs. evenly distributed on main couch
Certifications:	By UL to Standard 60601-1, CAN/CSA C22.3 No. 601-1