

# ENOCHS EXCEL 350

## General Examination Table

**PATIENT LOAD: 600 LBS.**  
evenly distributed on the table's main couch

**SPRING ASSISTED COUNTER BALANCE MECHANISM** for smooth easy adjustment of the back section - NO gas cylinder to replace

**ONE-PIECE UPHOLSTERED TOP** padded with 1.5" high density foam on a steel reinforced frame for maximum strength

**PAPER ROLL HOLDER**

Concealed when not in use, the **MULTI-POSITION STIRRUPS** offer length and lateral adjustment

**PAPER CUTTER**

**FIELD REVERSIBLE SIDE DRAWERS**

**EXTRA WIDE PULL-OUT LEG EXTENSION**

Durable 4" **DEEP MOULDED DRAWERS** with rounded edges and corners for easy cleaning

4" deep **IRRIGATION BASIN** moves independent of leg extension

Extra strength **STEEL REINFORCED SIDE GUIDE MOUNTING SYSTEM** for discipline non-binding movement

**FOOTEND STORAGE DRAWERS**

Easy to clean **SLIP-RESISTANT TREAD**

**PULL-OUT LOCKING STEP** engineered to accommodate the 600 lbs patient load

**FULLY WELDED STEEL BASE** construction – NO PLASTIC PARTS

### TABLE FEATURES:

- 2-YEAR Warranty
- Pull-Out Locking Step engineered to accommodate the 600 lbs patient load
- All Steel Base Construction – NO PLASTIC
- Spring Assisted counter balance allows smooth adjustment of couch back
- Easy to clean One-Piece Upholstered Top designed without cross-seams, decorative stitching or embossing
- Multi-position Stirrups with easy length and lateral adjustment
- 3.32 cubic feet of Storage in the five spacious Drawers
- Certified by UL to Standard 60601-1 CAN/CSA C22.3 No. 601-1
- ENOCHS Excel 350 replaces the Excel 71-01



The extra long pull-out foot step features a steel stabilizing leg designed to accommodate a 600 pound patient load capacity. The step easily extends out from the table and locks into position providing the patient with a stable surface when getting on and off the table. The step's all steel construction offers extraordinary strength and durability. The 17.5" long surface is covered in easy to clean slip-resistant tread.



New field reversible side drawers allows table to have the flexibility to be placed into any room configuration. Side drawers combined with the two footend drawers offer maximum storage capacity totaling 3.32 cubic feet. Nylatron® on the bottom drawer edges protect against wear and enhances the drawer's overall performance. An integrated stops eliminate accidental spillage of drawer contents.



The working end of the table has an exclusive ENOCHS design of a pullout padded leg extension that moves independent of the 4" deep irrigation basin. This separate basin design allows physician to easily transition from leg extension to basin and provides a fast easy clean up by the staff after the procedure. The Excel 350 also features multi-position stirrups that lock into position with the weight of the patient's foot.

# Specifications for the ENOCHS EXCEL 350 General Examination Table

<b>Warranty:</b>	2-YEAR
<b>Base Construction:</b>	Fully welded all steel base
<b>Back Lift Mechanism:</b>	Horizontal to 80° by spring assisted counter balance
<b>Couch Top:</b>	One-piece upholstered top, steel reinforced
<b>Perineal Recess:</b>	Not available on this model
<b>Pelvic Tilt:</b>	Not available on this model
<b>Leg Extension:</b>	Extends out 18" from end of seat section
<b>Irrigation Basin:</b>	4" Deep basin
<b>Electrical Outlets:</b>	Not available on this model
<b>Instrument Warmer: with ON/OFF Switch, Circuit Breaker &amp; Reset</b>	Not available on this model
<b>Stirrups:</b>	Multi-positional, lateral and length adjustment
<b>Storage:</b>	<b>3 Side Drawers:</b> 21.5" W x 15" D x 4" H <b>2 Footend Drawers:</b> 15.5" W x 15" D x 4" H <b>Total Storage Space:</b> 3.32 cubic feet
<b>Foot Step:</b>	Extra long surface engineered for 600 lb. patient load
<b>Paper Roll Holder:</b>	Accommodates 18" or 21" paper widths
<b>Accessories:</b>	knee crutches, drawer locks, stainless steel basin, restraint straps, IV pole, urology basin with drain and screen, footend accessory rails, caster base
<b>Mains Power Entry: with Circuit Breaker</b>	Not available on this model
<b>Duty Cycle:</b>	<b>Receptacles:</b> n/a <b>Instrument Warming Drawer:</b> n/a
<b>Ratings:</b>	<b>Table:</b> n/a <b>Receptacles:</b> n/a <b>Total Receptacle Load:</b> n/a
<b>Table Dimensions:</b>	<b>Width:</b> 27", <b>Height:</b> 34", <b>Length:</b> 55" (extends to 73")
<b>Table Weight:</b>	268 lbs.
<b>Patient Load:</b>	600 lbs. evenly distributed on main couch
<b>Certifications:</b>	By UL to Standard 60601-1, CAN/CSA C22.3 No. 601-1