

JD HONIGBERG INTERNATIONAL

X-RAY FILM PROCESSORS

X-Ray film processors models FI45/75/245 and Ecotank made for 110V and 60Hz usage ship from USA.

X-Ray film processors models FI45/75/245 and Ecotank made for 110V 50Hz and 220V 50/60Hz usage ship from Germany.

X-Ray film processor model L356 made for 110V and 60Hz usage ships from USA.

X-Ray film processor model L356 made for 220V and 60Hz usage ships from Japan.



FI-45

*The **FI-45** is the world's simplest automatic processor: simple to install, operate & maintain.*

*Plug & process is all you have to do! The revolutionary design of the **FI-45** reduces energy, water & chemical use, saving its owner money while delivering outstanding quality images.*

Features & Benefits

- Lightweight, tabletop processor
- 100 second processing speed, leading edge to leading edge
- Energy and water saving design limits run time to when film is actually being processed
- Anti-oxidation and anti-crystallization programs result in chemical savings and outstanding image quality
- All components are easily accessible for service
- Simple microprocessor control about the size of a credit card
- Ingenious shallow tank design reduces chemical oxidation resulting in cost & environmental savings
- Top film feed and top return for convenience
- New film transport rack design minimizes cleaning through reduced part count & ensures smooth film path
- Quality polyurethane & rubber rollers for outstanding image quality and no artifacts
- 2 year warranty
- Meets ISO, IEC, CE, TUV, CSA & UL quality standards
- Made in Germany

J D HONIGBERG INTERNATIONAL, INC

155 Pfingsten Road, Ste 150; Deerfield, IL 60015 USA

Tel: 847- 412 0200 / Fax: 847 – 412 0034 / E-Mail: medical@jdhintl.com / www.jdhmedical.com