

# IDI PRODUCT CATALOG

## IMAGING TABLES, ACCESSORIES & VIDEO INTEGRATION SYSTEMS

Quality products for surgical, interventional and diagnostic imaging

230V



# Table of Contents

## Imaging Tables and options

Aspect ISR G2 Table Accessories	3-6
Aspect 100-4T G2 Table Accessories	7-10
Aspect 100-4 Table Accessories	11-13
Aspect 100RTL Table Accessories	14-17
Aspect 100RT Table Accessories	18-21
Aspect 100T Table Accessories	22-25
Aspect 100E Table Accessories	26-29
Aspect 100UC <i>Plus</i> Urology Table	30-33
Aspect 100US Gyn Table	34-35

## Video Integration Systems and options

MDS	36-38
MDS D	39-42
ilex 32	43

## Comparison Matrices

Aspect Four-Way Float Top Tables Specifications & Accessory Matrix	44
Aspect Pain Management Table Specifications & Accessory Matrix	45
Aspect Urology & GYN Table Specifications & Accessory Matrix	46

## Customer Service

IDI Customer Service Contacts	Back page
-------------------------------	-----------

The part numbers listed in this catalog represent 230VAC tables with metric sized accessory rails. 115V systems are available. Part numbers for these can be provided upon request. The part numbers listed are the 230V systems.

**Aspect ISR G2 Table**  
**Catalog # A100-2200 or A100-2641**

**\$49,900.00**



ISR G2 A100-2200



ISR G2 A100-2641

**Aspect ISR G2 Table Description:**

**Mobile Imaging Table with Motorized Elevation, Trendelenburg Tilt, Iso Roll and Manual 4-way Tabletop Float. For Advanced Vascular Surgery and Interventional Procedures.**

- 230V is standard – includes metric sized accessory rails
- 115V is optional and must be specified at time of order.

**Aspect ISR G2 #A100-2200**

For Advanced Vascular Procedures; Includes short, recessed side rails 5" (127mm) in length at imaging end of table that facilitate use of removable headrest extension (option) and also facilitates attachment of catheter tray extension at narrow end of table. Two separate catheter tray extensions are available: One for the pedestal end of table and one for the imaging end of table.

Includes the following:

- X100-2492: Tabletop, 19.5/22" (495/559mm)/24" (610mm) X 84" (2134mm)
- A100-2170 Pendant hand control
- X100-2310 Tabletop Pad, 2" (51mm) thick.
- A270-0178: Panning hand control.
- C000-0328: Patient restraint straps.



**Aspect ISR G2 #A100-2641**

Contoured shape, 22" (559mm) wide at imaging end, 24" (610mm) wide at pedestal end, 24" (610mm) long. Does not include 5" (127mm) long rails at imaging end of top. Not compatible with optional headrest. Optional catheter tray attachment can be attached to pedestal end of table only.

Includes the following:







- X100-2426: Tabletop, 22" (559mm)/24" (610mm) x 84" (2134mm).
- A100-2170 Pendant hand control
- X100-2489 Tabletop Pad, 2" (51mm) thick.
- A270-0178: Panning hand control.
- C000-0328: Patient restraint straps.








## Aspect ISR G2 Accessories Included:

- Tabletop pad, 2" thick  
# X100-2310 for the X100-2492 table top or # X100-2489 for the X100-2426 table top.
- Patient Restraint Straps. # C000-0328
- Panning Hand Control, rail mounted. #A270-0173

## Aspect ISR G2 Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<b>Clamp-on Accessory Rail</b> 10mm x 25mm x 203mm clamp-on accessory rail.	A100-2379	\$495.00
	<b>Armboard</b> Quick-release pivoting arm board designed to mount to standard 10mm x 25mm EU accessory rails. It consists of a rail mounting armboard platform, pad with loop/hook fastener strip and two patient arm restraint straps. If the A310-062 armboard is going to be used on a table top edge without accessory rails, then a clamp-on accessory rail A100-2379 is required. L = 676 mm (26.63in.) W = 153 mm (6.00in.) Gauge = 44 mm (1.75in.)	A310-062	\$450.00
	<b>VAB Armboard</b> Carbon fiber, one-piece, extra wide armboard for vascular access procedures (body section slides under patient/table pad). Includes strap and pad.	A100-2244	\$1489.00
	<b>Arm Board</b> Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059	\$975.00
	<b>Headrest Extension - EU Rail</b> Articulating carbon fiber headrest, designed to mount to standard 10mm x 25mm EU accessory rails used on the imaging end of IDI mobile imaging table models ISR G2 (A100-2200) and 100-4T G2 (A100-2533). Consists of a carbon fiber platform and a pad with loop/hook fastener strips. Requires the table to have an X100-2492 style tabletop. It can also accommodate a flexible anesthesia "hoop" style screen holder, IDI part number C000-0965. L = 462 mm (18.19"), B= 588 mm (23.13"), Pad gauge=51 mm (2"), Adjustment angle= ±12°	A100-2505	\$4635.00
Catheter tray extensions are optional and sold separately; see part numbers A100-1769 (pedestal end catheter tray extension), , A100-2322 (imaging end catheter tray extension for use with table A100-1416) or A100-2387 (imaging end catheter tray extension for use with table A100-2598).			
	<b>Tabletop Catheter Tray Extension</b> 24.25"/616mm wide, 30"/762mm long. Mounts to pedestal end of table only. The maximum weight this extension can support is 20 LBS or 9 Kgs.	A100-1769	\$995.00
IDI Export Product Catalog	Images Diagnostics, Inc. D000-2032E		

## Aspect ISR G2 Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<p><b>Tabletop Catheter Tray Extension</b> 20 ¾"/527mm x 27"/686mm long. Mounts to imaging end of table only. Requires tabletop X100-2492. The maximum weight this extension can support is 20 LBS or 9 Kgs.</p>	A100-2322	\$995.00
	<p><b>Radiation Shield</b> Articulating, rail-mounted. Attaches to table side-rail at base end of table, either side, adjustable. Articulating joint allows shield to conform to table contour. Snap-on upper section above level of tabletop is removable. Designed to mount to standard 10mm x 25mm EU accessory rails used on the IDI mobile imaging table models ISR G2 (A100-2200), 100-4T G2 (A100-2533), 100-4 with 45.72cm/76.20cm (18in./30in.) wide tabletop (A100-2632) and 100-4 with 59.69cm/76.20cm (23.5in./30in.) wide tabletop (A100-2653) .</p>	A610-106	\$1995.00
	<p><b>Anesthesia Screen Holder</b> Flexible "hoop" style. Mounts directly to A100-2505 Headrest Extension. For use on accessory rails - a pair of rotating side rail clamps, part number C000-1065, are required. Rotating side rail clamps and clamp-on accessory rails are not required if table is ordered with the A100-2505 Head- rest Extension.</p>	C000-0965	\$389.00
	<p><b>Peripheral Extension</b> Carbon fiber extension, mounts to imaging end of table only. 22"/559mm wide,, 12"/305mm long. Maximum weight for this extension is 25 lbs or 11.4 Kgs. Requires the table to have an X100-2492 style tabletop..</p>	A100-2486	\$2795.00
	<p><b>Multi-position Arm Support</b> Attaches to IDI tables with accessory rails - or can be used with a clamp on accessory rail on tables or areas within a table that do not have accessory rails. Length: 23.5 inches (597mm).</p>	C000-1180	\$866.00

## Aspect ISR G2 Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<b>Pendant Hand Control</b> For Aspect ISR. Controls elevate and tilt functions only (Included with table).	A100-2170	\$450.00
	<b>Panning Hand Control</b> For Aspect ISR G2 series. (Included with table).	A100-0178	\$2150.00
	<b>ISR Table Foot Switch with Safety Hood.</b> Foot control for Aspect ISR G2 table. Foot control motions are: iso-roll, transverse travel, elevate/lower and Trendelenburg down/up. (Safety hood helps avoid accidental motorized table movement from foot switch).	A100-2071	\$1645.00
	<b>IDI Medical Grade Anti-Fatigue Mat</b> (22"/56cm x 66"/167.6cm) Heavy-duty, one-piece construction for high-traffic areas. Safe & non-toxic, contains no plasticizers and will not leach over time. Durable brushed texture.	C000-1181	\$357.00
	<b>Side-rail clamp</b> Rotating, with socket for attachment of accessories with 15.88mm (5/8") circular post. Requires A100-2379 clamp-on accessory rail.	C000-1065	\$395.00

## Shipping Crate Specifications

- Weight (approx.): 750 lbs. (340kg)
- Length: 96" (2438mm).
- Width: 42" (1067mm).
- Height: 47" (1194mm).

## Aspect 100-4T G2 Table Catalog # A100-2533 or A100-2594

\$36,900.00



100-4T G2 A100-2533



100-4T G2 A100-2594

### Aspect 100-4T G2 Table Description:

**Mobile Imaging Table with Motorized Elevation, Trendelenburg Tilt and manual 4-way tabletop float.**

- 230V is standard – includes metric sized accessory rails
- 115V is optional and must be specified at time of order.

### Aspect 100-4T G2 #A100-2533

For *Advanced Vascular* Procedures; Includes short, recessed side rails 5" (127mm) in length at imaging end of table that facilitate use of removable headrest extension (option) and also facilitates attachment of catheter tray extension at narrow end of table. Two separate catheter tray extensions are available: One for the pedestal end of table and one for the imaging end of table.

Includes the following:

- X100-2492: Tabletop 19.5"/495mm x 22"/559mm x 24"/610mm X 84"/2134mm
- X100-2310 Tabletop Pad, 2" (51mm) thick.
- A100-0179: Panning hand control.
- C000-0328: Patient restraint straps.



### Aspect 100-4T G2 #A100-2594

Contoured shape, 22" (559mm) wide at imaging end, 24" (610mm) wide at pedestal end, 24" (610mm) long. Does not include 5" (127mm) long rails at imaging end of top. Not compatible with optional headrest. Optional catheter tray attachment can be attached to pedestal end of table only.

Includes the following:

- X100-2426: Tabletop, 22" (559mm)/24" (610mm) x 84" (2134mm).
- X100-2489 Tabletop Pad, 2" (51mm) thick.
- A100-0179: Panning hand control
- \* C000-0328: Patient restraint straps



## Aspect 100-4T G2 Accessories Included:

- Tabletop pad, 2" thick  
# X100-2310 for the X100-2492 table top or # X100-2489 for the X100-2426 table top.
- Patient Restraint Straps. # C000-0328
- Panning Hand Control, rail mounted. #A100-0179






## Aspect 100-4T G2 Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<b>Clamp-on Accessory Rail</b> 10mm x 25mm x 203mm clamp-on accessory rail.	A100-2379	\$495.00
	<b>Armboard</b> Quick-release, rail mounted (requires Clamp-on Accessory Rail A100-2379, above)	A310-062	\$450.00
	<b>VAB Armboard</b> Carbon fiber, one-piece, extra wide armboard for vascular access procedures (body section slides under patient/table pad). Includes strap and pad.	A100-2244	\$1489.00
	<b>Arm Board</b> Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059	\$975.00
	<b>Headrest Extension - EU Rail</b> Articulating carbon fiber headrest, designed to mount to standard 10mm x 25mm EU accessory rails used on the imaging end of IDI mobile imaging table models ISR G2 (A100-2200) and 100-4T G2 (A100-2533). Consists of a carbon fiber plaoorm and a pad with loop/hook fastener strips. Requires the table to have an X100-2492 style tabletop. It can also accommodate a flexible anesthesia "hoop" style screen holder, IDI part number C000-0965. L = 462 mm (18.19"), B= 588 mm (23.13"), Pad gauge=51 mm (2"), Adjustment angle=±12°	A100-2505	\$4635.00






## Aspect 100-4T G2 Table Optional and Replacement Accessories:

Catheter tray extensions are optional and sold separately; see part numbers A100-1769 (pedestal end catheter tray extension), A100-2322 (imaging end catheter tray extension for use with table A100-2533) or A100-2387 (imaging end catheter tray extension for use with table A100-2594).

Photo	Description	Catalog Number	MSRP
	<b>Tabletop Catheter Tray Extension</b> 24.25"/616mm wide, 30"/762mm long. Mounts to pedestal end of table only. The maximum weight this extension can support is 20 LBS or 9 Kgs.	A100-1769	\$995.00
	<b>Tabletop Catheter Tray Extension</b> 20 3/4"/527mm x 27"/686mm long. Mounts to imaging end of table only. Requires tabletop X100-2492. The maximum weight this extension can support is 20 LBS or 9 Kgs.	A100-2322	\$995.00
	<b>Radiation Shield</b> Articulating, rail-mounted. Attaches to table side-rail at base end of table, either side, adjustable. Articulating joint allows shield to conform to table contour. Snap-on upper section above level of tabletop is removable.	A610-106	\$1995.00
	<b>Anesthesia Screen Holder</b> Flexible "hoop" style. Mounts directly to A100-2317 Headrest Extension. For use on accessory rails - a pair of rotating side rail clamps, part number C000-1065, are required. Rotating side rail clamps and clamp-on accessory rails are not required if table is ordered with the A100-2317 Headrest Extension.	C000-0965	\$389.00
	<b>Peripheral Extension</b> Carbon fiber extension, mounts to imaging end of table only. 22"/559mm wide, 12"/305mm long. Maximum weight for this extension is 25 lbs or 11.4 Kgs. - Requires tabletop X100-2492.	A100-2486	\$2795.00

## Aspect 100-4T G2 Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<p><b>Multi-position Arm Support</b> Attaches to IDI tables with accessory rails - or can be used with a clamp on accessory rail on tables or areas within a table that do not have accessory rails. Length: 23.5 inches (597mm).</p>	C000-1180	\$866.00
	<p><b>Side-rail clamp</b> Rotating, with socket for attachment of accessories with 15.88mm (5/8") circular post. Requires A100-2379 clamp-on accessory rail.</p>	C000-1065	\$395.00
	<p><b>IDI Medical Grade Anti-Fatigue Mat</b> (22"/56cm x 66"/167.6cm) Heavy-duty, one-piece construction for high-traffic areas. Safe &amp; non-toxic, contains no plasticizers and will not leach over time. Durable brushed texture.</p>	C000-1181	\$357.00
	<p><b>Foot Control for 1004T and 1004T G2</b> (must verify table serial # with IDI for compatibility)</p>	A100-2360	\$1475.00
	<p><b>Panning Hand Control, for 100-4T and G2 series. (Included with table).</b></p>	A270-0179	\$1950.00

## Shipping Crate Specifications

- Weight (approx.): 750 lbs. (340kg)
- Length: 96" (2438mm).
- Width: 42" (1067mm).
- Height: 47" (1194mm).

## Aspect 100-4 Table Catalog # A100-2653 and A100-2654

\$25,995.00



### Aspect 100-4 Description:

The Aspect 100-4 is designed for C-arm imaging in endovascular or other interventional procedures where a 4-way floating tabletop is required in a non-tilting table.

- 4-Way Manual Floating Tabletop
- Motorized Elevation
- 230V is standard – includes metric sized accessory rails
- 115V is optional and must be specified at time of order.

### Aspect 100-4 #A100-2653

Includes the following:

- X100-517: Tabletop, 18" (457mm)/30" (762mm) x 80" (2032mm) tabletop.
- X100-1780: Tabletop pad, 2" (51mm) thick, for 18" (457mm)/30" (762mm) x 80 (2032 mm) tabletop.
- A100-435: Tabletop width extender board (lateral arms support) for 18" (457mm) wide tabletop. /24" (610mm) X 84" (2134mm)
- A270-183 Panning hand control.
- C000-0328: Patient restraint straps.



### Aspect 100-4 #A100-2654

Includes the following:



- X100-982: Tabletop 23.5" (597mm)"/30" (711mm) x 80" (2032mm) contoured shape.
- X100-1767 Tabletop pad, 2" (51mm) thick, for 23.5" (597mm)"/30" (711mm) x 80" (2032mm) tabletop.
- A270-183 Panning hand control.
- C000-0328: Patient restraint straps.



## Aspect 100-4 Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<b>Tabletop Catheter Tray Extension</b> Mounts to pedestal end only.	A100-814	\$995.00
	<b>Armboard</b> Quick-release, rail mounted. Requires Clamp-on Accessory Rail A100-2379	A310-062	\$450.00
	<b>Clamp-on Accessory Rail</b> 10mm x 25mm x 203mm clamp-on accessory rail.	A100-2379	\$495.00
	<b>Arm Board</b> Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059	\$975.00
	<b>VAB Armboard</b> Carbon fiber, one-piece, extra wide armboard for vascular access procedures (body section slides under patient/table pad). Includes strap and pad.	A100-2244	\$1489.00
	<b>Radiation Shield</b> Articulating, rail-mounted. Attaches to table side-rail at base end of table, either side, adjustable. Articulating joint allows shield to conform to table contour. Snap-on upper section above level of tabletop is removable.	A610-106	\$1995.00
	<b>ID410 Extender Board</b> Tabletop Width Extender Board with pads: lateral arms support. Use with Aspect 100-4 tables with 18" or 22" wide tops.	18"/ 457mm: A100-435  22"/ 559mm: A100-780	18"/457mm: \$705.00  22"/559mm: \$705.00
	<b>Panning Hand Control for Aspect 100-4</b> (Included with table)	A270-183	\$1150.00

## Aspect 100-4 Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<p><b>Multi-position Arm Support</b> Attaches to IDI tables with accessory rails - or can be used with a clamp on accessory rail on tables or areas within a table that do not have accessory rails. Length: 23.5 inches (597mm).</p>	C000-1180	\$866.00
	<p><b>IDI Medical Grade Anti-Fatigue Mat</b> (22"/56cm x 66"/167.6cm) Heavy-duty, one-piece construction for high-traffic areas. Safe &amp; non-toxic, contains no plasticizers and will not leach over time. Durable brushed texture.</p>	C000-1181	\$357.00

## Aspect 100-4 Shipping Crate Specifications

- Weight (approx.): 635 lbs. (288kg)
- Length: 96" (2438mm)
- Width: 42" (1067mm)
- Height: 47" (1194mm)

## Aspect 100-RTL Mobile Imaging Table Catalog # A100-2637, A200-2647 & A100-2673

\$21,500.00



### Aspect RTL Description:

C-Arm table for pain management and other imaging procedures where motorized elevation, tilt, lateral roll and longitudinal motion are required.

- Motorized elevation
- Trendelenburg tilt
- Lateral roll
- Longitudinal motion
- 230V is standard
- 115V is optional and must be specified at time of order.

### Aspect RTL #A100-2637

Includes the following:

- X100-995: Rectangular tabletop 24" (610mm) x 80" (2032mm).
- X100-1742: Tabletop pad, 2" (51mm) thick, for 24" (610mm) x 80" (2032mm) top.
- K000-0197: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.

### Aspect RTL #A100-2647

Includes the following:

- X100-733: Tabletop with tapered end and facial cutout 22" (559mm) x 80" (2032mm).
- X100-1757: Tabletop pad, 2" (51mm) thick, for 22" (559mm) x 80" (2032mm) top with facial cutout.
- K000-0197: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.

## Aspect 100RTL #A100-2673

Includes the following:




- X100-1325: Rectangular tabletop 28" (711mm) x 80" (2032mm).
- X100-2083: Tabletop pad, 2" (51mm) thick, for 28" (711mm) x 80" (2032mm).
- K000-0197: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.



### Aspect 100RTL Table Optional and Replacement Accessories:




Photo	Description	Catalog Number	MSRP
	<b>Clamp-on Accessory Rail</b> 10mm x 25mm x 203mm clamp-on accessory rail.	A100-2379	\$495.00
	<b>Side-rail clamp</b> Rotating, with socket for attachment of accessories with 15.88mm (5/8") circular post. Requires A100-2379 clamp-on accessory rail.	C000-1065	\$395.00
	<b>Armboard</b> Quick-release, rail mounted. Requires Clamp-on Accessory Rail A100-2379	A310-062	\$450.00
	<b>Arm Board</b> Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059	\$975.00
	<b>Multi-position Arm Support</b> Attaches to IDI tables with accessory rails - or can be used with a clamp on accessory rail on tables or areas within a table that do not have accessory rails. Length: 23.5 inches (597mm).	C000-1180	\$866.00
	<b>Hand Control, 100RTL</b> (must verify table serial # with IDI for compatibility). Motions controlled: Elevation, trendelenburg, roll, slide. (Included with table).	K000-0197	\$305.00
	<b>Foot Control for 100RTL</b> Foot Control unit controls all table motions.	A100-1961	\$1535.00

## Aspect 100RTL Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	Scatter Guard radiation shield draping For 28" (711mm) x 80" (2032mm) wide Optional Tabletop X100-1325 28"x80" (71x203cm) Dimensions: 59" (150cm) long. 24" (61cm) draping below tabletop. Draping extends along both sides and head-- end of table. Supported by fabric section over tabletop. Not for use with arm board A310-062.	A610-108	\$1595.00
	<b>Radiation Shield</b> Fixed shield with integrated clamp system. For use with tables that do not have side rails. Clamps to table sides or at head-end of table. Snap-on upper section above level of tabletop is re- movable. Use with 100RTL, 100RT, 100T, 100E tables.	A610-080	\$1295.00
	<b>Scatter-Guard radiation shield draping</b> For for 24"/61cm rectangular tabletop X100-995. Dimensions: 59"/150cm long. 24"/61cm draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop <i>Note: Not for use with arm board A310-062.</i>	A610-096	\$1595.00
	<b>Scatter-Guard radiation shield draping</b> for 22"/56cm tabletop X100-733 models with 22"/56cm tapered-end tabletop with facial cutout. Dimensions: 59"/150cm long. 24"/61cm draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop. <i>Note: Not for use with arm board A310-062.</i>	A610-092	\$1595.00
	<ul style="list-style-type: none"> <li>• Facial Cushion, crescent shape.</li> <li>• Disposable covers for facial cushion, box of 50.</li> </ul>	Cushion: #C000-0597  Covers: #C000-0598	Cushion: \$70.00  Covers: \$40.00
	<b>IDI Medical Grade Anti-Fatigue Mat</b> (22"/56cm x 66"/167.6cm) Heavy-duty, one-piece construction for high- traffic areas. Safe & non-toxic, contains no plasticizers and will not leach over time. Durable brushed texture.	C000-1181	\$357.00



## Aspect 100RTL Tabletop options

	<p><b>Standard Tabletop X100-995</b>                  24" x 80" (61 x 203cm)                  Low attenuation, carbon fiber construction.                  Standard tops may be substituted with optional sizes/styles as shown below at no extra charge.</p>
	<p><b>Optional Tabletop X100-733 with facial cutout</b>                  22" x 80" (56 x 203cm). Low attenuation, carbon fiber construction. Optional tabletop available at no extra charge.                  Must specify top selection with order.</p>
	<p><b>Optional Tabletop X100-1325</b>                  28" x 80" (71 x 203cm)                  Wide for larger patients. Low attenuation, carbon fiber construction. Must specify top selection with order.</p>

## 100RTL Shipping Crate Specifications

- Weight (approx.): 605 lbs. (274kg)
- Length: 96" (2438mm).
- Width: 42" (1067mm).
- Height: 47" (1194mm)

## Aspect 100RT Mobile Imaging Table

Catalog # A100-2636, A100-2645 & A100-2661

\$17,500.00



### Aspect 100RT Description:

C-arm table for pain management and other imaging procedures where motorized elevation, tilt and lateral roll are required.

- Motorized Elevation
- Trendelenburgtilt
- Lateral Roll
- 230V is standard
- 115V is optional and must be specified at time of order.

### Aspect 100RT #A100-2636

Includes the following:

- X100-995: Rectangular tabletop 24" (610mm) x 80" (2032mm).
- X100-1742: Tabletop pad, 2" (51mm) thick, for 24" (610mm) x 80" (2032mm) top
- K000-0195: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.

### Aspect 100RT #A100-2645

Includes the following:

- X100-733: Tabletop with tapered end and facial cutout 22" (559mm) x 80" (2032mm).
- X100-1757: Tabletop pad, 2" (51mm) thick, for 22" (559mm) x 80" (2032mm) top with facial cutout.
- K000-0195: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.

## Aspect 100RT #A100-2661

Includes the following:





- X100-1325: Rectangular tabletop 28" (711mm) x 80" (2032mm).
- X100-2083: Tabletop pad, 2" (51mm) thick, for 28" (711mm) x 80" (2032mm).
- K000-0195: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.






## Aspect 100RT Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<b>Side-rail clamp</b> Rotating, with socket for attachment of accessories with 15.88mm (5/8") circular post. Requires A100-2379 clamp-on accessory rail.	C000-1065	\$395.00
	<b>Armboard</b> Quick-release, rail mounted. Requires Clamp-on Accessory Rail A100-2379	A310-062	\$450.00
	<b>Clamp-on Accessory Rail</b> 10mm x 25mm x 203mm clamp-on accessory rail.	A100-2379	\$495.00
	<b>Arm Board</b> Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059	\$975.00
	<b>Multi-position Arm Support</b> Attaches to IDI tables with accessory rails - or can be used with a clamp on accessory rail on tables or areas within a table that do not have accessory rails. Length: 23.5 inches (597mm).	C000-1180	\$866.00
	<b>Scatter Guard radiation shield draping</b> For 28"/71cm wide Optional Tabletop X100-1325 28"/ 71cm x 80"/ 203cm) Dimensions: 59"/150cm long. 24"/61cm draping below tabletop. Draping extends along both sides and head--end of table. Supported by fabric section over tabletop. <i>Not for use with arm board A310-062.</i>	A610-108	\$1595.00

## Aspect 100RT Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<b>Hand Control, 100RT</b> Must verify table serial # with IDI for compatibility. (Included with table).	K000-0195	\$305.00
	<b>Foot Control</b> for 100RT (must verify table serial # with IDI for compatibility).	A100-1960	\$1490.00
	<b>Radiation Shield</b> Fixed shield with integrated clamp system. For use with tables that do not have side rails. Clamps to table sides or at head-end of table. Snap-on upper section above level of tabletop is removable. Use with 100RTL, 100RT, 100T, 100E tables.	A610-080	\$1295.00
" data-bbox="131 448 276 528"/>	<b>Scatter-Guard radiation shield draping</b> For for 24"/61cm rectangular tabletop X100-995. Dimensions: 59"/150cm long. 24"/61cm draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop <i>Note: Not for use with arm board A310-062.</i>	A610-096	\$1595.00
" data-bbox="161 588 271 686"/>	<b>Scatter-Guard radiation shield draping</b> for 22"/56cm tabletop X100-733 models with 22"/56cm tapered-end tabletop with facial cutout. Dimensions: 59"/150cm long. 24"/61cm draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop. <i>Note: Not for use with arm board A310-062.</i>	A610-092	\$1595.00
	<ul style="list-style-type: none"> <li>• Facial Cushion, crescent shape.</li> <li>• Disposable covers for facial cushion, box of 50.</li> </ul>	Cushion: #C000-0597  Covers: #C000-0598	Cushion: \$70.00  Covers: \$40.00
	<b>IDI Medical Grade Anti-Fatigue Mat</b> (22"/56cm x 66"/167.6cm) Heavy-duty, one-piece construction for high-traffic areas. Safe & non-toxic, contains no plasticizers and will not leach over time. Durable brushed texture.	C000-1181	\$357.00

## Aspect 100RT Tabletop options

	<p><b>Standard Tabletop X100-995</b>                  24" x 80" (61 x 203cm)                  Low attenuation, carbon fiber construction.                  Standard tops may be substituted with optional sizes/styles as shown below at no extra charge.</p>
	<p><b>Optional Tabletop X100-733 with facial cutout</b>                  22" x 80" (56 x 203cm). Low attenuation, carbon fiber construction. Optional tabletop available at no extra charge.                  Must specify top selection with order.</p>
	<p><b>Optional Tabletop X100-1325</b>                  28" x 80" (71 x 203cm)                  Wide for larger patients. Low attenuation, carbon fiber construction. Must specify top selection with order.</p>

## 100RT Shipping Crate Specifications

- Weight (approx.): 571 lbs. (259kg.)
- Length: 96" (2438mm)
- Width: 42" (1067mm)
- Height: 47" (1194mm)

## Aspect 100T Mobile Imaging Table

Catalog # A100-2638, A100-2649 & A100-2659

\$15,500.00



### Aspect 100T Description:

C-arm table for pain management and other procedures where only motorized elevation and Trendelenburg tilt is required.

- Motorized Elevation
- Trendelenburg tilt
- 230V is standard
- 115V is optional and must be specified at time of order.

### Aspect 100T #A100-2638

Includes the following:

- X100-995: Rectangular tabletop 24" (610mm) x 80" (2032mm).
- X100-1742: Tabletop pad, 2" (51mm) thick, for 24" (610mm) x 80" (2032mm) top
- K000-0054: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.

### Aspect 100T #A100-2649

Includes the following:

- X100-733: Tabletop with tapered end and facial cutout 22" (559mm) x 80" (2032mm).
- X100-1757: Tabletop pad, 2" (51mm) thick, for 22" (559mm) x 80" (2032mm) top with facial cutout.
- K000-0054: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.

## Aspect 100T #A100-2659

Includes the following:





- X100-1325: Rectangular tabletop 28" (711mm) x 80" (2032mm).
- X100-2083: Tabletop pad, 2" (51mm) thick, for 28" (711mm) x 80" (2032mm).
- K000-0054: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.



## Aspect 100T Table Optional and Replacement Accessories:




Photo	Description	Catalog Number	MSRP
	<b>Clamp-on Accessory Rail</b> 10mm x 25mm x 203mm clamp-on accessory rail.	A100-2379	\$495.00
	<b>Side-rail clamp</b> Rotating, with socket for attachment of accessories with 15.88mm (5/8") circular post. Requires A100-2379 clamp-on accessory rail.	C000-1065	\$395.00
	<b>Arm Board</b> Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059	\$975.00
	<b>Multi-position Arm Support</b> Attaches to IDI tables with accessory rails - or can be used with a clamp on accessory rail on tables or areas within a table that do not have accessory rails. Length: 23.5 inches (597mm).	C000-1180	\$866.00
	<b>Scatter Guard radiation shield draping</b> For 28"/71cm wide Optional Tabletop X100-1325 28"/71cm x 80"/203cm) Dimensions: 59"/150cm long. 24"/61cm draping below tabletop. Draping extends along both sides and head--end of table. Supported by fabric section over tabletop. <i>Not for use with arm board A310-062.</i>	A610-108	\$1595.00
	<b>Hand Control</b> for 100T table. Elevate and tilt functions. (Included with table).	K000-0054	\$230.00

## Aspect 100T Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<p><b>Foot Control for 100T</b></p>	A100-500	\$995.00
	<p><b>Radiation Shield</b> Fixed shield with integrated clamp system. For use with tables that do not have side rails. Clamps to table sides or at head-end of table. Snap-on upper section above level of tabletop is removable. Use with 100RTL, 100RT, 100T, 100E tables.</p>	A610-080	\$1295.00
	<p><b>Scatter-Guard radiation shield draping</b> For for 24"/61cm rectangular tabletop X100-995. Dimensions: 59"/150cm long. 24"/61cm draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop <i>Note: Not for use with arm board A310-062.</i></p>	A610-096	\$1595.00
	<p><b>Scatter-Guard radiation shield draping</b> for 22"/56cm tabletop X100-733 models with 22"/56cm tapered-end tabletop with facial cutout. Dimensions: 59"/150cm long. 24"/61cm draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop. <i>Note: Not for use with arm board A310-062.</i></p>	A610-092	\$1595.00
	<ul style="list-style-type: none"> <li>• Facial Cushion, crescent shape.</li> <li>• Disposable covers for facial cushion, box of 50.</li> </ul>	<p>Cushion: #C000-0597 Covers: #C000-0598</p>	<p>Cushion: \$70.00 Covers: \$40.00</p>
	<p><b>IDI Medical Grade Anti-Fatigue Mat</b> (22"/56cm x 66"/167.6cm) Heavy-duty, one-piece construction for high-traffic areas. Safe &amp; non-toxic, contains no plasticizers and will not leach over time. Durable brushed texture.</p>	C000-1181	\$357.00



## Aspect 100T Tabletop options

	<p><b>Standard Tabletop X100-995</b>                  24" x 80" (61 x 203cm)                  Low attenuation, carbon fiber construction.                  Standard tops may be substituted with optional sizes/styles as shown below at no extra charge.</p>
	<p><b>Optional Tabletop X100-733 with facial cutout</b>                  22" x 80" (56 x 203cm). Low attenuation, carbon fiber construction. Optional tabletop available at no extra charge.                  Must specify top selection with order.</p>
	<p><b>Optional Tabletop X100-1325</b>                  28" x 80" (71 x 203cm)                  Wide for larger patients. Low attenuation, carbon fiber construction. Must specify top selection with order.</p>

## 100RT Shipping Crate Specifications

- Weight (approx.): 540 lbs./245kg
- Length: 96" (2438mm)
- Width: 42" (1067mm).
- Height: 47" (1194mm)

## Aspect 100E Mobile Imaging Table

Catalog # A100-2635, A100-2643 & A100-2657

\$12,995.00



### Aspect 100E Description:

C-arm table for pain management and other procedures where motorized elevation is the only movement required.

- 230V is standard
- 115V is optional and must be specified at time of order.

### Aspect 100E #A100-2635

Includes the following:

- X100-995: Rectangular tabletop 24" (610mm) x 80" (2032mm).
- X100-1742: Tabletop pad, 2" (51mm) thick, for 24" (610mm) x 80" (2032mm) top
- K000-0041: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.

### Aspect 100E #A100-2643

Includes the following:

- X100-733: Tabletop with tapered end and facial cutout 22" (559mm) x 80" (2032mm).
- X100-1757: Tabletop pad, 2" (51mm) thick, for 22" (559mm) x 80" (2032mm) top with facial cutout.
- K000-0041 Pendant hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.

## Aspect 100E #A100-2657




Includes the following:

- X100-1325: Rectangular tabletop 28" (711mm) x 80" (2032mm).
- X100-2083: Tabletop pad, 2" (51mm) thick, for 28" (711mm) x 80" (2032mm).
- K000-0041: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.




## Aspect 100E Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<b>Clamp-on Accessory Rail</b> 10mm x 25mm x 203mm clamp-on accessory rail.	A100-2379	\$495.00
	<b>Side-rail clamp</b> Rotating, with socket for attachment of accessories with 15.88mm (5/8") circular post. Requires A100-2379 clamp-on accessory rail.	C000-1065	\$395.00
	<b>Arm Board</b> Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059	\$975.00
	<b>Multi-position Arm Support</b> Attaches to IDI tables with accessory rails - or can be used with a clamp on accessory rail on tables or areas within a table that do not have accessory rails. Length: 23.5 inches (597mm).	C000-1180	\$866.00
	<b>Scatter Guard radiation shield draping</b> For 28"/71cm wide Optional Tabletop X100-1325 28"/ 71cm x 80"/ 203cm) Dimensions: 59"/150cm long., 24"/61cm draping below tabletop. Draping extends along both sides and head--end of table. Supported by fabric section over tabletop. <i>Not for use with arm board A310-062.</i>	A610-108	\$1595.00
	<b>Hand control for 100E.</b> (Included with table).	K000-0041	\$85.00

## Aspect 100E Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<p><b>Radiation Shield</b> Fixed shield with integrated clamp system. For use with tables that do not have side rails. Clamps to table sides or at head-end of table. Snap-on upper section above level of tabletop is removable. Use with 100RTL, 100RT, 100T, 100E tables.</p>	A610-080	\$1295.00
	<p><b>Scatter-Guard radiation shield draping</b> For for 24"/61cm rectangular tabletop X100-995. Dimensions: 59"/150cm long. 24"/61cm draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop <i>Note: Not for use with arm board A310-062.</i></p>	A610-096	\$1595.00
	<p><b>Scatter-Guard radiation shield draping</b> for 22"/56cm tabletop X100-733 models with 22"/56cm tapered-end tabletop with facial cutout. Dimensions: 59"/150cm long. 24"/61cm draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop. <i>Note: Not for use with arm board A310-062.</i></p>	A610-092	\$1595.00
	<p><b>Foot Control 100E</b></p>	A100-758	\$1535.00
	<ul style="list-style-type: none"> <li>• Facial Cushion, crescent shape.</li> <li>• Disposable covers for facial cushion, box of 50.</li> </ul>	<p>Cushion: #C000-0597</p> <p>Covers: #C000-0598</p>	<p>Cushion: \$70.00</p> <p>Covers: \$40.00</p>
	<p><b>IDI Medical Grade Anti-Fatigue Mat</b> (22"/56cm x 66"/167.6cm) Heavy-duty, one-piece construction for high-traffic areas. Safe &amp; non-toxic, contains no plasticizers and will not leach over time. Durable brushed texture.</p>	C000-1181	\$357.00

## Aspect 100E Tabletop options

	<p><b>Standard Tabletop X100-995</b>                  24" x 80" (61 x 203cm)                  Low attenuation, carbon fiber construction.                  Standard tops may be substituted with optional sizes/styles as shown below at no extra charge.</p>
	<p><b>Optional Tabletop X100-733 with facial cutout</b>                  22" x 80" (56 x 203cm). Low attenuation, carbon fiber construction. Optional tabletop available at no extra charge.                  Must specify top selection with order.</p>
	<p><b>Optional Tabletop X100-1325</b>                  28" x 80" (71 x 203cm)                  Wide for larger patients. Low attenuation, carbon fiber construction. Must specify top selection with order.</p>

## 100RT Shipping Crate Specifications

- Weight (approx.): 495 lbs. (225kg)
- Length: 96" (2438mm)
- Width: 42" (1067mm).
- Height: 47" (1194mm).

## Aspect 100-UCPLUS Urology/ Multi-purpose Catalog # A100-2639

\$37,250.00



Shown with optional leg holders (#A100-2247).



Shown with tabletop extension (A100-2675)

### Aspect 100UCPlus description:

Imaging table for urology and multi-purpose use. Urological accessories are easily removable so the table can be quickly configured for multi-purpose C-arm imaging or general patient positioning for non-urological procedures. A removable radiographic cassette tray provides KUB radiography compatibility with mobile x-ray machines, including both conventional film/grid systems and digital detectors.

- 230V is standard – includes metric sized accessory rails
- 115V is optional and must be specified at time of order.

### Aspect 100UCPlus #A100-2369





Includes the following:

- A100-2379 9.525mm x 25mm x 203.2mm clamp-on accessory rail
- #A100-2676 Stirrups Mounting Rail extension Kit
- X100-1818: tabletop pad, 28" (711mm) x 46" (1168mm) x 2" (51mm) thick
- A100-2675 Tabletop extension, 24" (610mm)W x 34" (864mm)L
- X100-1817: Tabletop extension pad, 24" (610mm) x 34" (864mm) x 2" (51mm) thick
- K000-0158: Pendant hand control
- A100-1540: Foot switch control
- Z100-1060: Removable x-ray shield flap, for perineal end of table
- C000-0571: Patient restraint straps
- A310-062: Arm board, rail mounted, with pad and strap
- C000-0595: Urology drain bag support hoop
- W100-1001: Radiographic cassette tray, removable
- C000-0593: Disposable drain bag & hose ass'y, box of 20
- C000-0492: Disposable covers for foot switch, box of 50
- C000-0631: disposable cover for x-ray shield flap (qty. of 1)
- C000-0645: disposable tabletop pad protector (qty. of 1)






## Aspect 100UCPlus Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<b>Lithotomy Leg holders (pair)</b> Boot-style Stirrups (pair), pneumatic-assist, with pair of side rail clamps (C000-1171). For patients weighing under 159kg.	A100-2678	\$5090.00
	<b>Heavy Duty Lithotomy Legholders (pair)</b> Heavy Duty, boot style, pneumatic assist with one-hand adjustment, with pads. For patients up to 525 lbs (238kkg). Includes pair of side-rail clamps.	A100-2679	\$7155.00
	<b>Physician Elbow Supports (pair)</b> Removable	C000-0594	\$1195.00
	<b>Armboard</b> Rail mounted, with pad and strap 1 arm board included with table.	A310-062	\$450.00
	<b>Split Leg Holders</b> For use with Aspect 100UC Plus Tables. Includes pads. Prone percutaneous, lateral or modified lateral nephrostomy. Includes two split leg holder clamps # C000-1217.	C000-1216	\$6400.00
	<b>Replacement Split Leg Holder Clamp</b> For use with split leg holders. # C000-1115	C000-1217	\$125.00
	<b>Stirrups Mounting Rail extension Kit</b> Included with table. For use with leg-holders (Leg-holders are sold separately).	#A100-2676	\$525.00

## Aspect 100UCPlus Optional and Replacement Accessories:






Photo	Description	Catalog Number	MSRP
	<b>Multi-position Arm Support</b> Attaches to IDI tables with accessory rails - or can be used with a clamp on accessory rail on tables or areas within a table that do not have accessory rails. Length: 23.5 inches (597mm).	C000-1180	\$866.00
	<b>Hand Control for 100UCPlus and 100UC.</b> Controls height up/down, trendelenburg/ revers, slide, cross table movements. (Included with table).	K000-0158	\$185.00
	<b>Foot Switch Control for 100UCPlus.</b> (Included with table).	A100-1540	\$1855.00
	<b>Side-rail clamp</b> Rotating, with socket for attachment of accessories with 15.88mm (5/8") circular post. Requires A100-2379 clamp-on accessory rail.	C000-1065	\$395.00

## Aspect 100UCPlus Consumables

Photo	Description	Catalog Number	MSRP
	<b>Vinyl Urology Table Covers</b> Box of 20	C000-0645	\$350.00
	<b>Uro Drain Bags with Hose</b> Disposable drain bag and hose assembly. Box of 20	C000-0593	\$265.00
	<b>Disposable Covers for Aspect 100UCPlus Foot Controls</b> 24"/610mm x 20"/508mm Box of 50	C000-0492	\$50.00
	<b>Uro Drain Support Hoop</b> Collapsible & removable drain bag support hoop.	C000-0595	\$350.00
	<b>3"/ 76mm Restraint Straps (Pair)</b> Fabric with velcro fasteners for Aspect 100UCPlus table.	C000-0571	\$110.00



## Aspect 100UCPlus Consumables

Photo	Description	Catalog Number	
	<p><b>Restraint Straps</b> One each 3" and 5" (76mm and 127mm) Fabric with velcro fasteners for Aspect 100UCPlus table.</p>	C000-0328	\$110.00
	<p><b>Uro Collection Unit</b> 5 gallon/18.9 liter capacity, Box of 10. Soft sided Vinyl container gently expands into a cube as it collects waste fluid when floor drain or suction is unavailable. Internal baffle system enables container to fill evenly as powder converts liquid to gel. Built in zipper and handles allow for storage, closure of drain bag and hose assembly for enhanced infection control during transport and disposal.</p>	C000-0612	\$275.00
	<p><b>Uro Collection Unit</b> 2.6 gal/10 liter capacity, Box of 10. Flexible plastic sided container expands into a cube as it collects waste fluid when floor drain or suction is unavailable. Internal baffle system enables container to fill evenly as powder converts liquid to gel. Integrated valve system allows user to shut off the flow of liquid from the drain bag in order to change collection units..</p>	C000-0905	\$199.00
	<p><b>Uro Collection Unit</b> 5.3 gal/20 liter capacity, Box of 10. Flexible plastic sided container expands into a cube as it collects waste fluid when floor drain or suction is unavailable. Internal baffle system enables container to fill evenly as powder converts liquid to gel. valve Flow valve system allows user to shut off the flow of liquid from the drain bag in order to change collection units.</p>	C000-0916	\$289.00
	<p><b>Removable X-ray Shield Flap</b> For perineal end of table, 0.5mm Pb equivalent.</p>	Z100-1060	\$170.00

## Aspect 100UCPlus Shipping Crate Specifications

- Weight (approx.): 815 lbs. (370kg)
- Length: 96" (2438mm)
- Width: 42" (1067mm).
- Height: 47" (1194mm).

**Aspect 100-US Gyn/Fertility Clinic Table**  
**Catalog # A100-2674**

**\$14,950.00**



**Aspect 100US description:**

GYN/Fertility Mobile Imaging Table. An economical C-arm table for hysterosalpingograms.




- 230VAC is standard.
- 115VAC is optional and must be specified at time of order.

**Aspect 100US #A100-2674**





Includes the following:

- A100-758 Foot control
- X100-1687: tabletop pad, 24" (610mm) x 46" (1168mm) x 1" (25 mm) thick.
- A100-2379: 10mm x 25mm x 203mm clamp-on accessory rail.
- C000-0328: Patient restraint straps.
- A100-2380: Leg support tabletop extension with pad.
- X100-2351: Leg support tabletop extension pad.
- A100-2248: Basic Knee Crutches [pair] with a pair of rotating side rail clamps[C000-1065] for patients weighing up to 400 lbs.(181 kg)

**Aspect 100US Optional and Replacement Accessories:**

Photo	Description	Catalog Number	MSRP
	<b>Armboard</b> Quick-release, rail mounted. Requires Clamp-on Accessory Rail A100-2379	A310-062	\$450.00
	<b>Clamp-on Accessory Rail</b> 10mm x 25mm x 203mm clamp-on accessory rail.	A100-2379	\$495.00
	<b>Hand control for 100US</b>	K000-0041	\$85.00

## Aspect 100US Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<b>Basic Knee Crutches</b> For patients up to 400 lbs. <b>(181 kg)</b> Basic pair with pads. Includes two C000-1065 rotating side rail clamps.	A100-2677	\$1595.00
	<b>Foot Control 100US</b> (Included with table.)	A100-758	\$650.00
	<b>Side-rail clamp</b> Rotating, with socket for attachment of accessories with 15.88mm ( $\frac{5}{8}$ " ) circular post. Requires A100-2379 clamp-on accessory rail.	C000-1065	\$395.00
	<b>Multi-position Arm Support</b> Attaches to IDI tables with accessory rails - or can be used with a clamp on accessory rail on tables or areas within a table that do not have accessory rails. Length: 23.5 inches (597mm).	C000-1180	\$866.00

## Aspect 100UCPlus Shipping Crate Specifications

Weight (approx.): 245 kg / 540 lbs.

Length: 244 cm / 96 in. Width: 107 cm.

42 in. Height: 119 cm / 47 in.



## MDS Description:

MDS Mobile Display System, Complete with High Definition Video Monitors, Router and tablet controller.

- 230V is standard
- 115V is optional and must be specified at time of order.

## MDS #A280-0335

Includes the following:

- Dual-Arm Mobile Display Stand with adjustable counterbalanced arms
- Adjustable shelves for placement of ancillary equipment (Qty. of 4)
- Bilateral input source cable manifold supporting DVI, SDI, VGA, S-Video, SOG and composite inputs
- Built-in isolation transformer and 8 outlet power strip for video and ancillary equipment AC power
- FSN 26" (660mm) HD LED Widescreen Multi-Modality Video Monitors (Qty of 2).
- FSN IPS1000 Video Router/Controller with 10.1" Tablet Controller and articulated support arm with lock
- Dual 8" (20cm) "Preview" LCD Reference Monitors, mounted to MDS stand center column
- Qty (1) External cable storage hook
- Qty (1) Kit - Accessory Basket Kit
- Protective covers with storage for 26" (660mm) FSN monitors (Qty of 2)
- Video Cable kit, including common cables/connectors needed for connection of C-arm and equipment that will be mounted onto the MDS Mobile Display System. Customer choice of one of the following cable kits is included with the MDS System.



**Customers must specify cable kit choice with order:**

- **A280-0325:** Cable kit for GE OEC 9800/9900 C-Arms and camera cable kit for systems mounted onto the MDS
- **A280-0328:** Cable kit for all DVI output C-Arms and camera cable kit for systems mounted onto the MDS
- **A280-0326:** Cable kit for Siemens VGA C-Arm with camera cable kit for systems mounted onto the MDS
- **A280-0494:** Cable kit for Siemens Cios C-Arm with Video Manager and camera cable kit for systems mounted on the MDS

## MDS Optional and Replacement Accessories:

### MDS Optional Cable Kits

Not Included but Recommended

<p><b>External Source Cable Kit</b> Add this Kit if the customer plans to add any external device that is not mounted directly onto the MDS stand other than the C Arm (i.e. ultra-sound, PACS, hemodynamics) Kit Includes: qty (1) each of 20 ft. (6.1m) S Video, (2) 15 ft. (4.6m) lockable HDMI, 20 ft. (6.1m) SVG cable and qty. (2) lockable DVI/HDMI adapters.</p>	A282-0003	\$470.00	
Photo	Description	Catalog Number	MSRP
	<p><b>Endo Shelf</b> Bilateral slide-out shelf, increases storage and shelf space For MDS or MDS D.</p>	A280-0415	\$1441.00
	<p><b>CO<sub>2</sub> gas brackets</b> Pair of holders for type E gas cylinders (for insufflator use) mounted to sides of MDS stand.</p>	A280-0401	\$475.00
	<p><b>Shelf with Nook</b> Kit includes Kit- Shelf Nook Kit includes: medical rail shelf with under mounted nook that will accommodate a DAU switcher, freeing up shelf space. Customer installation only. Nook size: 4"/102mm H x 19"/483mm W.</p>	A282-0008	\$1257.00
	<p><b>Medical Rail Shelf</b> Shelf Includes medical rails. For accessory attachments.</p>	A280-0130	\$850.00
	<p><b>Shelf with Keyboard</b> Kit includes medical rail shelf with sliding keyboard carriage. Suitable for most standard customer-supplied keyboards.</p>	A282-0005	\$1461.00
	<p><b>Accessory Basket Kit</b> For customer-supplied foot switch or other storage.</p>	A282-0007	\$674.00
	<p><b>Suction Canister Kit</b> Includes: 2 suction canister rail clamps and brackets suitable for 3000 ml canister.</p>	A282-0009	\$357.00

## MDS Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<b>Cable Hook</b> Cable storage J-hook for external cables.	A280-0570	\$265.00
	<b>IV Pole Kit</b> Includes two IV rail clamps suitable for mounting a 4' (1.2 m) 2-hook IV pole onto two medical rail equipped (A280-0130) shelves. Note: this option requires two shelves with medical rails.	A280-0508	\$633.00
	<b>MDS Monitor Protection Covers F</b> or 26"/660mm flat panel monitor. (2 covers are included with MDS System A280-0330)	X280-0262	\$169.00
	<b>MDS Articulated Arm</b> Kit includes multi-joint arm (18 lb/8kg max load), assembly including qty (1) rail clamp, and articulated arm with VESA 100mm x 100mm mount	A282-0016	\$1100.00
	<b>MDS 3rd Monitor with articulated arm</b> Includes 3rd monitor (19"/483mm FSN), articulated arm and cables.	A282-0006	\$4320.00

## MDS Installation

Installation is the responsibility of the dealer



## MDS D Description

The Image Diagnostics MDS D provides a lower cost solution for flexible monitor positioning of 26"/660mm HD medical grade monitors. Images are manually PIP'ed at the monitor and changed at the IPS-500 switching device. The MDS D is recommended for those facilities where image routing requirements remain fairly consistent from day to day.

- 230V is standard
- 115V is optional and must be specified at time of order.

## MDS D #A280-0585

Includes the following:

- Dual counter balanced arms
- Two 26" / 660mm high definition monitors
- Video switcher/scalar for switching source images
- (1) 19.75"/502mm) x 19.75" / 502mm equipment shelf for placement of ancillary equipment. Includes under shelf "nook" for the video switcher
- Dual articulating arms with positioning handles and low-force manual motion for easy horizontal or vertical positioning of monitors
- Self-leveling monitor arms provide 36"/914mm of fully counter-balanced vertical travel  
Arm radius extends horizontally to 74"/188cm
- Stand height (with arms fully retracted): 64.5"/164 cm
- Extended arm vertical height: 96"/244 cm
- Lockable double-wheeled 6"/152mm casters with cable pushers
- Front and rear mounted ergonomic transport and positioning handles
- Front cross mounted medical rail for accessory attachment
- Qty (1) accessory basket
- Qty (1) cable storage hook
- Qty (2) Padded, Protective Covers for monitors with storage compartment for 26"/660mm
- Built-in isolation transformer and 8-outlet AC power strip for connection of video components and ancillary equipment

## MDS D Cable Kits – [Included but Selection is Needed](#)




For connection of C-arm and Endoscopy Systems. Customer choice of one of the following cable kits is included with MDS D System A280-0330 at no extra charge. Customer must specify the one selected cable kit with MDS System order. Prices shown below apply only to orders for replacement or additional cable kits.

<b>MDS Cable kit for all high resolution SOG or composite output C-Arms</b> Video connection kit for all SOG or composite out C-Arms including GE OEC 9800, 9900 C-arms.	A280-0491	\$261.00
<b>MDS D Cable kit for all DVI Output C-arms</b> Suitable for GE OEC Elite, Ziehm RFD, Siemens Cios, Phillips Veradius C-Arms.	A280-0492	\$1650.00
<b>MDS D Cable kit for VGA C-arm including Siemens Arcadis</b> Video connection kit.	A280-0493	\$2098.00
<b>MDS D Cable kit for Siemens Cios with Video Manager C-arm</b> Video connection kit.	A280-0494	\$2878.00

## MDS D Optional Cable Kits

Not Included but Recommended

<b>External Source Cable Kit</b> Add this Kit if the customer plans to add any external device that is not mounted directly onto the MDS stand other than the C Arm (i.e. ultrasound, PACS, hemodynamics) Kit Includes: qty (1) each of 20 ft. (6.1m) S Video, (2) 15 ft. (4.6m) lockable HDMI, 20 ft. (6.1m) SVG cable and qty. (2) lockable DVI/HDMI adapters.	A282-0003	\$470.00
---	-----------	----------


Photo	Description	Catalog Number	MSRP
	<b>Endo Shelf</b> Includes medical rail shelf with two bilateral slide-out shelves for MDS or MDS D. Factory installation only.	A280-0415	\$1441.00
	<b>CO<sub>2</sub> Gas Brackets</b> Pair of holders for type E gas cylinders (for insufflator use) mounted to sides of MDS stand (customer installation).	A280-0401	\$475.00
	<b>Shelf with Nook</b> Kit includes Kit- Shelf Nook Kit includes: medical rail shelf with under mounted nook that will accommodate a DAU switcher, freeing up shelf space. Customer installation only . Nook size: 4"/102mm H x 19"/483mm W.	A282-0008	\$1257.00



## MDS D Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<b>Medical Rail Shelf</b> Shelf Includes medical rails. Add if equipment needs may overload the standard single shelf configuration.	A280-0130	\$850.00
	<b>Suction Canister Kit</b> Includes: 2 suction canister rail clamps and brackets suitable for 3000 ml canister.	A282-0009	\$357.00
	<b>Accessory Basket Kit</b> Wire basket suitable for accessory storage.	A282-0007	\$674.00
	<b>Shelf with Keyboard</b> Kit includes: (A280-0130) medical rail shelf with sliding keyboard carriage for customer-supplied keyboard.	A282-0005	\$1461.00
	<b>Cable Hook</b> Cable storage J-hook for external cables.	A280-0570	\$265.00
	<b>IV Pole Kit</b> Includes two IV rail clamps suitable for mounting a 4' (1.2 m) 2-hook IV pole onto two medical rail equipped (A280-0130) shelves. Note: this option requires two shelves with medical rails.	A280-0508	\$633.00
	<b>MDS D Monitor Protection Covers</b> For 26" (660mm) flat panel monitors. Two covers are included with MDS D system A280-0230.	X280-0262	\$169.00
	<b>MDS D Articulated Arm</b> Kit includes multi-joint arm (18 lb/8kg max load), assembly including qty (1) rail clamp, and articulated arm with VESA 100mm x 100mm mount	A282-0016	\$1100.00

## MDS D Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<p><b>MDS 3rd Monitor with articulated arm</b> Includes 3rd monitor (19" FSN), articulated arm and cables.</p>	A282-0006	\$4320.00

### MDS D Installation


Installation is the responsibility of the dealer

### MDS D Shipping Crate Specifications

- Weight (approx.): 825 lbs. (374kg)
- Length: 65 in. (1651mm)
- Width: 42 in. (1067mm).
- Height: 72 in. (1829mm)

## ilex 32 video monitor and stand Catalog # A280-0743

\$9800.00



### ilex 32 Description

The Ilex 32 monitor and stand provides a cost effective solution for those customers requiring an additional monitor or monitors with mirrored images in multiple locations.

- 230V is standard
- 115V is optional and must be specified at time of order.

### ilex 32 #A280-0743

Includes the following:

- 32"/81cm HD medical grade monitor with medical bios and dual inputs
- Single or split-screen images
- Heavy duty cart stand with tilt (Maximum height: 76"/193cm to top of monitor)
- Cable channel
- Accessory basket
- 15 foot power cord (4.6 m)
- Cable hook
- Qty. (1) C-Arm kit including a pair of input cables qty (2) (20 ft. (6.1m) BNC and qty. (2) 20ft. (6.1m) DVI cables)

### Ilex 32 Installation

Installation is the responsibility of the dealer

## Aspect Four-Way Float Table Specifications & Accessories

•=Standard feature ◊=Optional NA=Not available	100 - 4	100 - 4T G2	ISR G2	
<b>Elevate (Height with 2" pad)</b>	29.9 - 39.7" 75.9 - 100.8 cm	32 - 42" 81 - 107 cm	32 - 40" 81 - 102 cm	
<b>Trendelenburg Tilt</b>	NA	±12°	±12°	
<b>Lateral Tilt</b>	NA	NA	±12°	
<b>Longitudinal Tabletop Motion</b>	40" float 102 cm	32" float 81 cm	32" float 81 cm	
<b>Transverse Tabletop Motion</b>	10" float 25 cm	8" float 20 cm	8" float 20 cm	
<b>Standard Tabletop Size</b>	18/30" x 80" 46/76 x 203 cm	24"/22"/19.5" x 84" 60/55.8/49.5 x 213 cm	24"/22"/19.5" x 84" 60/55.8/49.5 x 213 cm	
<b>Optional Tabletop Size</b>	23.5-30" x 80" 59-76 x 203 cm	22"/24" x 84" 56 /61 x 213 cm	22"/24" x 84" 56 /61 x 213 cm	
<b>Imaging Area</b>	72" 182.9 cm	61" 155 cm	61" (155cm), 74" (188 cm) with optional peripheral extension	
<b>Patient Capacity</b>	450 lb./204 kg	450 lb./204 kg	450 lb./204 kg	
<b>Table Weight</b>	490 lb. 222 kg	605 lb. 274 kg	605 lb. 274 kg	
<b>Tabletop Pad</b>	•	•	•	
<b>Restraint Straps</b>	•	•	•	
<b>Pendant Hand Control</b>	NA	NA	•	
<b>4-way Panning Hand Control</b>	•	•	•	
<b>Foot Control</b>	NA	Custom order)	Custom order	
<b>Armboard, rail mounted</b>	NA	◊	◊	
<b>Armboard, carbon fiber</b>	◊	◊	◊	
<b>VAB Armboard</b>	◊ Use of the VAB arm board with the 100-4 is only recommended on the models with 23.5"the table top . Use is only possible with a patient's weight securing the VAB to the table.	◊	◊	
<b>Clamp-on Accessory Rail</b>	◊	◊	◊	
<b>Catheter Tray extension</b>	◊	◊	◊	
<b>Headrest extension</b>	NA	◊	◊	
<b>Table-Mounted Radiation Shield</b>	◊	◊	◊	
<b>Anesthesia Screen Holder</b>	NA	◊	◊	

## Aspect Pain Management Table Specifications & Accessories

•Standard feature ◊Optional NA Not available	100E	100T	100RT	100RTL
<b>Elevate (Height with 2" pad)</b>	28 - 39.8" 71 - 101 cm	30 - 41.8" 76 - 106 cm	32.4 - 44.2" 82 - 112 cm	31.2 - 41.2" 79.2 - 104.6 cm
<b>Trendelenburg Tilt</b>	NA	±15°	±15°	±15°
<b>Lateral Tilt</b>	NA	NA	20°	20°
<b>Longitudinal Tabletop Motion</b>	NA	NA	NA	10"
<b>Standard Tabletop Size</b>	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm	24" x 80" 61 - 203 cm
<b>Optional Tabletop Size</b>	22" x 80" 56 x 203 cm with facial cut-out	22" x 80" 56 x 203 cm with facial cut-out	22" x 80" 56 x 203 cm with facial cut-out	22" x 80" 56 - 203 cm
<b>Optional Tabletop Size</b>	28" x 80" 72 x 203 cm	28" x 80" 72 x 203 cm	28" x 80" 72 x 203 cm	28" X 80" 72 x 203 cm
<b>Imaging Area</b>	66" 167.6 cm	62" 157.5 cm	62" 157.5 cm	59" 149.9 cm
<b>Patient Capacity</b>	450 lb. 204 kg	450 lb. 204 kg	450 lb. 204 kg	450 lbs. 204 kg
<b>Table Weight</b>	290 lb. 131 kg	325 lb. 147 kg	425 lb. 192 kg	355 lbs. 161 kg
<b>Tabletop Pad (2" thick)</b>	•	•	•	•
<b>Restraint Straps</b>	•	•	•	•
<b>Facial Support Cushion</b>	•	•	•	•
<b>Pendant Hand Control</b>	•	•	•	•
<b>Foot Control</b>	◊	◊	◊	◊
<b>Armboard, rail mounted</b>	◊	◊	◊	◊
<b>Armboard, carbon fiber</b>	◊	◊	◊	◊
<b>VAB Armboard</b>	◊	◊	◊	◊
<b>Clamp-on Accessory Rail</b>	◊	◊	◊	◊
<b>Scatter-Guard Radiation Shield</b>	◊	◊	◊	◊

## Aspect Urology & GYN Table Specifications & Accessories

•Standard feature ◊Optional NA Not available	100 UC-Plus	100 US
Elevate (Height with pad)	33.3" – 43.1" (84.5 – 109.4 cm)	27" – 39" (68.5 – 99cm)
Trendelenburg Tilt	± 15°	NA
Transverse travel	± 3.5", 7" total (± 8.9 cm, 17.8 cm total)	NA
Longitudinal Tabletop Motion	10" (25.4 cm)	NA
Tabletop Size	28"W x 46"L (71 x 117 cm)	24"W x 46"L (61 x 117 cm)
Tabletop Extension Size	24"W x 34"L (61 x 86 cm)	NA
Imaging Area	Main tabletop: 25"W x 20"L ( 63.5 x 50.8 cm)	22"W x 32"L (56 x 81 cm)
Imaging Area with table extension	54" (137 cm)	NA
Patient Weight Capacity	Upper body on main table: 525lb. (238 kg) Upper body on extension: 400 lbs ( 181kg) Load capacity on extension only: 200 lbs (90.7kg)	450 lbs (204 kg)
Table Weight	490 lb. (222 kg)	375 lbs ((170kg)
Tabletop Pad (2" thick)	•	•
Tabletop extension with pad	•	Leg support tabletop extension
Hand control pendant & footswitch control	•	Footswitch only
Patient restraint straps	•	•
Removeable drain bag support hoop	•	NA
Disposable drain bag & hose assembly, box of 20	•	NA
Removable radiation shield flap for perineal end of table	•	NA
Armboard, rail-mounted with pad and strap	•	◊
Split Leg Stirrups	•	NA
Disposable covers for footswitch, box of 50	•	◊
Radiographic cassette tray	•	NA

**Diagnostics, Inc.**

Toll free:	877.304.5434
Main:	978.829.0009
Fax:	978.829.0027

