

IDI PRODUCT CATALOG

IMAGING TABLES, ACCESSORIES & VIDEO INTEGRATION SYSTEMS

Quality products for surgical, interventional and diagnostic imaging

115V



Table of Contents

Imaging Tables and options

Aspect ISR G2 Table Accessories	3-6
Aspect 100-4T G2 Table Accessories	7-10
Aspect 100-4 Table Accessories	11-13
Aspect 100RTL Table Accessories	14-17
Aspect 100RT Table Accessories	18-21
Aspect 100T Table Accessories	22-25
Aspect 100E Table Accessories	26-29
Aspect 100UC <i>Plus</i> Urology Table	30-33
Aspect 100US Gyn Table	34-35

Video Integration Systems and options

MDS	36-39
MDS D	40-43
ilex 32	44

Comparison Matrices

Aspect Four-Way Float Top Tables Specifications & Accessory Matrix	45
Aspect Urology & GYN Table Specifications & Accessory Matrix	46
Aspect Pain Management Table Specifications & Accessory Matrix	47

Customer Service

IDI Customer Service Contacts	Back page
-------------------------------	-----------

The part numbers listed in this catalog represent 115VAC tables with SAE sized accessory rails. 230V systems are available. Part numbers for these can be provided upon request. The part numbers listed in this catalog are for 115V systems.

Aspect ISR G2 Table
Catalog # A100-1902 or A100-2640

\$49,900.00



ISR G2 A100-1902



ISR G2 A100-2640

Aspect ISR G2 Table Description:

Mobile Imaging Table with Motorized Elevation, Trendelenburg Tilt, Iso Roll and Manual 4-way Tabletop Float. For Advanced Vascular Surgery and Interventional Procedures.

- 115VAC is standard.
- 230VAC is optional and must be specified at time of order.

Aspect ISR G2 #A100-1902

For Advanced Vascular Procedures; Includes short, recessed side rails 5" (127mm) in length at imaging end of table that facilitate use of removable headrest extension (option) and also facilitates attachment of catheter tray extension at narrow end of table. Two separate catheter tray extensions are available: One for the pedestal end of table and one for the imaging end of table.

Includes the following:

- A100-2170 Pendant hand control
- X100-2492: Tabletop 19.5"/495mm x 22"/559mm x 24"/610mm X 84"/2134mm
- X100-2310 Tabletop Pad, 2" (51mm) thick.
- A270-0163: Panning hand control.
- C000-0328: Patient restraint straps.



Aspect ISR G2 #A100-2640

Contoured shape, 22" (559mm) wide at imaging end, 24" (610mm) wide at pedestal end, 24" (610mm) long. Does not include 5" (127mm) long rails at imaging end of top. Not compatible with optional headrest. Optional catheter tray attachment can be attached to pedestal end of table only.

Includes the following:

- A100-2170 Pendant hand control
- X100-2426: Tabletop, 22" (559mm)/24" (610mm) x 84" (2134mm).
- X100-2489 Tabletop Pad, 2" (51mm) thick.
- A270-0163: Panning hand control.
- C000-0328: Patient restraint straps.







Aspect ISR G2 Accessories Included:

- Tabletop pad, 2" thick
X100-2310 for the X100-2492 table top or # X100-2489 for the X100-2426 table top.
- Patient Restraint Straps. # C000-0328
- Panning Hand Control, rail mounted. #A270-0173

Aspect ISR G2 Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	Clamp-on Accessory Rail Rail is 6" long x 3/8" thick x 1-1/8" high (152 x 10 x 29 mm)	A100-1007	\$495.00
	Armboard Quick-release, rail mounted (requires Clamp-on Accessory Rail A100-1007, above)	A310-056	\$450.00
	VAB Armboard Carbon fiber, one-piece, extra wide armboard for vascular access procedures (body section slides under patient/table pad). Includes strap and pad.	A100-2244	\$1489.00
	Arm Board Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059	\$975.00
	Headrest Extension Articulating carbon fiber headrest, Attaches to the imaging end of the table. Requires the table to have an X100-2492 style tabletop. 23-1/8" / 587mm wide x 12" / 305mm long.	A100-2317	\$4635.00
Catheter tray extensions are optional and sold separately; see part numbers A100-1769 (pedestal end catheter tray extension), , A100-2322 (imaging end catheter tray extension for use with table A100-1416) or A100-2387 (imaging end catheter tray extension for use with table A100-2598).			
	Tabletop Catheter Tray Extension 24.25"/616mm wide, 30"/762mm long. Mounts to pedestal end of table only. The maximum weight this extension can support is 20 LBS or 9 Kgs.	A100-1769	\$995.00
	Tabletop Catheter Tray Extension 20 3/4"/527mm x 27"/686mm long. Mounts to imaging end of table only. Requires tabletop X100-2492. The maximum weight this extension can support is 20 LBS or 9 Kgs.	A100-2322	\$995.00

Aspect ISR G2 Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<p>Radiation Shield Articulating, rail-mounted. Attaches to table side-rail at base end of table, either side, adjustable. Articulating joint allows shield to conform to table contour. Snap-on upper section above level of tabletop is removable.</p>	A610-051	\$1995.00
	<p>Anesthesia Screen Holder Flexible "hoop" style. Mounts directly to A100-2317 Headrest Extension. For use on accessory rails - a pair of rotating side rail clamps, part number C000-0746, are required. Rotating side rail clamps and clamp-on accessory rails are not required if table is ordered with the A100-2317 Headrest Extension.</p>	C000-0965	\$389.00
	<p>Peripheral Extension Carbon fiber extension, mounts to imaging end of table only. 22"/559mm wide, 12"/305mm long. Maximum weight for this extension is 25 lbs or 11.4 Kgs. Requires the table to have an X100-2492 style tabletop.</p>	A100-2486	\$2795.00
	<p>Multi-position Arm Support Attaches to IDI tables with accessory rails - or can be used with a clamp on accessory rail on tables or areas within a table that do not have accessory rails. Length: 23.5 inches (597mm).</p>	C000-1180	\$866.00

Aspect ISR G2 Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	Pendant Hand Control For Aspect ISR. Controls elevate and tilt functions only (Included with table).	A100-2170	\$450.00
	Panning Hand Control For Aspect ISR and G2 series. (Included with table).	A270-0163	\$2150.00
	ISR Table Foot Switch with Safety Hood. Foot control for Aspect ISR G2 table. Foot control motions are: iso-roll, transverse travel, elevate/lower and Trendelenburg down/up. (Safety hood helps avoid accidental motorized table movement from foot switch).	A100-2071	\$1645.00
	IDI Medical Grade Anti-Fatigue Mat (22"/559mm x 66"/1676mm) Heavy-duty, one-piece construction for high-traffic areas. Safe & non-toxic, contains no plasticizers and will not leach over time. Durable brushed texture.	C000-1181	\$357.00
	Side-rail clamp Rotating, with socket for attachment of accessories with 5/8 inch (16 mm) circular post. Requires A100-1007 clamp- on accessory rail.	C000-0746	\$395.00

Shipping Crate Specifications

- Weight (approx.): 750 lbs. (340kg)
- Length: 96" (2438mm).
- Width: 42" (1067mm).
- Height: 47" (1194mm).

Aspect 100-4T G2 Table Catalog # A100-1416 or A100-2598

\$36,900.00



100-4T G2 A100-1416



100-4T G2 A100-2598

Aspect 100-4T G2 Table Description:

Mobile Imaging Table with Motorized Elevation, Trendelenburg Tilt and manual 4-way tabletop float.

- 115VAC is standard.
- 230VAC is optional and must be specified at time of order.

Aspect 100-4T G2 # A100-1416

For *Advanced Vascular* Procedures; Includes short, recessed side rails 5" (127mm) in length at imaging end of table that facilitate use of removable headrest extension (option) and also facilitates attachment of catheter tray extension at narrow end of table. Two separate catheter tray extensions are available: One for the pedestal end of table and one for the imaging end of table.

Includes the following:

- X100-2492: Tabletop 19.5" (495mm)/22" (559mm)/ 24" (610mm) X 84" (2134mm)
- X100-2310 Tabletop Pad, 2" (51mm) thick.
- A270-0173: Panning hand control.
- C000-0328: Patient restraint straps.



Aspect 100-4T G2 # A100-2598

Contoured shape, 22" (559mm) wide at imaging end, 24" (610mm) wide at pedestal end, 24" (610mm) long. Does not include 5" (127mm) long rails at imaging end of top. Not compatible with optional headrest. Optional catheter tray attachment can be attached to pedestal end of table only.

Includes the following:

- X100-2426: Tabletop, 22" (559mm)/24" (610mm) x 84" (2134mm).
- X100-2489 Tabletop Pad, 2" (51mm) thick.
- A270-0173: Panning hand control
- * C000-0328: Patient restraint straps



Aspect 100-4T G2 Accessories Included:

- Tabletop pad, 2" thick
X100-2310 for the X100-2492 table top or # X100-2489 for the X100-2426 table top.
- Patient Restraint Straps. # C000-0328
- Panning Hand Control, rail mounted. #A270-0173

Aspect 100-4T G2 Table Optional and Replacement Accessories:


Photo	Description	Catalog Number	MSRP
	Clamp-on Accessory Rail Rail is 6" long x 3/8" thick x 1-1/8" high (152 x 10 x 29 mm)	A100-1007	\$495.00
	Armboard Quick-release, rail mounted (requires Clamp-on Accessory Rail A100-1007).	A310-056	\$450.00
	VAB Armboard Carbon fiber, one-piece, extra wide armboard for vascular access procedures (body section slides under patient/table pad). Includes strap and pad.	A100-2244	\$1489.00
	Arm Board Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059	\$975.00
	Headrest Extension Articulating carbon fiber headrest, Attaches to the imaging end of the table. Requires the table to have an X100-2492 style table-top. 23-1/8" / 587mm wide x 12" / 305mm long	A100-2317	\$4635.00

Aspect 100-4T G2 Table Optional and Replacement Accessories:

Catheter tray extensions are optional and sold separately; see part numbers A100-1769 (pedestal end catheter tray extension), A100-2322 (imaging end catheter tray extension for use with table A100-1416) or A100-2387 (imaging end catheter tray extension for use with table A100-2598).

Photo	Description	Catalog Number	MSRP
	<p>Tabletop Catheter Tray Extension 24.25"/616mm wide, 30"/762mm long. Mounts to pedestal end of table only. The maximum weight this extension can support is 20 LBS or 9 Kgs.</p>	A100-1769	\$995.00
	<p>Tabletop Catheter Tray Extension 20 3/4"/527mm x 27"/686mm long. Mounts to imaging end of table only. Requires tabletop X100-2492. The maximum weight this extension can support is 20 LBS or 9 Kgs.</p>	A100-2322	\$995.00
	<p>Radiation Shield Articulating, rail-mounted. Attaches to table side-rail at base end of table, either side, adjustable. Articulating joint allows shield to conform to table contour. Snap-on upper section above level of tabletop is removable.</p>	A610-051	\$1995.00
	<p>Anesthesia Screen Holder Flexible "hoop" style. Mounts directly to A100-2317 Headrest Extension. For use on accessory rails - a pair of rotating side rail clamps, part number C000-0746, are required. Rotating side rail clamps and clamp-on accessory rails are not required if table is ordered with the A100-2317 Headrest Extension.</p>	C000-0965	\$389.00
	<p>Peripheral Extension Carbon fiber extension, mounts to imaging end of table only. 22"/559mm wide, 12"/305mm long. Maximum weight for this extension is 25 lbs or 11.4 Kgs. - Requires tabletop X100-2492.</p>	A100-2486	\$2795.00

Aspect 100-4T G2 Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<p>Multi-position Arm Support Attaches to IDI tables with accessory rails - or can be used with a clamp on accessory rail on tables or areas within a table that do not have accessory rails. Length: 23.5 inches (597mm).</p>	C000-1180	\$866.00
	<p>Side-rail clamp Rotating, with socket for attachment of accessories with 5/8 inch (16mm) circular post. Requires A100-1007 clamp-on accessory rail.</p>	C000-0746	\$395.00
	<p>IDI Medical Grade Anti-Fatigue Mat (22"/559mm x 66"/1676mm) Heavy-duty, one-piece construction for high-traffic areas. Safe & non-toxic, contains no plasticizers and will not leach over time. Durable brushed texture.</p>	C000-1181	\$357.00
	<p>Foot Control for 1004T and 1004T G2 (must verify table serial # with IDI for compatibility)</p>	A100-2360	\$1475.00
	<p>Panning Hand Control, for 100-4T and G2 series. (Included with table).</p>	A270-0173	\$1950.00

Shipping Crate Specifications

- Weight (approx.): 750 lbs. (340kg)
- Length: 96" (2438mm).
- Width: 42" (1067mm).
- Height: 47" (1194mm).

Aspect 100-4 Table Catalog # A100-566 and A100-2655

\$25,995.00



Aspect 100-4 Description:

The Aspect 100-4 is designed for C-arm imaging in endovascular or other interventional procedures where a 4-way floating tabletop is required in a non-tilting table.

- 4-Way Manual Floating Tabletop
- Motorized Elevation
- 115VAC is standard.
- 230VAC is optional and must be specified at time of order.

Aspect 100-4 #A100-566

Includes the following:

- X100-517: Tabletop, 18" (457mm)/30" (762mm) x 80" (2032mm) tabletop.
- X100-1780: Tabletop pad, 2" (51mm) thick, for 18" (457mm)/30" (762mm) x 80 (2032 mm) tabletop.
- A100-435: Tabletop width extender board (lateral arms support) for 18" (457mm) wide tabletop. /24" (610mm) X 84" (2134mm)
- A270-126 Panning hand control.
- C000-0328: Patient restraint straps.



Aspect 100-4 #A100-2655

Includes the following:



- X100-982: Tabletop 23.5" (597mm)"/30" (762mm) x 80" (2032mm) contoured shape.
- X100-1767 Tabletop pad, 2" (51mm) thick, for 23.5" (597mm)"/30" (762mm) x 80" (2032mm) tabletop.
- A270-126 Panning hand control.
- C000-0328: Patient restraint straps.



Aspect 100-4 Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	Tabletop Catheter Tray Extension Mounts to pedestal end only.	A100-814	\$995.00
	Armboard Quick-release, rail mounted. Requires Clamp-on Accessory Rail A100-1007	A310-056	\$450.00
	Clamp-on Accessory Rail Rail is 6" long x 3/8" thick x 1-1/2" high (152 x 10 x 29 mm)	A100-1007	\$495.00
	Arm Board Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059	\$975.00
	VAB Armboard Carbon fiber, one-piece, extra wide arm-board for vascular access procedures (body section slides under patient/table pad). Includes strap and pad.	A100-2244	\$1489.00
	Radiation Shield Articulating, rail-mounted. Attaches to table side-rail at base end of table, either side, adjustable. Articulating joint allows shield to conform to table contour. Snap-on upper section above level of tabletop is removable.	A610-051	\$1995.00
	ID410 Extender Board Tabletop Width Extender Board with pads: lateral arms support. Use with Aspect 100-4 tables with 18" or 22" wide tops.	18": A100-435 22": A100-780	18": \$705.00 22": \$705.00
	Panning Hand Control for Aspect 100-4 (Included with table)	A270-126	\$1150.00

Aspect 100-4 Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<p>Multi-position Arm Support Attaches to IDI tables with accessory rails - or can be used with a clamp on accessory rail on tables or areas within a table that do not have accessory rails. Length: 23.5 inches (597mm).</p>	C000-1180	\$866.00
	<p>IDI Medical Grade Anti-Fatigue Mat (22"/559mm x 66"/1676mm) Heavy-duty, one-piece construction for high-traffic areas. Safe & non-toxic, contains no plasticizers and will not leach over time. Durable brushed texture.</p>	C000-1181	\$357.00

Aspect 100-4 Shipping Crate Specifications

- Weight (approx.): 635 lbs. (288kg)
- Length: 96" (2438mm)
- Width: 42" (1067mm)
- Height: 47" (1194mm)

Aspect 100-RTL Fluoroscopic Procedure Table Catalog # A100-1254, A200-2646 & A100-2672

\$21,500.00



Aspect RTL Description:

C-Arm table for aping management and other imaging procedures where motorized elevation, tilt, lateral roll and longitudinal motion are required.

- Motorized elevation
- Trendelenburg tilt
- Lateral roll
- Longitudinal motion
- 115VAC is standard.
- 230VAC is optional and must be specified at time of order.

Aspect RTL #A100-1254

Includes the following:

- X100-995: Rectangular tabletop 24" (610mm) x 80" (2032mm).
- X100-1742: Tabletop pad, 2" (51mm) thick, for 24" (610mm) x 80" (2032mm) top.
- K000-0197: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.

Aspect RTL #A100-2646

Includes the following:

- X100-733: Tabletop with tapered end and facial cutout 22" (559mm) x 80" (2032mm).
- X100-1757: Tabletop pad, 2" (51mm) thick, for 22" (559mm) x 80" (2032mm) top with facial cutout.
- K000-0197: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.

Aspect 100RTL #A100-2672

Includes the following:






- X100-1325: Rectangular tabletop 28" (711mm) x 80" (2032mm).
- X100-2083: Tabletop pad, 2" (51mm) thick, for 28" (711mm) x 80" (2032mm).
- K000-0197: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.






Aspect 100RTL Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	Side-rail clamp Rotating, with socket for attachment of accessories with 5/8 inch (16mm) circular post. Requires A100-1007 clamp-on accessory rail.	C000-0746	\$395.00
	Armboard Quick-release, rail mounted. Requires Clamp-on Accessory Rail A100-1007	A310-056	\$450.00
	Clamp-on Accessory Rail Rail is 6" long x 3/8" thick x 1-1/8" high (152 x 10 x 29 mm)	A100-1007	\$495.00
	Arm Board Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059	\$975.00
	Multi-position Arm Support Attaches to IDI tables with accessory rails - or can be used with a clamp on accessory rail on tables or areas within a table that do not have accessory rails. Length: 23.5 inches (597mm).	C000-1180	\$866.00
	Hand Control, 100RTL (must verify table serial # with IDI for compatibility). Motions controlled: Elevation, trendelenburg, roll, slide. (Included with table).	K000-0197	\$305.00
	Foot Control for 100RTL Foot Control unit controls all table motions.	A100-1961	\$1535.00

Aspect 100RTL Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	Scatter Guard radiation shield draping For 28" wide Optional Tabletop X100-1325 28"x80" (71x203cm) Dimensions: 59" long. 24" draping below tabletop. Draping extends along both sides and head--end of table. Sup- ported by fabric section over tabletop. <i>Not for use with arm board A310-056.</i>	A610-108	\$1595.00
	Radiation Shield Fixed shield with integrated clamp system. For use with tables that do not have side rails. Clamps to table sides or at head-end of table. Snap-on upper section above level of tabletop is re- movable. Use with 100RTL, 100RT, 100T, 100E tables.	A610-080	\$1295.00
	Scatter-Guard radiation shield draping For for 24" rectangular tabletop X100-995. Dimensions: 59" long. 24" draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop <i>Note: Not for use with arm board A310-056.</i>	A610-096	\$1595.00
	Scatter-Guard radiation shield draping for 22" tabletop X100-733 models with 22" tapered-end tabletop with facial cutout. Dimensions: 59" long. 24" draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop. <i>Note: Not for use with arm board A310-056.</i>	A610-092	\$1595.00
	<ul style="list-style-type: none"> • Facial Cushion, crescent shape. • Disposable covers for facial cushion, box of 50. 	Cushion: #C000-0597 Covers: #C000-0598	Cushion: \$70.00 Covers: \$40.00
	IDI Medical Grade Anti-Fatigue Mat (22"/559mm x 66"/1676mm) Heavy-duty, one-piece construction for high- traffic areas. Safe & non-toxic, contains no plasticizers and will not leach over time. Durable brushed texture.	C000-1181	\$357.00

Aspect 100RTL Tabletop options

	<p>Standard Tabletop X100-995 24" x 80" (61 x 203cm) Low attenuation, carbon fiber construction. Standard tops may be substituted with optional sizes/styles as shown below at no extra charge.</p>
	<p>Optional Tabletop X100-733 with facial cutout 22" x 80" (56 x 203cm). Low attenuation, carbon fiber construction. Optional tabletop available at no extra charge. Must specify top selection with order.</p>
	<p>Optional Tabletop X100-1325 28" x 80" (71 x 203cm) Wide for larger patients. Low attenuation, carbon fiber construction. Must specify top selection with order.</p>

100RTL Shipping Crate Specifications

- Weight (approx.): 605 lbs. (274kg)
- Length: 96" (2438mm).
- Width: 42" (1067mm).
- Height: 47" (1194mm)

Aspect 100RT Fluoroscopic Procedure Table

Catalog # A100-870, A100-2644 & A100-2660

\$17,500.00



Aspect 100RT Description:

C-arm table for pain management and other imaging procedures where motorized elevation, tilt and lateral roll are required.

- Motorized Elevation
- Trendelenburg tilt
- Lateral Roll
- 115VAC is standard.
- 230VAC is optional and must be specified at time of order.

Aspect 100RT #A100-870

Includes the following:

- X100-995: Rectangular tabletop 24" (610mm) x 80" (2032mm).
- X100-1742: Tabletop pad, 2" (51mm) thick, for 24" (610mm) x 80" (2032mm) top
- K000-0195: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.



Aspect 100RT #A100-2644

Includes the following:

- X100-733: Tabletop with tapered end and facial cutout 22" (559mm) x 80" (2032mm).
- X100-1757: Tabletop pad, 2" (51mm) thick, for 22" (559mm) x 80" (2032mm) top with facial cutout.
- K000-0195: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.



Aspect 100RT #A100-2660

Includes the following:

- X100-1325: Rectangular tabletop 28" (711mm) x 80" (2032mm).
- X100-2083: Tabletop pad, 2" (51mm) thick, for 28" (711mm) x 80" (2032mm).
- K000-0195: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.




Aspect 100RT Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	Side-rail clamp Rotating, with socket for attachment of accessories with 5/8 inch (16mm) circular post. Requires A100-1007 clamp-on accessory rail.	C000-0746	\$395.00
	Armboard Quick-release, rail mounted. Requires Clamp-on Accessory Rail A100-1007	A310-056	\$450.00
	Clamp-on Accessory Rail Rail is 6" long x 3/8" thick x 1-1/8" high (152 x 10 x 29 mm)	A100-1007	\$495.00
	Arm Board Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059	\$975.00
	Multi-position Arm Support Attaches to IDI tables with accessory rails - or can be used with a clamp on accessory rail on tables or areas within a table that do not have accessory rails. Length: 23.5 inches (597mm).	C000-1180	\$866.00
	Scatter Guard radiation shield draping For 28" wide Optional Tabletop X100-1325 28"x80" (71x203cm) Dimensions: 59" long. 24" draping below tabletop. Draping extends along both sides and head--end of table. Supported by fabric section over tabletop. <i>Not for use with arm board A310-056.</i>	A610-108	\$1595.00

Aspect 100RT Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	Hand Control, 100RT Must verify table serial # with IDI for compatibility. (Included with table).	K000-0195	\$305.00
	Foot Control for 100RT (must verify table serial # with IDI for compatibility).	A100-1960	\$1490.00
	Radiation Shield Fixed shield with integrated clamp system. For use with tables that do not have side rails. Clamps to table sides or at head-end of table. Snap-on upper section above level of tabletop is removable. Use with 100RTL, 100RT, 100T, 100E tables.	A610-080	\$1295.00
	Scatter-Guard radiation shield draping For for 24" rectangular tabletop X100-995. Dimensions: 59" long. 24" draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop <i>Note: Not for use with arm board A310-056.</i>	A610-096	\$1595.00
	Scatter-Guard radiation shield draping for 22" tabletop X100-733 models with 22" tapered-end tabletop with facial cutout. Dimensions: 59" long. 24" draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop. <i>Note: Not for use with arm board A310-056.</i>	A610-092	\$1595.00
	<ul style="list-style-type: none"> • Facial Cushion, crescent shape. • Disposable covers for facial cushion, box of 50. 	Cushion: #C000-0597 Covers: #C000-0598	Cushion: \$70.00 Covers: \$40.00
	IDI Medical Grade Anti-Fatigue Mat (22"/559mm x 66"/1676mm) Heavy-duty, one-piece construction for high-traffic areas. Safe & non-toxic, contains no plasticizers and will not leach over time. Durable brushed texture.	C000-1181	\$357.00

Aspect 100RT Tabletop options

	<p>Standard Tabletop X100-995 24" x 80" (61 x 203cm) Low attenuation, carbon fiber construction. Standard tops may be substituted with optional sizes/styles as shown below at no extra charge.</p>
	<p>Optional Tabletop X100-733 with facial cutout 22" x 80" (56 x 203cm). Low attenuation, carbon fiber construction. Optional tabletop available at no extra charge. Must specify top selection with order.</p>
	<p>Optional Tabletop X100-1325 28" x 80" (71 x 203cm) Wide for larger patients. Low attenuation, carbon fiber construction. Must specify top selection with order.</p>

100RT Shipping Crate Specifications

- Weight (approx.): 571 lbs. (259kg.)
- Length: 96" (2438mm)
- Width: 42" (1067mm)
- Height: 47" (1194mm)

Aspect 100T Fluoroscopic Procedure Table

Catalog # A100-658, A100-2648 & A100-2658

\$15,500.00



Aspect 100T Description:

C-arm table for pain management and other procedures where only motorized elevation and Trendelenburg tilt is required.

- Motorized Elevation
- Trendelburg tilt
- 115VAC is standard.
- 230VAC is optional and must be specified at time of order.

Aspect 100T #A100-658

Includes the following:

- X100-995: Rectangular tabletop 24" (610mm) x 80" (2032mm).
- X100-1742: Tabletop pad, 2" (51mm) thick, for 24" (610mm) (610mm) x 80" (2032mm) top
- K000-0054: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.

Aspect 100T #A100-2648

Includes the following:

- X100-733: Tabletop with tapered end and facial cutout 22" (559mm) x 80" (2032mm).
- X100-1757: Tabletop pad, 2" (51mm) thick, for 22" (559mm) x 80" (2032mm) top with facial cutout.
- K000-0054: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.

Aspect 100T #A100-2658

Includes the following:

- X100-1325: Rectangular tabletop 28" (711mm) x 80" (2032mm).
- X100-2083: Tabletop pad, 2" (51mm) thick, for 28" (711mm) x 80" (2032mm).
- K000-0054: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.




Aspect 100T Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	Clamp-on Accessory Rail Rail is 6" long x 3/8" thick x 1-1/8" high (152 x 10 x 29 mm)	A100-1007	\$495.00
	Side-rail clamp Rotating, with socket for attachment of accessories with 5/8 inch (16mm) circular post. Requires A100-1007 clamp-on accessory rail.	C000-0746	\$395.00
	Arm Board Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059	\$975.00
	Multi-position Arm Support Attaches to IDI tables with accessory rails - or can be used with a clamp on accessory rail on tables or areas within a table that do not have accessory rails. Length: 23.5 inches (597mm).	C000-1180	\$866.00
	Scatter Guard radiation shield draping For 28" wide Optional Tabletop X100-1325 28"x80" (71x203cm) Dimensions: 59" long. 24" draping below tabletop. Draping extends along both sides and head--end of table. Supported by fabric section over tabletop. <i>Not for use with arm board A310-056.</i>	A610-108	\$1595.00
	Hand Control for 100T table. Elevate and tilt functions. (Included with table).	K000-0054	\$230.00

Aspect 100T Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<p>Foot Control for 100T</p>	A100-500	\$995.00
	<p>Radiation Shield Fixed shield with integrated clamp system. For use with tables that do not have side rails. Clamps to table sides or at head-end of table. Snap-on upper section above level of tabletop is removable. Use with 100RTL, 100RT, 100T, 100E tables.</p>	A610-080	\$1295.00
	<p>Scatter-Guard radiation shield draping For for 24" rectangular tabletop X100-995. Dimensions: 59" long. 24" draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop <i>Note: Not for use with arm board A310-056.</i></p>	A610-096	\$1595.00
	<p>Scatter-Guard radiation shield draping for 22" tabletop X100-733 models with 22" tapered-end tabletop with facial cutout. Dimensions: 59" long. 24" draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop. <i>Note: Not for use with arm board A310-056.</i></p>	A610-092	\$1595.00
	<ul style="list-style-type: none"> • Facial Cushion, crescent shape. • Disposable covers for facial cushion, box of 50. 	<p>Cushion: #C000-0597</p> <p>Covers: #C000-0598</p>	<p>Cushion: \$70.00</p> <p>Covers: \$40.00</p>
	<p>IDI Medical Grade Anti-Fatigue Mat (22"/559mm x 66"/168cm) Heavy-duty, one-piece construction for high-traffic areas. Safe & non-toxic, contains no plasticizers and will not leach over time. Durable brushed texture.</p>	C000-1181	\$357.00

Aspect 100T Tabletop options

	<p>Standard Tabletop X100-995 24" x 80" (61 x 203cm) Low attenuation, carbon fiber construction. Standard tops may be substituted with optional sizes/styles as shown below at no extra charge.</p>
	<p>Optional Tabletop X100-733 with facial cutout 22" x 80" (56 x 203cm). Low attenuation, carbon fiber construction. Optional tabletop available at no extra charge. Must specify top selection with order.</p>
	<p>Optional Tabletop X100-1325 28" x 80" (71 x 203cm) Wide for larger patients. Low attenuation, carbon fiber construction. Must specify top selection with order.</p>

100T Shipping Crate Specifications

- Weight (approx.): 540 lbs. (245 kg)
- Length: 96" (2438mm)
- Width: 42" (1067mm).
- Height: 47" (1194mm)

Aspect 100E Fluoroscopic Procedure Table

Catalog # A100-402, A100-2642 & A100-2656

\$12,995.00



Aspect 100E Description:

C-arm table for pain management and other procedures where motorized elevation is the only movement required.

- 115VAC is standard.
- 230VAC is optional and must be specified at time of order.

Aspect 100E #A100-402

Includes the following:

- X100-995: Rectangular tabletop 24" (610mm) x 80" (2032mm).
- X100-1742: Tabletop pad, 2" (51mm) thick, for 24" (610mm) x 80" (2032mm) top
- K000-0041: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.



Aspect 100E #A100-2642

Includes the following:

- X100-733: Tabletop with tapered end and facial cutout 22" (559mm) x 80" (2032mm).
- X100-1757: Tabletop pad, 2" (51mm) thick, for 22" (559mm) x 80" (2032mm) top with facial cutout.
- K000-0041 Pendant hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.



Aspect 100E #A100-2656

Includes the following:

- X100-1325: Rectangular tabletop 28" (711mm) x 80" (2032mm).
- X100-2083: Tabletop pad, 2" (51mm) thick, for 28" (711mm) x 80" (2032mm).
- K000-0041: Hand control.
- C000-0328: Patient restraint straps.
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50.




Aspect 100E Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	Clamp-on Accessory Rail Rail is 6" long x 3/8" thick x 1-1/8" high (152 x 10 x 29 mm)	A100-1007	\$495.00
	Side-rail clamp Rotating, with socket for attachment of accessories with 5/8 inch (16mm) circular post. Requires A100-1007 clamp-on accessory rail.	C000-0746	\$395.00
	Arm Board Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059	\$975.00
	Multi-position Arm Support Attaches to IDI tables with accessory rails - or can be used with a clamp on accessory rail on tables or areas within a table that do not have accessory rails. Length: 23.5 inches (597 mm).	C000-1180	\$866.00
	Scatter Guard radiation shield draping For 28" wide Optional Tabletop X100-1325 28"x80" (71x203cm) Dimensions: 59" long. 24" draping below tabletop. Draping extends along both sides and head--end of table. Supported by fabric section over tabletop. <i>Not for use with arm board A310-056.</i>	A610-108	\$1595.00
	Hand control for 100E. (Included with table).	K000-0041	\$85.00
	Foot Control 100E	A100-758	\$650.00

Aspect 100E Table Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<p>Radiation Shield Fixed shield with integrated clamp system. For use with tables that do not have side rails. Clamps to table sides or at head-end of table. Snap-on upper section above level of tabletop is removable. Use with 100RTL, 100RT, 100T, 100E tables.</p>	A610-080	\$1295.00
	<p>Scatter-Guard radiation shield draping For for 24" rectangular tabletop X100-995. Dimensions: 59" long. 24" draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop <i>Note: Not for use with arm board A310-056.</i></p>	A610-096	\$1595.00
	<p>Scatter-Guard radiation shield draping for 22" tabletop X100-733 models with 22" tapered-end tabletop with facial cutout. Dimensions: 59" long. 24" draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop. <i>Note: Not for use with arm board A310-056.</i></p>	A610-092	\$1595.00
	<ul style="list-style-type: none"> • Facial Cushion, crescent shape. • Disposable covers for facial cushion, box of 50. 	<p>Cushion: #C000-0597</p> <p>Covers: #C000-0598</p>	<p>Cushion: \$70.00</p> <p>Covers: \$40.00</p>
	<p>IDI Medical Grade Anti-Fatigue Mat (22"/559mm x 66"/168cm) Heavy-duty, one-piece construction for high-traffic areas. Safe & non-toxic, contains no plasticizers and will not leach over time. Durable brushed texture.</p>	C000-1181	\$357.00

Aspect 100E Tabletop options

	<p>Standard Tabletop X100-995 24" x 80" (61 x 203cm) Low attenuation, carbon fiber construction. Standard tops may be substituted with optional sizes/styles as shown below at no extra charge.</p>
	<p>Optional Tabletop X100-733 with facial cutout 22" x 80" (56 x 203cm). Low attenuation, carbon fiber construction. Optional tabletop available at no extra charge. Must specify top selection with order.</p>
	<p>Optional Tabletop X100-1325 28" x 80" (71 x 203cm) Wide for larger patients. Low attenuation, carbon fiber construction. Must specify top selection with order.</p>

100E Shipping Crate Specifications

- Weight (approx.): 495 lbs. (225kg)
- Length: 96" (2438mm)
- Width: 42" (1067mm).
- Height: 47" (1194mm).

Aspect 100-UCPLUS Urology/ Multi-purpose Catalog # A100-1517

\$37,250.00



Shown with optional leg holders (#A100-2247).



Shown with tabletop extension (A100-1400)

Aspect 100UCPlus description:

Imaging table for urology and multi-purpose use. Urological accessories are easily removable so the table can be quickly configured for multi-purpose C-arm imaging or general patient positioning for non-urological procedures. A removable radiographic cassette tray provides KUB radiography compatibility with mobile x-ray machines, including both conventional film/grid systems and digital detectors.

- 115VAC is standard.
- 230VAC is optional and must be specified at time of order.

Aspect 100UCPlus #A100-1517





Includes the following:

- A100-1007 6"/152mm x 1.25"/32mm clamp-on accessory rail
- #A100-1876 Stirrups Mounting Rail extension Kit
- X100-1818: tabletop pad, 28" (711mm) x 46" (1168mm) x 2" (51mm) thick
- A100-1400 Tabletop extension, 24" (610mm)W x 34" (864mm)L
- X100-1817: Tabletop extension pad, 24" (610mm) x 34" (864mm) x 2" (51mm) thick
- K000-0158: Pendant hand control
- A100-1540: Foot switch control
- Z100-1060: Removable x-ray shield flap, for perineal end of table
- C000-0571: Patient restraint straps
- A310-062: Arm board, rail mounted, with pad and strap
- C000-0595: Urology drain bag support hoop
- W100-1001: Radiographic cassette tray, removable
- C000-0593: Disposable drain bag & hose ass'y, box of 20
- C000-0492: Disposable covers for foot switch, box of 50
- C000-0631: disposable cover for x-ray shield flap (qty. of 1)
- C000-0645: disposable tabletop pad protector (qty. of 1)





Aspect 100UCPlus Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<p>Lithotomy Legholders (pair) Boot style, pneumatic assist with one-hand adjustment, with pads. For patients up to 350 lbs (158kg). Includes pair of side-rail clamps.</p>	A100-2246	\$5090.00
	<p>Heavy Duty Lithotomy Legholders (pair) Heavy Duty, boot style, pneumatic assist with one-hand adjustment, with pads. For patients up to 525 lbs (238kg). Includes pair of side-rail clamps.</p>	A100-2247	\$7155.00
	<p>Physician Elbow Supports (pair) Removable</p>	C000-0594	\$1195.00
	<p>Armboard Rail mounted, with pad and strap 1 arm board included with table.</p>	A310-056	\$450.00
	<p>Split Leg Holders For use with Aspect 100UC Plus Tables. Includes pads. Prone percutaneous, lateral or modified lateral nephrostomy. Includes two split leg holder clamps # C000-1115.</p>	C000-1146	\$6400.00
	<p>Replacement Split Leg Holder Clamp For use with split leg holders. # C000-1146</p>	C000-1115	\$125.00
	<p>Stirrups Mounting Rail extension Kit Included with table. For use with leg-holders (Leg-holders are sold separately).</p>	#A100-1876	\$525.00




Aspect 100UCPlus Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	Multi-position Arm Support Attaches to IDI tables with accessory rails - or can be used with a clamp on accessory rail on tables or areas within a table that do not have accessory rails. Length: 23.5 inches (597mm).	C000-1180	\$866.00
	Hand Control for 100UCPlus and 100UC. Controls height up/down, trendelenburg/ revers, slide, cross table movements. (Included with table).	K000-0158	\$185.00
	Foot Switch Control for 100UCPlus. (Included with table).	A100-1540	\$1855.00
	Side-rail clamp Rotating, with socket for attachment of accessories with 5/8 inch. (16mm) circular post. Requires A100-1007 clamp- on accessory rail.	C000-0746	\$395.00

Aspect 100UCPlus Consumables

Photo	Description	Catalog Number	
	Vinyl Urology Table Covers Box of 20	C000-0645	\$350.00
	Uro Drain Bags with Hose Disposable drain bag and hose assembly. Box of 20	C000-0593	\$265.00
	Disposable Covers for Aspect 100UCPlus Foot Controls 24" x 20" Box of 50	C000-0492	\$50.00
	Uro Drain Support Hoop Collapsible & removable drain bag support hoop.	C000-0595	\$350.00

Aspect 100UCPlus Consumables

Photo	Description	Catalog Number	
	Uro Collection Unit 5 gallon/18.9 liter capacity, Box of 10. Soft sided Vinyl container gently expands into a cube as it collects waste fluid when floor drain or suction is unavailable. Internal baffle system enables container to fill evenly as powder converts liquid to gel. Built in zipper and handles allow for storage, closure of drain bag and hose assembly for enhanced infection control during transport and disposal.	C000-0612	\$275.00
 5.3 gal/20 liter	Uro Collection Unit 2.6 gal/10 liter capacity, Box of 10. Flexible plastic sided container expands into a cube as it collects waste fluid when floor drain or suction is unavailable. Internal baffle system enables container to fill evenly as powder converts liquid to gel. Integrated valve system allows user to shut off the flow of liquid from the drain bag in order to change collection units..	C000-0905	\$199.00
 2.6 gal/10 liter	Uro Collection Unit 5.3 gal/20 liter capacity, Box of 10. Flexible plastic sided container expands into a cube as it collects waste fluid when floor drain or suction is unavailable. Internal baffle system enables container to fill evenly as powder converts liquid to gel. valve Flow valve system allows user to shut off the flow of liquid from the drain bag in order to change collection units.	C000-0916	\$289.00

Aspect 100UCPlus Shipping Crate Specifications

- Weight (approx.): 815 lbs. (370kg)
- Length: 96" (2438mm)
- Width: 42" (1067mm).
- Height: 47" (1194mm).

Aspect 100-US Gyn/Fertility Clinic Table Catalog # A100-1479

\$14,950.00



Aspect 100US description:

GYN/Fertility Mobile Imaging Table. An economical C-arm table for hysterosalpingograms




- 115VAC is standard.
- 230VAC is optional and must be specified at time of order.

Aspect 100US #A100-1479




Includes the following:

- A100-758 Foot control
- X100-1687: tabletop pad, 24" (610mm) x 46" (1168mm) x 1" (25 mm) thick.
- A100-1007 6" (152mm) x 1.25" (32mm) clamp-on accessory rail.
- C000-0328: Patient restraint straps.
- A100-2380: Leg support tabletop extension with pad.
- X100-2351: Leg support tabletop extension pad.
- A100-2248: Basic Knee Crutches [pair] with a pair of rotating side rail clamps [C000-0746] for patients weighing up to 400 lbs (181kg)

Aspect 100US Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	Armboard Quick-release, rail mounted. Requires Clamp-on Accessory Rail A100-1007	A310-056	\$450.00
	Clamp-on Accessory Rail Rail is 6" long x 3/8" thick x 1-1/8" high.	A100-1007	\$495.00
	Hand control for 100US	K000-0041	\$85.00

Aspect 100US Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	Foot Control 100US (Included with table.)	A100-758	\$650.00
	Side-rail clamp Rotating, with socket for attachment of accessories with 5/8 inch (16mm) circular post. Requires A100-1007 clamp- on accessory rail.	C000-0746	\$395.00
	Multi-position Arm Support Attaches to IDI tables with accessory rails - or can be used with a clamp on accessory rail on tables or areas within a table that do not have accessory rails. Length: 23.5 inches (597mm).	C000-1180	\$866.00

Aspect 100US Shipping Crate Specifications

Weight (approx.): 245 kg / 540 lbs.
 Length: 244 cm / 96 in.
 Width: 107 cm. / 42"
 Height: 119 cm / 47 in.



MDS Description:

MDS Mobile Display System, Complete with High Definition Video Monitors, Router and tablet controller.

- 115VAC is standard.
- 230VAC is optional and must be specified at time of order.

MDS #A280-0330

Includes the following:

- Dual-Arm Mobile Display Stand with adjustable counterbalanced arms
- Adjustable shelves for placement of ancillary equipment (Qty. of 4)
- Bilateral input source cable manifold supporting DVI, SDI, VGA, S-Video, SOG and composite inputs
- Built-in isolation transformer and 8 outlet power strip for video and ancillary equipment AC power
- FSN 26" (660mm) HD LED Widescreen Multi-Modality Video Monitors (Qty of 2).
- FSN IPS1000 Video Router/Controller with 10.1" Tablet Controller and articulated support arm with lock
- Dual 8" (20cm) "Preview" LCD Reference Monitors, mounted to MDS stand center column
- Qty (1) External cable storage hook
- Qty (1) Kit - Accessory Basket Kit
- Protective covers with storage for 26" (660mm) FSN monitors (Qty of 2)
- Video Cable kit, including common cables/connectors needed for connection of C-arm and equipment that will be mounted onto the MDS Mobile Display System. Customer choice of one of the following cable kits is included with the MDS System.

Customers must specify cable kit choice with order:

- **A280-0325:** Cable kit for GE OEC 9800/9900 C-Arms and camera cable kit for systems mounted onto the MDS
- **A280-0328:** Cable kit for all DVI output C-Arms and camera cable kit for systems mounted onto the MDS
- **A280-0326:** Cable kit for Siemens VGA C-Arm with camera cable kit for systems mounted onto the MDS
- **A280-0494:** Cable kit for Siemens Cios C-Arm with Video Manager and camera cable kit for systems mounted on the MDS

MDS Optional and Replacement Accessories:

MDS Optional Cable Kits








Not Included but Recommended

External Source Cable Kit

Add this Kit if the customer plans to add any external device that is not mounted directly onto the MDS stand other than the C Arm (i.e. ultrasound, PACS, hemodynamics) Kit Includes: qty (1) each of 20 ft. (6.1m) S Video, (2) 15 ft. (4.6m) lockable HDMI, 20 ft. (6.1m) SVG cable and qty. (2) lockable DVI/HDMI adapters.

A282-0003

\$470.00

Photo	Description	Catalog Number	MSRP
	Endo Shelf Bilateral slide-out shelf, increases storage and shelf space For MDS or MDS D.	A280-0415	\$1441.00
	CO₂ gas brackets Pair of holders for type E gas cylinders (for insufflator use) mounted to sides of MDS stand.	A280-0401	\$475.00
	Shelf with Nook Kit includes Kit- Shelf Nook Kit includes: medical rail shelf with under mounted nook that will accommodate a DAU switcher, freeing up shelf space. Customer installation only. Nook size: 4" H x 19" W.	A282-0008	\$1257.00
	Medical Rail Shelf Shelf Includes medical rails. For accessory attachments.	A280-0130	\$850.00
	Shelf with Keyboard Kit includes medical rail shelf with sliding keyboard carriage. Suitable for most standard customer-supplied keyboards.	A282-0005	\$1461.00
	Accessory Basket Kit For customer-supplied foot switch or other storage.	A282-0007	\$674.00
	Suction Canister Kit Includes: 2 suction canister rail clamps and brackets suitable for 3000 ml canister.	A282-0009	\$357.00

MDS Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	<p>Cable Hook Cable storage J-hook for external cables.</p>	A280-0570	\$265.00
	<p>IV Pole Kit Includes two IV rail clamps suitable for mounting a 4' 2-hook IV pole onto two medical rail equipped (A280-0130) shelves. Note: this option requires two shelves with medical rails.</p>	A280-0508	\$633.00
	<p>MDS Monitor Protection Covers F or 26" (660mm) flat panel monitor. (2 covers are included with MDS System A280-0330)</p>	X280-0262	\$169.00
	<p>MDS Articulated Arm Kit includes multi-joint arm (18 lb. max load), assembly including qty (1) rail clamp, and articulated arm with VESA 100mm x 100mm mount</p>	A282-0016	\$1100.00
	<p>MDS 3rd Monitor with articulated arm Includes 3rd monitor (19" FSN), articulated arm and cables.</p>	A282-0006	\$4320.00

MDS Optional and Replacement Accessories:

Description	Catalog Number	MSRP
<p>MDS Installation options</p> <p>System installation by IDI, including two-day operator in-service training, is available for the USA and Canada. Dealers who choose to have IDI perform the technical installation only will receive a discount provided they supply two days of certified clinical training. Additional in-service training of more than two days is available on a billable basis only. Contact IDI regarding fees for additional training. Installation costs are not discountable and installation is only available in the USA and Canada.</p>		
<p>MDS Installation, technical integration and clinical training by IDI in the continental United States</p>	<p>Y-MDSINS-01</p>	<p>\$5000.00</p>
<p>MDS installation and integration only by IDI technician in Continental United States OEM agrees to supply 2-day certified clinical applications training to customer.</p>	<p>Y-MDSINS-04</p>	<p>\$2500.00</p>
<p>MDS Installation, integration and clinical training by IDI in Alaska or Hawaii 4-day installation, integration and clinical operation training by IDI.</p>	<p>Y-MDSINS-02</p>	<p>\$6000.00</p>
<p>MDS Installation, integration and clinical training by IDI in Canada 4-day installation, integration and clinical operation training by IDI.</p>	<p>Y-MDSINS-03</p>	<p>\$6000.00</p>

MDS D Video Integration System Catalog # A280-230

\$54,900.00



MDS D Description

The Image Diagnostics MDS D provides a lower cost solution for flexible monitor positioning of 26" HD medical grade monitors. Images are manually PIP'ed at the monitor and changed at the IPS-500 switching device. The MDS D is recommended for those facilities where image routing requirements remain fairly consistent from day to day.

- 115VAC is standard.
- 230VAC is optional and must be specified at time of order.

MDS D #A280-230

Includes the following:

- Dual counter balanced arms
- Two 26" (660mm) high definition monitors
- Video switcher/scalar for switching source images
- (1) 19.75" (502mm) x 19.75" (502mm) equipment shelf for placement of ancillary equipment. Includes under shelf "nook" for the video switcher
- Dual articulating arms with positioning handles and low-force manual motion for easy horizontal or vertical positioning of monitors
- Self-leveling monitor arms provide 36"/914mm of fully counter-balanced vertical travel
Arm radius extends horizontally to 74"
- Stand height (with arms fully retracted): 64.5"/164 cm
- Extended arm vertical height: 96"/244 cm
- Lockable double-wheeled 6"/152mm casters with cable pushers
- Front and rear mounted ergonomic transport and positioning handles
- Front cross mounted medical rail for accessory attachment
- Qty (1) accessory basket
- Qty (1) cable storage hook
- Qty (2) Padded, Protective Covers for monitors with storage compartment for 26"/660mm
- Built-in isolation transformer and 8-outlet AC power strip for connection of video components and ancillary equipment

MDS D Optional and Replacement Accessories:




MDS D Cable Kits – [Included but Selection is Needed](#)

For connection of C-arm and Endoscopy Systems. Customer choice of one of the following cable kits is included with MDS D System A280-0330 at no extra charge. Customer must specify the one selected cable kit with MDS System order. Prices shown below apply only to orders for replacement or additional cable kits.

MDS Cable kit for all high resolution SOG or composite output C-Arms Video connection kit for all SOG or composite out C-Arms including GE OEC 9800, 9900 C-arms.	A280-0491	\$261.00
MDS D Cable kit for all DVI Output C-arms Suitable for GE OEC Elite, Ziehm RFD, Siemens Cios, Phillips Veradius C-Arms.	A280-0492	\$1650.00
MDS D Cable kit for VGA C-arm including Siemens Arcadis Video connection kit.	A280-0493	\$2098.00
MDS D Cable kit for Siemens Cios with Video Manager C-arm Video connection kit.	A280-0494	\$2878.00

MDS D Optional Cable Kits


Not Included but Recommended

External Source Cable Kit Add this Kit if the customer plans to add any external device that is not mounted directly onto the MDS stand other than the C Arm (i.e. ultrasound, PACS, hemodynamics) Kit Includes: qty (1) each of 20 ft. (6.1m) S Video, (2) 15 ft. (4.6m) lockable HDMI, 20 ft. (6.1m) SVG cable and qty. (2) lockable DVI/HDMI adapters.	A282-0003	\$470.00	
Photo	Description	Catalog Number	MSRP
	Endo Shelf Includes medical rail shelf with two bilateral slide-out shelves for MDS or MDS D. Factory installation only.	A280-0415	\$1441.00
	CO₂ Gas Brackets Pair of holders for type E gas cylinders (for insufflator use) mounted to sides of MDS stand (customer installation).	A280-0401	\$475.00
	Shelf with Nook Kit includes Kit- Shelf Nook Kit includes: medical rail shelf with under mounted nook that will accommodate a DAU switcher, freeing up shelf space. Customer installation only . Nook size: 4" H x 19" W.	A282-0008	\$1257.00

MDS D Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	Medical Rail Shelf Shelf Includes medical rails. Add if equipment needs may overload the standard single shelf configuration.	A280-0130	\$850.00
	Suction Canister Kit Includes: 2 suction canister rail clamps and brackets suitable for 3000 ml canister.	A282-0009	\$357.00
	Accessory Basket Kit Wire basket suitable for accessory storage.	A282-0007	\$674.00
	Shelf with Keyboard Kit includes: (A280-0130) medical rail shelf with sliding keyboard carriage for customer-supplied keyboard.	A282-0005	\$1461.00
	Cable Hook Cable storage J-hook for external cables.	A280-0570	\$265.00
	IV Pole Kit Includes two IV rail clamps suitable for mounting a 4' (1.2 m) 2-hook IV pole onto two medical rail equipped (A280-0130) shelves. Note: this option requires two shelves with medical rails.	A280-0508	\$633.00
	MDS D Monitor Protection Covers For 26" (660mm) flat panel monitors. Two covers are included with MDS D system A280-0230.	X280-0262	\$169.00
	MDS D Articulated Arm Kit includes multi-joint arm (18 lb. max load), assembly including qty (1) rail clamp, and articulated arm with VESA 100mm x 100mm mount	A282-0016	\$1100.00

MDS D Optional and Replacement Accessories:

Photo	Description	Catalog Number	MSRP
	MDS 3rd Monitor with articulated arm Includes 3rd monitor (19" (483 mm) FSN), articulated arm and cables.	A282-0006	\$4320.00

MDS D Installation options

System installation by IDI, including two-day operator in-service training, is available for the USA and Canada. Dealers who choose to have IDI perform the technical installation only will receive a discount provided they supply two days of certified clinical training. Additional in-service training of more than two days is available on a billable basis only. Contact IDI regarding fees for additional training. Installation costs are not discountable and installation is only available in the USA and Canada.

MDS D Installation, technical integration and clinical training by IDI in the continental United States 3-day installation, integration and clinical operation training by IDI.	Y-MDSINS-01	\$5000.00
MDS D installation only by IDI technician in Continental United States OEM agrees to supply 2-day certified clinical applications training to customer.	Y-MDSINS-04	\$2500.00
MDS D Installation and clinical training by IDI in Alaska or Hawaii One day of installation, 1½ days of clinical operation training by IDI.	Y-MDSINS-02	\$6000.00
MDS D Installation, integration and clinical operation training by IDI. One day installation 1½ days of and clinical operation training by IDI.	Y-MDSINS-03	\$6000.00

MDS D Shipping Crate Specifications

- Weight (approx.): 825 lbs.(374kg)
- Length: 65 in.(1651mm)
- Width: 42 in.(1067mm).
- Height: 72 in.(1829mm)

ilex 32 video monitor and stand
Catalog # A280-0360

\$9800.00



ilex 32 Description

The Ilex 32 monitor and stand provides a cost effective solution for those customers requiring an additional monitor or monitors with mirrored images in multiple locations.

- 115VAC is standard.
- 230VAC is optional and must be specified at time of order.

ilex 32 #A280-0360

Includes the following:

- 32"/81cm HD medical grade monitor with medical bios and dual inputs
- Single or split-screen images
- Heavy duty cart stand with tilt (Maximum height: 76"/193cm to top of monitor)
- Cable channel
- Accessory basket
- 15 foot (4.6 m) power cord
- Cable hook
-

ilex 32 Installation options

The customer typically performs installation. Installation is not available outside of the United States.
 For Dealers who choose to have IDI perform the installation:

<p>ilex 32 Installation by IDI in the continental US Includes installation and 1 day of clinical applications training by IDI certified applications person</p>	<p>Y-ILEXINSO1</p>	<p>\$1800.00</p>
<p>ilex 32 Installation by IDI in Alaska or Hawaii Includes installation by IDI in Canada, Alaska, or Hawaii: includes one-day assembly and ½ day operator training</p>	<p>Y-ILEXINSO2</p>	<p>\$2500.00</p>

Aspect Four-Way Float Table Specifications & Accessories

•=Standard feature ◊=Optional NA=Not available	100 - 4	100 - 4T G2	ISR G2
Elevate (Height with 2" pad)	29.9 - 39.7" 75.9 - 100.8 cm	32 - 42" 81 - 107 cm	32 - 40" 81 - 102 cm
Trendelenburg Tilt	NA	±12°	±12°
Lateral Tilt	NA	NA	±12°
Longitudinal Tabletop Motion	40" float 102 cm	32" float 81 cm	32" float 81 cm
Transverse Tabletop Motion	10" float 25 cm	8" float 20 cm	8" float 20 cm
Standard Tabletop Size	18/30" x 80" 46/76 x 203 cm	24"/22"/19.5" x 84" 61/56/50 x 213 cm	24"/22"/19.5" x 84" 61/56/50 x 213 cm
Optional Tabletop Size	23.5-30" x 80" 60-76 x 203 cm	22"/24" x 84" 56/61 x 213 cm	22"/24" x 84" 56/61 x 213 cm
Imaging Area	72" 183 cm	61" 155 cm	61" (155cm), 74" (188 cm) with optional peripheral extension
Patient Capacity	450 lb./204 kg	450 lb./204 kg	450 lb./204 kg
Table Weight	490 lb. 222 kg	605 lb. 274 kg	605 lb. 274 kg
Tabletop Pad	•	•	•
Restraint Straps	•	•	•
Pendant Hand Control	NA	NA	•
4-way Panning Hand Control	•	•	•
Foot Control	NA	Custom order)	Custom order
Armboard, rail mounted	NA	◊	◊
Armboard, carbon fiber	◊	◊	◊
VAB Armboard	◊ Use of the VAB arm board with the 100-4 is only recommended on the models with 23.5" table tops . Use is only possible with a patient's weight securing the VAB to the table.	◊	◊
Clamp-on Accessory Rail	◊	◊	◊
Catheter Tray extension	◊	◊	◊
Headrest extension	NA	◊	◊
Table-Mounted Radiation Shield	◊	◊	◊
Anesthesia Screen Holder	NA	◊	◊

Aspect Urology & GYN Table Specifications & Accessories

•Standard feature ◊Optional NA Not available	100 UC-Plus	100 US
Elevate (Height with pad)	33.3" – 43.1" (85 – 109 cm)	27" – 39" (69 – 99cm)
Trendelenburg Tilt	± 15°	NA
Transverse travel	± 3.5", 7" total (± 9 cm, 18 cm total)	NA
Longitudinal Tabletop Motion	10" (25 cm)	NA
Tabletop Size	28"W x 46"L (71 x 117 cm)	24"W x 46"L (61 x 117 cm)
Tabletop Extension Size	24"W x 34"L (61 x 86 cm)	NA
Imaging Area	Main tabletop: 25"W x 20"L (64 x 51 cm)	22"W x 32"L (56 x 81 cm)
Imaging Area with table extension	54" (137 cm)	NA
Patient Weight Capacity	Upper body on main table: 525lb. (238 kg) Upper body on extension: 400 lbs (181kg) Load capacity on extension only: 200 lbs (91kg)	450 lbs (204 kg)
Table Weight	490 lb. (222 kg)	375 lbs ((170kg)
Tabletop Pad (2" thick)	•	•
Tabletop extension with pad	•	Leg support tabletop extension
Hand control pendant & footswitch control	•	Footswitch only
Patient restraint straps	•	•
Removeable drain bag support hoop	•	NA
Disposable drain bag & hose assembly, box of 20	•	NA
Removable radiation shield flap for perineal end of table	•	NA
Armboard, rail-mounted with pad and strap	•	◊
Split Leg Stirrups	•	NA
Disposable covers for footswitch, box of 50	•	◊
Radiographic cassede tray	•	NA

Aspect Pain Management Table Specifications & Accessories

•Standard feature ◊Optional NA Not available	100E	100T	100RT	100RTL
Elevate (Height with 2" pad)	28 - 39.8" 71 - 101 cm	30 - 41.8" 76 - 106 cm	32.4 - 44.2" 82 - 112 cm	31.2 - 41.2" 79 - 105 cm
Trendelenburg Tilt	NA	±15°	±15°	±15°
Lateral Tilt	NA	NA	20°	20°
Longitudinal Tabletop Motion	NA	NA	NA	10"
Standard Tabletop Size	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm	24" x 80" 61 - 203 cm
Optional Tabletop Size	22" x 80" 56 x 203 cm with facial cut-out	22" x 80" 56 x 203 cm with facial cut-out	22" x 80" 56 x 203 cm with facial cut-out	22" x 80" 56 - 203 cm
Optional Tabletop Size	28" x 80" 71 x 203 cm	28" x 80" 71 x 203 cm	28" x 80" 71 x 203 cm	28" X 80" 71 x 203 cm
Imaging Area	66" 168 cm	62" 157 cm	62" 157 cm	59" 150 cm
Patient Capacity	450 lb. 204 kg	450 lb. 204 kg	450 lb. 204 kg	450 lbs. 204 kg
Table Weight	290 lb. 132 kg	325 lb. 147 kg	425 lb. 193 kg	355 lbs. 161 kg
Tabletop Pad (2" thick)	•	•	•	•
Restraint Straps	•	•	•	•
Facial Support Cushion	•	•	•	•
Pendant Hand Control	•	•	•	•
Foot Control	◊	◊	◊	◊
Armboard, rail mounted	◊	◊	◊	◊
Armboard, carbon fiber	◊	◊	◊	◊
VAB Armboard	◊	◊	◊	◊
Clamp-on Accessory Rail	◊	◊	◊	◊
Scatter-Guard Radiation Shield	◊	◊	◊	◊

Image Diagnostics, Inc.

Since 1986, IDI has been a manufacturer of innovative products for surgical, interventional and diagnostic imaging.

Our product focus is in electro-mechanical solutions for fluoroscopic and digital radiographic applications.

IDI is committed to superior product design, value and customer support in everything we do.

Image Diagnostics, Inc.

310 Authority Drive
Fitchburg, Massachusetts
01420

Toll free: 877.304.5434

Main: 978.829.0009

Fax: 978.829.0027

ImageDiagnostics.com

Please direct all inquiries to:
sales@imagediagnostics.com

idi

Image Diagnostics Inc.



A Veteran-Owned Business