

MASS™ Medical Secure-A-Scope™ Storage systems

MASS™ Medical Storage Cabinets not only organize your scopes, they protect them as well. Specially designed for the length of endoscopes and colonoscopes, cushioned scope storage racks and padded back walls help eliminate costly repairs. Cabinets store 4, 8, 9, or 18 flexible Endoscopes.

Two filtered ambient air vents keep air moving in and out of the scope cabinet, even when the door is closed. For cleanliness and convenience, a removable drip tray comes standard in our scope storage cabinets.

Door Options - FlexLok™ tambour door, solid hinged doors or hinged doors with glass inserts.



Pass-through scope cabinet




Cabinet for 8 flexible Endoscopes



Cabinet for 9 Endoscopes



Cabinet especially designed for small-diameter scopes
(bronchoscopes, cytosopes, ureteroscopes, TEE probes)



Rotating scope holders make loading and removing scopes easy, safe and secure.



Autoclavable scope tray with liner



Scope Transport cart that is steam or machine washable

JD Honigberg International, Inc

650 Dundee Road, Northbrook, IL 60062 USA

Tel: (1) 847 412 0200 / Fax: (1) 847 412 0034 / E-Mail: medical@jdhintl.com / www.jdhmedical.com