

MASS

MEDICAL STORAGE



MASS PUTS YOU IN CONTROL

Contents

MASS™ Philosophy	2,3
MASS Metal™	
Exchange Carts	4, 5
Wall Storage Cabinets.....	6
Exchange Components	7
Exchange Cart Options	8
Carts for MRI	9
Endoscopy	
Secure-A-Scope™ Storage System....	10-12
Procedure Carts	13
Scope Transport Carts	13
Cardiac Cath Lab - EP Lab	
Wall Storage Cabinets.....	14
Wall Storage Cabinet Options.....	15
Mobile Cath Lab Storage.....	16
Procedure Carts	17
Surgery – Ambulatory Surgery	
MASS Metal™ Exchange Carts.....	18
Surgical Supplies "To Go".....	19
Additional MASS™ Models	20
Surgical Supply Wall Cabinets	21
Radiology – Angiography	
Interventional Supply Storage Cabinets....	22
Wall Storage Cabinets.....	23
Mobile Interventional Supplies.....	24
Radiology Supply Carts	25
Metal Cabinets and Carts.....	26
Radiology – CT	
CT Supply Storage	27
Design Worksheet	
MASS™ Small Series Carts	28, 29
Standard Cabinets.....	30
Metal Exchange Cabinets.....	31
Metal Exchange Carts	32
Typical Rooms.....	33
Accessories/Options	34, 35

Specifications are subject to change without notice.

MASS PUTS YOU IN CONTROL

SERVICES:

Storage Needs Analysis
Free Engineering CAD Design Service
Mobile and Architectural Solutions

MATERIALS:

Medical Grade Laminates
14 Gauge Powder Coated Metal
Stainless Steel
ABS Plastic

FEATURES:

Customizable Shelves, Bins,
Baskets, Drawers, Dividers
Pull Out and Rotating
Catheter Hooks (FIFO Glides)
Pull Out Suture Modules
Secure-A-Scope™ Easy Access
Endoscopy Holders
Custom Small Scope Holders
LEDs, UV Lighting, Timed Fans
Roll Top or Hinged Glass Doors
and much more.....

STATE OF THE ART INVENTORY MANAGEMENT:

Electronic Security
Tracking Software

SCOPE TRACKER™

INFECTION CONTROL COMPLIANCE:

SATISFY JOINT COMMISSION AND
MANUFACTURER STORAGE SPECS



ELECTROPHYSIOLOGY HYBRID OR CART

PUTS YOU IN CONTROL WITH STAINLESS STEEL

MASS Medical Storage now offers customized hospital-grade Stainless Steel products, favored by Infection Control and OR directors for the ease of cleaning and anti-rust properties of stainless steel materials. By customizing with MASS for large OR and Sterile Core builds, you can save money and precisely meet clinical storage and infection control standards without sacrificing the specialized componentry of other MASS Medical Storage products.

STAINLESS OR PORTABLE WORKSTATION



FIFO GLIDE™ CATHETER MANAGEMENT



EXTRA WIDE STAINLESS MOBILE STORAGE CARTS



With MASS Medical Storage, you will get the highest performance and strongest cabinets built today with our free engineering and architectural services to fit your space.

CALL
800.593.1900
www.massmedicalstorage.com

MASS METAL™ EXCHANGE CARTS

- Durable and light weight – Easy to move and use
- Premium 5" casters, front two locking for safety
- Interchangeable shelves, trays, and baskets to adapt to future needs
- High-impact rigid pvc wrap-around cart guard
- Four cart sizes to meet every need
- Heavy-duty, powder-coated 14 and 20 gauge steel construction
- No wasted space!

Design for
Today's
Storage Needs
Reconfigure for
Tomorrow's
Demands

MX6101DK
Full-Size Metal
Exchange Cart
with FlexLok™ Door and
optional Keyless Lock.
26"D x 45.12"W x 78.75"H





MX6101DK
Full Size Exchange Cart
with FlexLok™ Door and optional Keyless Lock.
26"D x 45.12"W x 78.75"H



MX6107
Half Size Exchange Cart
with FlexLok™ Door.
31"D x 28"W x 78.75"H



MX6102DK
Mid Size Exchange Cart
with FlexLok™ Door, with Stainless Steel Top
and optional Keyless Lock.
26"D x 45.12"W x 44"H



MX6108
Standard Size Exchange Cart
with FlexLok™ Door, with Stainless Steel
Top and optional Red Bumper.
26"D x 25.62"W x 44"H

MASS METAL™ **EXCHANGE CARTS**

Four cart sizes... Hundreds of options

- Mix and match interior components, trays, baskets, shelves and more!
- Easy to reconfigure for new packaging or procedures
- ABS Trays and Wire Baskets available in 3", 6", and 9" heights
- Components slide in and out easily on ABS side panels
- Easy to maneuver with lightweight construction and premium casters

**Visit Our
Web Site
MASSMedical
Storage.com**



MASS METAL™
Wall Storage Cabinets

MX4400D
Metal Exchange
Wall Storage
Cabinet
with FlexLok™ Door.
22"D x 41"W x
89.37"H



FIFOGLIDE™ CATHETER TRACKS
WITH ROTATING HOOKS



SLIDE-OUT SHELF



M4601
Upper Metal Wall
Cabinets
16"D x 41"W x 36"H



MX4401D
Lower Metal Exchange
Wall Cabinets
22"D x 41"W x 35.5"H



MX4402
Metal Exchange Half Cabinet
22"D x 21.62"W x 84.37"H



SLIDE-OUT SLOTSHELF™ WITH ADJUSTABLE DIVIDERS



ABS TRAYS WITH OPTIONAL DIVIDERS



WIRE BASKETS AVAILABLE IN THREE DIFFERENT SIZES



FIFOGLIDE™ CATHETER STORAGE MODULE




SLIDE-OUT SHELF



BINS AVAILABLE IN THREE CONVENIENT SIZES - 4", 6", 8" WIDTHS



SLIDE-OUT SUTURE MODULE



ELECTRONIC KEYLESS LOCK ADDS SECURITY AND CONVENIENCE

EXCHANGE COMPONENTS

- Select colored SlotShelf™ dividers and bins to coordinate interiors
- MASS™ modules make a difference! Baskets, trays, and shelves are 18.5" wide
- Interchangeable shelves, trays, and baskets - easily moved from carts to cabinets!

Contact Us Today for a QuikQuote™

800.593.1900



SIX INCH LOCKING METAL DRAWER



ID LABEL FOR BASKETS

EXCHANGE CART OPTIONS

Color code to improve organization and efficiency

- Color code carts for specific procedures
- Choose cart guard and door strip accent colors
- Select colored SlotShelf™ dividers and bins to coordinate interiors
- MASS™ modules make a difference! Baskets, trays, and shelves are 18.5" wide and easily moved from carts to cabinets!



SLIDE-OUT BIN MODULE



SLIDE-OUT SHELF



MX6108 RED



MX6108 BLUE



MX6108 BLACK



MX6108 GREEN Standard Size
Exchange Cart with FlexLok™ Door.
26"D x 25.62"W x 44"H



MX6108 GREEN
Pictured with closed FlexLok™ Door
and colored door strip.



**FRONT LOCKING PREMIUM
CASTER WHEELS**



SLIDE-OUT TRAYS WITH DIVIDERS



M3101HGMR
Full Size Modular Storage Cart
Hinged Doors with Glass Inserts.
26"D x 45.12"W x 79.25"H



M3101HGMR
Full Size Modular Storage Cart
Hinged Doors with Glass Inserts.
26"D x 45.12"W x 79.25"H

- Mix and match interior components to customize carts for specific supplies
- All exposed metal components and hardware non-ferrous
- Coil Cart
- Supply Cart
- Medication Cart

**Contact Us
Today for a
QuikQuote™
800.593.1900**

**If you don't
see what your
looking for,
call us!**



FIFOGLIDE™ CATHETER TRACKS.
Rotating stainless steel hooks.



**MRI CONDITIONAL CASTERS WITH
STAINLESS STEEL BEARINGS**




STAINLESS STEEL HINGES

Secure-A-Scope™ Storage System

- Revolutionary new design sets the gold standard for scope storage and accessibility
- Rotating scope holders make loading and removing scopes easy, safe, and secure
- Scope holders engineered to protect scope heads and support umbilicals
- Back wall of cabinet padded vinyl to prevent damage to distal tips
- Removable plastic drip tray
- Filtered ambient air vents
- Specialized holders designed for EUS, Enteroscopy, TEE

Ask to see a demonstration in your hospital today.

Call MASS™ or find your local rep on our website.



M1200LHG SECURE-A-SCOPE™ STORAGE SYSTEM

Locking Hinged Door with Glass Inserts.
Holds 9 Flexible Endoscopes.
13.62"D x 40.87"W x 88.62"H




M1800 SECURE-A-SCOPE™ STORAGE SYSTEM

with FlexLok™ Door.
Holds 18 Flexible Endoscopes.
18"D x 40.87"W x 95.25"H

Secure-A-Scope™ Exclusively From MASS™
Unique, Rotating Scope Holders
The Storage System To Protect All Your Scopes

- Cabinets to store 4, 8, 9, or 18 flexible Endoscopes
- Door Options - all models
 - FlexLok™ Door
 - Hinged Doors with Glass Inserts
 - Solid Hinged Doors



M1200
SECURE-A-SCOPE™
STORAGE SYSTEM
with FlexLok™ Door. Holds 9 Flexible Endoscopes.
12.87"D x 40.87"W x 94.12"H



M1202
SECURE-A-SCOPE™
STORAGE SYSTEM
with FlexLok™ Door. Holds 4 flexible Endoscopes.
12.87"D x 21.62"W x 94.12"H



M1802
SECURE-A-SCOPE™
STORAGE SYSTEM
with FlexLok™ Door. Holds 8 Flexible Endoscopes. 18"D x 21.62"W x 95.25"H



M2000LHG
Pass-through SECURE-A-SCOPE™
Scope Cabinet
with Locking Hinged Glass Doors. Holds 4 scopes on each side wall for capacity of 8 scopes. 23.12"D x 40.87"W x 93"H

- Cabinets to store 4/8 or 8/16 small diameter scopes such as Bronchoscopes, Cystoscopes, or Ureteroscopes



M1208
SECURE-A-SCOPE™
STORAGE SYSTEM
with FlexLok™ Door. Holds 9 small diameter scopes. 12.87"D x 40.87"W x 77"H



M1204
SECURE-A-SCOPE™
STORAGE SYSTEM
with FlexLok™ Door. Holds 4 small diameter scopes. 12.87"D x 21.62"W x 77"H



- Vertically and independently adjustable upper and lower scope supports
- Holds 4-8 small diameter scopes
- Hinged, 3-point locking, clear hinged door
- Wall mounted for loading and removing scopes easily, safely and securely
- Silicone inserts gently cushion and support any scopes. The inserts are easily, quickly cleaned or replaced
- Available for immediate shipment

ENDOSCOPY PROCEDURE CARTS



M3103
Full Size ERCP Cart
with FlexLok™ Door on upper
and lower sections.
26"D x 45.12"W x 71"H



M3104
Mid Size ERCP Cart
with FlexLok™ Door.
26"D x 45.12"W x 42"H



M3102ED
Dilator Storage Cart
with FlexLok™ Door.
26"D x 45.12"W x 42"H



M1201HG
**Small Diameter Wall
Mounted SECURE-A-SCOPE™
Scope Cabinet**
with hinged, 3-point locking, clear hinged
door. Holds 4-8 small scopes.
11.25"D x 19.125"W x 48.125"H, 51 lbs

SCOPE TRANSPORT CARTS

- Carts available with four or nine autoclavable trays (lids included)
- Heavy-duty handles for easier mobility

**This Cart is Steam or
Machine Washable!**



M3124
HDPE Scope Transport Cart
22.62"D x 27"W x 44"H



M6204SS
**Stainless Steel Scope
Transport Cart**
21"D x 24.4"W x 35.75"H



M3102
Mid-Size Modular Storage Cart
with FlexLok™ Door.
26"D x 45.12"W x 42"H



MX3106DK
Mid-Size Exchange Cart
with FlexLok™ Door, Keyless Lock, and Stainless
Steel Top and Rail. 26"D x 45.12"W x 42"H

PROCEDURE CARTS

- Stock one cart with all the supplies you need for a specific procedure
- The right size cart for every task
- Customized interiors and accessories
- Unlimited design options
 - Stent Cart
 - PTCA
 - Pacemaker Cart
 - Pediatric Supplies
 - Computer Cart
 - Emergency Cart
 - Equipment Cart
 - Bedside Recovery Cart



MX6101DHG
Full Size Metal Exchange Cart
Hinged Doors with Glass Inserts. Also pictured:
Suture Module, Locking six inch Metal Drawer,
and three FifoGlide™ Catheter Modules.
26"D x 45"W x 79.25"H



S6102
Stainless Procedure Cart
26"D x 45.12"W x 44"H



M2200HG

Wall Storage Cabinets

Hinged Doors with Glass Inserts, flat top pictured. 22.75"D x 40.87"W x 84"H

M1600HG

Wall Storage Cabinets

Hinged Doors with Glass Insert, flat top pictured. 16.75"D x 40.87"W x 84"H

WALL STORAGE CABINET OPTIONS

Door options:

- Hinged with Glass Inserts
- Solid Hinged Doors
- FlexLok™ Door (roll down)
- No Door



**M1602
Half Wall
Storage
Cabinets**

16"D x 21.62"W
x 89.5"H

**M2202
Half Wall
Storage
Cabinets**

22"D x 21.62"W x
91.12"H

Half Cabinet is 21.6" wide, but has over 17 cu. ft. of interior storage space.



**VERTICAL DIVIDER
MODULE FOR LONG
BOXES.**

**SLIDE-OUT
VERTICAL
DIVIDER
MODULE
FOR EASY
ACCESS OF
LONG BOXES.**

Slide-out Vertical
Divider Module



MOBILE CATH LAB STORAGE

- Create a totally mobile cath lab storage system
- More storage capacity than any other mobile storage space
- Fifth wheel provides greater support, easier rolling, and better directional control



M3101D
Full Size Divided Modular Storage Cart
with FlexLok™ Door. 26"D x 45.12"W x 79.25"H



M3101DHG
Full Size Divided Modular Storage Cart
Hinged Doors with Glass Inserts.
26"D x 45.12"W x 79.25"H



M3101
Full Size Modular Storage Cart
with FlexLok™ Door. 26"D x 45.12"W x 79.25"H



MX3101DK
Full Size Exchange Cart
with FlexLok™ Door. 26"D x 45.12"W x 79.25"H



M3101
Full Size Modular Storage Cart
with FlexLok™ Door. 26"D x 45.12"W x 79.25"H



MX3101DK
Full Size Exchange Cart
with FlexLok™ Door. 26"D x 45.12"W x 79.25"H



M3102
Mid-Size Modular Storage Cart
with FlexLok™ Door.
26"D x 45.12"W x 42.5"H



MX3106DK
Mid-Size Exchange Cart
with FlexLok™ Door, Keyless Lock, and Stainless
Steel Top and Rail. 26"D x 45.12"W x 43.5"H

PROCEDURE CARTS

- Stock one cart with all the supplies you need for a specific procedure
- The right size cart for every task
- Customized interiors and accessories
- Unlimited design options
 - Stent Cart
 - PTCA
 - Pacemaker Cart
 - Pediatric Supplies
 - Computer Cart
 - Emergency Cart
 - Equipment Cart
 - Bedside Recovery Cart



MX6101DHG
Full Size Metal Exchange Cart
Hinged Doors with Glass Inserts. Also pictured:
Suture Module, Locking six inch Metal Drawer,
and three FifoGlide™ Catheter Modules.
26"D x 45.12"W x 78"H



M3110
MASS™ Small Series Cart
22.62"D x 27"W x 44"H

MASS METAL™

EXCHANGE CARTS GET
THE JOB DONE



**SLIDE-OUT SLOTSHELF™
WITH ADJUSTABLE
DIVIDERS**



**ABS TRAYS WITH OPTIONAL
DIVIDERS**



**FIFOGLIDE™ CATHETER
STORAGE MODULE**



**ABS SIDE PANELS ALLOW
EASY EXCHANGE**



MX6101DHG
Full Size Metal Exchange Cart
Hinged Doors with Glass Inserts.
26"D x 45.12"W x 78"H



MX6107
Half Size Metal
Exchange Cart
with FlexLok™ Door.
31"D x 28"W x 78.75"H



MX6102DK
Mid Size Metal Exchange Cart
with FlexLok™ Door and optional Keyless Lock.
26"D x 45.12"W x 44"H



MX6108
Standard Size Metal
Exchange Cart
with FlexLok™ Door.
26"D x 25.62"W x 44"H

ARTHROSCOPY



ENT



GENERAL SURGERY



LIVER



NEURO



ORTHO



THORACIC



VASCULAR



SURGICAL SUPPLIES "TO GO"

- Trays, baskets, and shelves slide in and out quickly
- Create a totally mobile storage system
- Move supplies quickly and easily where needed
- Reduce inventory, save money
- 14 and 20 gauge powder-coated steel
- Lightweight and easy to transport
- ABS side panels permit easy module exchange
- Change is inevitable – choose the one storage system that adapts

GEN. SUPPLY



GI



MX6101DK
Full Size Metal
Exchange Cart
with FlexLok™ Door
and optional Keyless Lock.
26"D x 45.12"W x 78.75"H

MX6107
Half Size Metal Exchange Cart
with FlexLok™ Door.
31"D x 28"W x 78.75"H

ADDITIONAL MASS™ MODELS TO COMPLETE YOUR MOBILE STORAGE SYSTEMS

Portable • Flexible • Functional

- Move supplies and instruments room to room effortlessly
- Consolidate inventory and eliminate duplication



MX6101K
Full Size Metal Exchange Storage Cart
with FlexLok™ Door and optional Keyless Lock.
26"D x 45.12"W x 79.25"H



M3101DHG
Full Size Modular Storage Cart
Hinged Doors with Glass Inserts.
26"D x 45.12"W x 79.25"H



M3101
Full Size Modular Storage Cart
with FlexLok™ Door.
26"D x 45.12"W x 79.25"H



S6101HG
Full Size Stainless Catheter Cabinet
with FIFO Glides, and optional Keyless Lock.
22"D x 41"W x 91.25"H



MX3101DK
Full Size Exchange Cart
with FlexLok™ Door and optional Keyless Lock.
26"D x 45.12"W x 79.25"H



MX3101DK
Full Size Exchange Cart
with FlexLok™ Door and optional Keyless Lock.
26"D x 45.12"W x 79.25"H



MX4400DHG
Metal Exchange Wall Storage Cabinet
Hinged Doors with Glass Inserts.
22.75"D x 41"W x 84.37"H



MX4402
Metal Exchange Half Cabinet
with FlexLok™ Door. 22"D x 21.62"W x 84.37"H

SURGICAL SUPPLY WALL CABINETS

- Hinged Doors with Glass Inserts providing visibility for all inventory
- Four door options (see page 15)
- Modular System allowing easy, rapid utilization or changes in components
- Easily secured to wall
- Full and half widths available
- Sloped or flat tops



MX4400D
Metal Exchange
Wall Storage
Cabinet
with FlexLok™ Door.
22"D x 41"W x 84.37"H



INTERVENTIONAL SUPPLY STORAGE CABINETS

Door options:

- Hinged with Glass Inserts
- Solid Hinged Doors
- FlexLok™ Door (roll down)
- No Door
- More interior storage capacity than any other medical storage...over 34 cu. ft.
- Interior modules to store every size package with no wasted space
- Maximize narrow spaces
- Half cabinet – 21.6" wide



**M1600HG
Wall
Storage
Cabinets**
Hinged Doors
with Glass
Inserts.
16"D x
40.87"W x
84"H

**M2200HG
Wall
Storage
Cabinets**
Hinged Doors
with Glass
Inserts.
22"D x
40.87"W x
84"H



**M2200
Wall Storage Cabinets**
with FlexLok™ Door. 22"D x 41"W x 91.25"H

**M1600
Wall Storage Cabinets**
with FlexLok™ Door. 16" D x 41"W x 89.5"H



(Pictured—Two each of M1601 and M2201 with 82" continuous counter top)

M1601
Upper Wall Cabinets
16"D x 41"W x 37.08"H


M2201
Lower Wall Cabinets
22"D x 41"W x 37.5"H

WALL STORAGE CABINETS

- Creative designs to fit unique environments
- Organize your vast interventional supply inventory
- Complete your lab and maximize storage with upper and lower cabinets and MASS™ modular components

Contact Us Today for a QuikQuote™

800.593.1900



Hanger Bar

Vertical Divider Module

**VERTICAL DIVIDER
MODULE FOR
LONG BOXES**

Also available in a pullout/slideout module



**SLOTSHELF™ MODULES
WITH OPTIONAL
COLORED DIVIDERS**



**OPTIONAL COLORED
DECOR STRIP HANDLE**

**KEY-LOCKED FLEXLOK™
TAMBOUR DOOR ROLLS
UP INTO TOP OF
CABINET, SHOWN WITH
CLOSED DOOR**



**BINS AVAILABLE IN
THREE CONVENIENT
SIZES - 4", 6", 8"
WIDTHS**



**WIRE BASKETS
AVAILABLE IN NINE
DIFFERENT SIZES**



**EASY-OPEN DRAWERS
WITH ADJUSTABLE
DIVIDERS**

MOBILIZE INTERVENTIONAL SUPPLIES

- Create a totally mobile storage system
- More storage capacity than any other mobile storage system
- Door options:
 - Hinged with Glass Inserts
 - Solid Hinged Doors
 - FlexLok™ Door (roll down)
 - No Door
- Mix and match interior components for effective product storage utilization



M3101D
Full Size Divided Modular Storage Cart
with FlexLok™ Door. 26"D x 45.12"W x 79.25"H



MX6101K
Full Size Exchange Storage Cart
with FlexLok™ Door, and optional
Keyless Lock.
26"D x 45.12"W x 79.25"H



M3101DHG
Full Size Divided Modular Storage Cart
Hinged Doors with Glass Inserts.
26"D x 45.12"W x 79.25"H




M3101
Full Size Modular Storage Cart
with FlexLok™ Door. 26"D x 45.12"W x 79.25"H



**SLIDE-OUT
VERTICAL
DIVIDER
MODULE
FOR EASY
ACCESS OF
LONG BOXES**

Slide-out Vertical
Divider Module





Standard System
Full Width Shelves

Exchange System
Slide-out Components

MX3101DK
Full Size Exchange Cart
with FlexLok™ Door and Keyless Lock.
26"D x 45.12"W x 79.25"H



OPTIONAL KEYLESS LOCK



**ABS TRAY WITH
OPTIONAL DIVIDERS**

RADIOLOGY SUPPLY CARTS

- Unlimited design options
- Unique interior modular components
- Easy exchange of components on ABS Side Panels



MX3106DK
Mid-Size Exchange Cart
with FlexLok™ Door, Stainless Steel Top
and Rail and optional Keyless Lock.
26"D x 45.12"W x 42"H



M3102
Mid-Size Modular Storage Cart
with FlexLok™ Door.
26"D x 45.12"W x 42"H

METAL CABINETS AND CARTS FOR RADIOLOGY

- Modules interchangeable between cabinets and carts
- Durable 14 and 20 gauge steel construction
- Wire baskets, trays, shelves, and SlotShelf™ module with adjustable dividers



MX6102DK
Mid Size Metal Exchange Cart
with FlexLok™ Door and optional Keyless Lock.
26"D x 45.12"W x 44"H



MX6108
Standard Size Metal Exchange Cart
with FlexLok™ Door.
26"D x 25.62"W x 44"H



MX4400
Metal Wall Storage Cabinet
with FlexLok™ Door.
22"D x 41"W x 84.75"H



MX4400DHG
Metal Exchange Wall Storage Cabinet
Hinged Doors with Glass Inserts.
22"D x 41"W x 84.75"H



MX4402
Metal Exchange Half Wall Storage Cabinet
with FlexLok™ Door.
22"D x 21.6"W x 84.37"H



MX6107
Half Size Metal Exchange Cart
with FlexLok™ Door.
31"D x 28"W x 78.75"H



MX6101DHG
Full Size Metal Exchange Cart
Hinged Doors with Glass Inserts.
26"D x 45.12"W x 78.75"H



M1600
Wall Storage Cabinets
with FlexLok™
Door. 16"D x
40.87"W x 89.5"H

M2200
Wall Storage Cabinets
with FlexLok™
Door. 22"D x
40.87"W x 91.25"H

CT SUPPLY STORAGE

- Wall storage cabinets organize your room and eliminate clutter
- Custom cabinets without the custom price
- Mobile supply carts– the ultimate convenience
- Unlimited design possibilities
- MASS™ Small Series Carts (below) available in over 500 different configurations

Contact Us Today for a QuikQuote™

800.593.1900



M3110
MASS™ Small Series Cart
22.62"D x 27"W x 44"H



MX3106DK
Mid-Size Exchange Cart
with FlexLok™ Door, Stainless Steel Top and Rail and optional Keyless Lock.
26"D x 45.12"W x 43.5"H



M3101
Full Size Modular Storage Cart
with FlexLok™ Door.
26"D x 45.12"W x 79.25"H

ACCESSORIES

MASS Chrome™ Shelves

Pictured with Dividers



MASS Chrome™ Slanted Suture Shelves

Call for Available Dimensions



Cart Cover

Call for Available Dimensions



Medical Kanban Bins

Call for Available Dimensions



Transparent Bins



Standard Stackable Bins



Wide Bin System



Stackable Dividable Bins



MASS Chrome™ High Density Storage

Pictured with Upper Rail System

Call for Available Dimensions





**THREE BIN SIZES:
4", 6", AND 8" WIDE**



**EQUIPMENT ACCESS DOOR
(MASS™ SMALL SERIES
CART ONLY)**



IV POLE EXTENDS TO 36"



PROPRIETARY LOCK



**ELECTRONIC KEYLESS
LOCK ADDS SECURITY
AND CONVENIENCE**



**SLOTSHELF™ MODULES
WITH OPTIONAL COLORED
DIVIDERS**



**AUTOCCLAVABLE SCOPE
TRAY AND LINER**



**SIX INCH LOCKING METAL
DRAWER**



**METAL SHELF OR CLEAR
ACRYLIC DIVIDED TRAY
FOR ADJUSTABLE ABOVE
CART STORAGE OF SMALL
ITEMS**



**FIFOGLIDE™ CATHETER
STORAGE MODULE**



**FOLD-DOWN
WRITING SHELF
16 3/16" X 19 1/8"**



**GLOVE BOX HOLDER
TWO OR THREE
BOX HOLDERS**



LED LIGHTING



**CORD MANAGER
INDIVIDUALLY WRAP
POWER CORDS AND
LEADS. ELIMINATE
TANGLED CORDS**




**KEYLESS LOCK
OPTIONAL ON ALL
CARTS EXCEPT MASS™
SMALL SERIES CARTS**




**POWER STRIP
HOSPITAL-GRADE SIX-PLUG.
ACCESS VIA CORD PORTS IN
BACK OR SIDE WALL**


TYPICAL CARDIAC CATH LAB & INTERVENTIONAL RADIOLOGY



TYPICAL METAL EXCHANGE SURGERY CARTS



TYPICAL MOBILE SURGERY STORAGE DESIGNS





PUTS YOU IN CONTROL

MASS™ Medical Storage puts you in control of time, space & productivity throughout the GI Lab with a wide range of customizable solutions to store, transport and protect your Endoscopy equipment. We work with you to design and build an efficient storage system to best utilize & track inventory in various sizes & quantities including colonoscopes, EUS scopes, TEE probes, bronchoscopes, cystoscopes, laparoscopes and small or large-diameter scopes. We also offer a variety of transport & procedure carts ranging from standard, auto-clavable and motorized models.

Our systems satisfy the Joint Commission standards set by Infection Control and manufacturer requirements for storage and are the most highly compliant in the industry.

THE JOINT COMMISSION:

- 2 filtered ambient air vents
- Easy-to-clean removal drip tray
- No contact between scopes & drip tray
- Scopes separated from other scopes
- Scopes hang straight, not coiled
- Padded walls to protect distal tip



Our Secure-A-Scope cabinets feature a revolutionary new design through simple engineering, setting the standard for scope storage and accessibility.

SECURE A SCOPE:

- Simply Engineered
- Positive loading & unloading so user knows scopes are securely in place
- Highest hanging height in the industry
- Front holders rotate for easy access to rear scopes and even scope rotation
- High quality molded ABS holding system
- Contiguous holder with no welds, glues or cavities to prevent scope contamination

CALL

800.593.1900 for a quote

LARGE SCOPE



M1200LHG
Secure-A-Scope™
Storage System
 Locking Hinged Door with Glass Inserts.
 Holds 9 Flexible Endoscopes.

SMALL SCOPE



M1201HG
Small Diameter Wall Mounted
Secure-A-Scope™ Scope Cabinet
 with hinged, 3-point locking, clear hinged door.
 Holds 4-8 small scopes.



M6204SS
Stainless Steel
Scope Transport Cart



M3103
Full Size ERCP Cart
 with FlexLok™ Door on upper and lower sections.

TRANSPORT CARTS

PROCEDURE CARTS

SCOPE CABINET OPTIONS:

- Keyed, keyless or digital proximity lock
- LED lights for easy viewing & scope management
- Proprietary specialized EUS holder to securely store Ultrasound connector
- Powered Circulation Fan

MASS MEDICAL SERVICES:

- Storage Needs Analysis
- Free Engineering CAD Design Service
- Mobile & Architectural Solutions



Mix

&



MASS Medical Storage already puts you in control with customized services, features and a range of materials including hospital-grade stainless steel, medical-grade laminate, & powder-coated metal. Knowing budgets also come in a wide range, MASS now offers **Mix & Match** services.

With **Mix & Match** services, you can customize your storage needs to include the Infection Control benefits of stainless steel doors with the more economical brushed aluminum laminate cabinets which have the remarkable look of stainless steel. The **Mix & Match** option affords you the style, ease of cleaning, and performance of stainless steel while still being responsible to your bottom line.

At MASS Medical Storage, we'll do anything to keep you in control.

CALL
800.593.1900

MATCH

With MASS Medical Storage, you will get the highest performance and strongest cabinets built today with our free engineering and architectural services to fit your space.

www.massmedicalstorage.com



PUTS ARCHITECTS & EQUIPMENT PLANNERS IN CONTROL

MASS Medical Storage has worked in partnership with Architects and Equipment Planners for more than 15 years in evaluating inventory & performance needs to optimize space, budget and ergonomics. MASS is your design & engineer partner from common spec files to fully-customized storage systems.



MASS MEDICAL STORAGE PLANNING & DESIGN SERVICES:

- Valuable expertise in Infection Control and Joint Commission standards for medical grade storage in the OR, Hybrid OR, Cath Lab, Radiology, GI, ED, Ortho and all other acute intervention areas
- Design services at no additional cost
- Knowledgeable and accessible design & engineering professionals
- Efficient custom storage systems to maximize time, space & productivity
- Enhanced construction documents with custom CAD & Inventor product drawings
- Elevations

CALL
800.593.1900
for Elevations & Construction Documents



MEDICAL STORAGE

PREMIER

Tubing, Stopcocks

Contracted supplier of

AmerinetSM



Contract Holder



Minority Business
MANUFACTURER OF THE YEAR
2013

MedAssets[®]
Contracted Supplier



Minority Business
MANUFACTURER OF THE YEAR
2014



(800) 593-1900 • www.MASSMedicalStorage.com



MEDICAL STORAGE

7848 Barton Street
Lenexa, KS 66214

800-593-1900
Fax: 913-438-2275