

MED-BIN CARTS

SL36BIN5 (pictured)
Standard Line Med-Bin Cart with BEST® lock on cabinet, pull-out shelf. Includes 36 - 5" bins, Standard Package

SL28BIN5 (not shown)
Standard Line Med-Bin Cart with BEST® lock on cabinet, pull-out shelf. Includes 28 - 5" bins, Standard Package

SL54BIN3 (pictured)
Standard Line Med-Bin Cart with a BEST® lock on cabinet, pull-out shelf. Capacity of 54 - 3.5" bins, Standard Package

SL42BIN3 (not shown)
Standard Line Med-Bin Cart with a BEST® lock on cabinet, pull-out shelf. Capacity of 42 - 3.5" bins, Standard Package



Bin dividers & labels included with each cart

SL36BIN5XL (pictured)
Standard Line Med-Bin Cart with BEST® lock on cabinet, pull-out shelf. Includes 36 extra long 5.6" bins, Standard Package

SL24BIN5XL (pictured)
Standard Line Med-Bin Cart with BEST® lock on cabinet, pull-out shelf. Includes 24 extra long 5.6" bins, Standard Package

SL18BIN5XL (pictured)
Standard Line Med-Bin Cart with BEST® lock on cabinet, pull-out shelf. Includes 18 extra long 5.6" bins, Standard Package



Bin transfer cart options also available

Bin dividers & labels included with each cart



40442
Replacement
3.5" Bins
2.75"H x 3.5"W x 9.5"D

680533
Bin Label for
3.5" Bins
(150 Labels)

680513
10 Bin Dividers
for 3.5" Bins



40446
Replacement
5" Bins
2.75"H x 5"W x 9.5"D

680535
Bin Labels for
5" Bins
(100 Labels)

680515
10 Bin Dividers
for 5" Bins



40433
Replacement
5.6" Extra Long Bins
2.75"H x 5.6"W x 18"D

680537
Bin Labels for
5.6" Bins
(210 Labels)

40434
20 Bin Dividers
for 5.6" Bins

VL36BIN3 (pictured)
Value Line Med-Bin Cart with
key lock on cabinet, internal
narcotics box, pull-out shelf.
Includes 36 - 3.5" bins,
Standard Package




**Mechanical keyless entry
available on
Value Line Carts**



VL24BIN5 (pictured)
Value Line Med-Bin Cart with key lock on
cabinet, internal narcotics box, pull-out shelf.
Includes 24 - 5" bins, Standard Package

VL24MBIN5 (not shown)
Value Line Med-Bin Cart with keyless lock on
cabinet, internal narcotics box, pull-out shelf.
Includes 24 - 5" bins, Standard Package

VL36MBIN3 (not shown)
Value Line Med-Bin Cart with keyless lock on
cabinet, internal narcotics box, pull-out shelf.
Includes 36 - 3.5" bins, Standard Package



Bin Drawer shown with 3.5" bins

Harloff standard bin carts are available
with 5" and 3.5" bins. They extend
8.5" into the drawer and have a
small tab that fits under a
metal edge. There is 7" of
open space behind



VL12BIN5 (pictured)
Value Line Med-Bin Cart with
key lock on cabinet, internal
narcotics box, pull-out shelf. Includes
12 - 5" bins, Standard Package
(shown with a Specialty Package)

VL18BIN3 (pictured)
Value Line Med-Bin Cart with
key lock on cabinet, internal
narcotics box, pull-out shelf.
Includes 18 - 3.5" bins,
Standard Package

**Accessory packages available
for all cart styles**



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

**Standard Line Med-Bin Cart with 36 – 5” Bins Standard Package
 #SL36BIN5**

Specifications



Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Dual column medication cart with 36 – 5” wide medication bins. Dividers and labels are included, single width narcotics drawer
- BEST® brand lock on cabinet; Illinois lock on narcotics drawer
- Painted pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Full wrap around vinyl bumper, protecting the cart
- 5” ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish
- Offered in specialty package (#SL36BIN5SP) Included are: #68640-P1 drawer tray with adjustable dividers, #680551 eye drop tray, #AL680435 aluminum waste container and #684810 locking sharps container

Drawer Configuration:

Left Column Drawers:

- 1 ea – 4” multipurpose storage drawer
- 9 ea – 3.25” med bin drawers with 4 – 5” bins per drawer

Right Column Drawers:

- 1 ea – 4” Multipurpose Storage Drawer
- 1 ea – 9.75” Locking Narcotics Drawer
- 2 ea – 9.75” Multipurpose Storage Drawer

Thirteen available color choices:

Sand, Gloss Beige, Burgundy,
 Light Blue, Navy, Light Gray,
 Mauve, Red, Teal, Yellow,
 Hammer Tone Blue, Hammer Tone Gray
 and Hammer Tone Green

Assembled Cart Dimensions:

- 45.75” H x 43” W x 23.5” D
 (116.2 cm x 109.2 cm x 59.7 cm)
- Writing Surface: 44” (111.8 cm)
- Cart Weight: 220 lbs (100 kg)
- Approx Shipping Weight: 265 lbs (120.5 kg)



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Standard Line Med-Bin Cart with 54 – 3” Bins Standard Package #SL54BIN3

Specifications

Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Dual column medication cart with 54 – 3” wide medication bins. Dividers and labels are included, single width narcotics drawer
- BEST® brand lock on cabinet; Illinois lock on narcotics drawer
- Painted steel pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Full wrap around vinyl bumper, protecting the cart
- 5” ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish
- Offered in specialty package (#SL54BIN3SP) Included are: #68640-P1 drawer tray with adjustable dividers, #680551 eye drop tray, #AL680435 aluminum waste container and #684810 locking sharps container

Drawer Configuration:

Left Column Drawers:

- 1 ea – 4” multipurpose storage drawer
- 9 ea – 3.25” med bin drawers with 6 – 3” bins per drawer

Right Column Drawers:

- 1 ea – 4” Multipurpose Storage Drawer
- 1 ea – 9.75” Locking Narcotics Drawer
- 2 ea – 9.75” Multipurpose Storage Drawer

Assembled Cart Dimensions:

- 45.75” H x 43” W x 23.5” D
(116.2 cm x 109.2 cm x 59.7 cm)
- Writing Surface: 44” (111.8 cm)
- Cart Weight: 220 lbs (100 kg)
- Approx Shipping Weight: 265 lbs (120.5 kg)



Thirteen available color choices:

Sand, Gloss Beige, Burgundy,
Light Blue, Navy, Light Gray,
Mauve, Red, Teal, Yellow,
Hammer Tone Blue, Hammer Tone Gray
and Hammer Tone Green



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

**Standard Line Med-Bin Cart with 28 – 5” Bins, Standard Package
 #SL28BIN5**

Specifications

Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Dual column medication cart with 28 - 5” wide medication bins. Dividers and labels are included, single width narcotics drawer
- BEST® brand lock on cabinet; Illinois lock on narcotics drawer
- Painted pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Full wrap around vinyl bumper, protecting the cart
- 5” ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish
- Offered in specialty package (#SL28BIN5SP) Included are: #68640-P1 drawer tray with adjustable dividers, #680551 eye drop tray, #AL680435 aluminum waste container and #684810 locking sharps container



Drawer Configuration:

Left Column Drawers

- 1 ea – 4” multipurpose storage drawer
- 7 ea – 3.25” med bin drawers with 4 – 5” bins per drawer
- 1 ea – 6.5” multipurpose drawer with adjustable dividers

Right Column Drawers

- 1 ea – 4” multipurpose storage drawer
- 1 ea – 9.75” locking narcotics drawer
- 2 ea – 9.75” multipurpose storage drawer

Assembled Cart Dimensions:

- 45.75” H x 43” W x 23.5” D (116.2 cm x 109.2 cm x 59.7 cm)
- Writing surface: 44” (111.8 cm)
- Cart weight: 205 lbs (93.2 kg)
- Approx shipping weight: 250 lbs (113.7 kg)

Thirteen available color choices:

Sand, Gloss Beige, Burgundy, Light Blue, Navy, Light Gray, Mauve, Red, Teal, Yellow, Hammer Tone Blue, Hammer Tone Gray and Hammer Tone Green



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

**Standard Line Med-Bin Cart with 42 – 3.5” Bins, Standard Package
 #SL42BIN3**

Specifications



Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Dual column medication cart with 42 – 3.5” wide medication bins. Dividers and labels are included, single width narcotics drawer
- BEST® brand lock on cabinet; Illinois lock on narcotics drawer
- Painted pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Full wrap around vinyl bumper, protecting the cart
- 5” ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish
- Offered in specialty package (#SL42BIN3SP) Included are: #68640-P1 drawer tray with adjustable dividers, #680551 eye drop tray, #AL680435 aluminum waste container and #684810 locking sharps container

Drawer Configuration:

Left Column Drawers:

- 1 ea – 4” multipurpose storage drawer
- 7 ea – 3.25” med bin drawers with 6 – 3.5” bins per drawer
- 1 ea – 6.5” multipurpose drawer with adjustable dividers

Right Column Drawers:

- 1 ea – 4” multipurpose storage drawer
- 1 ea – 9.75” locking narcotics drawer
- 2 ea – 9.75” multipurpose storage drawer

Thirteen available color choices:

Sand, Gloss Beige, Burgundy,
 Light Blue, Navy, Light Gray,
 Mauve, Red, Teal, Yellow,
 Hammer Tone Blue, Hammer Tone Gray
 and Hammer Tone Green

Assembled Cart Dimensions:

- 45.75” H x 43” W x 23.5” D
 (116.2 cm x 109.2 cm x 59.7 cm)
- Writing surface: 44” (111.8 cm)
- Cart weight: 220 lbs (100 kg)
- Approx shipping weight: 265 lbs (120.5 kg)



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

Value Line Med Bin Cart with 24 – 5 inch Bins, Key Locking #VLT24BIN5

Harloff's value medication carts offer many of the same features our customers love in our large med carts at an economical price point. This economy med cart has become our fastest growing small med cart.



Color Shown is Hammertone Blue



Thirteen available color choices: (may be viewed at www.harloff.com)

- Sand, Gloss Beige, Burgundy, Light Blue, Yellow, Navy, Light Gray, Mauve, Red, Teal, Hammertone Blue, Hammer Tone Gray and Hammer Tone Green

Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Tall cabinet with key locking system with 4 keys
- Internal narcotics box in bottom 9" drawer
- Pull out side shelf
- Includes 20 Harloff 5" medication bins with dividers and labels
- Bumper frame with projected front casters and (4) corner bumpers
- 5" ball-bearing premium casters, all full swivel, - one with brake, one directional
- Painted steel top rail/push handle
- Removable plastic top
- Full extension ball-bearing drawer slides
- Cart is delivered fully assembled
- Pre-threaded accessory mounting holes
- Hard baked powder coat finish
- Certified "Latex Free"
- Twelve year manufacturer's warranty

Drawer Configuration:

- 1 - ea 3" multipurpose storage drawer
- 6 - ea 3" bin drawers including 4 - 5" bins each
- 1 - ea 9" drawer with 1 divider and locking narcotics box

Internal Drawer Dimensions:

- 23" W x 17" D (58.5 cm W x 43.2 cm D)

Assembled Cart Dimensions:

- 46.75" H x 29.5" W x 22.5" D (121.3 cm H x 75 cm W x 57.2 cm D)
- Writing surface: 43.25" (109.9 cm)
- Cart weight: 165 pounds (75 kg)
- Approx. shipping weight: 210 pounds (95.5kg)



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

Value Line Med Bin Cart with 36 – 3.5 inch Bins, Key Locking #VLT36BIN3

Harloff’s value medication carts offer many of the same features our customers love in our large med carts at an economical price point. This economy med cart has become our fastest growing small med cart.



Color Shown is Sand



Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Tall cabinet with key locking system with 4 keys
- Internal locking narcotics box in bottom 9” drawer
- Pull out side shelf
- Includes 36 Harloff 3.5” medication bins with dividers and labels
- Bumper frame with projected front casters and (4) corner bumpers
- 5” ball-bearing premium casters, all full swivel, - one with brake, one directional
- Painted steel top rail/push handle
- Removable plastic top
- Full extension ball-bearing drawer slides
- Cart is delivered fully assembled
- Pre-threaded accessory mounting holes
- Hard-baked powder coat finish
- Certified “Latex Free”
- Twelve year manufacturer’s warranty

Drawer Configuration:

- 1 – ea 3” multipurpose storage drawers
- 6 – ea 3” bin drawers including 6 - 3.5” bins each
- 1 – ea 9” drawer with 1 divider and locking narcotics box

Internal Drawer Dimensions:

- 23” W x 17” D (58.5 cm W x 43.2 cm D)

Assembled Cart Dimensions:

- 46.75” H x 29.5” W x 22.5” D (121.3 cm H x 75 cm W x 57.2 cm D)
- Writing surface: 43.25” (109.9 cm)
- Cart weight: 170 pounds (77.3 kg)
- Approx. shipping weight: 210 pounds (94.5 kg)

Thirteen available color choices:

(may be viewed at www.harloff.com)

- Beige, Burgundy, Light Blue, Dark Blue, Light Gray, Mauve, Red, Teal, Yellow, Hammertone Blue, Hammertone Gray, Hammertone Green and Sand



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

Value Line Med Bin Cart with 12 – 5 inch Bins, Key Locking #VLT12BIN5

Harloff’s value medication carts offer many of the same features our customers love in our large med carts at an economical price point. This economy med cart has become our fastest growing small med cart.

Specifications



Construction:

- Painted steel cabinet and drawers

Features:

- Tall cabinet with key lock with 4 keys
- Internal locking narcotics box in bottom 9” drawer
- Pull out side shelf
- Includes 12 Harloff 5” medication bins with dividers and labels
- Bumper frame with projected front casters and (4) corner bumpers
- 5” ball-bearing premium casters, all full swivel, - one with brake, one directional
- Painted steel top rail/push handle
- Removable plastic top
- Full extension ball-bearing drawer slides
- Cart is delivered fully assembled
- Pre-threaded accessory mounting holes
- Hard-baked powder coat finish
- Certified “Latex Free”
- Twelve year manufacturer’s warranty

Drawer Configuration:

- 1 - ea 3” multipurpose storage drawer
- 3 - ea 3” bin drawers including 4 - 5” bins each
- 1 - ea 9” drawer with 2 adjustable dividers
- 1 - ea 9” drawer with 1 divider and locking narcotics box

Internal Drawer Dimensions:

- 23” W x 17” D (58.5 cm W x 43.2 cm D)

Assembled Cart Dimensions:

- 46.75” H x 29.5” W x 22.5” D (121.3 cm H x 75 cm W x 57.2 cm D)
- Writing surface: 43.25” (109.9 cm)
- Cart weight: 150 pounds (68.2kg)
- Approx. shipping weight: 190 pounds (86.4 kg)

Thirteen available color choices:

(may be viewed at www.harloff.com)

- Sand, Gloss Beige, Burgundy, Light Blue, Yellow, Navy, Light Gray, Mauve, Red, Teal, Hammer Tone Blue, Hammer Tone Gray and Hammer Tone Green



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Value Line Med Bin Cart with 18 – 3.5 inch Bins, Key Locking #VLT18BIN3

Harloff's value medication carts offer many of the same features our customers love in our large med carts at an economical price point. This economy med cart has become our fastest growing small med cart.

Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Tall cabinet with key locking system with 4 keys
- Internal locking narcotics box in bottom 9" drawer
- Pull out side shelf
- Includes 18 Harloff 3.5" medication bins with dividers and labels
- Bumper frame with projected front casters and (4) corner bumpers
- 5" ball-bearing premium casters, all full swivel, - one with brake, one directional
- Painted steel top rail/push handle
- Removable plastic top
- Full extension ball-bearing drawer slides
- Cart is delivered fully assembled
- Pre-threaded accessory mounting holes
- Hard-baked powder coat finish
- Certified "Latex Free"
- Twelve year manufacturer's warranty

Drawer Configuration:

- 1 - ea 3" multipurpose storage drawer
- 3 - ea 3" bin drawers including 6 -3.5" bins each
- 1 - ea 9" drawer with 2 adjustable dividers
- 1 - ea 9" drawer with 1 divider and locking narcotics box

Internal Drawer Dimensions:

- 23" W x 17" D (58.5 cm W x 43.2 cm D)

Assembled Cart Dimensions:

- 46.75" H x 29.5" W x 22.5" D
(121.3 cm H x 75 cm W x 57.2 cm D)
- Writing surface: 43.25" (109.9 cm)
- Cart weight: 155 pounds (70.5 kg)
- Approx. shipping weight: 195 pounds (88.6 kg)



Color Shown is Sand



Thirteen available color choices:

(may be viewed at www.harloff.com)

- Sand, Gloss Beige, Burgundy, Light Blue, Yellow, Navy, Light Gray, Mauve, Red, Teal, Hammer Tone Blue, Hammer Tone Gray and Hammer Tone Green