

# Diaper & Specimen Scales

## Model PC-400



Constructed entirely of easy to clean stainless steel, the Model PC-400 is a compact battery powered scale. Extremely precise, the large LCD read-out displays weight in as fine as a 0.1 gram. Power is provided by either an internal rechargeable battery pack or the included AC adapter. The auto shut-off feature maximizes battery life.

Approvals:



### MAXIMUM CAPACITY

Model PC400-02 - 2 x 0.0002 lb (1000 x 0.1 g)

Model PC400-05 - 5 x 0.002 lb (2300 x 1 g)

Model PC400-10 - 10 x 0.005 lb (4500 x 2 g)

### ACCURACY

Model PC400-02 +/- 0.0002 lb (0.1 g)

Model PC400-05 +/- 0.002 lb (1 g)

Model PC400-10 +/- 0.005 lb (2 g)

### DIMENSIONS

Overall 8.75" L x 6.5" W x 2" H

Platter 6.5" x 6.5"

### DISPLAYED UNITS

lb, kg, g, oz, or lb:oz

### CONTROLS

Zero, On/Off, and Units

### OPTIONS

RS-232 Cable (P/N PCOPT05)

### POWER

115V AC Adapter, 6' Cord  
Rechargeable Battery, 10 hour life

### CONSTRUCTION

304 Stainless Steel

### WARRANTY

Scale - One Year / Battery - 90 Days

## Model PC-500



The PC-500 features IP66 water resistance to allow for thorough sanitizing. The 8.5" x 8.5" stainless steel platter can be easily removed and sanitized. Operates on standard alkaline batteries for 1,000 hours of continuous use. The auto shut-off feature will turn off the scale when it's not being used. An AC adapter is available as an option.

Approvals:



### MAXIMUM CAPACITY

Model PC500-05 - 5 x 0.002 lb (2200 x 1 g)

Model PC500-10 - 10 x 0.005 lb (4500 x 2 g)

Model PC500-25 - 25 x 0.01 lb (11,300 x 5 g)

### ACCURACY

Model PC500-05 +/- 0.002 lb (1 g)

Model PC500-10 +/- 0.005 lb (2 g)

Model PC500-25 +/- 0.01 lb (5 g)

### DIMENSIONS

Overall 12.75" L x 9.75" W x 5.26" H

Platter 8.5" x 8.5"

### DISPLAYED UNITS

lb, kg, g, oz, or lb:oz

### CONTROLS

Zero, On/Off, and Units

### OPTIONS

RS-232 Cable (P/N PC5OPT03)

### POWER

Six D Batteries (Included)  
or Optional AC Adapter

### CONSTRUCTION

Scale Platter - Stainless Steel  
Case - Plastic

### WARRANTY

One Year



**Doran Scales, Inc.** 1-800-264-4106  
1315 Paramount Parkway Fax: 630-879-0073  
Batavia, IL 60510 www.doranmedical.com