

4040-650-200

Power Hi-Lo & Power Back Exam Tables



Standard Features

- Power Hi-Lo and power backrest
- 650lb (294.84kg) patient weight capacity
- Removable upholstered top with seamless antimicrobial guard
- 3 year limited warranty
- Low access ADA compliant 18.5" (±0.5") ((46.9cm)(±1.27cm))
- Two function foot control with Hi-Lo and power back
- Large front storage drawer
- Side storage drawer - reversible
- Dual electro-mechanical actuator system for efficient speed
- 12 gauge steel extension leg pull-out to accommodate all patient heights
- Seamless upholstered leg cushion elevates 10° for patient comfort
- Paper roll holder (18" – 21") (45.7cm – 53.3cm)
- One touch, easy-glide stirrups with 4 lateral positions and friction lock length adjustment
- Pelvic tilt (7°)
- Magnetic hand control with OneTouch Patient Positioning™ – two programmable table positioning buttons, plus return-to-home
- Hospital grade receptacle 110V with child safety cover (patient left)
- Stainless steel drain pan (1.5" deep) (3.81cm)

Optional Accessories

- **014** – Leg extension booster pad to elevate leg section (*field installable*)
- **146** – OneTouch Wheelbase™ for infection control, ADA compliant
- **251** – Bierhoff knee crutches (*pair/field installable*)
- **469** – Welch Allyn Light bracket capabilities (NOT actual bracket) (*specify patient right/left*)
- **822** – Stainless steel drain pan with spout (*factory installed only*)
- **870** – Drawer warmer 110V (*factory installed only*)
- **4040-650-200F** – Fire retardant foam (CAL133)

Table Specifications

- **Height:**
 - Min. 18.5" (± .5") (46.9cm) (±1.27cm) at seat
 - Max. 37" at seat (93.9cm)
- **Width:** 27" (68.6cm)
- **Length:** 57" – 76" with foot section extended (144.7cm - 193cm)
- **Leg rest:** 21" W x 13" D (53.3cm W x 33cm D)
- **Front Storage Drawer:** 19" W x 12" D (48.26cm W x 30.48cm D)
- **Rear Storage Drawer:** 12.5" W x 24.5" D (31.75cm W x 62.23cm D)
- **Electrical:** 120 VAC, 60HZ
- **Stirrups:** Four lateral positions, Infinitely adjustable

Colors



4040-650-300

Power Hi-Lo & Power Back Exam Tables



Standard Features

- Power Hi-Lo with power backrest
- 650lb (294.84kg) patient weight capacity
- Removable upholstered top with seamless antimicrobial guard
- OneTouch Wheelbase™ for infection control
- 3 year limited warranty
- Low access ADA compliant 18.5" (±0.5") ((46.9cm)(±1.27cm))
- Two function foot control with Hi-Lo and power back
- Large front storage drawer
- Side storage drawer - reversible
- Dual electro-mechanical actuator system for efficient speed
- 12 gauge steel extension leg pull-out to accommodate all patient heights
- Seamless upholstered leg cushion elevates 10° for patient comfort
- Paper roll holder (18" – 21") (45.7cm – 53.3cm)
- One touch, easy-glide stirrups with 4 lateral positions and friction lock length adjustment
- Pelvic tilt (7°)
- Magnetic hand control with OneTouch Patient Positioning™ – two programmable table positioning buttons, plus return-to-home
- Hospital grade receptacle 110V with child safety cover (patient left)
- Drawer warmer 110V
- Stainless steel drain pan (1.5" deep) (3.81cm)

Optional Accessories

- **014** – Leg extension booster pad to elevate leg section
- **251** – Bierhoff knee crutches (*pair/field installable*)
- **822** – Stainless steel drain pan with spout (*factory installed only*)
- **4040-650-300F** – Fire retardant foam (CAL133)

Table Specifications

- **Height:**
 - Min. 18.5" (± .5") (46.9cm) (±1.27cm) at seat
 - Max. 37" at seat (93.9cm)
- **Width:** 27" (68.6cm)
- **Length:** 57" – 76" with foot section extended (144.7cm - 193cm)
- **Leg rest:** 21" W x 13" D (53.3cm W x 33cm D)
- **Front Storage Drawer:** 19" W x 12" D (48.26cm W x 30.48cm D)
- **Rear Storage Drawer:** 12.5" W x 24.5" D (31.75cm W x 62.23cm D)
- **Electrical:** 120 VAC, 60HZ
- **Stirrups:** Four lateral positions, Infinitely adjustable

Colors



4070-650-200

Power Hi-Lo & Manual Back Exam Tables



Standard Features

- Power Hi-Lo with air-spring assisted manual back adjustment
- 650lb (294.84kg) patient weight capacity
- Removable upholstered top with seamless antimicrobial guard
- 3 year limited warranty
- Low access ADA compliant 18.5" (±0.5") ((46.9cm)(±1.27cm))
- One function foot control with Hi-Lo
- Large front storage drawer
- Side storage drawer - reversible
- Dual electro-mechanical actuator system for efficient speed
- 12 gauge steel extension leg pull-out to accommodate all patient heights
- Seamless upholstered leg cushion elevates 10° for patient comfort
- Paper roll holder (18" – 21") (45.7cm – 53.3cm)
- One touch, easy-glide stirrups with 4 lateral positions and friction lock length adjustment
- Pelvic tilt (7°)
- Magnetic hand control with OneTouch Patient Positioning™ – two programmable table positioning buttons, plus return-to-home
- Hospital grade receptacle 110V with child safety cover (patient left)
- Stainless steel drain pan (1.5" deep) (3.81cm)

Optional Accessories

- **014** – Leg extension booster pad to elevate leg section (*field installable*)
- **146** – OneTouch Wheelbase™ for infection control, ADA compliant
- **251** – Bierhoff knee crutches (*pair/field installable*)
- **469** – Welch Allyn Light bracket capabilities (NOT actual bracket) (*specify patient right/left*)
- **822** – Stainless steel drain pan with spout (*factory installed only*)
- **870** – Drawer warmer 110V (*factory installed only*)
- **4070-650-200F** – Fire retardant foam (CAL133)

Table Specifications

- **Height:**
Min. 18.5" (± .5") (46.9cm) (±1.27cm) at seat
Max. 37" at seat (93.9cm)
- **Width:** 27" (68.6cm)
- **Length:** 57" – 76" with foot section extended (144.7cm - 193cm)
- **Leg rest:** 21" W x 13" D (53.3cm W x 33cm D)
- **Front Storage Drawer:** 19" W x 12" D (48.26cm W x 30.48cm D)
- **Rear Storage Drawer:** 12.5" W x 24.5" D (31.75cm W x 62.23cm D)
- **Electrical:** 120 VAC, 60HZ
- **Stirrups:** Four lateral positions, Infinitely adjustable

Colors



4070-650-300

Power Hi-Lo & Manual Back Exam Tables



Standard Features

- Power Hi-Lo with air-spring assisted manual back adjustment
- 650lb (294.84kg) patient weight capacity
- Removable upholstered top with seamless antimicrobial guard
- OneTouch Wheelbase™ for infection control
- 3 year limited warranty
- Low access ADA compliant 18.5" (±0.5") ((46.9cm)(±1.27cm))
- One function foot control with Hi-Lo
- Large front storage drawer
- Side storage drawer - reversible
- Dual electro-mechanical actuator system for efficient speed
- 12 gauge steel extension leg pull-out to accommodate all patient heights
- Seamless upholstered leg cushion elevates 10° for patient comfort
- Paper roll holder (18" – 21") (45.7cm – 53.3cm)
- One touch, easy-glide stirrups with 4 lateral positions and friction lock length adjustment
- Pelvic tilt (7°)
- Magnetic hand control with OneTouch Patient Positioning™ – two programmable table positioning buttons, plus return-to-home
- Hospital grade receptacle 110V with child safety cover (patient left)
- Drawer warmer 110V
- Stainless steel drain pan (1.5" deep) (3.81cm)

Optional Accessories

- **014** – Leg extension booster pad to elevate leg section
- **251** – Bierhoff knee crutches (*pair/field installable*)
- **822** – Stainless steel drain pan with spout (*factory installed only*)
- **4070-650-300F** – Fire retardant foam (CAL133)

Table Specifications

- **Height:**
 - Min. 18.5" (± .5") (46.9cm) (±1.27cm) at seat
 - Max. 37" at seat (93.9cm)
- **Width:** 27" (68.6cm)
- **Length:** 57" – 76" with foot section extended (144.7cm - 193cm)
- **Leg rest:** 21" W x 13" D (53.3cm W x 33cm D)
- **Front Storage Drawer:** 19" W x 12" D (48.26cm W x 30.48cm D)
- **Rear Storage Drawer:** 12.5" W x 24.5" D (31.75cm W x 62.23cm D)
- **Electrical:** 120 VAC, 60HZ
- **Stirrups:** Four lateral positions, Infinitely adjustable

Colors



4010

4000 Series Power Procedure / Exam Tables



Four Function Power Procedure Table



Standard Features

- Ultra comfort seamless top with antimicrobial guard
- ADA compliant
- 19" (48.3cm) low access seat height
- Power Hi-Lo
- Power back
- Power tilt
- Power leg section
- Four function hand control (hi-lo, back, leg rest, and seat tilt)
- Articulating headrest
- Stainless steel leg pullout
- Removable stainless steel drain pan (3.5qt./3.3L)
- Two position upholstered kneeler/footrest
- One touch, easy-glide stirrups with 3 lateral positions and friction lock length adjustment
- Concealed paper roll holder (18"-21") (45.7cm-53.3cm)
- Trendelenburg position
- Electro-mechanical actuator system
- 600lb (272.2kg) patient weight capacity
- 3 year limited warranty

Optional Accessories

- 035 – Arm board assembly (specify patient right or left/factory installed)
- 041 – Patient assist armrests (pair/factory installed)
- 045 – Accessory rail (pair/specify head or seat section/factory installed)
- 065 – Vision screen
- 143 – Easy-Glide locking casters (set of four/ factory installed)
- 251 – Bierhoff knee crutches (pair)
- 375 – Safety rails full assembly (pair/factory installed)
- 411 – Four function programmable foot control (simultaneous use of foot and hand controls)
- 810 – Restraint straps assembly (pair/factory installed)
- 822 – Stainless steel drain pan with spout (factory installed)
- 852 – Hospital grade receptacle
- 4010F – Fire retardant foam

Table Specifications

- **Height:** Min. 19" (48.3cm) at seat
Max.35" (88.9cm) at seat
- **Width:** 27" (68.6cm)
- **Length:** 71.5" (181.6cm) with foot section extended; 85" (215.9cm) with foot section and headrest extended
- **Leg rest:** 16.25" (41.3cm) W x 9" (22.9cm) L
- **Head rest:** 15" (38.1cm) W x 8.5" (21.6cm) H
- **Electrical:** 120 VAC, 60 HZ or 220/240VAC, 60HZ
- **Stirrups:** Infinitely adjustable

Colors



4011

4000 Series Power Procedure / Exam Tables



Four Function Power Procedure Table



*Photo shown with optional accessory #046 – Accessory rail

Standard Features

- Ultra comfort seamless top with antimicrobial guard
- ADA compliant
- 19" (48.3 cm) low access seat height
- Power Hi-Lo
- Power back
- Power tilt
- Power leg section
- Four function programmable foot control (hi-lo, back, leg rest, and seat tilt) with three memory positions
- Articulating headrest
- Stainless steel leg pullout
- Removable stainless steel drain pan (3.5qt./3.3L)
- Two position upholstered kneeler/footrest
- One touch, easy-glide stirrups with 3 lateral positions and friction lock length adjustment
- Concealed paper roll holder (18"-21") (45.7cm-53.3cm)
- Trendelenburg position
- Electro-mechanical actuator system
- 600lb (272.2kg) patient weight capacity
- 3 year limited warranty

Optional Accessories

- 035 – Arm board assembly (specify patient right or left/factory installed)
- 041 – Patient assist armrests (pair/factory installed)
- 046 – Accessory rail (pair/specify head or seat section/factory installed)
- 065 – Vision screen
- 143 – Easy-Glide locking casters (set of four/ factory installed)
- 251 – Bierhoff knee crutches (pair)
- 375 – Safety rails full assembly (pair/factory installed)
- 415 – Four function hand control (simultaneous use of foot and hand controls)
- 810 – Restraint straps assembly (pair/factory installed)
- 822 – Stainless steel drain pan with spout (factory installed)
- 852 – Hospital grade receptacle
- 4011F – Fire retardant foam

Table Specifications

- **Height:** Min. 19" (48.3cm) at seat
Max. 35" (88.9cm) at seat
- **Width:** 27" (68.6cm)
- **Length:** 71.5" (181.6cm) with foot section extended; 85" (215.9cm) with foot section and headrest extended
- **Leg rest:** 16.25"(41.3cm) W x 9" (22.9cm) L
- **Head rest:** 15" (38.1cm) W x 8.5" (21.6cm) H
- **Electrical:** 120 VAC, 60 HZ or 220/240VAC, 60HZ
- **Stirrups:** Infinitely adjustable

Colors



5020

5000 Series Power Procedure Tables



Power Procedure Table



Standard Features

- Ultra comfort seamless top with antimicrobial guard
- 24" (61cm) low access seat high
- Power Hi-Lo
- Power backrest
- Power tilt
- Power leg section
- Four function hand control
- Articulating headrest
- Stainless steel leg pullout
- Removable stainless steel drain pan (1.9 Quart / 1.8L)
- Two position upholstered kneeler / footrest
- One touch, easy-glide stirrups with 3 lateral positions and friction lock length adjustment
- Concealed paper roll holder (18" – 21") (45.7cm – 53.3cm)
- Trendelenburg position
- Electro-mechanical actuator system
- 400lb (181.4kg) weight capacity
- 3 year limited warranty

Optional Accessories

- **035** – Armboard assembly (*specify patient right or left/factory installed*)
- **041** – Patient assist armrests (*factory installed*)
- **046** – Accessory rail pair (*specify head or seat section/ factory installed*)
- **065** – Vision screen
- **143** – Easy-Glide locking casters (*set of four, factory installed*)
- **251** – Bierhoff knee crutches (*pair*)
- **375** – Safety rails full assembly (*pair*)
- **410** – Four function programmable foot control (*substitution for hand control/60-*)
- **411** – Four function programmable foot control (*simultaneous use of foot and hand controls*)
- **810** – Restraint straps (*factory installed*)
- **822** – Stainless steel drain pan with spout
- **852** – Hospital grade receptacle 110V (*factory installed*)
- **855** – Hospital grade receptacle 220V (*factory installed*)
- **856** – 220/240 Volt (*factory installed*)
- **5020F** – Fire retardant foam

Table Specifications

- **Height:** Min. 24" (61cm) at seat
Max. 36" (91.4cm) at seat
- **Width:** 26.5" (67.3cm)
- **Length:** 71.5" (181.6cm) with foot section extended, 85" (215.9cm) with foot section and head rest extended
- **Leg Rest:** 16.25" (41.3cm) W x 9" (22.9cm) H
- **Head Rest:** 15" (38.1cm) W x 8.5" (21.6cm) H
- **Electrical:** 120 VAC, 60 HZ or 220/240VAC, 60HZ
- **Stirrups:** Infinitely adjustable

Colors



5005

5000 Series Bariatric Treatment Tables



Two Function Bariatric Power Table



Standard Features

- Ultra comfort seamless top with antimicrobial guard and built in pillow
- ADA Compliant with 19" (48.3cm) low access
- All-steel heavy duty construction
- Power Hi-Lo
- Two function hand control
- Concealed paper roll holder (18" – 21") (45.7cm – 53.3cm)
- Large, full extension storage drawer on roller bearing drawer guides
- 10" (25.4cm) leg section
- Electro-mechanical actuator system
- 800lb (363kg) weight capacity
- 3 year limited warranty

Optional Accessories

- **035** – Armboard assembly (*specify patient right or left/factory installed*)
- **046** – Accessory rails (*pair/specify head or seat section/factory installed*)
- **450** – Two function foot control (*substitution for hand control/64-*)
- **451** – Simultaneous two function foot and hand control (*69-*)
- **802** – Paper cutter band
- **810** – Restraint straps assembly
- **856** – 220/240 Volt (*factory installed*)
- **5005F** – Fire retardant foam

Table Specifications

- **Height:** Min. 19" (48.3cm) at seat
Max. 35" (88.9cm) at seat
- **Width:** 30" (76.2cm)
- **Length:** 70" (177.8cm) to 80" (203.2cm) with leg section extended
- **Storage drawer:** 16" D x 19" W x 4.5" H (40.6cm D x 48.3cm W x 11.4cm H)
- **Electrical:** 120 VAC, 60 HZ or 220/240 VAC, 60HZ

Colors



5050

5000 Series Bariatric Treatment Tables



Bariatric Table with Power Back



Standard Features

- Ultra comfort seamless top with antimicrobial guard and built in pillow
- All-steel heavy duty construction
- One function hand control
- Concealed paper roll holder (18" – 21") (45.7cm – 53.3cm)
- Large, full extension storage drawer on roller bearing drawer guides
- Fully articulating 10" (25.4cm) leg section
- Reinforced pull-out footstep with non-slip pad
- Electro-mechanical actuator system
- 800lb (363kg) weight capacity
- 3 year limited warranty

Optional Accessories

- **430** – One function foot control (*substitution for hand control/64-*)
- **431** – Simultaneous one function foot control and hand control (*64-*)
- **856** – 220/240 Volt (*factory installed*)
- **5050F** – Fire retardant foam

Table Specifications

- **Height:** 31" (78.7cm) at seat section
- **Width:** 30" (76.2cm)
- **Length:** 70" (177.8cm) to 80" (203.2cm) with leg section extended
- **Storage drawer:** 16" D x 27" W x 4.5" H (40.6cm D x 68.6cm W x 11.4cm H)
- **Foot step:** 13.3" L x 17.5" W x 7.8" H (40.6cm L x 68.6cm W x 11.4cm H)
- **Electrical:** 120 VAC, 60 HZ or 220/240 VAC, 60HZ

Colors



5015

5000 Series Power Procedure Tables



Podiatry Chair with Three Function Programmable Foot Control



Standard Features

- Ultra comfort seamless top with antimicrobial guard
- Power Hi-Lo
- Power backrest
- Power tilt
- Three function memory position foot control
- 24" (61cm) low access seat height
- Adjustable, sliding, removable stainless steel debris tray (5.6 Quart / 5.3L)
- Adjustable leg section
- Clear vinyl foot section cover
- Electro-mechanical actuator system
- 375lb (170.1kg) weight capacity
- 3 year limited warranty

Optional Accessories

- **143** – Easy-glide locking casters (*set of 4/ factory installed*)
- **457** – Simultaneous non-programmable hand control
- **852** – Hospital grade receptacle 110V (*factory installed*)
- **855** – Hospital grade receptacle 220V (*factory installed*)
- **856** – 220/240 Volt (*factory installed*)
- **961** – Mayo tray assembly
- **963** – Light post assembly, 35 watt halogen
- **964** – Vision screen with accessory mounting bracket
- **969** – Complete accessory package (*includes 961, 963, 964*)
- **5015F** – Fire retardant foam

Table Specifications

- **Height:** Min. 24" (61cm) at seat
Max. 36" (91.4cm) at seat
- **Width:** 26.5" (67.3cm)
- **Length upright:** 50" to 59" (127cm to 149.9cm)
- **Length flat:** 74" to 83" (188cm to 210.8cm)
- **Foot section:** 21" W x 15" L (53.3cm x 38.1cm)
- **Debris Tray:** 9" W x 15" L x 2.5"D (22.9cm W x 38.1cm L x 6.4cm D)
- **Electrical:** 120 VAC, 60 HZ or 220/240VAC, 60HZ

Colors



5016

5000 Series Power Procedure Tables



Podiatry Chair with Three Function Programmable Foot Control



Standard Features

- Ultra comfort seamless top with antimicrobial guard
- Power Hi-Lo
- Power backrest
- Power tilt
- Three function memory position foot control
- 19" (48.3cm) low access seat height
- Adjustable, sliding, removable stainless steel debris tray (5.6 Quart / 5.3L)
- Adjustable leg section
- Clear vinyl foot section cover
- Electro-mechanical actuator system
- 600lb (272.2kg) weight capacity
- 3 year limited warranty

Optional Accessories

- **143** – Easy-glide locking casters (*set of 4/ factory installed*)
- **457** – Simultaneous non-programmable hand control
- **852** – Hospital grade receptacle 110V (*factory installed*)
- **855** – Hospital grade receptacle 220V (*factory installed*)
- **856** – 220/240 Volt (*factory installed*)
- **961** – Mayo tray assembly
- **963** – Light post assembly, 35 watt halogen
- **964** – Vision screen with accessory mounting bracket
- **969** – Complete accessory package (*includes 961, 963, 964*)
- **5016F** – Fire retardant foam

Table Specifications

- **Height:** Min. 19" (48.3cm) at seat
Max. 35" (89cm) at seat
- **Width:** 26.5" (67.3cm)
- **Length upright:** 50" to 59" (127cm to 149.9cm)
- **Length flat:** 74" to 83" (188cm to 210.8cm)
- **Foot section:** 21" W x 15" L (53.3cm x 38.1cm)
- **Debris Tray:** 9" W x 15" L x 2.5"D (22.9cm W x 38.1cm L x 6.4cm D)
- **Electrical:** 120 VAC, 60 HZ or 220/240VAC, 60HZ

Colors



8677

8000 Series Power Phlebotomy Chairs



Phlebotomy Chair with Hi-Lo and Power Back



Standard Features

- Ultra Comfort Seamless Top with antimicrobial guard
- 24" low access height (61cm)
- Power Hi-Lo
- Power backrest
- Two function hand control
- Self-positioning armrests
- Self-positioning leg section
- Reclines to 30° from horizontal position
- Electro-mechanical actuator system
- Available in all UMF Medical upholstery color options
- 375lb (170.1kg) patient weight capacity
- 3 year limited warranty

Optional Accessories

- **032** – Optional backrest with articulating headrest (*factory installed*)
- **033** – Fold-away table (13.5" W x 10" L x .75" D) (34.3cm W x 25.4cm L x 1.9cm H)
- **143** – Easy-Glide locking casters (*set of 4/factory installed*)
- **420** – Two function foot control (*substitution for hand control*)
- **421** – Simultaneous two function foot and hand controls
- **856** – 220/240 Volt (*factory installed*)
- **8677F** – Fire retardant foam

Table Specifications

- **Seat Height:**
Min. 24" (61cm) at seat
Max. 36" (91.4cm) at seat
- **Overall Width:** 31.5" (80cm)
- **Seat Width:** 18" (45.7cm)
- **Length (fully extended):** 67" (170.2cm)
- **Arm Length:** 23" (58.4cm)
- **Electrical:** 120 VAC, 60 HZ or 220/240 VAC, 60HZ
- **Backrest Width:** 17" Max / 12" Min. (43.18cm Max/ 30.48cm Min.)

Colors



8678

8000 Series Power Phlebotomy Chairs



Phlebotomy Chair with Hi-Lo and Power Back



Standard Features

- Ultra Comfort Seamless Top with antimicrobial guard
- 19" (48.3cm) ADA compliant low access height
- Power Hi-Lo
- Power backrest
- Two function hand control with OneTouch Patient Positioning™
- Self-positioning armrests
- Self-positioning leg section
- Reclines to 30° from horizontal position
- Electro-mechanical actuator system
- Available in all UMF Medical upholstery color options
- 375lb (170.1kg) patient weight capacity
- 3 year limited warranty

Optional Accessories

- **032** – Optional backrest with articulating headrest (*factory installed*)
- **033** – Fold-away table (13.5" W x 10" L x .75" D) (34.3cm W x 25.4cm L x 1.9cm H)
- **143** – Easy-Glide locking casters (*set of 4/factory installed*)
- **421** – Simultaneous two function foot and hand controls
- **856** – 220/240 Volt (*factory installed*)
- **8678F** – Fire retardant foam

Table Specifications

- **Seat Height:**
Min. 19" (48.3cm) at seat
Max. 36" (91.4cm) at seat
- **Overall Width:** 31.5" (80cm)
- **Seat Width:** 18" (45.7cm)
- **Length (fully extended):** 67" (170.2cm)
- **Arm Length:** 23" (58.4cm)
- **Electrical:** 120 VAC, 60 HZ or 220/240 VAC, 60HZ
- **Backrest Width:** 17" Max / 12" Min. (43.18cm Max/ 30.48cm Min.)

Colors



5240

5200 Signature Series Manual Exam Tables



Manual Exam Table with Four Storage Drawers



Standard Features

- 500lb capacity (226.8kg)
- Removable ultra-comfort seamless top
- 5 year limited warranty
- Galvanized steel base, 18 gauge, powder coated
- Easy-reach articulating back adjustment
- One-touch, easy glide stirrups with 4 lateral positions and friction lock length adjustment
- Slide out leg section with removable upholstered leg rest
- Removable 9"x9" plastic drip pan (22.9cm x 22.9cm)
- Pull-out all steel non slip step 7" high (17.8cm)
- 21" Paper roll holder (53.3cm)
- 2 front drawers on full extension roller bearing glides (18.25"x15.5"x4" each or 46.4cm x 39.4cm x 10.2cm)
- 2 reversible side storage drawers on full extension roller bearing glides (18.25"x15.5"x4" each or 46.4cm x 39.4cm x 10.2cm)

Optional Accessories

- **010** – Leg extension booster pad
- **037** – Armboard assembly (*specify patient right or left side / factory installed*)
- **045** – Accessory rail (*pair / specify head or seat section / factory installed*)
- **251** – Bierhoff knee crutches (*pair*)
- **350** – Pelvic tilt
- **365** – I.V. pole (*factory installed*)
- **375** – Safety rails assembly (*pair / factory installed*)
- **463** – Light bracket capabilities (DOES NOT INCLUDE ACTUAL BRACKET)
- **701** – Removable standard premium top
- **701F** – Removable standard premium top with CAL-133 foam
- **810** – Restraint straps assembly (*factory installed*)
- **820** – Stainless steel drain pan (*factory installed*)
- **822** – Stainless steel drain pan with spout (*factory installed*)
- **849** – Drawer warmer 220V (*factory installed*)
- **851** – Drawer warmer 110V (*factory installed*)
- **852** – Hospital grade receptacle 110V with child safe cover (*factory installed / patient left side*)
- **855** – Hospital grade receptacle 220V with child safe cover (*factory installed / patient left side*)
- **5240F** – Fire retardant foam
- **5240B** – Base only
- SilverSan antimicrobial powder coated base

Table Specifications

- **Height:** 31.8" (80.8cm)
- **Length:** 58.5" – 77.5" (148.6cm – 196.9cm)
- **Width:** 27" (68.6cm)
- **Stirrups:** Infinitely adjustable

Colors



5240-145

5200 Signature Series Manual Exam Tables Exam Tables



Manual Exam Table with Side Step and Four Storage Drawers



Standard Features

- 500lb capacity (226.8kg)
- Removable ultra-comfort seamless top
- 5 year limited warranty
- Galvanized steel base, 18 gauge, powder coated
- Easy-reach articulating back adjustment
- One touch, easy-glide stirrups with 4 lateral positions and friction lock length adjustment
- Slide out leg section with removable upholstered leg rest
- Removable 9" by 9" plastic drip pan (22.9cm x 22.9cm)
- Pull out steel non-slip side step 7" high (17.8cm)
- 21" paper roll holder (53.3cm)
- 2 front drawers on full extension roller bearing glides (18.25"x15.5"x4" each or 46.4cm x 39.4cm x 10.2cm)
- 2 reversible side storage drawers on full extension roller bearing glides (18.25"x15.5"x4" each or 46.4cm x 39.4cm x 10.2cm)

Optional Accessories

- 010 – Leg extension booster pad
- 307 – Armboard assembly (*specify patient left or right side / factory installed*)
- 045 – Accessory rail (*pair / specify head or seat section / factory installed*)
- 251 – Bierhoff knee crutches (*pair*)
- 350 – Pelvic tilt
- 365 – I.V. pole (*factory installed*)
- 375 – Safety rails assembly (*pair / factory installed*)
- 463 – Light bracket capabilities (DOES NOT INCLUDE ACTUAL BRACKET)
- 701 – Removable standard premium top
- 701F – Removable standard premium top with CAL-133 foam
- 810 – Restraint straps assembly (*factory installed*)
- 820 – Stainless steel drain pan (*factory installed*)
- 822 – Stainless steel drain pan with spout (*factory installed*)
- 849 – Drawer warmer 220V (*factory installed*)
- 851 – Drawer warmer 110V (*factory installed*)
- 852 – Hospital grade receptacle 110V with child safe cover (*factory installed / patient left side*)
- 855 – Hospital grade receptacle 220V with child safe cover (*factory installed / patient left side*)
- 5240-145F – Fire retardant foam
- 5240-145B – Base only
- SilverSan antimicrobial powder coated base

Table Specifications

- **Height:** 31.8" (80.8cm)
- **Length:** 58.5" – 77.5" (148.6cm – 196.9cm)
- **Width:** 27" (68.6cm)
- **Stirrups:** Infinitely adjustable

Colors



5250

5200 Signature Series Manual Exam Tables



Manual Exam Table with Five Drawers



Standard Features

- 500lb capacity (226.8kg)
- Removable ultra-comfort seamless top
- 5 year limited warranty
- Galvanized 18 gauge, powder coated steel base
- Easy-reach articulating back adjustment
- One touch, easy-glide stirrups with 4 lateral positions and friction lock length adjustment
- Slide out leg section with removable upholstered leg rest
- Removable 9"x9" plastic drip pan (22.9cm x 22.9cm)
- Pull out all steel non-slip step 7" high (17.8cm)
- 21" Paper roll holder (53.3cm)
- 2 front drawers on full extension roller bearing glides (18.25"x15.5"x4" each or 46.4cm x 39.4cm x 10.2cm)
- 3 reversible side storage drawers on full extension roller bearing glides (18.25"x15.5"x4" each or 46.4cm x 39.4cm x 10.2cm)

Optional Accessories

- **010** – Leg extension booster pad
- **037** – Armboard assembly (*specify patient right or left side / factory installed*)
- **045** – Accessory rail (*pair / specify head or seat section / factory installed*)
- **251** – Bierhoff knee crutches (*pair*)
- **350** – Pelvic tilt
- **365** – I.V. pole (*factory installed*)
- **375** – Safety rails assembly (*pair / factory installed*)
- **463** – Light Bracket Capabilities (DOES NOT INCLUDE ACTUAL BRACKET)
- **701**– Removable standard premium top
- **701F** – Removable standard premium top with CAL-133 foam
- **810** – Restraint straps assembly (*factory installed*)
- **820** – Stainless steel drain pan (*factory installed*)
- **822** – Stainless steel drain pan with spout (*factory installed*)
- **849** – Drawer warmer 220V (*factory installed*)
- **851** – Drawer warmer 110V (*factory installed*)
- **852** – Hospital grade receptacle 110V with child safe cover (*factory installed / patient left side*)
- **855** – Hospital grade receptacle 220V with child safe cover (*factory installed / patient left side*)
- **5250F** – Fire retardant foam
- **5250B** – Base only
- SilverSan antimicrobial powder coated base

Table Specifications

- **Height:** 31.8" (80.8cm)
- **Length:** 57.5" - 74" (146.1cm - 188cm)
- **Width:** 27" (68.6cm)

Colors



5900

5900 Series Pediatric Tables



Pediatric Examination Table



Standard Features

- Ultra comfort seamless top with antimicrobial guard
- Heavy duty, all-metal storage drawers
- Heavy duty roller bearing drawer glides
- All metal cabinet doors with commercial grade hinges
- Padded front rail for infant safety
- 18" (45.7cm) paper roll holder
- Digital scale operated by 6 D-cell batteries (not included) with auto shut-off feature
- SECA® stainless steel measuring rule with paddles
- Weighs in both pounds/ounces and kilograms
- Tare weight capability
- 3 year limited capability

Optional Accessories

- **091** – KG scale only
- **092** – Pounds scale only
- **365** – I.V. pole for 5200 series tables (*factory installed only*)
- **502** – Double door lock (*factory installed only*)
- **505** – Single drawer lock (*factory installed only*)

Table Specifications

- **Overall height:** 44.5" (113cm)
- **Work surface height:** 35" (88.9cm)
- **Width:** 48" (121.9cm)
- **Depth:** 25" (63.5cm)
- **Drawers:** 16.5" W x 17" D x 4" H (42cm x 43cm x 10cm)
- **Overall door opening:** 40.25" W x 17.25" H (102cm x 44cm)
- **Weight capacity:** 50lb (22.7kg)

Colors



5585

5500 Series Treatment Table



H-Brace Treatment Table with Two Drawers



Standard Features

- 350lb (158.8kg) weight capacity
- Ultra comfort seamless top with antimicrobial guard
- 5 year limited warranty
- All-welded steel construction
- Available in all upholstery color options
- Two all metal full extension storage drawers
- Leveling glides

Optional Accessories

- **800** – Paper cutter band
- **804** – Paper rod holder
- **5585F** – Fire retardant foam

Table Specifications

- **Height:** 31.5" (80cm)
- **Length:** 70" (177.8cm)
- **Width:** 26.5" (67.3cm)
- **Storage drawers:** 27" W x 16" D x 4.2" H (67cm x 41cm x 11cm)

Colors



5588

5500 Series Treatment Table



H-Brace Treatment Table



Standard Features

- 350lb (158.8kg) weight capacity
- Ultra comfort seamless top with antimicrobial guard
- 5 year limited warranty
- All-welded steel construction
- Available in all upholstery color options
- Leveling glides

Optional Accessories

- **800** – Paper cutter band
- **804** – Paper rod holder
- **5585F** – Fire retardant foam

Table Specifications

- **Height:** 31.5" (80cm)
- **Length:** 70" (177.8cm)
- **Width:** 26.5" (67.3cm)

Colors



5570

5500 Series Treatment Tables



Treatment Table with Adjustable Backrest



*Shown with optional accessory 015 – Foot section and pad.

Standard Features

- 500lb (226.8kg) weight capacity
- Ultra comfort seamless top with antimicrobial guard
- Manual, Adjustable Backrest
- Available in all upholstery color options
- All-welded steel construction
- Leveling glides

Optional Accessories

- **015** – Foot Section and pad
- **800** – Paper cutter band
- **804** – Paper rod holder
- **826** – Wood laminate storage shelf
- **5570F**– Fire retardant foam

Table Specifications

- **Height:** 31.5" (80cm)
- **Length:** 70"- 81" with leg section (177.8cm - 205.74cm)
- **Width:** 26.5" (67.3cm)

Colors



8612

8000 Series Phlebotomy / ENT Chairs



ENT Chair with Foot Operated Pump



Standard Features

- 250lb (113.4kg) patient weight capacity
- Ultra comfort seamless top with antimicrobial guard
- 3 year limited warranty
- Sturdy steel frame
- Available in all UMF Medical upholstery options
- 360° rotation with swivel lock
- 10° elevation from horizontal position
- Manual foot pedal height adjustment
- Manually adjustable backrest
- Manually adjustable footrest
- Hinged arms

Optional Accessories

- **031** – Optional backrest with articulating headrest

Table Specifications

- **Min. Seat Height:** 23" (58.4cm)
- **Max. Seat Height:** 38" (96.5cm)
- **Overall Width:** 31.5" (80cm)
- **Seat Width:** 18" (45.72cm)
- **Arm Length:** 23" (58.4cm)
- **Length:** 71" (180.3cm) fully extended

Colors



8675

8000 Series Phlebotomy / ENT Chairs



Phlebotomy Chair with Foot Operated Pump



Standard Features

- 250lb (113.4kg) patient weight capacity
- Ultra comfort seamless top with antimicrobial guard
- 3 year limited warranty
- Sturdy steel frame
- Available in all UMF Medical upholstery options
- 360 degree rotation with swivel lock
- 10° elevation from horizontal position
- Manual foot pedal height adjustment
- Interlinked backrest and leg section for reclining patient to a horizontal position
- Hinged arms

Table Specifications

- **Length:** 67" (170.2cm) fully extended
- **Min. Seat Height:** 23" (58.4cm)
- **Max. Seat Height:** 38" (96.5cm)
- **Overall Width:** 31.5" (80cm)
- **Seat Width:** 18" (45.72cm)
- **Arm Length:** 23" (58.4cm)

Colors



8690

8600 Series Blood Draw Chairs



Bariatric and Parent with Child Phlebotomy Chair



*Photo shown with Side tray accessory #120

Standard Features

- Ultra comfort seamless upholstery with antimicrobial guard
- Accommodates bariatric patients and parent with child
- Available in all upholstery color options
- Padded, flip-up, height adjustable arm board
- Heavy duty, all-welded steel frame
- 750lb (340.2kg) patient weight capacity
- Adjustable leg levelers
- 3 year limited warranty

Optional Accessories

- 120 – Side tray

Table Specifications

- **Width:** 48.5" (123cm)
- **Depth:** 27.75" (70cm)
- **Height:** 39.875" (101cm)

Seat dimensions:

- **Height:** 21.25" (54cm)
- **Width:** 26.625" (68cm)
- **Between arm boards:** 30" (76cm)

Colors

