



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* [www.harloff.com](http://www.harloff.com)

## Classic Line, 2 Drawer with Lower Storage Compartment, Treatment Cart, Standard Package #6200

### Specifications



#### Construction:

- Sturdy, uni-body steel construction

#### Features:

- Tall cabinet
- 21" Storage compartment, locking doors and adjustable shelf
- Key lock with 4 keys
- 5" (12.7cm) ball-bearing premium casters, all full swivel one with brake, one tracking
- Full wrap-around padded bumper on separate steel frame
- Aluminum push handle / side rail
- Ball-bearing full extension drawer slides
- Cart is delivered fully assembled
- Removable plastic top
- Pre-drilled threaded holes for side mounted accessories
- Latex free
- Twelve year warranty

#### Drawer Configuration:

- (1) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- Storage Compartment: 21-inch (53.3 cm) high opening, with adjustable internal shelf

#### Internal Drawer Dimensions:

- 20" W x 17" D (50.8 cm W x 43.2 cm D)

#### Assembled Cart Dimensions

- 41.5" H x 32" W x 22" D (105.4 cm H x 81.3 cm W x 55.9 cm D)
- Cart Net Weight: 115 pounds (52.2 kg)
- Shipping dimensions: 39"D x 28"W x 48"H
- Approximate shipping weight: 130 pounds (59 kg)

#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* [www.harloff.com](http://www.harloff.com)

**Classic Line, Six Drawer Workstation (Tall) with Key Lock,  
Standard Package  
#6450**

**Specifications**



**Color Shown is Beige.  
Additional Colors Available**

**Construction:**

- Sturdy, uni-body steel construction

**Features:**

- Key lock with two keys
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Ball-bearing full extension drawer slides
- 5" ball-bearing casters, all full swivel - one directional, one with brake
- Removable plastic top
- Pre-drilled threaded holes for side mounted accessories
- Two aluminum side rails / push handles
- Stabilizing pontoon bumpers
- Twelve year warranty on manufactured parts
- Latex free
- Durable, hard-baked powder coat finish

**Drawer Configuration:**

- (4) – 3"
- (1) – 6"
- (1) – 12"

**Internal Drawer Dimensions:**

- 23" W x 17" D (58.42 cm W x 43.18 cm D)

**Assembled Cart Dimensions:**

- 42.5" H x 32" W x 22" D  
(107.95 cm H x 81.3 cm W x 55.9 cm D)
- Cart net weight: 129 lbs (58.5 kg)
- Approximate shipping weight: 170 lbs (77.1 kg)
- Approximate shipping dims: 53" H x 33" W x 26" D  
(134.6 cm H x 83.8 cm W x 66 cm D)



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## 6 Drawer Workstation (Tall) with Basic Electronic Pushbutton Lock, Standard Package #7450E

### Specifications



Color Shown is Beige.  
Additional Colors Available

#### Construction:

- Sturdy, uni-body steel construction

#### Features:

- Basic electronic pushbutton lock
- Key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Ball-bearing full extension drawer slides
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Removable plastic top
- Pre-drilled threaded holes for side mounted accessories
- Two aluminum side rails / push handles
- Stabilizing pontoon bumpers
- Latex free
- Durable, hard-baked powder coat finish
- Twelve year warranty on manufactured parts

#### Drawer Configuration:

- (4) – 3" (7.6 cm)
- (1) – 6" (15.2 cm)
- (1) – 12" (30.5 cm)

#### Internal Drawer Dimensions:

- 23" W x 17" D (58.42 cm W x 43.18 cm D)

#### Assembled Cart Dimensions:

- 42.5" H x 32" W x 22" D  
(108 cm H x 81.3 cm W x 55.9 cm D)
- Cart net weight: 131 lbs (59.4 kg)
- Approximate shipping weight: 168 lbs (76.2 kg)
- Approximate shipping dims: 53" H x 33" W x 26" D  
(134.6 cm H x 83.8 cm W x 66 cm D)



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* [www.harloff.com](http://www.harloff.com)

## Long Term Care Treatment Cart – Standard Package #6220

### Specifications

#### Construction:

- Painted steel

#### Features:

- Tall cabinet
- Dual flush-folding full height utility doors with slam lock feature
- Corner bumpers
- 5" (12.7cm) sealed ball-bearing premium casters, all full swivel – one with brake, one directional
- 20 standard size storage bins – 5.5" W x 9.5" D x 2.5" H (14 cm x 24.1 cm x 6.4 cm) located on 5 fixed shelves (4 bins per shelf)
- One roll-out shelf at bottom with deep storage tub
- Tub dims: 22.38" W x 17.38" D x 8" H (56.9 cm x 44.2 cm x 20.4 cm)
- Removable plastic top
- Cart is delivered fully assembled
- Certified "Latex Free"
- Hard-baked powder coat finish
- Five year warranty



#### Thirteen available color choices:

Sand, Gloss Beige, Burgundy,  
Light Blue, Navy, Light Gray,  
Mauve, Red, Teal, Yellow,  
Hammer Tone Blue, Hammer Tone Gray  
and Hammer Tone Green

#### Assembled Cart Dimensions:

- 41.5" H x 32" W x 22" D  
(105.4 cm x 81.3 cm x 55.9 cm)
- Cart Net Weight: 147 lbs (66.9 kg)
- Approx. Shipping Weight: 182 lbs (82.8 kg)



650 Ford Street  
Colorado Springs, CO 80915  
Toll: (800)433-4064 \* Ph: (719)637-0300  
Fax: (719)597-8273 \* www.harloff.com

## Bin Treatment Cart with Locking Doors, Specialty Package #6221



### Specifications

#### Construction:

- Sturdy, painted steel construction

#### Features:

- Tall cabinet
- Dual, flush-folding, full height utility doors with slam lock feature
- Corner bumpers
- 5" (12.7cm) ball-bearing premium casters, all full swivel – one with brake, one directional
- 16 standard size storage bins – 5.6" W x 12" L x 2.5" H (14.2 cm x 30.5 cm x 6.4 cm) located on 4 fixed shelves. (4 bins per shelf)
- 4 long size storage bins – 5.6" W x 18" L x 2.5" H (14.2 cm x 45.7 cm x 6.4cm) located on 1 fixed shelf
- One roll-out shelf at bottom with deep storage tub
- Tub dimensions: 22.38" W x 17.38" D x 8" H (56.9 cm x 44.2 cm x 20.4 cm)
- Removable plastic top
- Cart is delivered fully assembled (sharps and waste containers are customer installed)
- Latex free
- Hard-baked powder coat finish
- Twelve year warranty

#### Drawer Configuration:

- 4 ea, shelves with 4 standard bins (16 total of bin 40432: 5.6" W x 12"L)
- 1 shelf with 4 long bins (4 total of bin 40433: 5.6" W x 18" L)
- Roll out bottom shelf with storage tub

#### Assembled Cart Dimensions

- 41.5" H x 28.5 W x 19.9" D (105.4 cm x 72.4 cm x 50.5 cm)
- Cart net weight: 147 lbs (66.9 kg)
- Approximate shipping weight: 182 lbs (82.8 kg)

#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

#### Accessories Included:

- Stainless steel top rail / push handle
- 3 Gallon sharps container and bracket
- 2 Gallon waste container

#### Internal Drawer Dimensions:

- 20" W x 17" D (50.8 cm W x 43.2 cm D)





650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## V-Series, Two Drawer Procedure Cart with Locking Storage Compartment #V30-2K

Harloff's V-Series Carts are value-priced for use in a variety of hospital departments. This drawer configuration is often a popular choice for treatment, procedure and isolation applications.

### Specifications

#### Construction:

- Steel construction

#### Features:

- Key locking drawers and storage compartment; includes four keys
- Aluminum side rail / push handle
- Tall cabinet
- 21" storage compartment with locking doors and adjustable shelf
- Stabilizing bumper frame with extended front casters and (4) corner bumpers
- 5" (12.7 cm) ball-bearing casters, all full swivel - two with brake
- Pre-drilled threaded holes for V-Series side mounted accessories
- Ball-bearing full extension drawer slides
- Twelve year manufacturer's warranty
- Ships with wheels installed and ready to use
- V-Series accessories are customer installed
- Latex free
- Durable, hard-baked powder coat finish

#### Internal Drawer Dimensions:

- 23" W x 16.75" D (58.4 cm W x 42.5 cm D)

#### Assembled Cart Dimensions

- 40" H x 29.5" W x 22" D (101.6cm H x 74.9 W x 55.9 D)
- Cart weight: 90 lbs (40.8 kg)
- Approx shipping weight: 113 pounds (51.2 kg)

#### Drawer Configuration:

- (1) – 3" (7.6 cm)
  - (1) – 6" (15.2 cm)
  - Storage compartment – 21" (53.3 cm)
- Integrated Pull Handle on all drawers



#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## V-Series, Six Drawer Procedure Cart with Key Lock #V30-6K

Harloff's V-Series Carts are value-priced for use in a variety of hospital departments. This drawer configuration is often a popular choice for treatment, procedure and isolation applications.

### Specifications

#### Construction:

- Steel construction

#### Features:

- Key lock with two keys
- Aluminum side rail / push handle
- 30" of vertical drawer space
- Stabilizing bumper frame with extended front casters and (4) corner bumpers
- 5" (12.7 cm) ball-bearing casters, all full swivel - two with brake
- Pre-drilled threaded holes for V-Series side mounted accessories
- Ball-bearing full extension drawer slides
- Twelve year manufacturer's warranty
- Ships with wheels installed and ready to use
- V-Series accessories are customer installed
- Latex free
- Durable, hard-baked powder coat finish

#### Internal Drawer Dimensions:

- 23" W x 16.75" D (58.4 cm W x 42.5 cm D)

#### Assembled Cart Dimensions

- 40" H x 29.5" W x 22" D  
(101.6 cm H x 74.9 cm W x 55.9 cm D)
- Cart weight: 118 pounds (53.5 kg)
- Approximate shipping weight: 141 pounds (63.9 kg)

#### Drawer Configuration:

- (4) – 3" (7.6 cm)
  - (1) – 6" (15.2 cm)
  - (1) – 12" (30.5 cm)
- Integrated Pull Handle on all drawers



#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* [www.harloff.com](http://www.harloff.com)

## V-Series, Five Drawer Procedure Cart with Key Lock #V24-5K

Harloff's V-Series Carts are value-priced for use in a variety of hospital departments. This drawer configuration is often a popular choice for treatment, procedure and isolation applications.

### Specifications

#### Construction:

- Steel construction

#### Features:

- Key lock with two keys
- Aluminum side rail / push handle
- 24" of vertical drawer space
- Stabilizing bumper frame with extended front casters and (4) corner bumpers
- 5" (12.7 cm) ball-bearing casters, all full swivel - two with brake
- Pre-drilled threaded holes for V-Series side mounted accessories
- Ball-bearing full extension drawer slides
- Twelve year manufacturer's warranty
- Ships with wheels installed and ready to use
- V-Series accessories are customer installed
- Latex free
- Durable, hard-baked powder coat finish

#### Internal Drawer Dimensions:

- 23" W x 16.75" D (58.4 cm W x 42.5 cm D)

#### Assembled Cart Dimensions

- 34" H x 29.5" W x 22" D  
(86.4 cm H x 74.9 cm W x 55.9 cm D)
- Cart weight: 101 pounds (45.8 kg)
- Approximate shipping weight: 123 pounds (55.8 kg)

#### Drawer Configuration:

- (3) – 3" (7.6 cm)
  - (1) – 6" (15.2 cm)
  - (1) – 9" (22.9 cm)
- Integrated Pull Handle on all drawers



#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black





650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* [www.harloff.com](http://www.harloff.com)

## Three Drawer Mini24 Instrument Cart (Tall Cabinet) with 3" Casters, Key Lock #4746K

Harloff's Mini24 Line offers a 24" wide cabinet with more internal storage space than our standard Mini Line. Select the cart best suited for your needs from our multiple cabinet heights, drawer configurations, and locking system choices



Color Shown is Beige,  
Additional Colors Available.

### Shipping Dimensions:

- 57" H x 33" W x 26" D  
(144.8 cm H x 83.8 cm W x 66 cm D)
- Approximate shipping weight: 119 pounds (54.1 kg)

## Specifications

### Construction:

- Painted steel cabinet and drawers

### Features:

- Space-saver design allows supply storage in areas with space limitations
- Cabinet with 24 vertical inches (61 cm) drawer space
- Key lock with two keys
- 3" (7.6 cm) single-wheel stem caster set, full swivel – two with brakes
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year manufacturer's warranty
- Cart requires caster installation and accessory assembly after delivery
- Durable hard-baked powder coat finish

### Internal Drawer Dimensions:

- 18.5" W x 16.75" D (47 cm x 42.5cm)

### Drawer Configuration:

- (1) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- (1) - 9 inch (22.86)
- 11 inches open space with non-adjustable bottom shelf

### Assembled Cart Dimensions:

- 43.25" H x 26.5" W x 21.5" D  
(109.855 cm H x 67.31 cm W x 54.61 cm D)
- Cart net weight: 80 pounds (36.4 kg)



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Three Drawer Mini24 Instrument Cart with 3" Casters, Basic Electronic Lock #4736E

Harloff's Mini24 Line offers a 24" wide cabinet with more internal storage space than our standard Mini Line. Select the cart best suited for your needs from our multiple cabinet heights, drawer configurations, and locking system choices.



Color Shown is Light Blue,  
Additional Colors Available

### Shipping Dimensions:

- 57" H x 33" W x 26" D  
(144.8 cm H x 83.8 cm W x 66 cm D)
- Approximate shipping weight: 113 pounds (51.4 kg)

### Specifications

#### Construction:

- Painted steel cabinet and drawers

#### Features:

- Space-saver design allows supply storage in areas with space limitations
- Cabinet has 24 vertical inches (61 cm) drawer space
- Basic electronic pushbutton lock with keylock override
- 3" (7.6 cm) single-wheel stem caster set, full swivel – two with brakes
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year manufacturer's warranty
- Cart requires caster installation and accessory assembly after delivery
- Durable hard-baked powder coat finish available in seventeen colors

#### Internal Drawer Dimensions:

- 18.5" W x 16.75" D (47 cm x 42.5cm)

#### Assembled cart dimensions:

- 43.25" H x 26.5" W x 21.5" D  
(109.855 cm H x 67.31 cm W x 54.61 cm D)
- Cart net weight: 74 pounds (33.6 kg)

#### Drawer Configuration:

- (2) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- 17 inches open space with non-adjustable bottom shelf



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* [www.harloff.com](http://www.harloff.com)

## Three Drawer Mini24 Instrument Cart with 3" Casters, Key Lock #4736K

Harloff's Mini24 Line offers a 24" wide cabinet with more internal storage space than our standard Mini Line. Select the cart best suited for your needs from our multiple cabinet heights, drawer configurations, and locking system choices



Color Shown is Light Blue,  
Additional Colors Available.

### Shipping Dimensions:

- 57" H x 33" W x 26" D  
(144.8 cm H x 83.8 cm W x 66 cm D)
- Approximate shipping weight: 113 pounds (51.4 kg)

### Specifications

#### Construction:

- Painted steel cabinet and drawers

#### Features:

- Space-saver design allows supply storage in areas with space limitations
- Cabinet with 24 vertical inches (61 cm) drawer space
- Key lock with two keys
- 3" (7.6 cm) single-wheel stem caster set, full swivel – two with brakes
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year manufacturer's warranty
- Cart requires caster installation and accessory assembly after delivery
- Durable hard-baked powder coat finish

#### Internal Drawer Dimensions:

- 18.5" W x 16.75" D (47 cm x 42.5cm)

#### Drawer Configuration:

- (2) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- 17 inches open space with non-adjustable bottom shelf

#### Assembled Cart Dimensions:

- 43.25" H x 26.5" W x 21.5" D  
(109.855 cm H x 67.31 cm W x 54.61 cm D)
- Cart net weight: 74 pounds (33.6 kg)



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Six Drawer Mini Line Anesthesia Cart with Key Lock, Standard Package #3156K

### Specifications



#### Construction:

- Painted steel cabinet and drawers

#### Features:

- Space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Cabinet with 30 vertical inches (76.2 cm) drawer space
- Key lock with two keys
- 3" (7.6 cm) single-wheel stem caster set, full swivel – two with brakes
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Cart requires caster installation and accessory assembly after delivery
- Durable hard-baked powder coat finish

#### Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (2) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

#### Internal Drawer Dimensions:

- 12.6" W x 16.75" D (32 cm x 42.5cm)

#### Assembled Cart Dimensions

- 40.75" H x 18" W x 18"D  
(103.5 cm H x 45.7 cm W x 45.7 cm D)
- Cart net weight: 87 pounds (39.5 kg)
- Approximate shipping weight: 107 pounds (48.5 kg)

#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Five Drawer Mini Line Anesthesia Cart with Key Lock, Standard Package #3145K

### Specifications

#### Construction:

- Painted steel cabinet and drawers

#### Features:

- Space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Cabinet with 24 vertical inches (61 cm) drawer space
- Key lock with two keys
- 3" (7.6 cm) single-wheel stem caster set, full swivel – two with brakes
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Cart requires caster installation and accessory assembly after delivery
- Durable hard-baked powder coat finish

#### Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

#### Internal Drawer Dimensions:

- 12.6" W x 16.75" D (32 cm x 42.5cm)

#### Assembled Cart Dimensions

- 34.4" H x 18" W x 18"D  
(87.4 cm H x 45.7 cm W x 45.7 cm D)
- Cart net weight: 74 pounds (33.6 kg)
- Approximate shipping weight: 94 pounds (42.5 kg)



#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black





650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## **Six Drawer Mini Line Anesthesia Cart with Extended Stabilizing Bumper 5" Wheels, with Key Lock #3256K**

Harloff's Mini Line of carts offer a narrow design with reduced footprint. Suitable for hospitals, clinics and doctor offices with both space and budget constraints

### **Specifications**

#### **Construction:**

- Painted steel cabinet and drawers

#### **Features:**

- Space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Cabinet with 30 vertical inches (76.2 cm) drawer space
- Key lock with two keys
- Bumper system extends wheel base and stabilizes cart
- 5" ball-bearing full swivel casters – one brake, one directional
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish

#### **Drawer Configuration:**

- (3) - 3-inch (7.6 cm)
- (2) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

#### **Internal Drawer Dimensions:**

- 12.6" W x 16.75" D (32cm x 42.5cm)

#### **Assembled Cart Dimensions:**

- 43.75" H x 19.5" W x 25" D  
(110.5 cm H x 49.5 cm W x 63.5 cm D)
- Cart net weight: 97 pounds (44 kg)
- Approximate shipping weight: 107 pounds (48.5 kg)



#### **Sixteen available color choices:**

##### **\*Antimicrobial Finish Available**

Sand, Gloss Beige, Burgundy, Light Blue, Navy, Light Gray, Mauve, Red, Teal, Yellow, Hammer Tone Blue, Hammer Tone Gray, Hammer Tone Green, Hammer Tone Red, Hammer Tone Black and White



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Five Drawer Mini Line Anesthesia Cart with Extended Stabilizing Bumper and 5" Wheels, with Key Lock #3245K

Harloff's Mini Line of carts offer a narrow design with reduced footprint. Suitable for hospitals, clinics and doctor offices with both space and budget constraints

### Specifications

#### Construction:

- Painted steel cabinet and drawers

#### Features:

- Space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Cabinet with 24 vertical inches (61 cm) drawer space
- Key lock with two keys
- Bumper system extends wheel base and stabilizes cart
- 5" ball-bearing full swivel casters – one brake, one directional
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish

#### Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

#### Internal Drawer Dimensions:

- 12.6" W x 16.75" D (32cm x 42.5cm)

#### Assembled Cart Dimensions:

- 37" H x 19.5" W x 25" D  
(94.0 cm H x 49.5 cm W x 63.5 cm D)
- Cart net weight: 84 pounds (38.1 kg)
- Approximate shipping weight: 104 pounds (47.1 kg)



#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Sand, Gloss Beige, Burgundy, Light Blue, Navy, Light Gray, Mauve, Red, Teal, Yellow, Hammer Tone Blue, Hammer Tone Gray, Hammer Tone Green, Hammer Tone Red, Hammer Tone Black and White



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* [www.harloff.com](http://www.harloff.com)

## ETC Line, Five Drawer, Treatment Cart #ETC-5



### Specifications

#### Construction:

- Sturdy, uni-body steel construction

#### Features:

- Low height make these carts ideal for bedside or under a counter
- Key lock with two keys
- 3" stem casters, two with brake
- Roller drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Hard-baked powder coat finish

#### Drawer Configuration:

- (3) – 2.75" (7.0 cm)
- (2) – 5.75" (14.6 cm)

#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

#### Assembled Cart Dimensions

- 27.25" H x 23.38" W x 15.13" D  
(69.2 cm H x 59.4 W x 38.4 D)
- Cart weight: 75 lbs (34.1 kg)
- Approx. Shipping Weight: 81 pounds (36.7 kg)

\*\* PLEASE NOTE THIS CART IS SHORTER THAN STANDARD COUNTER HEIGHTS



650 Ford Street  
Colorado Springs, CO 80915  
Toll: (800)433-4064 \* Ph: (719)637-0300  
Fax: (719)597-8273 \* [www.harloff.com](http://www.harloff.com)

## ETC Line, Four Drawer, Treatment Cart #ETC-4



### Specifications

#### Construction:

- Sturdy, uni-body steel construction

#### Features:

- Low height make these carts ideal for bedside or under a counter
- Key lock with two keys
- 3" stem casters, two with brake
- Roller drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Hard-baked powder coat finish

#### Drawer Configuration:

- (1) – 2.75" (6.98 cm)
- (3) – 5.75" (14.6 cm)

#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

#### Assembled Cart Dimensions

- 27.25" H x 23.38" W x 15.13" D  
(69.2 cm H x 59.4 W x 38.4 D)
- Cart weight: 48 lbs (21.8 kg)
- Approx. Shipping Weight: 63 pounds (28.7 kg)

\*\* PLEASE NOTE THIS CART IS SHORTER THAN STANDARD COUNTER HEIGHTS



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## ETC Line, Three Drawers with Lower Compartment, Treatment Cart #ETC-3



### Specifications

#### Construction:

- Sturdy, uni-body steel construction

#### Features:

- Low height make these carts ideal for bedside or under a counter
- Key lock with two keys
- 3" stem casters, two with brake
- Roller drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Hard-baked powder coat finish

#### Drawer Configuration:

- (3) – 2.75" (6.98 cm)
- (1) – 12" (30.5 cm) Storage compartment

#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

#### Assembled Cart Dimensions

- 27.25" H x 23.38" W x 15.13" D  
(69.2 cm H x 59.4 W x 38.4 D)
- Cart weight: 48 lbs (21.8 kg)
- Approx. Shipping Weight: 63 pounds (28.7 kg)

\*\* PLEASE NOTE THIS CART IS SHORTER THAN STANDARD COUNTER HEIGHTS





650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Five Drawer Mini Line Phlebotomy Cart Package with Key Lock #3245KPB



### Features

- Wrap-around Stabilizing Bumper
- 5-inch (12.7 cm) premium plate casters
- Cabinet has 24 vertical inches (60.96 cm) drawer space
- Key-locking system with 2 keys
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish available in 16 colors

### Accessories Included:

- 2-tier raised back rail system (MINI-RBR)
- Extended Push Handle
- Single Glove Box Holder (AL2362)
- Aluminum Waste Container with Lid (AL680436)
- Two 6-bin tilt bin organizer (683805)

### Other Mini Line Phlebotomy Cart Options:

- Locking Bemis ® Sharps Holder with 5 quart Container (684801)
- Upgrade to 3 Box Glove Holder (680422UPG)
- Chart Holder (680417)

### Specifications

#### Construction:

- Painted steel cabinet and drawers

#### Internal Drawer Dimensions:

- 15" W x 16.75" D (38.1cm x 42.5cm)

#### Assembled Cart Dimensions

- 52.3" H x 27.5" W x 21.5" D  
(133 cm H x 70 cm W x 54.6 cm D)
- Height to Work Surface – 33.8"
- Cart net weight: 108 lbs (49.1 kg)
- Approximate shipping weight: 133 pounds (28.2 kg)

#### Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Five Drawer Mini Line Phlebotomy Cart Package with Breakaway Lock #3245BPB



### Features

- Wrap-around Stabilizing Bumper
- 5-inch (12.7 cm) premium plate casters
- Cabinet has 24 vertical inches (60.96 cm) drawer space
- Breakaway Lock with 50 numbered seals
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish available in 16 colors

### Accessories Included:

- 2-tier raised back rail system (MINI-RBR)
- Extended Push Handle
- Single Glove Box Holder (AL2362)
- Aluminum Waste Container with Lid (AL680436)
- Two 6-bin tilt bin organizer (683805)

### Other Mini Line Phlebotomy Cart Options:

- Locking Bemis ® Sharps Holder with 5 quart Container (684801)
- Upgrade to 3 Box Glove Holder (680422UPG)
- Chart Holder (680417)

### Specifications

#### Construction:

- Painted steel cabinet and drawers

#### Internal Drawer Dimensions:

- 15" W x 16.75" D (38.1cm x 42.5cm)

#### Assembled Cart Dimensions

- 52.3" H x 27.5" W x 21.5" D  
(133 cm H x 70 cm W x 54.6 cm D)
- Height to Work Surface – 33.8"
- Cart net weight: 108 lbs (49.1 kg)
- Approximate shipping weight: 133 pounds (28.2 kg)

#### Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Five Drawer Mini Line Phlebotomy Cart Package with Turn Latch #3245LPB



### Features

- Wrap-around Stabilizing Bumper
- 5-inch (12.7 cm) premium plate casters
- Cabinet has 24 vertical inches (60.96 cm) drawer space
- Turn Latch to secure drawers while moving cart
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish available in 16 colors

### Accessories Included:

- 2-tier raised back rail system (MINI-RBR)
- Extended Push Handle
- Single Glove Box Holder (AL2362)
- Aluminum Waste Container with Lid (AL680436)
- Two 6-bin tilt bin organizer (683805)

### Other Mini Line Phlebotomy Cart Options:

- Locking Bemis ® Sharps Holder with 5 quart Container (684801)
- Upgrade to 3 Box Glove Holder (680422UPG)
- Chart Holder (680417)

### Specifications

#### Construction:

- Painted steel cabinet and drawers

#### Internal Drawer Dimensions:

- 15" W x 16.75" D (38.1cm x 42.5cm)

#### Assembled Cart Dimensions

- 52.3" H x 27.5" W x 21.5" D  
(133 cm H x 70 cm W x 54.6 cm D)
- Height to Work Surface – 33.8"
- Cart net weight: 108 lbs (49.1 kg)
- Approximate shipping weight: 133 pounds (28.2 kg)

#### Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Six Drawer Mini Line Phlebotomy Cart Package with Key Lock #3256KPB



### Features

- Wrap-around Stabilizing Bumper
- 5-inch (12.7 cm) premium plate casters
- Cabinet has 30 vertical inches (76.2 cm) drawer space
- Key-locking system with 2 keys
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish available in 16 colors

### Accessories Included:

- 2-tier raised back rail system (MINI-RBR)
- Extended Push Handle
- Single Glove Box Holder (AL2362)
- Aluminum Waste Container with Lid (AL680436)
- Two 6-bin tilt bin organizer (683805)

### Other Mini Line Phlebotomy Cart Options:

- Locking Bemis ® Sharps Holder with 5 quart Container (684801)
- Upgrade to 3 Box Glove Holder (680422UPG)
- Chart Holder (680417)

### Specifications

#### Construction:

- Painted steel cabinet and drawers

#### Internal Drawer Dimensions:

- 15" W x 17" D (38.1 cm W x 43.2 cm D)

#### Assembled Cart Dimensions

- 58.5" H x 27.5" W x 21.5" D  
(148.6 cm H x 70 cm W x 54.6 cm D)
- Height to Work Surface – 40"
- Cart net weight: 128 lbs (58.2 kg)
- Approximate shipping weight: 153 pounds (69.6 kg)

#### Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (2) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Six Drawer Mini Line Phlebotomy Cart Package with Breakaway Lock #3256BPB



### Features

- Wrap-around Stabilizing Bumper
- 5-inch (12.7 cm) premium plate casters
- Cabinet has 30 vertical inches (76.2 cm) drawer space
- Breakaway Lock with 50 numbered seals
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish available in 16 colors

### Accessories Included:

- 2-tier raised back rail system (MINI-RBR)
- Extended Push Handle
- Single Glove Box Holder (AL2362)
- Aluminum Waste Container with Lid (AL680436)
- Two 6-bin tilt bin organizer (683805)

### Other Mini Line Phlebotomy Cart Options:

- Locking Bemis ® Sharps Holder with 5 quart Container (684801)
- Upgrade to 3 Box Glove Holder (680422UPG)
- Chart Holder (680417)

### Specifications

#### Construction:

- Painted steel cabinet and drawers

#### Internal Drawer Dimensions:

- 15" W x 16.75" D (38.1cm x 42.5cm)

#### Assembled Cart Dimensions

- 58.5" H x 27.5" W x 21.5" D  
(148.6 cm H x 70 cm W x 54.6 cm D)
- Height to Work Surface – 40"
- Cart net weight: 128 lbs (58.2 kg)
- Approximate shipping weight: 153 pounds (69.6 kg)

#### Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (2) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black





650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Six Drawer Mini Line Phlebotomy Cart Package with Turn Latch #3256LPB



### Features

- Wrap-around Stabilizing Bumper
- 5-inch (12.7 cm) premium plate casters
- Cabinet has 30 vertical inches (76.2 cm) drawer space
- Turn Latch to secure drawers while moving cart
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish available in 16 colors

### Accessories Included:

- 2-tier raised back rail system (MINI-RBR)
- Extended Push Handle
- Single Glove Box Holder (AL2362)
- Aluminum Waste Container with Lid (AL680436)
- Two 6-bin tilt bin organizers (683805)

### Other Mini Line Phlebotomy Cart Options:

- Locking Bemis ® Sharps Holder with 5 quart Container (684801)
- Upgrade to 3 Box Glove Holder (680422)
- Chart Holder (680417)

### Specifications

#### Construction:

- Painted steel cabinet and drawers

#### Internal Drawer Dimensions:

- 15" W x 17" D (38.1 cm W x 43.2 cm D)

#### Assembled Cart Dimensions

- 58.5" H x 27.5" W x 21.5" D  
(148.6 cm H x 70 cm W x 54.6 cm D)
- Height to Work Surface – 40"
- Cart net weight: 128 lbs (58.2 kg)
- Approximate shipping weight: 153 pounds (69.6 kg)

#### Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (2) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Four Drawer Mini Line Phlebotomy Cart Package with Key Lock #3254KPB



### Features

- Wrap-around, stabilizing bumper
- 5-inch (12.7 cm) premium plate casters
- Cabinet with 30 vertical inches (76.2 cm) drawer space
- Key-locking system with 2 keys
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish available in 16 colors

### Accessories Included:

- 2-tier raised back rail system (MINI-RBR)
- Extended Push Handle
- Single Glove Box Holder (AL2362)
- Aluminum Waste Container with Lid (AL680436)
- Two 6-bin tilt bin organizer (683805)

### Other Mini Line Phlebotomy Cart Options:

- Locking Bemis ® Sharps Holder with 5 quart Container (684801)
- Upgrade to 3 Box Glove Holder (680422UPG)
- Chart Holder (680417)

### Specifications

#### Construction:

- Painted steel cabinet and drawers

#### Internal Drawer Dimensions:

- 15" W x 17" D (38.1 cm W x 43.2 cm D)

#### Assembled Cart Dimensions

- 58.5" H x 27.5" W x 21.5" D  
(148.6 cm H x 70 cm W x 54.6 cm D)
- Height to Work Surface – 40"
- Cart net weight: 128 lbs (58.2 kg)
- Approximate shipping weight: 153 pounds (69.6 kg)

#### Drawer Configuration:

- (1) - 3-inch (7.6 cm)
- (3) - 9-inch (22.9 cm)

#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Four Drawer Mini Line Phlebotomy Cart Package with Breakaway Lock #3254BPB



### Features

- Wrap-around Stabilizing Bumper
- 5-inch (12.7 cm) premium plate casters
- Cabinet has 30 vertical inches (76.2 cm) drawer space
- Breakaway Lock with 50 numbered seals
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish available in 16 colors

### Accessories Included:

- 2-tier raised back rail system (MINI-RBR)
- Extended Push Handle
- Single Glove Box Holder (AL2362)
- Aluminum Waste Container with Lid (AL680436)
- Two 6-bin tilt bin organizer (683805)

### Other Mini Line Phlebotomy Cart Options:

- Locking Bemis ® Sharps Holder with 5 quart Container (684801)
- Upgrade to 3 Box Glove Holder (680422UPG)
- Chart Holder (680417)

### Specifications

#### Construction:

- Painted steel cabinet and drawers

#### Internal Drawer Dimensions:

- 15" W x 16.75" D (38.1cm x 42.5cm)

#### Assembled Cart Dimensions

- 58.5" H x 27.5" W x 21.5" D  
(148.6 cm H x 70 cm W x 54.6 cm D)
- Height to Work Surface – 40"
- Cart net weight: 128 lbs (58.2 kg)
- Approximate shipping weight: 153 pounds (69.6 kg)

#### Drawer Configuration:

- (1) - 3-inch (7.6 cm)
- (3) - 9-inch (22.9 cm)

#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Four Drawer Mini Line Phlebotomy Cart Package with Turn Latch #3254LPB



### Features

- Wrap-around, stabilizing bumper
- 5-inch (12.7 cm) premium plate casters
- Cabinet with 30 vertical inches (76.2 cm) drawer space
- Turn Latch to secure drawers while moving cart
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish available in 16 colors

### Accessories Included:

- 2-tier raised back rail system (MINI-RBR)
- Extended Push Handle
- Aluminum Waste Container with Lid (AL680436)
- Single Glove Box Holder (AL2362)
- Two 6-bin tilt bin organizer (683805)

### Other Mini Line Phlebotomy Cart Options:

- Locking Bemis ® Sharps Holder with 5 quart Container (684801)
- Upgrade to 3 Box Glove Holder (680422UPG)
- Chart Holder (680417)

### Specifications

#### Drawer Configuration:

- (1) - 3-inch (7.6 cm)
- (3) - 9-inch (22.9 cm)

#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

#### Construction:

- Painted steel cabinet and drawers

#### Internal Drawer Dimensions:

- 15" W x 17" D (38.1 cm W x 43.2 cm D)

#### Assembled Cart Dimensions

- 58.5" H x 27.5" W x 21.5" D  
(148.6 cm H x 70 cm W x 54.6 cm D)
- Height to Work Surface – 40"
- Cart net weight: 128 lbs (58.2 kg)
- Approximate shipping weight: 153 pounds (69.6 kg)