

UNIT-DOSE BOX CARTS

AL420BOX (pictured)
OptimAL Line Unit-Dose Box Cart with lighter all aluminum construction. Includes narcotics drawer, pull-out shelf and BEST® lock on cabinet. Capacity of 420 boxes, Standard Package

SL420BOX (not shown)
Standard Line Unit-Dose Box Cart with BEST® lock on cabinet, pull-out shelf. Capacity of 420 boxes, Standard Package



AL294BOX (not shown)
OptimAL Line Unit-Dose Box Cart with lighter aluminum construction. Includes narcotics drawer, pull-out shelf and BEST® lock on cabinet. Capacity of 294 boxes, Standard Package

SL294BOX (pictured)
Standard Line Unit-Dose Box Cart with BEST® lock on cabinet, pull-out shelf. Capacity of 294 boxes, Standard Package

AL290BOXSC (not shown)
OptimAL Line Unit-Dose Box Cart with lighter all aluminum construction. Includes narcotics box in bottom drawer, pull-out shelf and BEST® lock on cabinet. Capacity of 290 boxes, Standard Package

SL290BOXSC (pictured)
Standard Line Unit-Dose Box Cart with BEST® lock on cabinet, pull-out shelf and Narcotics box in bottom drawer. Capacity of 290 boxes, Standard Package



680510E - Value Line Unit-Dose Box Divider Set
680553 - Standard and OptimAL Line Unit-Dose Box Divider Set for 20.25" wide drawers
680554 - Standard and OptimAL Line Unit-Dose Box Divider Set for 27.5" wide drawers

VL216BOX (pictured)
Value Line Unit-Dose Box Cart with key lock, internal narcotics box, pull-out shelf. Capacity of 216 boxes, Standard Package



Mechanical keyless entry available on Value Line Carts



Economy Accessory Packages



MCAP3 - Standard Line and OptimAL Line Accessory Package for 20.5" wide drawers

MCAP4 - Standard Line and OptimAL Line Accessory Package for 27.5" wide drawers

VLAP - Value Line Accessory Package for all Value Line Carts

STANDARD, VALUE & OPTIMAL LINE

VLOPT575 (pictured)
Value Line OPTI-PAK®
Compatible Cart with key lock
on cabinet internal narcotics box,
pull-out shelf. Capacity of 575
packs, Standard Package

SLOPT745 (pictured)
Standard Line OPTI-PAK® Cart
with BEST® lock on cabinet,
pull-out shelf. Capacity of 745
packs, Standard Package



OPTI-PAK Drawer

AL2397 - OPTI-PAK® Set, 12 trays /
60 bins for use in SLOPT carts (2 trays /
10 bins per drawer)

VL2397 - OPTI-PAK® Set, 48 Bins,
12 trays, 2 trays per drawer



WL294BOX-CM



WL108BOX-SO

WL108BOX-SO (pictured)
Wood look mobile medication
furniture for Unit-Dose boxes.
Shown in Solara Oak.
Includes internal narcotics box,
drawer dividers and Best® lock on
cabinet. Capacity of 108 boxes.

WL108BOX-CM (not shown)
Wood Look Unit-Dose
Box Cart in Cherry Mahogany.

WL294BOX-CM (pictured) Wood look
mobile medication furniture for Unit-Dose boxes.
Includes internal narcotics box, drawer dividers,
pull-out shelf and Best® lock on cabinet.
Capacity of 294 boxes, Standard Package
The WL294BOX is available in Cherry Mahogany only

Special order cabinet colors are available with minimum
quantity cabinet purchase and extended lead time.
Call Harloff for details



Cherry Mahogany
(CM)



Solara Oak
(SO)



Gray Granite



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Standard Line 294 Capacity Unit Dose Medication Box Cart, Standard Package #SL294BOX

Specifications



Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Dual column medication cart with 294 medication box capacity, single width narcotics drawer
- BEST® brand lock on cabinet; Illinois lock on narcotics drawer
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Pull-out side shelf
- Full wrap around vinyl bumper, protecting the cart
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish
- Offered in specialty package (#SL294BOXSP) Included are:
#68640-P1 drawer tray with adjustable dividers, #680551 eye drop tray, #AL680435 aluminum waste container and #684810 locking sharps container

Drawer Configuration:

Left Column Drawers:

- 1 ea – 4" multipurpose storage drawer
- 6 ea – 3.25" med box drawers with internal dividers, capacity 42 boxes per drawer
- 1 ea – 9.75" multipurpose storage drawer with adjustable dividers

Right Column Drawers:

- 1 ea – 4" multipurpose storage drawer
- 1 ea – 9.75" locking narcotics drawer
- 2 ea – 9.75" multipurpose storage drawer



Thirteen available color choices:

Sand, Gloss Beige, Burgundy,
Light Blue, Navy, Light Gray,
Mauve, Red, Teal, Yellow,
Hammer Tone Blue, Hammer Tone Gray
and Hammer Tone Green

Assembled Cart Dimensions:

- 45.75" H x 43" W x 23.5" D
(116.2 cm x 109.2 cm x 59.7 cm)
- Writing surface: 44" H (111.8 cm)
- Cart weight: 210 lbs (95.5 kg)
- Approx shipping weight: 255 lbs (116 kg)



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

OptimAL All Aluminum 294 Capacity Unit Dose Medication Cart, Standard Package #AL294BOX

Specifications

Construction:

- Painted aluminum cabinet and drawers

Features:

- Dual column medication cart with 294 unit dose box capacity, single width narcotics drawer
- BEST® brand lock on cabinet and Illinois lock on narcotics drawer
- Painted pull-out shelf
- Ball-bearing, full extension, self-closing drawer slides
- Plastic replaceable top with integrated push handles
- Full wrap around vinyl bumper, protecting the cart
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish
- Offered in specialty package (#AL294BOXSP) which includes #68640-P1 drawer tray with adjustable dividers, #680551 eye drop tray, #AL680436 aluminum waste container with lid and #684810 locking sharps container

Drawer Configuration:

Left Column Drawers:

- 1 ea – 4" multipurpose storage drawer
- 6 ea – 3.25" med box drawers with internal dividers, capacity 42 boxes per drawer
- 1 ea – 9.75" multipurpose storage drawer with adjustable dividers

Right Column Drawers:

- 1 ea – 4" multipurpose storage drawer
- 1 ea – 9.75" locking narcotics drawer
- 2 ea – 9.75" multipurpose storage drawer

Assembled Cart Dimensions:

- 45.75" H x 43" W x 23.5" D
(116.2 cm x 109.2 cm x 59.7 cm)
- Writing surface: 44" H (111.8 cm)
- Cart weight: 162 lbs (73.6 kg)
- Approx shipping weight: 197 lbs (89.5 kg)

Thirteen available color choices:

Sand, Gloss Beige, Burgundy,
Light Blue, Navy, Light Gray,
Mauve, Red, Teal, Yellow,
Hammer Tone Blue, Hammer Tone Gray
and Hammer Tone Green





650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Standard Line 294 Capacity Unit Dose Medication Box Cart, Standard Package #SL420BOX



Specifications

Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Dual column medication cart, single width narcotics drawer
- Capacity of 420 unit-dose boxes, 42 per drawer
- BEST® brand lock on cabinet; Illinois lock on narcotics drawer
- Painted steel pull-out side shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Full wrap around vinyl bumper, protecting the cart
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish

Offered in specialty package (#SL420BOXSP) which includes: BEST® lock on narcotics drawer, #68640-P1 drawer tray with adjustable dividers, #680551 eye drop tray, AL680435 aluminum waste container and #684810 locking sharps container

Thirteen available color choices:

Beige (-BG), Burgundy (-BU), Light Blue (-LB), Dark Blue (-DB), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammertone Dark Blue (-HTDB), Hammertone Grey (-HTGY), Hammertone Green (-HTGN) and Sand (-SD).

Assembled Cart Dimensions

- 45.75" H x 43" W x 23.5" D
(116.2 cm x 109.2 cm x 59.7 cm)
- Writing Surface: 44" (111.8 cm)
- Cart weight: 230 lbs (104.6 kg)
- Approx shipping weight: 275 lbs (125 kg)

Drawer configuration:

Left Column Drawers:

- 1 ea – 4" multipurpose storage drawer
- 9 ea – 3.25" med box drawers with internal divider set, 42 box capacity per drawer

Right Column Drawers:

- 1 ea – 4" multipurpose storage drawer
- 1 ea – 9.75" locking narcotics drawer
- 2 ea – 9.75" multipurpose storage drawer



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

OptimAL All Aluminum 420 Capacity Unit Dose Medication Cart, Standard Package #AL420BOX

Specifications

Construction:

- Painted aluminum cabinet and drawers

Features:

- Dual column medication cart with 420 unit dose box capacity, single width narcotics drawer
- BEST® brand lock on cabinet and Illinois lock on narcotics drawer
- Painted pull-out shelf
- Ball-bearing, full extension, self-closing drawer slides
- Plastic replaceable top with integrated push handles
- Full wrap around vinyl bumper, protecting the cart
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish
- Offered in specialty package (#AL420BOXSP) which includes #68640-P1 drawer tray with adjustable dividers, #680551 eye drop tray, AL680436 aluminum waste container with lid and #684810 locking sharps container

Drawer Configuration:

Left Column Drawers:

- 1 ea – 4" multipurpose storage drawer
- 9 ea – 3.25" med box drawers with internal dividers, capacity 42 boxes per drawer

Right Column Drawers:

- 1 ea – 4" multipurpose storage drawer
- 1 ea – 9.75" locking narcotics drawer
- 2 ea – 9.75" multipurpose storage drawer

Thirteen available color choices:

Sand, Gloss Beige, Burgundy,
Light Blue, Navy, Light Gray,
Mauve, Red, Teal, Yellow,
Hammer Tone Blue, Hammer Tone Gray
and Hammer Tone Green

Assembled Cart Dimensions:

- 45.75" H x 43" W x 23.5" D
(116.2 cm x 109.2 cm x 59.7 cm)
- Writing surface: 44" H (111.8 cm)
- Cart weight: 172 lbs (78.1 kg)
- Approx shipping weight: 207 lbs (94 kg)





650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Value Line Medication Cart with Key Lock

VLT216BOX

Harloff's value medication carts offer many of the same features our customers love in our large med carts at an economical price point. This economy med cart has become our fastest growing small med cart.



Color Shown is Light Gray



Thirteen available color choices:

(may be viewed at www.harloff.com)

• Beige, Burgundy, Light Blue, Dark Blue, Light Gray, Mauve, Red, Teal, Yellow, Hammertone Blue, Hammertone Gray, Hammertone Green and Sand

Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Tall cabinet with key locking system with 4 keys
- Internal locking narcotics box in bottom 9" drawer
- Capacity of up to 216 unit dose boxes or 60 sections for strip packaging
- Pull out Side shelf
- Bumper frame with projected front casters and (4) corner bumpers
- 5" ball-bearing premium casters, all full swivel – one with brake, one directional
- Painted steel top rail/push handle
- Removable plastic top
- Full extension ball-bearing drawer slides
- Cart is delivered fully assembled
- Pre-threaded accessory mounting holes
- Hard-baked powder coat finish
- Certified "Latex Free"
- Twelve year manufacturer's warranty

Drawer Configuration:

- 1 - ea 3" multipurpose storage drawer
- 6 - ea 3" drawers with internal unit dose dividers
- 1 - ea 9" drawer with 1 divider and locking narcotics box

Internal Drawer Dimensions:

- 23" W x 17" D (58.5 cm W x 43.2 cm D)

Assembled Cart Dimensions:

- 46.75" H x 29.5" W x 22.5" D (121.3 cm H x 75 cm W x 57.2 cm D)
- Writing surface: 43.25" (109.9 cm)
- Cart weight: 170 pounds (77.3 kg)
- Approx. shipping weight: 210 pounds (95.5 kg)