

SS8096

SS8000 Series Carts and Tables



Stainless Steel Mobile Utility Table with Full Perimeter Guard Rails



Standard Features

- Stainless steel construction
- Table top and shelf top are bent and hemmed for safety and strength
- Four sided guard rails on top and shelf
- 3" (7.6cm) swivel casters
- 3 year limited warranty

Optional Accessories

- 133 – Swivel casters locking 3" (7.6cm) (Set of 4)

Specifications

- **Height:** 34" (86.4cm)
- **Width:** 20" (50.8cm)
- **Depth:** 16" (40.6cm)



SS8153

SS8000 Series Carts and Tables



Stainless Steel Utility Table with One Drawer, One Shelf, and Three-sided Guard Rail



Standard Features

- Welded stainless steel construction
- Stainless steel drawer
- Three-sided guard rail on top
- Convenient storage shelf
- 2" (5.1cm) swivel casters
- 3 year limited warranty

Optional Accessories

- **117** – Swivel casters locking 2" (set of 4)
- **165** – Stainless Steel push handle
- **480** – Rubber crutch tips in place of casters (table height will be 2" (5.1cm)) lower)
- **505** – Single drawer lock (*factory installed*)

Specifications

- **Height:** 35.3" (89.5cm)
- **Width:** 20" (50.8cm)
- **Depth:** 16" (40.6cm)




SS8000

SS8000 Series Carts and Tables



Stainless Steel Instrument Table 24" x 72" (61cm x 152cm) with Shelf



Standard Features

- Welded, all stainless steel construction
- 16 gauge legs, 1.25" (3.18cm) outside diameter
- Edges hemmed for added rigidity and safety
- Top and shelf stiffeners for additional strength
- 2" (5.1cm) ball bearing swivel casters
- 3 year limited warranty

Optional Accessories

- **117** – Swivel casters locking 2" (5.1cm) (Set of 4)
- **480** – Rubber crutch tips in place of casters (table will be 2" (5.1cm)) lower

Specifications

- **Height:** 34" (86.4cm)
- **Width:** 72" (182.9cm)
- **Depth:** 24" (61cm)
- **Distance between shelves:** 19" (48.3cm)




SS8004

SS8000 Series Carts and Tables



Stainless Steel Instrument Table 24" x 60" (61cm x 152cm) with Shelf



Standard Features

- Welded, all stainless steel construction
- 16 gauge legs, 1.25" (3.18cm) outside diameter
- Edges hemmed for added rigidity and safety
- Top and shelf stiffeners for additional strength
- 2" (5.1cm) ball bearing swivel casters
- 3 year limited warranty

Optional Accessories

- **117** – Swivel casters locking 2" (5.1cm) (Set of 4)
- **480** – Rubber crutch tips in place of casters (table will be 2" (5.1cm)) lower

Specifications

- **Height:** 34" (86.4cm)
- **Width:** 60" (152.4cm)
- **Depth:** 24" (61cm)
- **Distance between shelves:** 19" (48.3cm)



SS8008

SS8000 Series Carts and Tables



Stainless Steel Instrument Table 24" x 48" (61cm x 122cm) with Shelf

Standard Features

- Welded, all stainless steel construction
- 16 gauge legs, 1.25" (3.18cm) outside diameter
- Edges hemmed for added rigidity and safety
- Top and shelf stiffeners for additional strength
- 2" (5.1cm) ball bearing swivel casters
- 3 year limited warranty



Optional Accessories

- **117** – Swivel casters locking 2" (5.1cm) (Set of 4)
- **479** – Rubber crutch tips in place of casters (table will be 2" (5.1cm)) lower

Specifications

- **Height:** 34" (86.4cm)
- **Width:** 48" (121.9cm)
- **Depth:** 24" (61cm)
- **Distance between shelves:** 19" (48.3cm)



SS8010

SS8000 Series Carts and Tables



Stainless Steel Instrument Table 20" x 36" (51cm x 91cm) with Shelf



Standard Features

- Welded, all stainless steel construction
- 16 gauge legs, 1.25" (3.18cm) outside diameter
- Edges hemmed for added rigidity and safety
- Top and shelf stiffeners for additional strength
- 2" (5.1cm) ball bearing swivel casters
- 3 year limited warranty

Optional Accessories

- **117** – Swivel casters locking 2" (5.1cm) (Set of 4)
- **479** – Rubber crutch tips in place of casters (table will be 2" (5.1cm)) lower

Specifications

- **Height:** 34" (86.4cm)
- **Width:** 36" (91.4cm)
- **Depth:** 20" (50.8cm)
- **Distance between shelves:** 19" (48.3cm)



SS8012

SS8000 Series Carts and Tables



Stainless Steel Instrument Table 18" x 33" (46cm x 84cm) with Shelf



Standard Features

- Welded, all stainless steel construction
- 16 gauge legs, 1.25" (3.18cm) outside diameter
- Edges hemmed for added rigidity and safety
- Top and shelf stiffeners for additional strength
- 2" (5.1cm) ball bearing swivel casters
- 3 year limited warranty

Optional Accessories

- **117** – Swivel casters locking 2" (5.1cm) (Set of 4)
- **479** – Rubber crutch tips in place of casters (table will be 2" (5.1cm)) lower

Specifications

- **Height:** 34" (86.4cm)
- **Width:** 33" (83.8cm)
- **Depth:** 18" (45.7cm)
- **Distance between shelves:** 19" (48.3cm)



SS8014

SS8000 Series Carts and Tables



Stainless Steel Instrument Table 20" x 48" (51cm x 122cm) with Shelf



Standard Features

- Welded, all stainless steel construction
- 16 gauge legs, 1.25" (3.18cm) outside diameter
- Edges hemmed for added rigidity and safety
- Top and shelf stiffeners for additional strength
- 2" (5.1cm) ball bearing swivel casters
- 3 year limited warranty

Optional Accessories

- **117** – Swivel casters locking 2" (5.1cm) (Set of 4)
- **479** – Rubber crutch tips in place of casters (table will be 2" (5.1cm)) lower

Specifications

- **Height:** 34" (86.4cm)
- **Width:** 48" (121.9cm)
- **Depth:** 20" (50.8cm)
- **Distance between shelves:** 19" (48.3cm)



SS8016

SS8000 Series Carts and Tables



Stainless Steel Instrument Table 16" x 30" (41cm x 76cm) with Shelf



Standard Features

- Welded, all stainless steel construction
- 16 gauge legs, 1.25" (3.18cm) outside diameter
- Edges hemmed for added rigidity and safety
- Top and shelf stiffeners for additional strength
- 2" (5.1cm) ball bearing swivel casters
- 3 year limited warranty

Optional Accessories

- **117** – Swivel casters locking 2" (5.1cm) (Set of 4)
- **479** – Rubber crutch tips in place of casters (table will be 2" (5.1cm)) lower

Specifications

- **Height:** 34" (86.4cm)
- **Width:** 30" (76.2cm)
- **Depth:** 16" (40.6cm)
- **Distance between shelves:** 19" (48.3cm)



SS8018

SS8000 Series Carts and Tables



Stainless Steel Instrument Table 16" x 20" (40.6cm x 51cm) with Shelf



Standard Features

- Welded, all stainless steel construction
- 16 gauge legs, 1.25" (3.18cm) outside diameter
- Edges hemmed for added rigidity and safety
- Top and shelf stiffeners for additional strength
- 2" (5.1cm) ball bearing swivel casters
- 3 year limited warranty

Optional Accessories

- **117** – Swivel casters locking 2" (5.1cm) (Set of 4)
- **165** – Stainless Steel push handle
- **479** – Rubber crutch tips in place of casters (table will be 2" (5.1cm) lower)

Specifications

- **Height:** 34" (86.4cm)
- **Width:** 20" (50.8cm)
- **Depth:** 16" (40.6cm)
- **Distance between shelves:** 19" (48.3cm)



SS8360

SS8300 Series Stands and Receptacles



Stainless Steel Twin Basin Stand with Shelf



Standard Features

- Heavy duty stainless steel tubular framework
- Two 13" (33cm) basins, 7 Quart (6.6L) capacity
- Rubber bumpers inside basin holder
- Single shelf below basins
- 2" (5.1cm) ball bearing swivel casters
- 3 year limited warranty

Optional Accessories

- 117 – Swivel casters locking 2" (5.1cm) (set of 4)

Specifications

- **Height:** 32" (81.3cm)
- **Basin diameter:** 13" (33cm)
- **Basin capacity (ea.):** 7 Quart (6.6L)



SS8365

SS8300 Series Stands and Receptacles



Stainless Steel Single Basin Stand with Shelf



Standard Features

- Heavy duty stainless steel tubular framework
- 13" (33cm) basin, 7 Quart (6.6L) capacity
- Rubber bumpers inside basin holder
- Single shelf below basin
- 2" (5.1cm) ball bearing swivel casters
- 3 year limited warranty

Optional Accessories

- 117 – Swivel casters locking 2" (5.1cm) (set of 4)

Specifications

- **Height:** 32" (81.3cm)
- **Basin diameter:** 13" (33cm)
- **Basin capacity:** 7 Quart (6.6L)

