

ENOCHS

EXAMINING ROOM FURNITURE

BLOOD DRAWING CHAIR



Specifications for Blood Drawing Chair and Side Chair



	201 BLOOD DRAWING CHAIR Extra Wide Design	200 SIDE CHAIR
Warranty:	4 year	4 year
Base Construction:	fully welded, steel rod base	fully welded, steel rod base
Base Finish:	chrome	chrome
Seat & Back Finish:	contoured seat and back smooth surface, polypropylene	contoured seat and back smooth surface, polypropylene
Armboards:	two padded armboards, covered in durable Naughahyde®, each armboard: 6"W x 19"L	N/A
Armboard Adjustment:	lateral: full height: 25" to 38"	N/A
Side Cabinet:	Stainless steel cabinet with one drawer, designed to mount under armboard on left or right-hand side depending on room configuration	N/A
Overall Dimensions:	43" W x 22.75" D x 33.25" H	19.75" W x 22.75" D x 33.25" H
Seating Area Dimensions:	27" W x 17" D x 17.5" H	18.25" W x 17" D x 17.5" H
Cabinet Dimensions:	12" W x 10.5" D x 2.5" H	N/A
Drawer Dimensions:	11" W x 9.5" D x 1.75" H	N/A
Patient Load on Seat:	225 lbs.	225 lbs.
Product Weight:	56 lbs.	18 lbs.
Shipping:	assembly required	no assembly required
Color Selection:	Black	Black