



Vivo

FEATURES

- Ideal for long-term care facilities, including retirement and nursing homes.
- Easy grip arm rests to help aid users out and into the seated position.
- The curved arm provides a more warm and homely aesthetic compared to a straight institutional styled arm rest.
- The backrest tilt and rising of the footrest allows for a completely horizontal resting position.
- The chair is standard with large dual wheel casters, two of them locking.
- Adjustable, removable headrest is fastened with Velcro, for easy movement and cleaning.
- Manual gas assisted recliner armchair with independent movements of the backrest and leg rest, by two gas springs.
- Mechanism is made of steel and powder coat painted black.
- The arms and frames are made of beech wood, in clear coat finish.
- All upholstered pieces are double stitched for a more durable finish and pleasing aesthetic finish.
- Upholstered and padded with CFC free, high-density and high-resiliency polyurethane foam.
- The seat and back are both webbed with Pirelli webbing for extra added comfort.
- Metal tubular steel push handle, powder coated in silver.

OPTIONS

- Closed, fully upholstered armrest.
- Foot-rest.
- Height adjustable armrest which fold away.
- Fire retardant padding.



8480F
Fixed Arms



8480A
Adjustable Arms

Max Height:	50"	1270 mm	50"	1270 mm
Max Width:	25.5"	650 mm	25.5"	650 mm
Max Depth:	28.5"	725 mm	28.5"	725 mm
Arm Height:	27"	685 mm	21-27"	535-685 mm
Back Height:	31"	790 mm	31"	790 mm
Seat Height:	21"	535 mm	21"	535 mm
Seat Width:	19"	485 mm	19"	485 mm
Seat Depth:	21"	535 mm	21"	535 mm
Wt/Chair:	80 lbs.	36.3 kg	80 lbs.	36.3 kg
Volume:	26 ft ³	0.74 m ³	26 ft ³	0.74 m ³
COM Fabric:	4 yd	4.5 m	4 yd	4.5 m
Per Carton:	1		1	

Vivo

