



mayo stands



foot-pedal operated

- Convenient foot-pedal operated lock adjusts height from 37" to 62½"
- All welded stainless steel construction
- Base designed to move under low clearance equipment with maximum stability
- 2" dual wheel swivel casters



P-1065-SS

- Removable extra large 20" x 25" stainless steel tray for easy cleaning
- 50 lbs. weight capacity evenly distributed



P-1066-SS

- Removable 16¼" x 21¼" stainless steel tray for easy cleaning
- 50 lbs. weight capacity evenly distributed



P-1068-SS

- Removable 12⁵/₈" x 19¹/₈" stainless steel tray for easy cleaning
- 50 lbs. weight capacity evenly distributed



P-1069-SS

- Removable large 16¼" x 21¼" stainless steel tray for easy cleaning
- Four 2" dual wheel swivel casters for increased mobility
- 50 lbs. weight capacity evenly distributed

hand operated

- Adjust height from 37" to 61½" - friction lock
- All welded stainless steel construction
- Base designed to move under low clearance equipment with maximum stability
- 2" dual wheel swivel casters



P-1066-A-SS

- Removable 16¼" x 21¼" stainless steel tray for easy cleaning
- 50 lbs. weight capacity evenly distributed



P-1068-A-SS

- Removable 12⁵/₈" x 19¹/₈" stainless steel tray for easy cleaning
- 50 lbs. weight capacity evenly distributed



P-1069-A-SS

- Removable large 16¼" x 21¼" stainless steel tray for easy cleaning
- 50 lbs. weight capacity evenly distributed



P-66

- Removable 19¹/₈"L x 12⁵/₈"W stainless steel tray on powder-coated frame
- Adjusts height from 32" to 53"
- Two 2" double ball bearing casters
- 50 lbs. weight capacity evenly distributed



P-65 Medical Tray Unit

- Removable 19¹/₈"L x 12⁵/₈"W stainless steel tray on chrome plated frame
- Two 2" double ball bearing casters
- 40 lbs. weight capacity evenly distributed

out of the **box** ready to use

