## Manufacturers of Quality X-Ray Protective Products

**Since 1958** 







**SHIELDING INTERNATIONAL™** has been a leading manufacturer of x-ray protective garments and accessories for the medical supplies industry for more than fifty years.

Since 1958, we have provided the most flexible, comfortable and longest wearing x-ray protective garments on the market.



- **1-7** Frontal Aprons
- 8-9 Vest Skirts
- **10-11** Wrap Around Aprons
- **12** Buckle Only Apparel
- 13 Half Aprons Demi Aprons
- 14-15 Dental
- 16-17 Pediatrics
- 18-19 Gonad Shields Thyroid Collars & Accessories
- 20-21 Gloves & Mitts
- 22-25 Eyewear

28-29 Porta-Shields, Curtains and Blockers

26-27 Apron Racks & Hangers

- 30-31 CT Shields Positioning Sandbags • Table & Panel Scatter Panel Shields
- 32-33 Mobile Barriers Windows & Frames
- 34-35 Sizing & Protection Options
- 36-39 Fabrics and Embroidery Options
- 40-41 Custom Options Orders & Shipping • Apron Care • Warranty



#### SHIELDING INTERNATIONAL<sup>™</sup> is an ISO 9001: CE, and FDA registered company.



		000
<b>Standard Protection</b> RL Regular Lead <sup>1</sup>	Benefits Most Economical 100% Recyclable	<b>Protection Value</b> Provides 0.5mm Pb Equivalent Protection @ 100 kVp 96% Attenuative
<b>Optional Protection</b> TL Lite Ply <sup>1</sup>	Benefits 20% Lighter than Regular Lead 100% Recyclable	<b>Protection Value</b> Provides 0.5mm Pb Equivalent Protection @ 100 kVp 95% Attenuative
<b>Optional Protection</b> NL Non Lead <sup>2</sup>	<b>Benefits</b> 25% Lighter than Regular Lead Lead Free & Environmentally Friendly	<b>Protection Value</b> Provides 0.5mm Pb Equivalent Protection @ 100 kVp 97% Attenuative

Α



A) Standard Coat Apron—Model SQR, ST, SV Three Closure Options: (T) Ties, (QR) Quick-Release Buckle, (V) Hook and Loop

For sizing, protection, fabrics and custom options refer to pages 34-40.

Product Conformance: In compliance with <sup>1</sup> ASTM test method F2547-06 and <sup>2</sup> IEC 61331-1:2014 core protection.

# Our Most Popular And Economical

#### **B) Economy Coat Apron—Model E** Tie Closures

• Available in our Selection of Solid Color Vinyl (V), Nylon (NY) or Diamond Ripstop (DRS) See fabrics on pages 37-38.

1

### Be Back Smart...




Features elastic back panels for added comfort and support, with hook and loop fastening flaps that wrap around waist for an adjustable fit. Shoulder pads included.



 Standard Protection: 0.5mm Pb Equivalency
 For Sizing, Protection, Fabrics and Custom Options, refer to Pages 34-40



D



### E) BRX Back Relief Apron — Model BRX

(Pictured with BRX Belt Option) The BRX Belt is a brace belt with a unique bionic system designed to help with lower back stress or strain by providing incredible support and comfort. The pulley system cinches the belt. The system and cords are backed by a lifetime warranty. The BRX belt can be adjusted for various body positions and is comfortable when worn outside our x-ray apparel (may also be used with other apron styles).\*

#### D) Back Relief Apron—Model BR

Back Relief Apron transfers the weight from the shoulders to the hips with the aid of a 6" wide elastic BR belt or for added support order Model BRX with brace belt. Adjustable back shoulder strap keeps apron snug. Shoulder pads included.\*

\* Indicate proper belt size from size chart when ordering.BR and BRX belt may also be ordered separately.

### **Frontal Aprons**

#### BRX Belt (E)

Belt Model No.	Size	Fits Waist
BRX-S	Small	25-30" (63.5-76.2cm)
BRX-M	Medium	30-35″ (76.2-88.9cm)
BRX-L	Large	35-40" (88.9-101.6cm)
BRX-XL	XLarge	40-45" (101.6-114.3cm)
BRX-XXL	XXLarge	45-50" (114.3-127cm)
BRX-XXXL	XXXLarge	50-55" (127-139.7cm)

#### BR Belt (D)

Belt Model No.	Size	Fits Waist
BR-S	Small	28-38" (71.1-96.5cm) 36-48" (91.4-121.9cm) 44-58" (111.8-147.3cm)
BR-M	Medium	36-48" (91.4-121.9cm)
BR-L	Large	44-58" (111.8-147.3cm)



- Standard Protection: 0.5mm Pb Equivalency
- For Sizing, Protection and Custom Options,
- refer to Pages 34-40



Ε

D



## **Suited for Special Procedures**



#### F) Surgical Drop Away Apron—Model SDA

Surgical Drop Away Apron offers top of the shoulder hook and loop fastening flaps for an adjustable fit. Features elastic back panels for added comfort and support. For quick removal, apron releases at shoulder and waist. Shoulder pads included.

Standard Protection: 0.5mm Pb Equivalency • For Sizing, Protection, Fabrics and Custom Options, refer to Pages 34-40

#### G) Surgical Drop Off Apron —Model SDO

Features back criss-cross flaps, with front hook and loop closure that assure easy removal from both over or under the surgical gown without disturbing the sterile field. Shoulder pads included.







• Standard Protection: 0.5mm Pb Equivalency • For Sizing, Protection, Fabrics and Custom Options, refer to Pages 34-40

4 www.shieldingintl.com



### **Frontal Aprons**

#### H) Three Quarter Wrap Apron— Model TQW

Three Quarter Wrap Apron, especially designed for the fluoro field, features unique side wrap protection with adjustable elastic back panels and hook and loop fastening flaps that wrap around waist for added support and comfortable fit. Apron also includes adjustable padded shoulders.

• Optional QR Buckle Closure available— Model TQWQR

#### I) Additional Side Wrap

Apron Width 18-20 - 4" each side 22-24 - 5" each side 27 - 6" each side

## For a Versatile Fit...

#### J) Versa-Adjustable—Model VA

Features two hook and loop fastening flaps at the waist, which aid in relieving weight off the shoulders. Its comfortable fit is attained by varying the position of the flaps. Shoulder pads included. Comes standard with hook and loop closure.

• Available with optional QR Buckle Closure as shown. Model VAQR

#### K) Adjust a Fit Apron—Model AF

Features two hook and loop fastening flaps at the waist, which aid in relieving weight off the shoulders and includes adjustable padded shoulders. Its comfortable fit is attained by varying the position of the flaps. Shoulder pads included.

Available with optional QR Buckle Closure. Model AFQR

Shown with optional QR **Buckle** Closure



 Standard Protection: 0.5mm Pb Equivalency • For Sizing, Protection, Fabrics and Custom Options, refer to Pages 34-40

### L) Fast Wrap—Model FW

Fast Wrap Apron offers both comfort and ease of wearing with single hook and loop fastener that wraps around the front. Adjustable back shoulder strap keeps apron snug. Shoulder pads included.







#### M) Uro-Apron—Model U

Designed for Urology, this apron features a unique tri-flap design for added comfort. Two hook and loop fastening flaps at the waist aid in relieving weight off the shoulders. Shoulder pads included.



• Standard Protection: 0.5mm Pb Equivalency • For Sizing, Protection, Fabrics and Custom Options, refer to Pages 34-40



## Two Pieces to Lighten the Load



Designed with two fully overlapping front panels, this vest and skirt combo reduces the weight by approximately 15% over our traditional and Adjust-A-Fit Vest Skirt Combos. Vest features hook and loop closure and shoulder pads for added comfort. For skirt description, refer to Model SO below.

**B) Skirt—Model SO** (when ordered separately) Skirt wraps around and overlaps in the front, offering extra protection. Fits securely on hips with extra wide hook and loop fastener and 1" quick release\* buckle for added support. Comes standard with large side pocket and hanging loops. \*2" quick release buckle available upon request.

Interested in Non-Slip Skirt Backing? See Custom Options page 40.





Adjust-A-Fit Vest Features wrap around protection and adjustable padded shoulders from the arm hole to the top of the shoulder, enabling the user to adjust the vest to a more secure and comfortable fit. Because of the versatility in fit, it accommodates a wide range of wearers. For skirt description, refer to page 8.



For Sizing, Protection, Fabrics and Custom Options, refer to Pages 34-40
Lighter weight protection options are recommended for all Vest Skirt models.



#### D) Vest Skirt Combo—Model VS

Proud holders of the original design and patent. Our standard Vest features all around protection with overlap in the front, shoulder pads and three panel hook and loop closure.

#### E) Flex Back Vest Skirt--Model FBVS

Flex Back Vest features wrap around protection with flex back support and adjustable hook and loop flaps at the waist and shoulder for added comfort and support, especially during long procedures. Shoulder pads included. For skirt description, refer to page 8.



### **Wrap Around Aprons**

В



#### A) Wrap Around Apron—Model WA

Wrap Around Apron is designed for front and back protection. This apron features extra wide shoulders, distributing weight over a greater area to reduce fatigue and three panel hook and loop fastener. Shoulder pads included.

#### B) Wrap Around Back Relief Apron—Model WABR

Wrap Around Back Relief Apron offering the same features of our Wrap Around Apron, transfers the weight from the shoulders to the hips with the aid of a 6" wide elastic belt. Comes standard with BR Belt, which may be ordered to fit various sized individuals. Indicate belt size from size chart when ordering. Shoulder pads included. BR Belt may be ordered separately. (Model BR)

#### Adjustable BR Belt Size Chart

Size	Fits Waist
	28-38" (7.1-96.5cm)
BR-M	36-48" (91.4-121.9cm)
BR-L	44-58" (111.8-147.3cm)

\* Lighter Weight Protection Options

### Looking for Extra **Back Support?**

Order with optional **BRX Belt**—Model WABRX For more info, refer to page 3

### C) Wrap Around Full Overlap Apron — Model WAFO

С

Designed with two fully overlapping front panels, this apron reduces the weight by approximately 15% over our traditional Wrap Around model. Features extra wide shoulders, hook and loop fastener and shoulder pads for added comfort.



For Sizing, Protection, Fabrics and Custom Options, refer to Pages 34-40 Lighter weight protection options are recommended on all Wrap Around Models

Recommended





#### D) Wrap Around Flex Back—Model WAFB

Features dual wrap elastic flex back with hook and loop side to front flaps while providing front and back protection for an adjustable and secure all around fit. Fits a wide range of sizes. Shoulder pads included.

#### (Not Pictured) Wrap Around Full Overlap Back Relief Apron -**Model WAFOBR**

Offering the same features of our Wrap Around Full Overlap Apron, this apron is approximately 15% lighter than our traditional Wrap Around Back Relief model and transfers the weight from the shoulders to the hips with the aid of a 6" wide elastic belt. Comes standard with BR Belt (see opposite page). Indicate belt size from size chart when ordering. BR Belt (Model BR) may also be ordered separately.

Looking for Extra Back Support? Order with optional BRX Belt – Model WAFOBRX (See page 3 for more information)





No Velcro<sup>®</sup>
 No Hook & Loop
 Buckles Only!

#### A) QR Vest Skirt Combo—Model QRVS

Quick-Release Buckle Vest Skirt features 2" quick release buckles for closure at the chest and waist. Also includes 1" buckle closure at the shoulder for secure fastening. Shoulder pads included.

**QR Skirt**—**Model QRSO** (If ordered separately) Quick-Release Buckle Skirt features a 1" quick release\* buckle closure on the inside and overlaps and closes with a 2" quick release buckle on the outside. Includes a large side pocket and hanging loops. \*2" quick release buckle available upon request.

Interested in non-slip backing? See custom options page 40

### **Buckle Aprons**

#### B) Wrap Around Quick-Release Apron—Model WAQR

Wrap Around Quick-Release Apron features extra wide shoulders, distributing weight over a greater area to reduce fatigue. Closing with 2" quick release buckles at the chest and waist and a 1" quick release buckle at the shoulder. Shoulder pads included.



• For Sizing, Protection, Fabrics and Custom Options, refer to pages 34-40 • Lighter weight protection options are recommended for both models.

Α

В

 Standard Protection: 0.5mm Pb Equivalency Royal Blue Nylon in stock for immediate shipment • For Sizing, Protection, Fabrics and Custom Options, refer to Pages 34-40

С



Half Apron is designed for front or back protection and is available with either Tie, Quick-Release Belt or Hook and Loop closure.

- Specify Closure next to P/N: (T) Ties, (V) Hook & Loop, (QR) Buckle
- In a hurry? In stock for immediate shipment: Royal Blue Nylon.

Model No.	Size: W x L
801	24x24"(61x61cm)
802	24x20" (61x50.8cm)
803	18x16″ (45.7x40.6cm)
804	10x10" (25.4x25.4cm)



### **Half Aprons**

### D) Demi Aprons

Demi Aprons feature convenient waist ring, protecting a range of sizes securely without ties or fasteners. The soft PVC covered spring steel ring holds the apron in place. Aprons may also be ordered with optional Quick-Release Buckle closure.

- In a hurry? In stock for immediate shipment: Royal Blue
- Nylon with waist ring. • Rack may be ordered separately. Model 830R

Waist Ring	QR Buckle	Size: W x L	Waist Measurements
830S	830QR-S	8.5x10" (21.6x25.4cm)	10-20" (25.4-50.8cm)
830M	830QR-M	12x12" (30.5x30.5cm)	21-30" (53.3-76.2cm)
830L	830QR-L	16x14" (40.6x35.6cm)	31-40" (78.7-101.6cm)
830XL	830QR-XL	20x16" (50.8x40.6cm)	41-50" (104.1-127cm)
830	830QR	Set of 4 w/ Rack	

830R Rack Only 14" x 9" x 3.5" (35.6cm x 22.9cm x 8.9cm)

#### E) Demi Rack — Model 830R

Sturdy metal rack is designed to hold set of four demi aprons using minimum wall space.



A) Panoramic Cape—Model 310 Panoramic Cape drapes over shoulders, shielding upper torso, front and back. Fastens in back with hook and loop closure. Overall: Width 23"(58.4cm) Length 25"(63.5cm)

#### B) Panoramic Poncho—Model 300PP

Panoramic Poncho features full front and partial back radiologic protection during the panoramic dental exam. It also features hook and loop shoulder closures for quick application. **300PP: (Adult)** Front 27" (68.6 cm) Back 16" (40.6 cm) Width 24" (61 cm)

**300PPC: (Child)** Front: 12" (30.5 cm) Back 6" (15.2 cm) Width 23" (58.4 cm)



Full Pano Apron provides excellent front and back protection. Features hook and loop shoulder closure and ties for quick application.

Front/Back: Width 23" (58.4 cm) Length 26.5" (67.3 cm)

• Protection, Fabrics and Custom Options, refer to Pages 35-40



D

 Standard Protection: 0.3mm Pb Equivalency or select an Alternate Protection
 Protection, Fabrics and Custom Options, refer to Pages 35-40 G



## **Dental Protection**

#### D) Kling-Kuver<sup>™</sup>--Model 300KK

Kling-Kuver<sup>™</sup> Dental Patient Apron is designed for complete patient protection of the neck, throat, chest, shoulders and gonadal area. Single hook and loop collar closure with non slip back makes this apron fast and easy to use. Includes hanging loops.

**300KK: (Adult)** Width 23.5" (59.7cm) Length 28" (71.1cm) **30KKC: (Child)** Width 20" (50.8cm) Length 20.5" (52cm)

#### E) Dental Tech Apron—Model 301

Dental Technologist Apron provides excellent dental staff radiologic protection. Ties provide for fast adjustment and a secure fit. Shown standing. One Size: Width 24" (61cm) Length 37" (94cm)

#### F) Dental Patient Apron – Model 300 & 300C

Dental Patient Apron is easy to apply. Its non-slip backing keep apron firmly in position by actually "gripping" the patients clothing. No ties or fasteners involved. Available in adult and children sizes. Includes hanging loops.

**300: (Adult)** Width 23" (58.4 cm) Length 26" (66 cm) **300C: (Child)** Width 19" (48.3 cm) Length 21" (53.3 cm)

### **Dental Aprons**

### G) Dental Thyroid Collar—Model 311NSL

Dental Thyroid Collar provides neck and thyroid protection during dental radiography. Fastens in back with hook and loop closure.

Model No.	Length	Bib
		9x3.75" (22.9x9.5cm)
311NSL	25.5″ (64.8 cm)	10x5.5" (25.4x14cm)





## Protection for Pediatrics

Α

### A) Pediatric Coat Apron (Model P-ST)

Pediatric Coat Apron features frontal protection with tie closure for easy on and off. Available in our wide selection of fabrics. Includes shoulder pads.

Model No.	Size	Measures: W x L
P-ST-S	Small	12x20" (30.5x50.8cm)
P-ST-M	Medium	14x22" (35.6x53.3cm)
P-ST-L	Large	16x24" (40.6x61cm)

### B) Pediatric Half Apron (Model P-HA)

Pediatric Half Apron with convenient hook and loop fastener. Available in our wide selection of fabrics.

Model No.	Size	Measures: W x L
P-HA-S	Small	10x10" (25.4x25.4cm)
P-HA-M	Medium	12x10" (30.5x25.4cm)
P-HA-L	Large	14x12" (35.6x25.4cm)

### C) Pediatric Spinal Stole (Model P-SS)

Pediatric Spinal Stole features lead-free back for total spine viewing and non-slip underside. Fully adjustable. Available in our wide selection of fabrics. One size.

 Model No.
 Size: L x W

 P-SS
 Small
 14x5.25" (35.6x13.3cm)

Standard Protection: 0.5mm Pb Equivalency
For Protection, Fabrics and Custom Options, refer to Pages 35-40



Standard Protection:

0.5mm Pb Equivalency

or select an alternate

protection

### D) Pediatric Thyroid Collar (Model P-TC)

Pediatric Thyroid Collar features hook and loop rear closure. Available in our wide selection of fabrics or in Royal Blue Nylon for immediate shipment.

Model No.	Length	Bib
P-TC-S	15" (38.1cm)	3″ (7.6cm)
P-TC-L	18" (45.7cm)	4" (10.2cm)



### E) Ped Pal<sup>™</sup>

Pediatric shield protects children from radiation. Features soft n' plush non-slip microfiber backing for added comfort.

- Includes Hanging Loops for Safe & Easy Storage
- Shield Measures 16"W x18"L
- Choice of Three Children Friendly Designs

Model No.	Description
PP-LBUG-RL	Ped Pal - Lady Bug
PP-KITTY-RL	Ped Pal - Kitty
PP-BEAR-RL	Ped Pal - Bear



Ped Pals™ Put Pediatric Patients at ease.

> • Standard Protection: 0.5mm Pb Equivalency or select an alternate protection









### A) Gonadal & Ovarian Belt System—Model 822BBV

Features four attachable gonad shields in assorted sizes that interchange to one hook and loop belt. Set of 4 includes 2 small, 1 medium and 1 large shield.

#### D) MP Shield – (Model MPS)

Multi Purpose Shield may be used for shielding patients or equipment. Nonslip backing keeps shield in place. May be ordered in standard RL, TL Lite Ply, or \*NL Non Lead protection options. Available in three sizes.



0.5mm Pb Equivalency Protection Type	Small: 18x24" (45.7x61cm)	Medium: 24x24″ (61x61cm)	Large: 24x36″ (61x91.4cm)
RL Reg Lead	MPS-S-RL	MPS-M-RL	MPS-L-RL
TL Lite Ply	MPS-S-TL	MPS-M-TL	MPS-L-TL
*NL Non Lead	MPS-S-NL	MPS-M-NL	MPS-L-NL

\*NL Non Lead recommended for (non x-ray) protection purposes from scatter radiation due to seed implants. The metals contained in this item are safe for standard disposal.

### Standard Protection: 0.5mm Pb Equivalency Styles A+B in stock for immediate shipment in Royal Blue Nylon or Light Blue Vinyl.

"Model 806AT"

"Model 807AT"

#### B) Gonadal Shields—Model 806

Gonad shields are available with or without optional Ties, Hook and Loop closure or QR Buckle. If closure is not specified, plain shields will be ordered. Available in three sizes.

• Specify Closure next to P/N: (T) Ties, (V) Hook & Loop, (QR) Buckle

Model No.	Size: W x L
806S	3x5" (7.6x12.7cm)
806M	6x7" (15.2x17.8cm)
806L	7x9" (17.8x22.9cm)
806AT	Set of 3

### C) Ova Shields—Model 807

Specify Closure next to P/N: (T) Ties, (V) Hook & Loop, (QR) Buckle

Model No.	Size: W x L
807S	3x4" (7.6x10.1cm)
807M	4.5x5.5" (11.4x13.9cm)
807L	7x6"(17.8x15.2cm)
807AT	Set of 3



#### J) Lead Sleeve—Model 808

Provides arm protection while allowing complete hand mobility. Hook and loop closure for fast and easy securing. Covered in Light Blue Vinyl.

Model No.	Length:	Diameter:
808S	19" (48.3cm)	6-12" (15.2-30.5cm)
808M	21″ (53.3cm)	7-13" (17.8-33cm)
808L	23" (58.4cm)	8-16" (20.3-40.6cm)

#### (For Apron Sleeves see Custom Options page 40)





#### **Comfort Fit Thyroid Collar**— Model 641CF

Features flexible protective design, but without the top binding which makes collar more comfortable to wear. Rear hook and loop closure.

Model No.	Length	Bib
641CFS	24" (61cm)	9x3.75" (22.9x9.5cm)
641CFL	26.5" (67.3cm)	9x5.5" (22.9x14cm)

#### Standard Thyroid Collar— Model 600

Standard Thyroid Collar is designed for maximum thyroid protection. Includes binding and rear hook and loop closure.

হৈছে	Model No.	Length	Bib
	600S	18.5″ (47cm)	10x6.5″ (25.4x16.5cn
	600L	22" (55.9cm)	10x6.5″ (25.4x16.5cn

#### Ultra Thyroid Collar—Model 621U

Ultra Thyroid Collar is our softest and most flexible, offering the best of comfort without the binding. Rear hook and loop closure.

Model No.	Length	Bib
621US	24" (61cm)	6.75x5.75" (17.1x1
621UL	26.5″ (67.3cm)	6.75x5.75" (17.1x1

## **Apron Collar Options**



### All Items on both pages 0.5mm Pb Equivalency

Accessories

### F) Radiation Safety Cap— Model RSC Adjustable hook and loop tab in back

allows for custom fit. Size/diameter is measured with tape around head, just above the ears. Available in three sizes.

#### Model No. Size 21"/53.3cm RSC-S RSC-M 22"/55.9cm 23.5″ / 59.7cm

## G) Spinal Stole—Model 810

Spinal Stole, designed for use in spinal radiograph procedures, features nonslip backing and lead-free back neck area providing total spine viewing. Fully adjustable.

In a Hurry? In stock for immediate shipment: Royal **Blue Nylon** 

#### Model Size

- 810 Adult: 25"L x 9"W (63.5x22.9cm)
- 810C Child: 20"L x 5.25"W (50.8x13.3cm)







G







14.6cm) (14.6cm)

#### Standard Protection: 0.5mm Pb Equivalency

- For Protection and Fabric Options, refer to Pages 35-38
- In a hurry? In stock for immediate shipment in Royal Blue Nylon, except Model 600 in Light **Blue Vinyl**

#### New Style Thyroid Collar— Model 611NS

New Style Thyroid Collar offers a trim, yet shielded and comfortable design. Includes binding and rear hook and loop closure.

Model No.	Length	Bib
611NSS	22.5" (57.1cm)	9x3.75" (22.9x9.5cm)
611NSL	25.5″ (64.8cm)	8.5″ x 5.5″ (21.6x14cm)

#### Quick-Release Collar—Model 631QR

Quick-Release Collar offers the same look and comfort as our model 621U, but without the hook and loop closure. Rear buckle closure.

Model No.	Length	Bib
631QRS	24" (61cm)	6.5x5.5" (16.5x14cm)
631QRL	25.5″ (64.8cm)	6.5x5.5″ (16.5x14cm)

н

### H) Tethered Collar—Model -T

Collar tethered to apron on right side by a ring for easy removal. Available on all models.

#### I) Attached Collar—Model 600A

Sewn directly to any frontal apron.

J) Shin Guards — Model SG

For long bone (tibia) protection, Shin Guards feature hook and loop closure for easy adjustment and comfort. Shin Guards are wider at the top tapering down to the foot. Available in two sizes or may be custom ordered. Sold in pairs.

/lodel	Shin Length
G-S	13.5″ (34.3 cm)
G-L	17" (43.2 cm)



## **Thyroid Collars**

## **Special Procedure Gloves**



## Attenuator-X

		Mo	odel No.	Size
Radiation Protection Gl	oves	AX	-6.5	6 1⁄2
A) Attenuator-X—Model AX			-7	7
			-7.5	7 1⁄2
Lead Free		AX	-8	8
<ul> <li>Natural Rubber Latex</li> <li>Powder Free</li> </ul>		AX	-8.5	8 1⁄2
<ul> <li>Powder Free</li> <li>Rolled Cuff for Ease of Donning and Removal</li> </ul>		AX	-9	9
<ul> <li>Exceptional Tactile Sensitivity</li> </ul>	Protect	tion Data	:*	
and Flexibility <ul> <li>High Attenuation</li> </ul>	60kVp	80kVp	100kVp 46%	120kVp
0.35mm Thickness	60%	53%	46%	41%

#### **Sterile Attenuative Gloves**

Model No.

IBG-6.5

For use during interventional procedures, diagnostic heart catheterization, coronary angioplasty, angiocardiography, pain management, orthopedics and fluoroscopy procedures.

\* Test results in accordance with ASTM Test Method F2547-06

Size

6 1/2

7 1/2

8

8 1⁄2

9

### **B)** Radiation **Attenuating Surgical** Glove—Model IBG

Glove—Model IBG	IBG-7
Lead Free	IBG-7.5
<ul><li>Non Latex</li><li>Powder Free</li></ul>	IBG-8
High Attenuation	IBG-8.5
<ul> <li>High Tactile Sensitivity</li> <li>0.18mm Thickness</li> </ul>	IBG-9

#### **Protection Data: \***

60kVp	80kVp	100kVp	120kVp
56.5%	47.1%	41.6%	37.7%

#### **C) Attenuating Specialty** Gloves—Model ASG

• Lead Fr	ee			Model No.	Size	37
Non Latex			ASG-6.5	6 ½		
<ul> <li>Excellent Flexibility, Dexterity and Touch Sensitivity</li> </ul>			ASG-7	7		
0.23mm Thickness				ASG-7.5	7 ½	
Protection Data: *				ASG-7.5	8	
60kVp	80kVp	100kVp	120Kvp	ASG-8.5	8 ½	
45%	36%	29%	24%	ASG-9	9	

\*Based on mean attenuation and combined results of fingers and cuff.

#### D) Leather Lead Gloves—Model 200L

This 15" one-piece, dip-molded lead glove features durable goat skin leather cover and foam liner. Available in pairs or single gloves.



www.shieldingintl.com

#### E) Vinyl Lead Gloves—Model 100V

Seamless Lead Vinyl Gloves are a 15" one-piece, dipmolded glove with supple vinyl cover and foam liner. Available in pairs or single gloves. Available with removable, washable Blue Sheath Cotton Cover-Model 100S.



#### Standard Protection: 0.5mm Pb Equivalency

## **Mitts and Mittens**

#### F) Open Palm Mitts—Model 1040PM

Open Palm Vet Mitt shields the back of the hand and features hook and loop wrist strap closure and elastic loops to secure fingers for quick removal. Covered in Dark Blue Vinyl.



#### Standard Protection: 0.5mm Pb Equivalency on all non sterile gloves, mitts & mittens.



#### H) Lead Mittens—Model 101MIT

Lead Mittens are 15" one-piece, dipmolded mitten featuring soft and supple black goat skin leather or blue vinyl cover with foam liner. Provides full hand and palm protection.

Model No.	Color
101MIT	Black Leather
101MIT-BV	Blue Vinyl





Model No. 101SLMIT 101SLMIT-BV

Color Dark Blue Burgundy

Model No.

100V-DKB

100V-BUR

С

J) Optional Removable Glove Liner

Both glove styles now available with removable washable foam liner. Add (-R) to Model No. when ordering. (example 100V-DKB-R)

#### G) Vet Mitts—Model 104VM

Designed for vets as a cost effective solution, the Vet Mitt features wrist slit allowing fingers to exit for freedom of movement during radiography. Contour cut to fit the arm and covered in Royal Blue Nylon. Custom coverings optional.




## Gloves

Color Black Leather Blue Vinyl

#### I) Slit Lead Mittens-Model 101SLMIT

Slit Mitten features same black leather cover or blue vinyl with foam liner and comes with base of the finger slit allowing for dexterity and hand protection.





• Front Lenses: 0.75mm Pb Equivalency Optical Lead Glass Side Shields: 0.35mm Pb Equivalency Lead Foil



White 53UFO-W

Silver 95-S

Black 95-BL

Lead Foil

BI—Bifocal

Orange

77-0

Yellov

77-Y

Black Tortoise White

Wt. (approx.) Lens Type Plano 92 grams Plano 92 grams Plano 92 grams


Black 96-B

С C) Ultra Wrap – Model 96 Red 95-R Silver 96-S Unifit Wrap Around Frame Wrap Around Design for Added Coverage

- Rubberized Nose Pads
- Secure Fit Nose Bridge
- Medium to Larger Frame Size
- Front Lenses: 0.75mm Pb Equivalency Approximate Wt. 78 grams

### E) Astro Flex—Model 70

Astro Flex eyewear features full field panoramic view with adjustable temple length and frame angle for a uniform fit. Includes side shield protection.








#### F) Max Fit—Model 97

- Front Lenses: 0.75mm Pb Equivalency Optical Lead Glass



with uni-fit ear and nose pieces for comfortable wear. Fits well for those with larger facial features. Includes side shield



#### D) Hornet—Model 77

• Approximate Wt. 71 grams

Optional RX: SV, BI

B) Sport Wrap – Model 95


Rubberized Temple Bars & Nose Pads

• Front Lenses: 0.75mm Pb Equivalency

• Light Weight Unifit Frame

R

Hornet eyewear's popular wrap-around design provides excellent coverage and features foam brow bar for comfortable flexibility and rubber temple tips and nosepads for non-slip wear. Fitted side shields increase scatter radiation protection. Includes side shield protection.





22 www.shieldingintl.com

#### Eyewear manufactured and inspected by board certified licensed optician. 0.75mm Pb equivalent optical lead glass - 98% attenuation up to 100 kVp. \*RX lenses meet minimum 0.75mm Pb equivalency (+/-.08) at 100 kVp.

Black	Red	Silver	Purple/Grey	Blue/Yellow	Wt. (Aprox.)
98UFB	98UFR	98UFS	98UFPG	98UFBY	73 grams



#### \* Rx Lenses meet minimum .075mm Pb Equivalency (+/- .08) at 100k Vp

#### H) Wrap Around—Model 53

Wrap Around eyewear features comfortable molded wrap styling in durable tortoise nylon frame with uni-fit bridge and leaded side shields for added protection. Provides excellent peripheral visibility. Available Rx\* lenses.

Titanium Oval eyewear features ultra lightweight composition and adjustable non-slip silicone nose pads for maximum long-wearing

comfort. Fitted side shields increase scatter radiation protection.

• Optional Rx\*: SV – Single Vision, BI – Bifocal and TRI – Trifocal

Model No.	Lens Type	Wt. (approx.)
53	Plano	90 grams

• Front Lenses: 0.75mm Pb Equivalency

• Side Shields: 0.50mm Pb Equivalency

• Optional Rx\*: SV – Single Vision, BI – Bifocal and TRI – Trifocal



Model No.	Lens Type	Wt. (approx.)
57	Plano	68 grams

K) Titanium Oval—Model 55

• Front Lenses: 0.75mm Pb Equivalency • Side Shields: 0.35mm Pb Equivalency Lead Foil

Available with Rx\* lenses.

#### J)Rad-Tech—Model 57

Rad-Tech eyewear features lightweight design with stainless metal frame, wrap temple and non-slip adjustable nose pads for stylish, yet comfortable fit. Fitted side shields increase scatter radiation protection. Available with Rx\* lenses.

- Front Lenses: 0.75mm Pb Equivalency
- Side Shields: 0.35mm Pb Equivalency Lead Foil • Optional Rx\*: SV – Single Vision, BI – Bifocal and TRI – Trifocal



Model No.	Lens Type	Wt. (approx.)
55	Plano	60 grams



#### L) Standard Metal—Model 50

Standard Metal eyewear is constructed of optical steel base metal with gold electro-plate and adjustable non-slip silicon nose pads for maximum comfort. Fitted side shields increase scatter radiation protection. Available with Rx\* lenses.

- Front Lenses: 0.75mm Pb Equivalency
- Side Shields: 0.35mm Pb Equivalent Lead Foil
- Optional Rx\*: SV Single Vision, BI Bifocal and TRI Trifocal

#### Model No. Lens Type Wt. (approx.)

50 Plano 69 grams



### **Designer Eyewear Offers Protection, Comfort & Style**



Model No.	Color	Lens Type	Wt. (approx.)	
82-GB	Gloss Black	Plano	70 grams	
82-MB	Matte	Plano		
82-T	Tortoise	Plano		
-			-	

#### N) Wiley X Brick – Model 83

- Square Lens Wrap Around Frame
- Rubberized Arms for Added Grip
- Foam Facial Seal for Splash Protection
- Front Lenses: 0.75mm Pb Equivalency
- Wt. 75 grams




#### Eyewear manufactured and inspected by board certified licensed optician. 0.75mm Pb equivalent optical lead glass - 98% attenuation up to 100 kVp. \*RX lenses meet minimum 0.75mm Pb equivalency (+/-.08) at 100 kVp.

#### M) Wiley X-PRO – Model 82

Matte Black Gloss Black

Curved frame provides contoured coverage around eyes. Features: non-slip rubberized temples and nose piece for added comfort.

- Front Lenses: 0.75mm Pb Equivalency
- Optional Rx\*: SV Single Vision Only Tortoise 82-T



	Model No.	Color	Lens Type	Wt. (approx.)
	85-BL	Black	Plano	63 grams
Res .	85-R	Red	Plano	

### O) CHIC – Model 85

Offers 'Chic' look and functionality. Features comfortable head-hugging temple design. Includes side shield protection.

- Front Lenses: 0.75mm Pb Equivalency
- Side Shields: 0.35mm Pb Equivalency Lead Foil
- Optional Rx\*: SV Single Vision, BI Bifocal

#### P) JACKIE-O – Model 86

Features sunglass-style big lenses and Star Décor with head hugging temple bars. Includes side shield protection.

Black 86-BL

Blue 86-B

- Front Lenses: 0.75mm Pb Equivalency
- Side Shields: 0.35mm Pb Equivalency Lead Foil
- Optional Rx\*: SV Single Vision, BI Bifocal

Model No.	Color	Lens Type	Wt. (approx.)
86-B	Blue w/Black	Plano	63 grams
86-BL	Black w/White	Plano	







#### D) Combination Apron and Glove Holder—Model 820


Combination Apron and Glove Holder offers quick, convenient storage for an apron and a pair of gloves. Powder coat and chrome steel construction. Dims. 27.75"W x 17"H (70.5 cm x 43.2 cm) Wt. 8.5 lbs (3.9 kg)

E) Multi-Apron Glove Rack—Model 821

Multi-Apron Glove Rack is ideal for storing multiple pairs of gloves and aprons. Also has convenient bar to hold thyroid collars. Mounts quickly and easily to any wall and requires a space only 19" wide. Powder coated heavy steel construction.

E

Dims. 19"W x 20"H (48.3 cm x 50.8 cm) Wt. 10 lbs (4.5 kg)



F) Wall Mount Swing Arm Apron Rack—Model 811

by-side for a 10 hanger wall console unit. Right Folding **811R** Dims. 27.5"W x 13"H (69.7 cm x 33 cm) Left Folding 811L Dims. 27.5"W x 13"H (69.7 cm x 33 cm)



**813-2** (Holds 1 Apron) Dims. 5" x 4" (12.7 cm x 10.2 cm) Wt. 1.5 lbs (680 g)

### A) Valet Rack—Model 816

Valet Rack can hold up to 300 lbs of aprons and can fit into any budget. Equipped with heavy duty casters and 10 closed loop hangers, this unit is extremely mobile, allowing aprons, thyroid shield, etc. to be transported easily to workstations. Additional hangers Model 816H may be purchased separately. Dims. 36.5"W x 22.5"D x 53"H (92.7 cm x 57.2 cm x 134.6 cm) Wt. 36 lbs (16.3 kg)

#### B) Porta-Rak—Model 817

Porta-Rak is a totally mobile unit for up to 8 aprons and designed for where space is at a premium. Sturdy aluminum five leg base with nylon casters uses only 24" of floor space. Dims. 24"W x 50"H (61 cm x 127 cm) Wt. 29 lbs (13.2 kg)

#### C) Mobile Swing Arm Apron Rack – Model 811M


Features steel construction with white powder coating and 10 individual chrome shoulder shape swing arms for easy apron access and hanging. Large locking casters for easy movement. Holds 10 aprons.

Dims: 22"W x 36"D x 53"H (55.9 cm x 91.4 cm x 134.6 cm) Wt. 48 lbs (21.8 kg)

Width 21" (53.3cm) Height 10" (25.4cm)



### **Porta-Shields**



Model No.	Shield Size: W x L	
1001	24x24" (61 x 61cm)	
1001-36	24x36" (61 x 91.4cm)	
1001-48	24x48" (61 x 121.9cm)	



Blue Nylon.



Covered in Royal Blue Nylon.

A) Porta-Shield—Model 1001

• In a hurry? In stock for immediate shipment!

B) Graduated Porta-Shield—Model 1002

Porta-Shield is ideal for patient protection during diagnostic

procedures and features adjustable height from 36" to 60" (91.4 cm to 152.4 cm) with sturdy five leg base for stable

mobility. Shield available in three sizes and covered in Royal

Features seven leaf shield with graduated effect, it fits easily

and comfortably over the patient and is perfect for cath lab use. Shield width is 25" (63.5 cm) and length of the longest panel is 28" (71.1 cm), decreasing in increments of 2" (5 cm).

Adjustable from 40" to 64" (101.6 cm to 162.6 cm) in height.

Seven leaf shield is ideal for general radiology, as well as fluoro applications. Shield width is 25" (63.5 cm) and length is 28" (71.1 cm). Adjustable from 40" to 64" (101.6 cm to 162.6 cm) in

C) Non Graduated Porta-Shield—Model 1003





#### E) Drape Lead Curtain—Model DLC

Drape Lead Curtains, ceiling mounted with optional curtain track, are an ideal and economical solution for temporary or permanent set-ups requiring radiation protection where a fixed, protective wall is not practical. Curtain sections easily hang and fasten together to cover large areas.

\* Maximum width per single panel is 40" (101.6 cm) \*\* Allow 2" (5 cm) overlap for multiple curtain panels.

Curtain Type	0.3mm Pb Equiv.	0.5mm Pb Equiv.	1.0mm Pb Equiv.
	Model No.	Model No.	Model No.
With Panels** Without Panels* Optional Curtain Track	DLC-3P DLC-3 DLC-CT	DLC-5P DLC-5	DLC-10P DLC-10







### **Vinyl Blockers**

**Coated Lead Vinyl Blockers** feature a smooth sanitary, non-absorbent surface for easy cleaning and are stocked in a range of sizes.



Available in 0.5mm and 1.0mm Pb lead equivalency. Color: Beige

#### • All sizes in stock for immediate shipment.

Size	Pb 0.5mm	Pb 1.0mm
4x10" (10.2x25.4cm)	SB1.5	SB1
5x8" (12.7x20.3cm)	SB2.5	SB2
5x12"(12.7x30.5cm)	SB3.5	SB3
6x10" (15.2x25.4cm)	SB4.5	SB4
5.5x14" (14x35.6cm)	SB5.5	SB5
7x11″(17.8x28cm)	SB6.5	SB6
7x17" (17.8x43.2cm)	SB7.5	SB7
8.5x14" (21.6x35.6cm)	SB8.5	SB8
8.5x17" (21.6x43.2cm)	SB9.5	SB9
10x12" (25.4x30.5cm)	SB10.5	SB10
11x14" (28x35.6cm)	SB11.5	SB11



Shielding's X-Curtains are manufactured specific to the needs of our customers' exclusive cabinet x-ray system engineering, color and protection (lead Pb equivalency) level requirements.

\* Standard X-Curtain (Model: **XCURTAIN**) features durable dual sided PVC vinyl outer covering impregnated with leaded vinyl.

\* Shielding's Specialty Curtain (Model: XCURTAIN/SC) features low coefficient outer coating with mesh scrim reinforcer impregnated with leaded vinyl, making this product preferable to applications that are processing lighter weight items (i.e., purses, wallets, jackets, etc.).

#### Features:

Lead Equivalencies: 0.25, 0.35, 0.5, 0.7 and 1.0mm, depending on level of protection and type of X-Curtain. Alternate equivalencies may be manufactured, but require a minimum. Outer Covering: Model **XCURTAIN** available in black, charcoal grey or beige. Model **XCURTAIN/SC** in black only.

Dimensions: Available in sheet size (48" x 96" / 121.9 cm x 243.8 cm) or custom die cut to your engineering specification.

D



#### D) Adjustable Mobile Shield—Model 1901

Adjustable Mobile Shield is easily portable and ideal for shielding patients in such procedures as chest, skull and spine radiography. Includes beige anti-microbial weblon covering. Adjustable from 24" to 48" (61 cm to 121.9 cm) in height. Shield Width 24" (61 cm). Chrome frame.

\*Standard Protection .05mm Pb Equivalency



### Lead Vinyl Sheeting

Coated Lead Vinyl Sheeting drapes readily for on the job applications wherever flexible lead sheeting is required and is offered in a variety of lead equivalencies. Sold in sheet size (48" x 96" / 121.9 cm x 243.8 cm) or cut to specifications. Cutting charge may apply. Color: Beige

Model No.	Pb Equivalency
901	0.25mm
902	0.5mm
903	1.0mm

(Black available for technical applications, subject to full sheet or volume ordering.)



### **Shields & Blockers**





### **CT** Shields

## **Immobilization Sand Bags and Table Scatter Shields**









A) CT Breast Shields offer up to a 57% reduction of radiation exposure during CT scan without significant changes in image quality, including artifacts and beam hardening effects. Shield consists of a thin (1mm) piece of bismuth impregnated synthetic rubber covering that mounts firmly on an offset quarter inch foam base. Designed for multiple uses. Available in Adult and Pediatric sizes.

	Model No.	Adult Size	Quantity	
	CT-BS-M	Medium	1 each	
	CT-BS-L	Large	1 each	
	CT-BS-COV	Disposable Cover	10 pack	
	Model No.	Pediatric Size	Quantity	
CT-PBS-N Nec		Neonate	10 pack	

CI-PBS-N	Neonate	10 pack
CT-BPS-1	1 Year Old	8 pack
CT-BPS-5	5 Year Old	6 pack
CT-BPS-10	10 Year Old	4 pack

B) CT Eye Shields, designed for use in head and facial CT scans, mold over the nasal bridge with the eye shield portions laying flat over the patient's closed eyes allowing for about 50% dose savings to the eye lens. Shield consists of a single piece thin (1mm) layer of bismuth impregnated synthetic rubber with adhesive back strip to ensure a tight seal with minimal air gaps. Intended as a single use product to prevent against potential biological contamination. Model No. CT-ES (20 pack)

**C) CT Thyroid Shields**, when used for the soft tissues of the neck, offer a 60% dose savings to the radiosensitive thyroid gland during CT procedures. Shield consists of a thin (1mm) piece of bismuth impregnated synthetic rubber and features a curved notch in the center to ensure proper positioning. Lateral bow ties are also attached to the interior flaps of the shield to accommodate larger necks or larger thyroid glands that extend beyond the clavicle or sternum. Adhesive back holds the shield firmly in place as patient breathes and swallows. Model No. CT-TS (10 pack)



#### A) Cervical Sandbag Set -Model CSS1012

10" x 12" (25.4 cm x 30.5 cm) 10 lbs (4.5 kg)





Α

R

**C)** General Sandbags come with grommets. If grommets are not desired, omit (-G) from the Model No. (i.e. GS612).

#### Table Scatter Shield—Model TSS-Custom

Custom table scatter shield features a dual or single sided saddle back design, allowing a secure fit to any style imaging table. The overlapping paneled shields are secured to a padded surface that simply drape over the tabletop, beneath the mattress. Single sided shields include self adhesive hook and loop for securing to table.

- Pleated lead panels enable staff and c-arm mobility
- Customized to specific table dimensions, panel lengths, or to include table end flaps
- Constructed of 0.5mm Pb equivalent panels
- Select from our wide range of solid color nylons, or our durable anti-microbial healthcare weblon coverings.



### B) Angiographic Sandbag Set—

7" x 12" (17.8 cm x 30.5 cm) 7 lbs (3.2 kg)



GS716-G



Model ASS712-G





#### Scatter Panel Shields – Model SPS-Custom

Scatter Panel Shield replacements for various OEM fluoro table shield assemblies.

- Customizable to Specific Panel Size
- Constructed of 0.5mm Pb Equivalent Panels
- Select from Our Durable Anti-Microbial Healthcare Weblon Coverings.

# 7x16″ (17.8x40.6cm) m) 7.5 lbs (3.4) kg)





### **Mobile Barriers**

#### **Regular and Jumbo Clear Pb Barriers**

Clear Pb Acrylic Mobile X-Ray Barriers feature optically clear and shatter resistant lead acrylic windows. Available in Regular and Jumbo size, as well as a variety of sizes for shielding 1, 2 or more people. Opaque lower panels finished in beige. Attractive, durable, easy to assemble and maneuver.

Barrier Specifications	Regular	Jumbo
<b>Model No.</b> Overall Size	<b>56-603</b> 75″H x 30″W (190.5x76.2cm)	<b>56-604</b> 75″ H x 30″ W (190.5x76.2cm)
Window Panel Size	24″ H x 30″W (61x76.2cm)	60″ H x 30″ W (152.4x76.2cm)
Lead Equivalency	0.5mm	0.5mm
Opaque Panel Size	48″ H x 30″ W (121.9x76.2cm)	12″ H x 30″ W (30.5x76.2cm)
Lead Equivalency	Min 1/32" (0.8mm)	Min 1/32″ (0.8mm)
Est. Shipping Weight	114 lbs / 51.7 kg	114 lbs / 51.7 kg

In a hurry? Above models in stock for immediate shipment. Subject to stock on hand.

## Leaded Windows and Frames

Leaded Glass Windows are made of high quality optical grade material, offering excellent radiation protection and superb visual clarity. Available in custom sizes up to 48" x 96". \*(LG = Lead Glass). For cost effectiveness, refer to standard sizes in chart at bottom of page.

Lead Equivalent	Thickness	X-Ray Peak Voltage
1.6 – 2.0mm (1/16" – 5/16")	7 – 9mm	150 kVp

Leaded Acrylic Windows feature 0.3mm or 0.5mm lead equivalency, as well as alternate equivalencies up to 2.0mm. Refer to chart for stock sizes. \*(LA = Lead Acrylic). Custom sizes also available.

Telescopic Lead Frames are manufactured for ease of installation in rough opening after wallboard is installed. Constructed from 16 gauge "JETCOAT" (HGDS) steel for rust resistance with paintable surface and fully lined with 1/16" lead lining. Adjusts for walls 4-1/4" (10.8 cm) to 7" (17.8 cm) thick. Fits into rough opening of 1" larger than glass size. \*(LF = Lead Frame).

> Lead Glass/Acrylic and Frames Sold Separately\* \*Crating Charges Apply

Size	Lead Glass/Acrylic* Model No.	Lead Frame Model No.	Size	Lead Glass/Acrylic* Model No.	Lead Frame Model No.
8x10″/20.3x25.4cm	LG810	LF810	24x36" / 61.0x91.4cm	LG2436/LA2436	LF2436
10x10″/25.4x25.4cm	LG1010	LF1010	24x48"/ 61.0x121.9cm	LG2448/LA2448	LF2448
10x12″/25.4x30.5cm	LG1012/LA1012	LF1012	30x30" / 76.2x76.2cm	LG3030	LF3030
12x12" / 30.5x30.5cm	LG1212/LA1212	LF1212	30x36" / 76.2x91.4cm	LG3036/LA3036	LF3036
12x16″ / 30.5x40.6cm	LG1216/LA1216	LF1216	32x40"/81.3x101.6cm	LG3240/LA3240	LF3240
12x18″/30.5x45.7cm	LG1218/LA1218	LF1218	36x36" / 91.4x91.4cm	LG3636/LA3636	LF3636
12x20″/ 30.5x50.8cm	LG1220	LF1220	36x48" / 91.4x121.9cm	LG3648/LA3648	LF3648
12x24″/30.5x61.0cm	LG1224/LA1224	LF1224	36x60"/91.4x152.4cm	LG3660/LA3660	LF3660
14x18" / 35.6x45.7cm	LG1418	LF1418	36x72" / 91.4x182.9cm	LG3672/LA3672	LF3672
16x20″/40.6x50.8cm	LG1620	LF1620	36x84" / 91.4x213.4cm	LG3684/LA3684	LF3684
16x24″/40.6x61.0cm	LG1624/LA1624	LF1624	36x96" / 91.4x243.8cm	LG3696	LF3696
18x24″ / 45.7x61.0cm	LG1824/LA1824	LF1824	48x48" / 121.9x121.9cm	LG4848/LA4848	LF4848
18x26″ / 45.7x66.0cm	LG1826	LF1826	48x60"/121.9x152.4cm	LG4860/LA4860	LF4860
20x24" / 50.8x61.0cm	LG2024	LF2024	48x72" / 121.9x182.9cm	LG4872/LA4872	LF4872
24x24" / 61.0x61.0cm	LG2424/LA2424	LF2424	48x84" / 121.9x213.4cm	LG4884/LA4884	LF4884
24x30" / 61.0x76.2cm	LG2430/LA2430	LF2430	48x96″/121.9x243.8cm	LG4896/LA4896	LF4896

#### **Mobile Screens**

Mobile floor screens feature 1.0mm Pb equivalent virgin lead lining, framed with heavy metal channeling, heavy-duty stainless steel legs and locking casters to ensure stability. Lead glass window has a minimum 2.0mm lead equivalent protection and is positioned with center of view 5'3" (160cm) up from finished floor. Plastic laminate skins covered in beige or grey. Screens and windows available in additional sizes.

Barrier Dimensions	8″ x 10″ Window 20.2 cm x 25.4 cm	10″ x 12″ Window 25.4 cm x 30.5 cm	12″ x 12″ Window 30.5 cm x 30.5 cm
24x72" (61x182.9cm)	MS2472-8	MS2475-10	MS2472-12
24x82" (61x208.3cm)	MS2482-8	MS2482-10	MS2482-12
30x72" (76.2x182.9cm)	MS3072-8	MS3072-10	MS3072-12
30x82" (76.2x208.3cm)	MS3082-8	MS3082-10	MS3082-12
36x72" (91.4x182.9cm)	MS3672-8	MS3672-10	MS3672-12
36x84" (91.4x213.4cm)	MS3684-8	MS3684-10	MS3684-12

Note: A crating charge applies to all orders.



## **Windows & Frames**

## **Ordering an apron is easy!**

Select your Model No. from pages 1-12 Use loose measurements to find the correct apron width and length from the charts below. Choose your protection from charts on facing page.

		l Aprons		Three C	Quarter Wrap <i>I</i>	Aprons
Example Part # for Frontal and TQW Aprons	Apron Width	Chest		Apron Width	Chest	Wrap Width
Standard Coat Apron with Ties – Reg. Lead Model - Width X Length - Protection		<b>30" – 34"</b> 76.2 cm – 86.4 cm		18″ 45.72 cm	<b>30″ – 34″</b> 76.2 cm – 86.4 cm	<b>4</b> ″ 10.2 cm
ST - 24 X 36 - RL	<b>20″</b> 50.8 cm	<b>34" – 38"</b> 86.4 cm – 96.5 cm		20″ 50.8 cm	34" – 38" 86.4 cm – 96.5 cm	<b>4″</b> 10.2 cm
Example Part # for Wrap Around Apron Wrap Around – Small Female – TL Lite Ply	<b>22″</b> 55.88 cm	38" – 42" 96.5 cm – 106.7 cm		22″ 55.88 cm	38" – 42" 96.5 cm – 106.7 cm	<b>5″</b> 12.7 cm
Model - Size X Length - Protection	24″ 60.96 cm	<b>42" – 46"</b> 106.7 cm – 116.8 cm		24″ 60.96 cm	42" – 46" 106.7 cm – 116.8 cm	5″ 12.7 cm
WA - SF X 36 - TL	27″ 68.58 cm	<b>47" – 52"</b> 119.4 cm – 132.1 cm		<b>27</b> ″ 68.58 cm	47" – 52" 119.4 cm – 132.1 cm	<b>6″</b> 15.2 cm
Unisex Wrap Around Aprons	Male Wrap Ar	ound Aprons		Female	Wrap Around	Aprons
Size Chest Waist Hips Siz	ze Chest	Waist	Hips	Size Che	st Wai	st Hips

			up to				up to				up to
Small	34" – 37"	27" – 31"	<b>39″</b>	Small	35″ – 38″	28" – 32"	<b>40″</b>	Small	34" – 35"	26" – 30"	<b>38″</b>
(S)	86.4 cm – 94 cm	68.6 cm – 78.7 cm	99.1 cm	(SM)	88.9 cm – 96.5 cm	71.1 cm – 81.3 cm	101.6 cm	(SF)	86.4 cm – 88.9 cm	66.4 cm – 76.2 cm	96.5 cm
Medium	38" – 41"	31" -35"	<b>43″</b>	Medium	<b>39″ – 42″</b>	<b>32″ -36″</b>	<b>44″</b>	Medium	36" – 37"	30" – 34"	<b>40″</b>
(M)	96.5 cm – 104.1 cm	78.7 cm – 88.9 cm	109.2 cm	(MM)	99.1 cm – 106.7 cm	81.3 cm – 91.4 cm	111.8 cm	(MF)	91.4 cm – 94 cm	76.2 cm – 86.36 cm	101.6 cm
Large	42" – 45"	35" – 39"	<b>47″</b>	Large	<b>43" – 46"</b>	<b>36" – 40"</b>	<b>48″</b>	Large	38" – 40"	34" – 38"	<b>43″</b>
(L)	106.7 cm – 114.3 cm	88.9 cm – 99.1 cm	119.4 cm	(LM)	109.2 cm – 116.8 cm	91.4 cm – 101.6 cm	129.5 cm	(LF)	96.5 cm – 101.8 cm	86.36 cm – 96.52 cm	109.2 cm
X-Large	<b>46" – 49"</b>	<b>39″ – 43″</b>	<b>51″</b>	X-Large	<b>47" – 50"</b>	40" – 44"	<b>52″</b>	X-Large	<b>41" – 44"</b>	<b>38" – 42"</b>	<b>47″</b>
(XL)	116.8 cm – 124.5 cm	99.1 cm – 109.2 cm	129.5 cm	(XLM)	119.4 cm – 127 cm	101.6 cm – 111.8 cm	132.1 cm	(XLF)	104.1 cm – 111.8 cm	96.52 cm – 106.68 cm	119.4 cm
XX-Large	<b>50" – 53"</b>	43" – 47"	<b>55″</b>	XX-Large	51" – 54"	<b>44" – 48"</b>	<b>56″</b>	XX-Large	<b>45" – 48"</b>	<b>42" – 46"</b>	<b>51″</b>
(XXL)	127 cm -134.6 cm	109.2 cm -119.4 cm	139.7 cm	(XXLM)	129.5 cm -137.2 cm	111.8 cm -121.9 cm	142.2 cm	(XXLF)	114.3 cm –121.9 cm	106.68 cm -116.84 cm	129.5 cm
	Your Approximate Height         4'11" – 5'3"         5'3" – 5'5"         5'5" – 5'8"         5'8" – 5'11"         5'11" – 6'2"           149.9cm – 160 cm         160 cm – 165.1 cm         165.1 cm – 172.7 cm         172.7 cm – 180.3 cm         180.3 cm – 187.9 cm										

Tour Approximate neight	149.9cm – 160 cm	160 cm – 165.1 cm	165.1 cm
Choose Length*	<b>30″</b> 76.2 cm	<b>32″</b> 81.3 cm	34″ 86.4 cm
*Wrap Around Aprons are 2" she	orter in back		

Sizing

Custom made (non-stock) aprons are non-returnable.

#### **Unisex Vest Skirt**

Size	Chest	Vest Length	Waist	Hips	Skirt Lengt
Small	34" – 37"	18″	27" – 31"	<b>39″</b>	21″
(S)	86.4 cm – 94 cm	45.7 cm	68.6 cm – 78.7 cm	99.1 cm	53.3 cm
Medium	38" – 41"	19″	31" – 35"	<b>43″</b>	22″
(M)	96.5 cm – 104.1 cm	48.3 cm	78.7 cm – 88.9 cm	109.2 cm	55.9 cm
Large	42" – 45"	<b>20″</b>	35" – 39"	<b>47″</b>	22″
(L)	106.7 cm – 114.3 cm	50.8 cm	88.9 cm – 99.1 cm	119.4 cm	55.9 cm
X-Large	<b>46" – 49"</b>	21″	<b>39″ – 43″</b>	51″	22″
(XL)	116.8 cm – 124.5 cm	53.3 cm	99.1 cm – 109.2 cm	129.5 cm	55.9 cm
XX-Large	<b>50" – 53"</b>	22″	<b>43″ – 47″</b>	<b>55″</b>	23″
(XXL)	127 cm -134.6 cm	55.9 cm	109.2 cm -119.4 cm	139.7 cm	58.4 cm

#### **Male Vest Skirt**

Size	Chest	Vest Length	Waist	Hips up to	Skirt Length
Small	35" – 38"	18″	28" – 32"	<b>40″</b>	21″
(SM)	88.9 cm – 96.5 cm	45.7 cm	71.1 cm – 81.3 cm	101.6 cm	53.3 cm
Medium	<b>39″ – 42″</b>	<b>19″</b>	32" -36"	<b>44″</b>	22″
(MM)	99.1 cm – 106.7 cm	48.3 cm	81.3 cm – 91.4 cm	111.8 cm	55.9 cm
Large	43" – 46"	<b>20″</b>	<b>36" – 40"</b>	<b>48″</b>	22″
(LM)	109.2 cm – 116.8 cm	50.8 cm	91.4 cm – 101.6 cm	129.5 cm	55.9 cm
X-Large	47" – 50"	21″	<b>40" – 44"</b>	<b>52″</b>	22″
(XLM)	119.4 cm – 127 cm	53.3 cm	101.6 cm – 111.8 cm	132.1 cm	55.9 cm
XX-Large	51" – 54"	<b>22″</b>	<b>44″ – 48″″</b>	<b>56″</b>	23″
(XXLM)	129.5 cm -137.2 cm	55.9 cm	111.8 cm -121.9 cm	142.2 cm	58.4 cm

	Example Part # for Vest Skirt Combo - TL Lite Ply							
		Vest Skirt -		Small Male -		TL Lite Ply		
	Mode		Model - Size -		Protection			
Female V	est Skirt	VS	-		SM -		TL	
Size	Ch	est	Vest Length		Waist	Hips up to	Skirt Length	
Small (SF)	34" – 35" 86.4 cm – 88.9 cm		17″ 43.2 cm	26" – 30" 66.04 cm – 76.2 cm		38″ 96.5 cm	<b>20″</b> 50.8 cm	
Medium (MF)	36" – 37" 91.4 cm – 94 cm		18″ 45.7 cm	<b>30″-34″</b> 76.2 cm – 86.36 cm		<b>40″</b> 101.6 cm	20″ 50.8 cm	
Large (LF)	38" – 40" 96.5 cm – 101.8 cm		19″ 48.3 cm	86.3	34" – 38" 6 cm – 96.52 cm	<b>43″</b> 109.2 cm	21″ 53.3 cm	
X-Large (XLF)	41" – 44" 104.1 cm – 111.8 cm		20″ 50.8 cm	<b>38″ – 42″</b> 96.52 cm – 106.68 cm		<b>47″</b> 119.4 cm	21″ 53.3 cm	
XX-Large (XXLF)	45" – 48" 114.3 cm –121.9 cm		21″ 53.3 cm	42" – 46" 106.68 cm -116.84 cm		<b>51″</b> 129.5 cm	<b>22″</b> 55.9 cm	
Size	Armhole Unisex Male Armhole Female							
Small		9.5″	24.1 cn	n	8″	20	.3cm	


36″

91.4 cm

38″

96.5 cm

Size	Armhole	Unisex Male	Armhole F	emale
Small	9.5″	24.1 cm	8″	20.3cm
Medium	10.5″	26.7 cm	8.5″	21.6 cm
Large	11″	27.9 cm	9″	22.9 cm
X-Large	11.5″	29.2 cm	9.5″	24.1 cm
XX-Large	12″	30.5 cm	10″	25.4 cm



VESTS-SKIRTS-WRAP		RL - REGULAR LEAD <sup>1</sup>							TL – TL Lite Ply <sup>1</sup> NL – Non Lead <sup>2</sup>				
Projection when ordering.			Front		Front Overlap*		Back		Front		Front Overlap*		Back
		t /Wrap rt:	0.5mm l 0.3mm l		1.0mm P 0.6mm P	•	0.3mm Pb 0.3mm Pb	•	0.5mm Pb 0.25mm Pt	•	1.0mm Pb 0.5mm Pb	•	0.25mm Pb Eq. 0.25mm Pb Eq.
standard.			RLL					TLL		Ν	LL		
			Front		Front Overlap*		Back		Front		Front Overlap*		Back
*Overlap does not apply to Flex Back Vest or Wrap Around Flex Back Apron		t /Wrap rt:	0.35mm 0.25mm		0.7mm P 0.5mm P		0.25mm Pl 0.25mm Pl	•	0.35mm Pk 0.25mm Pk	•	0.7mm Pb 0.5mm Pb	•	0.25mm Pb Eq. 0.25mm Pb Eq.
			RL - ALT					TL-ALT		NL	NL-ALT		
			Front Front Overlap*		erlap*	Back		Front		Front Overlap*		Back	
		t /Wrap rt:	0.25mm Pb eq 0.5mm Pb eq 0.25mm Pb eq 0.5mm Pb eq		0.25mm Pb eq 0.25mm Pb eq		0.7mm Pb eq 0.5mm Pb eq		0.35mm Pb eq 0.25mm Pb eq		0.7mm Pb eq 0.5mm Pb eq		
FULL OVERLAP STYLE-		RL - RE	GULAR LE	R LEAD <sup>1</sup> TL			TL – TL	L Lite Ply <sup>1</sup>		NL – Non Lead <sup>2</sup>			
VESTS & WRAP AROUND	F	Front	Front		Overlap* Back		Front		Front		t Overlap* Back		
Add Proper Initials for Optional Protection when ordering, otherwise Regular Lead is standard.		0.25mm	n Pb eq	5.0mm	Pb eq	0.25mr	n Pb eq	0.25mr	n Pb eq	0.5mm	Pb eq	0.25mr	n Pb Eq.
			Alt	ernate fro	ont or bac	k protect	on also available in 0.3mm, 0.35mm & 0.5mm Pb equivalencies.						

<sup>1</sup> Denotes Lead Equivalency of RL, RLL & RL-ALT Regular Lead and TL, TLL & TL-ALT Lite Ply @ 100 kVp. <sup>2</sup> Denotes Lead Equivalency of NL, NLL & NL-ALT Non Lead @ 100 kVp.

NOTE: For Proper protection in your working environment please contact a radiation physicist regarding adequate apron Pb equivalency

Contact us regarding prices on alternate protection options. Discounts may apply when reducing Pb equivalency.

## **Prints** Perfect for your Personality and Profession

DESIGNER FABRICS (DF)





Bermuda (DF)



Butterflies (DF)

Patriotic (DF)

Kaleidoscope (DF)



Owls (DF)

American Heart (DF)

Route 66 (DF)

Chocolate Floral (DF)





Aztec Sun (DF)

Red & White Stars (DF)

Tiki Blue (DF)





Lollipop Floral (DF)

c

Č,

Royal Paws (DF)









Hearts of Color (DF)

Splatter (DF)

1

Dog Paws (DF)

Green Surf (DF)

Denim (DF)





Dog Bones (DF)



Compass (DF)







Wild Cat Pink (DF)







Red Plaid (DF)



**DESIGNER FABRICS** (DF) (Continued)

Blue Squares (DF)



DIAMOND **RIPSTOP** (DRS)





Red (DRS)

Purple (DRS)

Royal (DRS)





Lavender (DRS)

Cinnamon (DRS)

Fushia (DRS)



Fabric and binding colors are approximate representations of actual colors/patterns available. Economy aprons (Model E) are available in Vinyl (V), Nylon (NY) or Diamond Ripstop (DRS) only. Diamond Ripstop (DRS) and Designer Fabrics (DF) feature a black nylon inside cover and black binding. Weblon (W) aprons feature a black or corresponding inside nylon cover. Black binding is standard on most apparel, except for royal binding on royal nylon and light blue vinyl, unless otherwise indicated from binding selection at time of order. Two tone coverings subject to up charge. To avoid snagging or tearing of fabric, all hook and loop fasteners must be secured properly.



Elephants (DF)







Red Flames (DF)



Desert Camo (DF)













Storm Camo (DF)

\*More prints available at shieldingintl.com

www.shieldingintl.com

Zebra (DF)



Athletes (DF)







Gray Paisley (DF)





Farm Animals (DF)

Faux Leopard (DF)











Pink Camo (DF)









Gray Squares (DF)



Black (DRS)



Slate Blue (DRS)



Green Plaid (DF)



Charcoal (DRS)



Marine Blue (DRS)



Hot Pink Plaid (DF)



Silver (DRS)



Kelly Green (DRS)



Tan (DRS)



Brown (DRS)

### **Fabrics & Colors**



Hot Pink (NY)





Mauve (NY)









Black (NY)

Purple (NY)



Navy (NY)









Lime Green (NY)

Tan (NY)

Coffee (NY)



Charcoal (NY)

Silver (NY)





## **Other Embroidery Options**



**Bottom Corner** 

Collar Bib



Shielding Intl.

Garamond

Shielding Intl.

Block

Shielding Intl.

Script

Shielding Intl.

Micro Scan

Shielding Intl.

Impress

Can't find your letter style?

Call for assistance.

location for one price!

		broideri	
Fo	or additional design	ns (under 10,000 stite	ches) see Dakota (
	Boots	Bone & Paw	Butterfly
		- 	
_	Football Helmet	Golf Crest	Sports
	×	A CONTRACTOR	The
	Swirl Sun	Tropical Fish	Palm Trees
_	Heart Bear	Heart	Bandaged Hear
			Å
	Skull and Bones	X-Ray Man	Skeleton
	X	and the second s	
_	Neuro	Neurology	Cardiology
_	Bone & Joint	Spine	Caducus
	RI	XIRAX	Sath



### igns

Collectibles website: www.dakotacollectibles.com

## **Custom Apron Options**

#### **Pocket Styles**



Patch Pocket



Slit Pocket



ID Pocket



Sleeves

**Optional Sleeve—Model 812** For upper arm protection. Features inner fastener for girth adjustment. Permanently sewn onto apron or select a hook and loop or snap fastener for easy attachment or removal.

Model Length – Measure from Shoulder to Elbow No.



Non-Slip Skirt Backing

# **Axillary Wings**



Model AW Permanently sewn onto apron, Axilliary Wings provide added protection without sacrificing mobility.

## **Ordering Information**

#### **Fast and Friendly Service**

referral to a select dealer nearest you.



Call Us! 1-800-292-2247 or 541-475-7211 Monday – Friday: 7am – 3:30pm PST



Email Us! info@shieldingintl.com We'll Reply ASAP

#### **Manufacturing and Delivery Time:**

Stock **SHIELDING INTERNATIONAL™** products ordered with select stock coverings and closures will usually be shipped within 24-48 hours of order receipt. These items include Economy & Standard Aprons, Half & Demi Aprons, Thyroid Collars, Gonad Shields, Gloves, Apron Racks, Mobile Shields, CT Shields, Non Rx Glasses and other select items. We will make every effort to ship all other custom made orders, with exception of high volume and special orders (ie., prescription glasses, lead or acrylic glass or frames and barriers) within 10 to 12 working days of confirmed order.

this service.

#### Shipping:

Orders are normally shipped with a leading major transportation carrier via ground service, unless otherwise specified, and ship F.O.B. Madras, Oregon or from a predetermined special order location(s).

claim is filed due to shipping damage.

#### **Apron Care & Guidelines:**

To maintain safe performance, as well as keep your x-ray apparel looking good for many years to come, SHIELDING **INTERNATIONAL™** strongly recommends that each garment be thoroughly inspected upon receipt and at regular intervals, as well as properly stored when not in use. Each piece of SHIELDING INTERNATIONAL<sup>™</sup> apparel is accompanied by an Apron Care & Guideline sheet for guidance in inspection and proper care of our products.

- on a SHIELDING heavy duty chrome hanger (Item# 200H) or equivalent, not on a hook.
- loop closed.

Failure to properly care for your apron may result in harm to the apron, as well as void its warranty.

**Guarantee & Warranty Disclaimer** SHIELDING INTERNATIONAL<sup>™</sup> warrants all Regular Lead and TL Lite Ply x-ray protective aprons from material defects and workmanship for a period of (2) years from date of purchase. (Non Lead aprons are limited to a 1 year warranty.) This guarantee is limited to the repair or replacement of the garment in case of manufacturer defect. The guarantee does not cover normal wear and tear, failure to properly care for the product, or any other

responsibility. **Apron Recycling Program:** 

# of your old lead aprons.

www.shieldingintl.com

and shoulder strap.

Apron Bag—Model AB

Conveniently holds one apron.

wrap around a sturdy cardboard

tube, while keeping it safe and secure. Includes carrying handle


(Size: 8"W x 29"L / 20.3cm x 73.7cm)

Tube shape design allows apron to





8" to 16" (20.3cm to 40.6cm) 812



### **Additional Apron Options**

- Attached Thyroid Collar (as shown on page 19)
- Tethered Thyroid Collar
- (as shown on page 19)
- Tu-Tone Covering (Specify Inside & Outside Fabric Choice)
- Removable Name Tag
- Additional 'QR' Quick-Release Buckle
- Skirt Suspenders



(Specify 1" or 2") • Built-in Vest Belt

#### To Order: Orders may be placed through a preferred or nearby SHIELDING INTERNATIONAL™ dealer. Feel free to contact us for a



Visit Us on the Web! www.shieldingintl.com



Fax Us! 541-475-6628 All Day, Every Day

Mail Us! Shielding International 2150 NW Andrews Drive Madras, OR 97741 USA

• Need it in a rush? SHIELDING INTERNATIONAL<sup>™</sup> offers a 3-day RUSH SERVICE for 15% on total cost of order or a 5 day RUSH SERVICE for 10% of the total cost. (min. \$10.00 rush fee per order). Contact us regarding large guantity rush orders to verify

• IMPORTANT: Be sure to inspect all merchandise at time of delivery, report any damage immediately, note any discrepancies on the Bill of Lading (if applicable) and RETAIN ALL PACKAGING MATERIAL (whether damage is visible or concealed) in the event a

• STORAGE: Don't fold, crease, drape or sit down in your SHIELDING garments. We strongly recommend that our aprons be hung

• Hook and loop closures must be secured properly to avoid snagging or tearing of fabric. Always store apron with hook and

Take advantage of SHIELDING INTERNATIONAL™'s free apron recycling program for the return/recycle of your retired lead aprons! Simply contact us for a return authorization to arrange for the recycling of your old aprons. All you pay is the freight and SHIELDING INTERNATIONAL<sup>™</sup> will assure the safe and environmentally friendly disposal



2150 NW Andrews Dr. Madras, OR 97741-0120 Toll Free (800) 292-2247 Tel: (541) 475-7211 Fax (541) 475-6628 Email info@shieldingintl.com www.shieldingintl.com

Distributed By: