



100US

Gyn/Fertility Mobile Imaging Table
Aspect™ Series table line



Ideal C-Arm imaging solution for Hysterosalpingograms

- Carbon fiber tabletop with low x-ray attenuation.
- Detachable boarding extension.
- Unobstructed fluoroscopy imaging area: 21"W x 32"L (53 x 81 cm).
- Motorized tabletop elevation: 27" – 38" (66 – 96,5 cm).
- 24"W x 46"L tabletop (61 x 117 cm).
- Patient weight capacity: 450 lbs (204 kgs).
- Standard US table side-rail: compatible with many lithotomy leg holders and clamps.
- 120 VAC operation with battery back-up operation (230VAC optional).
- Easy table transport with 5" (12,7 cm) locking casters.
- Table weight = 375 lbs (170 kgs)
- Warranty: One year labor / three years parts*



Accessories included with package A100-1479

- Basic knee crutches with two side rail clamps
- Tabletop pad
- Foot control
- Patient restraint straps

Optional accessories available:

- Armboard
- Table hand control

idi
Image Diagnostics, Inc.
Safety, Simplicity, Quality

IDI IMAGING TABLES

IDI Interventional / Endovascular Tables

IDI's Aspect line of interventional and endovascular tables deliver extended positioning and imaging capabilities.



IDI Urology Tables

For over 25 years IDI has manufactured the market leading mobile table for Urology.



IDI Spinal Imaging Pain Management Tables

C-arm tables for GI and pain management procedures. Tables with tilt, roll and longitudinal tabletop movements are available.



IDI VIDEO INTEGRATION SYSTEMS

ilex55 Video Integration System

The first mobile large field 4K UHD display solution with built in integration. A State-of-the-Art Display Solution for Vascular, GI and Neurosurgery.



MDS Video Integration System

Advanced Video Integration for multi-modality image visualization in Urology, GI & Vascular surgery.



ilex32 Monitor Stand

Portable image display when you just need a bigger monitor. Available in both HD and UHD-4K versions



Image Diagnostics, Inc.

310 Authority Drive Fitchburg, MA, 01420 USA Tel: 978-829-0009 Fax: 978-829-0027
www.imagediagnostics.com *Please reference IDI Warranty Document D000-2070 L100-3115A

idi
Image Diagnostics, Inc.
Safety, Simplicity, Quality