

SSCOR Duet

A Powerful Substitute for Wall Vacuum in an Emergency



The **SSCOR Duet®** is designed for the hospital crash cart, patient transport, surgi-center backup, long term care facilities and more.



Easily Attached to Walls or Carts
With the optional **SSCOR Duet** retention bracket.



Powerful and Portable
With adjustable vacuum.



Universal Canister Holder
Enables you to use the canister of your choice for standardized inventory control.

SSCOR

Complete information at www.sscor.com.

The **SSCOR** Duet battery care system will automatically charge and monitor itself.

Leave it plugged in all the time. The battery cannot be overcharged.



SSCOR Duet with Retention Bracket - Model 2314B



SSCOR Duet Flat Back - Model 2314

The **SSCOR** Duet® contains a battery care system that is designed to gently hold the battery at optimal charge, without harmful overcharging.

While plugged into AC power, the internal battery control system will:

- Bring the battery to full charge, illuminating the green light on the side of the unit
- Drop back to “float” when the battery reaches full charge, eliminating damage or destruction of the battery caused by overcharging

The battery care system monitors the battery during operation, informing the operator of anticipated running time and will disconnect a depleted battery to protect the battery from destructive deep discharge.

The **SSCOR** Duet’s automatic capacity check will warn you when the battery begins to lose capacity. (Eventually all batteries lose capacity and need to be replaced.)

SSCOR will replace a **SSCOR** Duet battery that does not operate to specifications within three years of a **SSCOR** Duet purchase.

Specifications:

Performance:	>30lpm airflow at the end of a patient tube; <-50mmHg - >-550mmHg
Dimensions:	17" x 9" x 5.25"
Weight:	10.05 lbs. (Model 2314), 10.65 lbs. (Model 2314B)
Charger:	Internal automatic charger offers fast charge up to full and then gentle float to maintain charge. LED's indicate charging status and operation mode.
Power:	12v DC Sealed Lead Acid Battery, 100 - 240V AC Input
Run Time:	45 minutes +/- 10% when battery fully charged at full capacity

Standard items included with product:

Canister, 6' patient tube, HI-D® suction tip, AC power cord and Quickstart Guide

SSCOR

SSCOR, Inc.
11064 Randall St. Sun Valley, CA 91352
+1 818 504 4054
Lit-2314-3-B

JD HONIGBERG INTERNATIONAL INC.
Deerfield, IL 60015
E-mail: medical@jdhintl.com
www.jdhmedical.com
U.S. Patent No. D270948



SG-80-SS

utility table with shelf, drawer and rail

overall size

20"W x 16"D x 34"H

drawer size

15¼"W x 12½"L x 3"H (I.D)



SG-80-A-SS

utility table with shelf, two drawers and rail

overall size

20"W x 16"D x 34"H

drawer size

15¼"W x 12½"L x 3"H (I.D)



SG-80-B-SS

no rail

utility table with shelf and drawer

overall size

20"W x 16"D x 34"H

drawer size

15¼"W x 12½"D x 3"H (I.D)



SG-80-C-SS

no rail

utility table with shelf and two drawers

overall size

20"W x 16"D x 34"H

drawer size

15¼"W x 12½"D x 3"H (I.D)



SG-81-SS

utility table with shelf and rail

overall size

20"W x 16"D x 34"H

features

- Convenient and easy to use
- Stainless steel, Steril-Gard construction eliminates full sub-tops
- Improves the control of undesirable substances
- Full extending drawer hangers (except SG-81-SS)
- Sound-proofed reinforcing channels for quiet use and added strength
- 3" double ball bearing, rubber casters for silent mobility
- 1½" high 3-sided top rail

options

- Part No. 46-B Drawer lock with key
- Part No. 06 3" casters with brakes

features

- Stainless steel Steril-Gard construction eliminates full sub-tops where undesirable substances can collect
- Easy to clean
- Lightweight but strong
- 3" rubber, double ball bearing casters for silent mobility
- Soundproofed reinforcing channels for quiet use

Mid-Size Tables with Shelf

	Dimensions
SG-82-SS	20"W x 16"D x 34"H
SG-84-SS	30"W x 16"D x 34"H
SG-86-SS	33"W x 18"D x 34"H
SG-88-SS	36"W x 20"D x 34"H



SG-86-SS

Mid-Size Tables without Shelf

	Dimensions
SG-83-SS	20"W x 16"D x 34"H



SG-87-SS

options

- 3 or 4 sided top and shelf rails available for all tables.
- Stainless steel drawers, or tote boxes may be mounted under top.
- Table height can be modified

Tote Boxes

Impact-resistant, high-density polyethylene.

- CDS-TB-5" (17¾"L x 14½"W x 5⅞"D)
- CDS-TB-10" (17¾"L x 14½"W x 10"D)

Tote Box Hangers


- SG-TBH-1 16"D SG Tables
- SG-TBH-2 18"D SG Tables
- SG-TBH-3 20"D SG Tables
- SG-TBH-4 24"D SG Tables

Glides/ Rubber Feet/ Casters

- #15 Glides
- #13 Rubber Feet
- #06 Locking Casters

IV Pole for SG Tables

- P-7202-1 available on lower shelf models only.



Contact Info: **JD Honigberg International, Inc.** Email: medical@jdhintl.com



4000 SE Columbia Way • Vancouver, Washington 98661



- Convenient foot-pedal operated lock adjusts height from 37" to 62½"
- All welded stainless steel construction
- Base designed to move under low clearance equipment with maximum stability
- 2" dual wheel swivel casters



P-1065-S

- Removable extra large 20" x 25" stainless steel tray for easy cleaning
- 50 lbs. weight capacity evenly distributed



P-1066-SS


- Removable 16¼" x 21¼" stainless steel tray for easy cleaning
- 50 lbs. weight capacity evenly distributed



- Removable large 16¼" x 21¼" stainless steel tray for easy cleaning
- Four 2" dual wheel swivel casters for increased mobility
- 50 lbs. weight capacity evenly distributed

Contact info: **JD Honigberg International, Inc.**

Email: medical@jdhintl.com





Pedigo's newly designed lidded hampers feature the Quiet-Close™ lid dampener for improved safety and added control of airborne bacteria and odors.

Linen hampers are constructed of heavy gauge steel and feature a triangular shape for efficient use of space. Double ball bearing 3" swivel casters make for easy mobility. Strong, heavy gauge shelf keeps bag from touching floor and simplifies bag removal.

**Dimensions**

32"H with a 57¼" perimeter

- Stainless steel
- Includes rolling bumpers
- Foot operated
- Accommodates 18" bags

Pedigo's Quiet-Close™ hydraulic lid dampener slows the closure of the lid, reducing the spread of bacteria and protects fingers and hands from being pinched.

- Foot operated
- Eliminates slamming
- Quiet operation
- Tested for 50,000 cycles



18" Bag Hamper Opening




warranty

Contact Info: **JD Honigberg International, Inc.** Email: medical@jdhintl.com

800.822.3501 • FAX 360.696.1700 • www.pedigo-usa.com



4000 SE Columbia Way • Vancouver, Washington 98661



05/16 MM-70

- Heavy, tip-resistant, all welded design
- 2" swivel dual-wheel casters for easy mobility
- 21¼" diameter base on all stands
- Tube: O.D. = 1"
- Easy to clean surfaces

- Stainless steel
- Hand operated, friction lock height adjustment
- Height adjustable from 68½" to 117½"
- Five leg base
- 70 lbs. overall weight capacity
- 10 lbs. per hook
(recommended 30 lbs. max. evenly distributed on top hooks for stability)



- Stainless steel
- Quick and easy one-hand foot pedal height adjustment
- Height adjustable from 68½" to 117½"
- 70 lbs. overall weight capacity
- 10 lbs. per hook
(recommended 30 lbs. max. evenly distributed on top hooks for stability)




Contact info: **JD Honigberg International, Inc.**

Email: medical@jdhintl.com

800.822.3501 • FAX 360.696.1700 • www.pedigo-usa.com



4000 SE Columbia Way • Vancouver, Washington 98661



chrome series

weight capacity rating of 135 kgs.

dual wheel locking casters option

- 66 Locking Casters - (activated upon standing)
- 66-A Locking Casters - (activated upon sitting)
- 01-B Adaptor (required for both 66 & 66-A)



- **P-36**
 - 14" round cushioned seat
 - 2" double ball bearing casters
 - Height adjustment (screw shaft) of 17, 7" to 26 1/4"

foot-operated

weight capacity rating of 135 kgs.

locking casters option

- 66 Locking Casters - (activated upon standing)
- 66-A Locking Casters - (activated upon sitting)



- **P-528-GS**
 - 14 1/2" round fully padded seat
 - 23 1/2" diameter aluminum base
 - 2" dual wheel casters
 - Height adjustment (pneumatic) of 19 1/4" to 24 1/2"

five series

weight capacity rating of 135 kgs.

locking casters option

- 66 Locking Casters - (activated upon standing)
- 66-A Locking Casters - (activated upon sitting)



- **P-521-GS**
 - 14½" round fully padded seat with back
 - 23½" diameter aluminum base
 - 2" dual wheel casters
 - Height adjustment (pneumatic) of 19¼" to 23¾"
 - Adjustable backrest



- **P-526-GS**
 - 14½" round round fully padded seat
 - 23½" diameter aluminum base
 - 2" dual wheel casters
 - Height adjustment (pneumatic) of 19¼" to 23¾"

Medical PPE Storage Carts

Infection Control PPE Cart for Coronavirus Disease.


Medical PPE Storage Cart allows critical COVID-19 PPE inventory for infection control to be stored in dedicated isolation wards or quickly travel to temporary isolation units.

Infectious disease outbreaks, such as Coronavirus Disease 2019, may strain central supply and limit your department's ability to respond.

Standard and Transmission-Based Precautions for Infection Control of Isolation patients require proper PPE usage. Most notably is that new PPE must be donned with each entry into a patient room¹

The MASS™ Medical PPE Storage Cart has enough capacity to store large quantities of Coronavirus PPE including: Goggles or disposable full-face shields, NIOSH-approved N95 respirators, Gowns, and Nonsterile gloves²

1. The Centers for Disease Control (CDC). (2019). 2007 Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Healthcare Settings [PDF file]. Retrieved from <https://www.cdc.gov/infectioncontrol/pdf/guidelines/isolation-guidelines-H.pdf>
2. The Centers for Disease Control (CDC). (2020). COVID-19 Personal Protective Equipment (PPE) for Healthcare Personnel [JPG file]. Retrieved from <https://www.cdc.gov/coronavirus/2019-ncov/images/COVID-19-PPE.jpg>



Infection Control PPE Cart for Coronavirus Disease
66cm D x 115cm W x 198cm H

Materials: Powder Coated Metal. Door: Hinged Glass; Keyed Locking.

(2-3 days Lead Time)

JD HONIGBERG INTERNATIONAL, INC

Deerfield, IL 60015 USA

Email: medical@jdhintl.com / Tel : (1) 847-412-0200

Point of Use Cart

Modular Surge Supply Storage

Point of Use Carts provide your facility with scalable storage to manage medical surge events. Temporary surge sites are an atypical environment with unique storage demands. With our transparent Lucent™ High Density Modular Storage supply trays, you can customize each supply tray quickly and easily. Keep your supplies organized, streamlined, and easily accessible.

- Easily adjustable dividers for modifying the size of individual supply compartments
 - Transparent trays provide visibility of stored supplies, increasing efficiency and productivity
 - Easy-to-use labeling system
- Highly maneuverable with lightweight construction and premium casters.

Material

- Antimicrobial HDPE
 - Continuous protection against pathogenic microorganism

Storage System

- Exchange panels
 - Effortlessly reconfigure interior components
- Lucent™ High Density Modular Storage supply trays
 - High visibility of stored supplies
 - Increase efficiency and productivity

Mobility

- Premium casters
 - Highly maneuverable
 - Durable construction



Modular Surge Supply Storage
65cm D x 47cm W x 33.84cm H

(2-3 days Lead Time)

JD HONIGBERG INTERNATIONAL, INC

Deerfield, IL 60015 USA

Email: medical@jdhintl.com / Tel : (1) 847-412-0200



Expand Capacity to Treat COVID-19 Patients

With the rise of the COVID-19 pandemic, Healthcare systems are rapidly being stretched and some overwhelmed as the demand for critical respiratory care expands beyond intended capacity. Elective procedures are being postponed to free needed capacity, and in locations where the crisis is intense, caregivers are having to turn away some infected patients.

Our treatment recliners are engineered for 24/7 use, and can help you respond to the need for more capacity to evaluate, triage, or treat patients. Champion is increasing production on a model that is ideal to relieve the pressure from this critical situation to get it to you as soon as possible. This manual recliner is made with buck-suede vinyl and can ship in 2-3 days.



CLINICAL BENEFITS:

- Durable steel frame for **24/7 use**.
- Can be used for **patient transport**.
- Swing-away arms for **lateral patient transfer**.
- Wall-saver mechanism, **saving floor space**.
- Open design is **easy to clean and disinfect**.
- Flat-body Trendelenburg for acute episodes.
- Two, fold-away side tables.
- Cost-effective capital purchase that can be absorbed into different departments within the healthcare system after the crisis is over.

Now more than ever, we are fully committed to our customers, and to whatever positive role we can have in halting the COVID-19 pandemic.

JD HONIGBERG INTERNATIONAL, INC.

Deerfield IL 60015 USA – Tel: (1) 847-412-0200 – Email: medical@jdhintl.com


Rapid Response Cart

Available to purchase for facilities responding to COVID-19 issues. This cart is offered during natural disasters and critical health concerns to provide a secure, mobile cart to meet immediate demands for storage.

* RAPID RESPONSE SHIPPING *

This cart is offered in the **color Sand** and **ships within 48-72 hours** from receipt of your order. For quantities over five, please contact us.

Dimensions: 34" H x 29.5" W x 22" D (86.4 cm H x 74.9 cm W x 55.9 cm D)



RRC344SD

Emergency Carts

Color: Red

Lead Time:

5 Business Days

Quantity:

5 or less



6400Q

Dimensions:

45.75" H x 32" W x 22" D
(116.2 cm H x 81.3 cm W x 55.9 cm D)

Dimensions:

38.5" H x 32" W x 22" D
(97.8 cm H x 81.3 cm W x 55.9 cm D)



6301Q

Dimensions:

44.5" H x 38" W x 22" D
(113 cm H x 96.5 cm W x 55.9 cm D)



6401Q

Infection Control and Procedure Carts

Color: Yellow

Lead Time:

5 Business Days

Quantity:

5 or less



3154KQ

Dimensions:

40.75" H x 19.25" W x 19.75" D
(103.5 cm H x 48.9 cm W x 50.2 cm D)

Color: Yellow

Lead Time:

6 Business Days

Quantity:

5 or less



6450YL

Dimensions:

42.5" H x 32" W x 22" D (107.95 cm H x 81.3 cm W x 55.9 cm D)

Color: Sand **Quantity:** 5 or less

Lead Time: 5 Business Days



6450SD

Dimensions:

42.5" H x 32" W x 22" D (107.95 cm H x 81.3 cm W x 55.9 cm D)

COVID-19 WARMING SOLUTIONS

All items ship within 2-3 weeks

ivNow Fluid Warmers

- Can provide from 25 to 75 liters of warmed fluid during a 12-hour period.
- Easy-to-Use: Place the IV bag in the pod to activate the unit. Fluid is automatically warmed up to 104°F (40°C) in 30 minutes or less.
- Capacity: 1-3 pods (holds up to 3-liter bag size)
- Configurations: Countertop, mobile equipment pole stand, or wall.
- Available via UPS shipment.



ivNow-1



ivNow-2



ivNow-3
with equipment pole

Fluid Warmers - Cleared for IV and irrigation fluids



EC250L

- Capacity: 12 1-liter bottles or 16 1-liter bags
- 56 cm H; 47 cm W; 60 cm D
- Available via UPS shipment
- Countertop mode
- Stackable



EC350L

- Capacity: 24 1-liter bottles/ bags
- 71 cm H; 47 cm W; 60 cm D
- Floor or countertop model
- Stackable

Blanket Warmers



EC250

- Capacity: 4-6 blankets
- 56 cm H; 47 cm W; 60 cm D
- Available via UPS shipment
- Countertop mode
- Stackable



EC350

- Capacity: 10-12 blankets
- 71 cm H; 47 cm W; 50 cm D
- Floor or countertop model
- Stackable

Contact Information

JD Honigberg International, Inc.
Deerfield, IL 60015 USA
medical@jdhintl.com
Tel: (1) 847-412-0200



SHIELDING^{INC.}

INTERNATIONAL

SPECIFICATIONS SHEET

Frontal Protection Apron model ECONOMY (E)



- Apron is available with 2,5cm tie strap closures only.
- Available in our selection of solid color nylon, vinyl and our wide selection of diamond ripstop fabrics! Optional protection or customization not available.
- Available for immediate shipment in Royal Blue Nylon (61cm width x 81, 86, 91 & 97cm lengths).
- Fabrics limited to Nylon, Vinyl or Diamond Ripstop only.
- Standard Protection: 0.5mm Pb Equivalency
- Available with Lead, T Lite Lead or Non Lead protective material.

Two Pieces to Lighten the Load



A) Full Overlap Vest Skirt—Model FOVS
Designed with two fully overlapping front panels, this vest and skirt combo reduces the weight by approximately 15% over our traditional and Adjust-A-Fit Vest Skirt Combos. Vest features hook and loop closure and shoulder pads for added comfort. For skirt description, refer to Model SO below.

B) Skirt—Model SO *(when ordered separately)*
Skirt wraps around and overlaps in the front, offering extra protection. Fits securely on hips with extra wide hook and loop fastener and 1" quick release* buckle for added support. Comes standard with large side pocket and hanging loops. *2" quick release buckle available upon request.

Interested in Non-Slip Skirt Backing?
See Custom Options page 40.



C) Adjust-A-Fit Vest Skirt—Model AFVS
Adjust-A-Fit Vest Features wrap around protection and adjustable padded shoulders from the arm hole to the top of the shoulder, enabling the user to adjust the vest to a more secure and comfortable fit. Because of the versatility in fit, it accommodates a wide range of wearers. For skirt description, refer to page 8.



- For Sizing, Protection, Fabrics and Custom Options, refer to Pages 34-40
- Lighter weight protection options are recommended for all Vest Skirt models.



D) Vest Skirt Combo—Model VS
Proud holders of the original design and patent. Our standard Vest features all around protection with overlap in the front, shoulder pads and three panel hook and loop closure.

E) Flex Back Vest Skirt—Model FBVS
Flex Back Vest features wrap around protection with flex back support and adjustable hook and loop flaps at the waist and shoulder for added comfort and support, especially during long procedures. Shoulder pads included. For skirt description, refer to page 8.

Vests & Skirts



Wrap Around Aprons

A



A) Wrap Around Apron—Model WA

Wrap Around Apron is designed for front and back protection. This apron features extra wide shoulders, distributing weight over a greater area to reduce fatigue and three panel hook and loop fastener. Shoulder pads included.

B) Wrap Around Back Relief Apron—Model WABR

Wrap Around Back Relief Apron offering the same features of our Wrap Around Apron, transfers the weight from the shoulders to the hips with the aid of a 6" wide elastic belt. Comes standard with BR Belt, which may be ordered to fit various sized individuals. Indicate belt size from size chart when ordering. Shoulder pads included. BR Belt may be ordered separately. (Model BR)

Adjustable BR Belt Size Chart

Size	Fits Waist
BR-S	28-38" (7.1-96.5cm)
BR-M	36-48" (91.4-121.9cm)
BR-L	44-58" (111.8-147.3cm)

* Lighter Weight Protection Options Recommended

Looking for Extra Back Support?

Order with optional
BRX Belt—Model WABRX
For more info, refer to page 3



C) Wrap Around Full Overlap Apron — Model WAFO

Designed with two fully overlapping front panels, this apron reduces the weight by approximately 15% over our traditional Wrap Around model. Features extra wide shoulders, hook and loop fastener and shoulder pads for added comfort.



- For Sizing, Protection, Fabrics and Custom Options, refer to Pages 34-40
- Lighter weight protection options are recommended on all Wrap Around Models

D) Wrap Around Flex Back—Model WAFB

Features dual wrap elastic flex back with hook and loop side to front flaps while providing front and back protection for an adjustable and secure all around fit. Fits a wide range of sizes. Shoulder pads included.

(Not Pictured) Wrap Around Full Overlap Back Relief Apron — Model WAFOBR

Offering the same features of our Wrap Around Full Overlap Apron, this apron is approximately 15% lighter than our traditional Wrap Around Back Relief model and transfers the weight from the shoulders to the hips with the aid of a 6" wide elastic belt. Comes standard with BR Belt (see opposite page). Indicate belt size from size chart when ordering. BR Belt (Model BR) may also be ordered separately.

Looking for Extra Back Support? Order with optional BRX Belt — Model WAFOBRX (See page 3 for more information)

ThermoSafe®

DiagnosticShipper

95 kPa Bags

Specimen transport bags are liquid tight and tamper evident.

The patented adhesive closure and continuous seal meet or exceed standards of all major air carriers and DOT.

Polyethylene construction and adhesive closure provides added security during handling by personnel.

Model	Inside Dimensions [in.]	Case Pack
951	6 x 9 $\frac{3}{8}$	250
952	9 x 12	250
953	9 x 16	250
954	13 $\frac{1}{4}$ x 17 $\frac{7}{8}$	100
955	10 $\frac{1}{4}$ x 14	250
*958	13 $\frac{1}{4}$ x 17 $\frac{7}{8}$	100

* Bag without biohazard logo.



ThermoSafe®

InfectiousShipper

Meets U.S. Department of Transportation (DOT) and IATA requirements for transporting infectious substances (UN 2814 and UN 2900).

4 models– accommodate shipments from 4 to 6 tubes/vials.

Pre-assembled packaging includes everything needed for shipping infectious agents in accordance with DOT Regulation 49 CFR Parts 171 to 178 and IATA Packing Instructions 602.

*Note: Compliance with UN/DOT regulations governing shipments of infectious substances in these approved shippers is contingent on the use of all supplied components. The Shipper is ultimately responsible for determining the suitability of packaging for use with contents. Consult appropriate regulations. **




Dry Ice Overpack

670

Substance Shippers protect frozen samples for over four days. The EPS foam container and lid accommodate either Model 630 or Model 660B shipper, plus 25 lb [11 kg] of dry ice. Foam insert holds in place while dry ice sublimates. Includes one Class 6 Infectious Substance label and one Class 9 Miscellaneous Dangerous Goods label (for dry ice), and shipping carton.

Model	Vessel/ Dimensions [in.]	Internal Cushioning	Capacity	Outer Dimensions LxWxH [in.]	Supplied Labeling	Additional Features
630	Aluminum can 3 x 6	Foam insert	Four 16 x 100 mm tubes	5 ½ x 5 ½ x 11 ¼	Infectious substance	Includes 4 x 10 ml blood collection tubes
660B		Absorbent material, flexible pressure vessel & bubble wrap material	Six 16 x 100 mm tubes maximum volume 60 ml	5 ½ x 5 ½ x 11 ¼	Infectious substance	Flexible pressure vessel accepts various sized tubes
670	N/A	EPS Insulated Shipper	Overpack for Models 630 and 660B	14 x 14 x 14 ⅝	Infectious substance & miscellaneous dangerous goods	Foam insert holds up to 11.3 kg [25 lb] of dry ice in place
INF2000	0.5 liter transparent container	Compartmentalized bubble pouch	Six tubes	4 x 4 x 6 ¾	Infectious substance & diagnostic	Liquid-tight container withstands up to 95kPa pressure
INF3000	1.0 liter transparent container	Compartmentalized bubble pouch	Up to 50 ml	6 ¼ x 6 ¼ x 7 ½	Infectious substance	Liquid-tight container withstands up to 95kPa pressure; vessel is autoclavable for re-use



ThermoSafe®
InfectiousShipper



Infectious Shipper

630

ThermoSafe shippers include EPS foam container, flexible foam insert, clear plastic press-lock bag, secondary packaging and absorbent material for full capacity of blood, corrugated carton, and step-by-step packing instructions.



Infectious Shipper

INF3000

Infectious shippers include pressure vessel with stop-lock lid, compartmentalized bubble pouch, 50 ml absorbent sheet, corrugated outer carton with required markings, shipper's declarations, and easy-to-follow instructions.

JD HONIGBERG INTERNATIONAL, INC
Deerfield, IL 60015 USA

Email: medical@jdhintl.com / www.jdhmedical.com

Envoy®

EFFICIENT SEASONAL SOLUTIONS

Numerous configurations available to meet a variety of needs

The ThermoSafe Envoy series features a vast selection of temperature ranges, utilizing standardized materials with the lowest packaging cost available. These efficient and easy-to-assemble shippers are ideal for shipping diagnostic samples.

- Ideal for 1-2 day shipments
- Designed to our real-world ISCSilver® ambient temperature profile
- ZERO BENCH-TIME® refrigerant design means no waiting which greatly simplifies and accelerates the assembly process
- Four container sizes and modular designs give you the flexibility to change refrigerants to meet your temperature criteria
- FEWER components, FASTER assembly, INCREASED operational efficiencies
- Molded containers utilize natural convection to reduce weight and overall shipper size
- Manufactured in the U.S.



SPECIFICATIONS:

- Temperature Ranges: 2-8°C, 15-25°C, < -20°C
- Duration:
 - Refrigerated and CRT: Minimum of 2 days
 - Frozen: Minimum of 4 days
- Packout Configuration: Seasonal
- Ambient Profile: ISCSilver
- Insulation: Expanded Polystyrene (EPS)

AVAILABLE MODELS

	PAYLOAD CAPACITY (liters)	WEIGHT lbs. [kg.]		PAYLOAD ID in. [mm.]			SHIPPER OD in. [mm.]			PART NUMBER	
		Summer	Winter	L	W	H	L	W	H	Summer	Winter
REFRIGERATED 2-8°C	2.5	17.4 [7.9]	12.4 [5.6]	9.01 [229]	7.99 [203]	2.13 [54]	16.7 [425]	15.75 [400]	16.26 [413]	E3R2S	E3R2W
	6.4	23.4 [10.6]	20.8 [9.4]	8.86 [225]	7.87 [200]	5.63 [143]	16.7 [425]	15.75 [400]	21.26 [540]	E6R2S	E6R2W
	12.1	36.1 [16.4]	28.9 [13.1]	11.14 [283]	9.6 [244]	6.89 [175]	20.51 [521]	19.01 [483]	20.75 [527]	E12R2S	E12R2W
	24.1	63.3 [28.7]	43.0 [19.9]	19.76 [502]	11.26 [286]	6.61 [168]	29.02 [737]	20.51 [521]	21.26 [540]	E24R2S	E24R2W
CRT 15-25°C	4.5	14.2 [6.4]	10.0 [4.5]	7.99 [203]	7.01 [178]	4.88 [124]	16.7 [424]	15.75 [400]	16.26 [413]	E5C2S	E5C2W
	8.1	19.2 [8.7]	14.2 [6.4]	7.99 [203]	7.01 [178]	8.87 [225]	16.7 [424]	15.75 [400]	21.26 [540]	E8C2S	E8C2W
	18.2	28.8 [13.1]	19.4 [8.8]	11.14 [283]	9.61 [244]	10.47 [266]	20.51 [521]	19.01 [483]	20.75 [527]	E18C2S	E18C2W
	39.6	49.8 [22.6]	32.8 [14.9]	19.76 [502]	11.26 [286]	10.87 [276]	29.02 [737]	20.51 [521]	21.26 [540]	E40C2S	E40C2W
FROZEN -20°C	2.5	23.8 [10.8]		9.01 [229]	7.99 [203]	2.13 [54]	16.7 [424]	15.75 [400]	16.26 [413]	E3F4	
	5.6	31.8 [14.4]		9.01 [229]	7.99 [203]	4.76 [121]	16.7 [424]	15.75 [400]	21.26 [540]	E6F4	
	12.1	49.2 [22.3]		11.14 [283]	9.6 [244]	6.89 [175]	20.51 [521]	19.01 [483]	20.75 [527]	E12F4	
	24.1	88.9 [40.3]		19.76 [502]	11.26 [286]	6.61 [168]	29.02 [737]	20.51 [521]	21.26 [540]	E24F4	

Please visit thermosafe.com/envoy for more detailed information.

ABOUT THERMOSAFE

Sonoco ThermoSafe, a business unit of Sonoco (NYSE: SON), is the leading global provider of temperature assurance packaging for the safe and efficient transport of pharmaceuticals, biologics, vaccines and other temperature-sensitive products. Sonoco ThermoSafe shipping solutions mitigate risk for our customers and ensure product efficacy throughout the extremes of a supply chain. For more information about Sonoco ThermoSafe's complete line of packaging solutions, please visit the company website.



COVID-19 (Coronavirus) Response

CDC GUIDELINES

USED FOR THE CARE OF ALL PATIENTS

Standard Precautions

Applies to blood, all body fluids, non-intact skin and mucous membranes. Utilize transmission-based precautions (contact, droplet) in addition to standard precautions as indicated.

- WASH HANDS FOR 20 SECONDS WITH SOAP AND WATER OR USE AN ALCOHOL-BASED HAND GEL**
 - After touching blood, body fluids, secretions, excretions, and contaminated equipment
 - After gloves are removed
 - Before and after patient contact
- WEAR GLOVES**
 - When touching blood, body fluids, secretions, excretions, mucous membranes, non-intact skin and contaminated items/equipment
- WEAR A MASK, EYE PROTECTION OR FACE SHIELD**
 - To protect mucous membranes of the eyes, nose and mouth during procedures likely to generate splashes/sprays of blood, body fluids, secretions and excretions
- WEAR A GOWN**
 - During procedures likely to generate splashes or sprays of blood, body fluids, secretions or excretions
- CLEANING AND DISINFECTION**
 - Clean common use equipment before removing from room
- AVOID/PREVENT INJURY:**
 - Use sharps with safety features
 - Sharps containers should be leak-proof and puncture resistant
 - Minimize needles in connections
 - Never reuse used needles
 - Do not bend, break or remove used needles
 - Place used needles/sharps in appropriate sharps containers
 - Never attempt to remove sharps from containers
 - Do not cut, crush, handle contact lenses, supply connectors or tip bath in any patient care area
 - Staff with respiratory lesions or weeping dermatitis should not provide direct patient care until the condition has resolved

www.BowmanDispensers.com • 800-962-4660
©2020 BOWMAN® Manufacturing Company, Inc. (1-800-962-4660) (BMT-HPD/CDC) Guidelines, 2007



TO PREVENT THE SPREAD OF INFECTION

Contact Precautions

- Perform hand hygiene before entering and prior to exiting the room.**
- Put on gloves before room entry. Discard gloves before exiting the room.**
- Put on gown before room entry. Discard gown before exiting the room. Do not wear the same gown and gloves for the care of more than one person.**
- Use dedicated or disposable equipment. Clean and disinfect reusable equipment before leaving room and using it on another person.**
- Contact Infection Prevention if you have questions.**

www.BowmanDispensers.com • 800-962-4660
©2020 BOWMAN® Manufacturing Company, Inc. (1-800-962-4660) (BMT-HPD/CDC) Guidelines, 2007



TO PREVENT THE SPREAD OF INFECTION

Airborne Precautions

- Perform hand hygiene before entering and prior to exiting the room.**
- Keep door closed at all times AIR. Negative air pressure recommended.**
- Put on a fit-tested N-95 or higher level respirator before entering room. Remove respirator after exiting the room and closing the door.**
- Contact Infection Prevention if you have questions.**

www.BowmanDispensers.com • 800-962-4660
©2020 BOWMAN® Manufacturing Company, Inc. (1-800-962-4660) (BMT-HPD/CDC) Guidelines, 2007



EYE PROTECTION/GOGGLES

PUBLIC AREA CONTAINMENT

AVOID INFECTION: Cover your cough and keep hands clean

BD111-0012



BD211-0012






KS101-0029



COVID-19 (Coronavirus) Response

Phase II - Patient Care Area Focus

CDC GUIDELINES



EYE PROTECTION/GOGGLES

PATIENT CARE AREA CONTAINMENT

ISOLATION SUPPLIES DISPENSERS



BD637-0012



BD670-0012



BD607-0012



BD602-0012



BD615-0012
with additional MP-047

PATIENT CARE AREA CONTAINMENT cont.



BD602-0012

- ✓ Gowns
- ✓ Face Masks
- ✓ Gloves
- ✓ Eye Protection



PD003-0111

- ✓ Eye Protection



FP-017

- ✓ Eye Protection or
- ✓ N-95 Masks



BK102-0111

- ✓ N-95 Masks
- ✓ Eye Protection/Goggles



FM300-0212

- ✓ N-95 Masks

ENVIRONMENTAL CONTAINMENT



ED-097

- ✓ Hand Gel
- ✓ Surface Wipes
- ✓ Tissue or Gloves

All consumable items are for demonstration purposes only and are not available for purchase through BOWMAN®

Distributed by:



JD Honigberg International, Inc.
Deerfield IL 60015 USA
Email: medical@jdhintl.com

www.JDHintl.com · 847-412-0200








PRODUCT VIDEOS:
[BowmanDispensers.com/
Product-Videos](http://BowmanDispensers.com/Product-Videos)

**MADE
IN THE
USA**

SPECIFICATIONS

Model Number	5400 ^N	5400 ^W
Weight Capacity	500 lbs. (227 kg)	500 lbs. (227 kg)
Overall Length	83" (211 cm)	83" (211 cm)
Overall Width	29.5" (75 cm)	33.5" (85 cm)
Height Range		
Low	23" (58 cm)	23" (58 cm)
High	32.75" (83 cm)	32.75" (83 cm)
Litter Positioning		
Fowler Backrest	0°-90°	0°-90°
Trendelenburg	+12°	+12°
Reverse Trendelenburg	-7°	-7°
Patient Surface	77" x 27" (196 x 69 cm)	77" x 31" (196 x 77 cm)
Side Rails	15" x 63 ⁵ / ₈ " (38 x 162 cm)	15" x 63 ⁵ / ₈ " (38 x 162 cm)
Casters	8" (20 cm)	8" (20 cm)
Warranty	1 year parts and labor	1 year parts and labor

OPTIONS

Accessory	N	W	Accessory	N	W	Accessory	N	W
X-Ray Board			IV Pole Options:			Additional Options:		
Full	750203	750204	Fixed Height	2101	same	Digital Scale	750214	same
Partial	750205	750206	Adjustable / Removable ..	5984001	same	Bed Exit Monitor	750216	same
Mattresses:			Adjustable / Attached.....	750501	750502	Extra Release Pedal - Hydraulic for stretchers with TDS	750105	same
3" Standard (black).....	5828002	5831002	Fold Down	59175001	same	iS6™ Instant Steer 6 Wheel	750103	same
3" Premium (blue).....	5850002	5852002	IV Pole Caddy	59163001	same	Extra Release Pedal - Hydraulic for stretchers with iS6™ ..	750106	same
4" Standard (black).....	5829002	5832002	Push/Pull Bar Options:			Storage Bag Hook.....	750215	same
4" Premium (blue).....	5851002	5853002	Push/Pull Bar (<i>Removable</i>) ..	5985001	5986001	Paper Roll Holder	5966001	5966003
4" Pressure Relief	5833002	5834002	Bar Pad	5704001	5705001	Quick Release O ₂ Holder	750213	same
5" Premium (blue) [†]	5855002	5856002	Chart Holder	5967001	same	O ₂ Tank Holder (Vertical)..	5981001	same
[†] Special order only.			Monitor Tray Board	59115001	59116001	Head/Foot Extension	5996001	same
Restraint Straps Options:			(includes chart holder 5967001)			Head/Foot 3" Pad	5708001	same
Velcro (pair).....	5903001	same	Knee Flex			Stirrups Set.....	59100001	same
Buckle (pair).....	5904001	same	Hand Crank	750201	750202	Color-Coded Bumper Options:		
Side Rail Options:			Hydraulic	750208	750209	Red 	750401	same
CuVerro® Rails	750303	same	 Leg Lift	 Knee Flex	 Trauma	Teal 	750402	same
Rail Pad (pair)	5716002	same						

Contact info: JD Honigberg International, Inc



800.822.3501 • FAX 360.696.1700 • www.pedigo-usa.com

Email: medical@jdhintl.com

4000 SE Columbia Way • Vancouver, Washington 98661

02/15 MM-80

stretcher series

5400

STANDARD FEATURES

- 500 lbs. weight capacity
- True Direction Steering (TDS)
- Side rails - full length, single release
- Center column design with precision ball bearings
- Single Trendelenburg assist cylinder (positive/negative - pneumatic)
- Pneumatic assist fowler 400 N (90#) force
- Dual side-mounted hydraulic foot pedals
- Brake & Steer activated foot pedals
- Four corner IV pole receptacles
- Four restraint strap positions
- Full PVC perimeter bumper
- 8" easy mobility casters
- Latex-free



stretcher series

5400

universal procedure stretcher




stretcher series

5400

The 5400 Series universal transport stretcher is versatile, yet economical. It's the right stretcher for ED, OR, PACU, or any department where cost, efficiency and quality are required.



- 1 RAIL SYSTEM:**
Easy to wipe clean and disinfect, full length tuck-away side rails allow for easier patient transfer by creating a zero transfer gap. Rails feature a release latch centrally located on either side of the litter top under the patient surface.



- 2 PNEUMATIC ASSIST FOWLER:**
Provides additional assistance for the safe and easy raising and lowering of the patient.





- 3 HYDRAULIC HEIGHT ADJUSTMENT FOOT PEDALS:**
Hydraulic Foot Pedals are located on each side of the stretcher within the patient care zone. Color-coded pedals and function graphics are quickly and easily understood.



- 4 CENTER COLUMN SUPPORT:**
Sturdy steel column with single hydraulic lift cylinder design uses sixteen precision ball-bearings to smoothly raise and lower the patient.



- 5 TRUE DIRECTION STEERING (TDS)**
For improved stretcher control, engage the True Direction Steering system. The steer caster located at the patient's right foot end of the stretcher will lock once the steering pedal is applied.



PATIENT LIFTS:
The 5400 Series Stretcher has 5⁷/₈" H of clearance between floor and cowlings to easily accommodate and utilize patient lifts.


stretcher options



FOLD DOWN IV POLE
Available as an option, the fold down IV pole tucks out of the way when not needed and remains attached to the stretcher when in use.



IV CADDY
Quickly and easily secure an IV pole to a stretcher using Pedigo's IV Caddy. Allows one caregiver to manage the safe transport of the patient, stretcher and IV equipment. Stores under litter top when not in use.



CuVerro®
Antimicrobial Copper Surfaces

Now as an added protection in the ongoing battle with bacteria, Pedigo offers the option of side rails made from **CuVerro® bactericidal copper alloy**.

Copper is an intrinsic bactericidal element and **CuVerro®** is the ONLY solid touch surface metal registered by the EPA to continuously kill bacteria¹ that cause infections and pose a risk to human health. Lab tests show more than 99.9% of bacteria¹ associated with disease and infection, including MRSA¹, are killed within two hours when in contact with copper. Bacteria will continue to be killed 24 hours a day, week after week, between regular cleanings.

¹ CuVerro® alloys have proven effectiveness on MRSA (Methicillin-Resistant Staphylococcus aureus), Staphylococcus aureus, Enterobacter aerogenes, Pseudomonas aeruginosa, and Escherichia coli O157:H7 and VRE (Vancomycin-Resistant Enterococcus). CuVerro is a registered trademark of GBC Metals LLC, d/b/a Olin Brass.

✓ **LATEX-FREE**



ZERO TRANSFER GAP
Tuck-Away side rails allows the caregiver to position the 5400 Series Stretchers flush with the bed for a zero transfer gap and easier patient transfer. No expensive or cumbersome transfer devices are required.

stretcher series

P-172

The P-172 Series general transport stretcher provides an economical patient transport solution. This non-hydraulic stretcher series features a 1³/₈" heavy-gauge O.D. tubular frame, continuous rubber bumper, long-life, double ball bearing, 8" swivel casters for quiet movement, and convenient hand holds on each corner.



P-172-C shown with rolling bumpers

P-172


- Four combination swivel lock & brake lock casters mounted on all corners for maximum control.



P-172 CONFIGURATIONS


	Swivel Caster	Brake & Steer Caster	P-170-1 Shelf	P-170-2 Rails	P-170-3 Fowler	P-170-4 IV Pole	P-170-5 Strap	P-170-6 Bumpers	P-170-8 Mattress	P-170-9 O2 Tank Holder
P-172		4								
P-172-C		4	•	•	•	•			•	•

stretcher options




MATTRESS PAD

2¹/₂" thick black soft-vinyl mattress pad measures 26"W x 77"L. Antimicrobial.




FOWLER HEAD REST

Stainless steel construction with an adjustable range from 0° - 45°.



IV POLE

Two hook pole includes holder and storage tube.



LOWER SHELF

Stainless steel storage shelf.

SPECIFICATIONS

Model Number	P-172
Weight Capacity	350 lbs. (159 kg)
Overall Length	78" (198 cm)
Overall Width	28" (71 cm)
Fixed Height	33¼" (85 cm)
Optional Litter Positioning Fowler Backrest	3°, 12°, 26°, 38°, 44°
Patient Surface	77" x 26" (196 x 66 cm)
Optional Side Rails	12" x 43" (30 x 109 cm)
Casters	8" (20 cm)
Warranty	1 year parts and labor

STANDARD FEATURES

- 350 lbs. weight capacity
- Designed as a cost efficient transport stretcher
- Four Brake & Steer Casters
- Heavy-gauge, 1⅜" O.D. tubular frame
- Fixed litter top height
- Full perimeter bumper
- Convenient hand holds on each corner
- Long-life polyurethane, double ball bearing, 8" swivel casters for quiet movement

✓ **LATEX-FREE**

OPTIONS

Accessory	Accessory	Accessory
Storage Shelf: Shelf P-170-1	IV Pole: Fixed Height P-170-4	Mattress: Pad (2½" thick black) P-170-8
Side Rail: Stainless Steel P-170-2	Restraint Strap: Velcro..... P-170-5	O2 Holder: E Size Tank Holder P-170-9 <i>(requires storage shelf P-170-1)</i>
Fowler: Stainless Steel P-170-3	Rolling Bumpers: 7½" round rubber P170-6	Bed Exit Monitor 750216

stretcher options



O₂ HOLDER

Bracket for E-size oxygen tank. Requires lower shelf.




RESTRAINT STRAP

Blue Velcro restraint strap.



CASTERS

8" long life, double ball bearing casters.



BUMPERS

7½" rolling rubber bumpers



Product Guide Covering Disinfectants & Germicides

DISINFECTANT SPRAY FOR HOSPITAL USE

CLO08 / CLO09 - FLORAL CLO10 / CLO11 - LEMON

DISINFECTANT SPRAY FOR HEALTH CARE USE

CLO14 / CLO15 - LEMON CLO18 - COUNTRY FRESH

CLO71 - AROMATIC (50 STATE COMPLIANT)

CLO75 - LAVENDER (50 STATE COMPLIANT) **NEW**

CLO81 - LEMON (50 STATE COMPLIANT)

- *HIV-1 (AIDS Virus), Methicillin Resistant Staphylococcus aureus (MRSA), Vancomycin Resistant Enterococcus faecalis (VRE), *Polio Virus Types 1 and 2, *Influenza A2/Hong Kong Virus, Mycobacterium Tuberculosis (TB) and *Herpes Simplex Virus Types 1 and 2 On Hard, Nonporous, Inanimate Environmental Surfaces
- Kills Pandemic 2009 H1N1 Influenza A Virus* [formerly called "swine flu"]
- Kills Most Germs - Controls Mold and Mildew - Eliminates Odors
- Bactericidal - Staphylocidal - Tuberculocidal
Fungicidal - Pseudomonocidal - *Virucidal

Disinfects: Disinfects on such bacterial and viral carrying surfaces as toilet seats, telephones, garbage and refuse cans, urinals, athletic equipment, light switches, bed pans, empty hampers and other out of the way places. Disinfects hard, non-porous surface of shoes and slippers.

Controls Mold and Mildew on mattresses, awnings, shower curtains, leather goods. Controls mildew and mold growth on painted or finished walls, ceilings and woodwork in closets, basements, attics, storage areas, summer cottages, damp and musty corners.

Deodorizes: Reduces annoying odors. Use to control odors in hospitals and nursing homes. In private homes, use in kitchen, bathroom, sick rooms, basements. Effectively suppresses odors caused by mold and mildew.

DISINFECTANT BATHROOM CLEANER CL866 - FLORAL

- Disinfects Influenza A2/Japan Virus* and Herpes Simplex Virus* Types 1 and 2 on Hard, Inanimate, Non-porous Environmental Surfaces
- Kills Pandemic 2009 H1N1 Influenza A Virus* [formerly called "swine flu"]
- ANTIBACTERIAL - Kills Common Kitchen Germs
- Bactericidal - *Virucidal - Fungicidal
- Cleans and Deodorizes
- Spray and Wipe Shining Clean

GERMICIDAL CLEANER

CL873 - FLORAL CL876 - COUNTRY FRESH

This foaming all-purpose cleaner cleans, disinfects and deodorizes in one easy step. Spray it on and wipe it off. Cleans and disinfects without scouring or scratching surfaces. Kills most common kitchen germs. Will disinfect Influenza A2/Japan Virus and Herpes Simplex Virus Types 1 and 2 on hard, inanimate, nonporous environmental surfaces. This is an EPA registered product that saves time by performing two operations in one step.

SURFACE SANITIZING WIPES CL913

- No-Rinse Food Contact Sanitizer Wipe
- Kills 99.9% of bacteria in just 30 seconds**
- **Staphylococcus aureus (Staph), Klebsiella pneumoniae
- Antibacterial formula sanitizes food and non-food contact surfaces

GENERAL PURPOSE DISINFECTANT WIPES CL915

- Effectively cleans, disinfects and deodorizes any hard, non-porous, inanimate environmental surface.
- Kills 99.9% of bacteria* in 15 seconds.
- Kills Pandemic 2009 H1N1 INFLUENZA A VIRUS (Formerly Called 'Swine Flu')

CL989 BROAD SPECTRUM GERMICIDAL

& DISINFECTANT WIPE 180 wipes **NEW**

- Kills 99.9% of Bacteria in Just 15 Seconds.
- Eliminates Bacteria In Just 4 Minutes • Deodorizes
- Fast Drying; Leaves No Dull or Sticky Residue
- Lemon Fragrance • Bleach Free
- Kills Pandemic 2009 H1N1 Influenza A Virus (formerly called swine flu)



JD Honigberg International, Inc. | Deerfield, IL 60015 USA | medical@jdhintl.com | www.jdhmedical.com



Emergency

Premier All Aluminum

6 Drawer Emergency Cart: 30"
Single Color (Red or Beige), Breakaway Lock.
Model: PAR-30



Standard features:

- Side shelf- Plastic top
- Push handles (2)
- Accessory panels (3)
- Locking casters (2)
- Swivel caster
- Tracking caster
- Breakaway locking bar (length of bar is equal to Cart size)
- Stabilizing frame with bumper
- Durable aluminum mounting tracks
- Carts require some assembly

Overall dimensions:

45 1/2"H x 25"D x 32 3/4"W / 115.6
115.6cm x 63.5 cm x 83.2 cm
(including bumper/handles/casters)

Drawer Size:

17 1/16"D x 22 1/4"W
43.3cm x 56.5cm

Drawer sizes:

3"H(3), 6"H(2), 9"H (1)
7.6cm (3), 15.2cm (2), 22.9cm (1)

Weight: 115 lbs / 52.2 kgs

Model PAR-30 - Red (pictured)

Color: Red

LATEXfree



Note: Carts may require some assembly.



Specialty

Standard Steel

4 Drawer Isolation Cart with Package

Model: AKL-Y-400IP



Model AKL-Y-400IP (pictured)

Color: Yellow

Standard features:

- AKL-Y-4 Four drawer cart
- PSR-1 Heavy Duty Railing
- APT-2 Plastic top
- AUH-1 Utility Hooks
- ACR-1 Coat Racks

Drawer Size:

17 1/16" D x 29 1/4"W

43.6cm Dx 74.3cm W

Drawer Configuration:

3"H (1), 9"H (3)

7.6cm H (1), 22.9cm H (3)

Weight: 189 lbs / 85.7 kgs

Note: Carts may require some assembly.

LATEXfree

	SB — SLATE BLUE
	SFG — SEA FOAM GREEN
	L — LAVENDER
	M — MAIVE
	Y — YELLOW
	CB — TAUPE
	T — TEAL
	B — BEIGE
	AB — CRASH CART BLUE
	CR — CRANBERRY
	BG — BLUE GREY
	HG — HUNTER GREEN
	AA — DEEP BLUE
	AR — RED

J.D. Honigberg International, Inc. 155 Pfingsten Rd, Suite 150, Deerfield, IL 60015 USA

Tel: 847-412-0200 • Fax: 847-412-0034 • medical@jdhintl.com

SAVE II™

SIMPLIFIED AUTOMATED VENTILATOR

The SAVE II is not just the markets smallest compressor driven ventilator, it is the only one that simplifies automated ventilation using a height chart to dial in a targeted tidal volume.

Supplemental O2 and Patient circuits

Replacement circuit is inexpensive.

Height presets enable minimally trained provider ability to deliver targeted therapy.

Confirm button prevents accidental setting changes.

Manual trigger enables rescuer to deliver an appropriate tidal volume every time regardless of stress levels/skill (BLS CPR Mode).

Rechargeable battery delivers air to the patient for 10+ hours, running and charging simultaneously when plugged in.



Heart icon flashes 100 times per minute to help rescuer time compressions (BLS CPR Mode)

Alarm dashboard lights up alerting provider of any adverse conditions

Battlefield

- » Point of Injury
- » CASEVAC
(Air & Ground Ambulances)
- » Battalion Aid Station
- » Backup for Forward Surgical Teams
- » Nation Building Missions

Hospitals

- » Intra Hospital Transport
- » Crash Cart
- » Surge Ventilator
- » Mass Casualty
- » Pan Flu
- » Backup Ventilator

Emergency Services

- » Advanced Life Support CPR
- » Basic Life Support CPR
- » Inter Hospital Transport
- » Respiratory Failure

JD Honigberg International | 155 N Pfungsten Rd, Ste 150, Deerfield, IL 60015, USA

Tel: (1) 847-412-0200 | www.jdhmedical.com | medical@jdhintl.com